

BIOMOLECULE CHARACTERIZATION & PURIFICATION — AT — EXTREMELY AFFORDABLE — PRICES —

AGGREGATES • ADC • mAb • BIOSIMILARS

**NEW
UHPLC Columns**



1.8 μ m SEC-X150 and X300

Replace Waters[®] BEH 1.7 μ m SEC columns to:

- Save Money
- Minimize Adsorption
- Maximize Separation Range



Welcome to Affordable, Higher Resolution Size Exclusion Chromatography

With Yarra™ SEC Columns you can:

- **Save Money** with extremely affordable prices
- **Achieve Better Results** through larger exclusion ranges and higher efficiencies
- **Get Enhanced Recovery** using more inert Yarra particles and Bio-Inert hardware
- **Gain Time** with faster more productive HPLC/UHPLC runs on the new 1.8 μm SEC-X150 and X300
- **Feel at Ease** knowing you have an unmatched product guarantee

Better Protein/Peptide SEC Starts With Yarra

Affordable Pricing and Unmatched Product Guarantee	4
Dependable Results	5

NEW 1.8µm

with Bio-Inert
Hardware

Exceptional Separation Power (Insulin / mAbs)	6-7
Expanded Separation Window	8
Better Recovery with Greater Inertness	9
HPLC/UHPLC System Versatility	10
Incredible Performance Upgrade	11
Excellent Lifetime and Stability	12

3µm

Simple Upgrade for 5µm SEC Methods	11
Excellent Lifetime and Stability	12

5µm

High Performance Preparative SEC	13
--	----

Select the Right Column for Your Application	14
Ordering Information.....	15

Affordable Pricing and Unmatched Product Guarantee

Yarra SEC columns are highly affordable and are accompanied by a product guarantee that ensures performance satisfaction, making it extremely easy to upgrade current SEC methods!

guarantee

If Yarra™ analytical columns do not provide you with at least an equivalent separation as any other GFC column of similar porosity, type, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Yarra vs. Waters® BEH SEC



Yarra 300 x 7.8 mm		VS.	Waters* 1.7 µm 150 x 4.6 mm	
SEC-X150	SEC-X300		BEH125 SEC	BEH200 SEC
1.8	1.8	Particle Size (µm)	1.7	1.7
150	300	Pore Size (Å)	125	200
1K-450K	10K-700K	MW Range in Native Conditions (Da)	1K-80K	10K-450K
>30,000	>30,000	Efficiency (plates/column)	>30,000	>30,000

*Waters specifications taken from Waters website.

Yarra vs. Tosoh Bioscience® TSKgel®



Yarra 300 x 7.8 mm			VS.	TSKgel® 300 x 7.8 mm		
SEC-2000	SEC-3000	SEC-4000		G2000SWxl	G3000SWxl	G4000SWxl
3	3	3	Particle Size (µm)	5	5	8
145	290	500	Pore Size (Å)	125	250	450
1K-300K	5K-700K	15K-1,500K	MW Range in Native Conditions (Da)	5K-150K	10K-500K	20K-7,000K
48,000	48,000	38,000	Efficiency (minimum theoretical plates)	20,000	20,000	16,000

*All TSKgel specifications were taken from Tosoh Bioscience 2004-5 Laboratory Products Catalog

Comparative separations may not be representative of all applications.


Expert Manufacturing for Dependable Performance and Reliability



Reproducible GFC columns require extreme detail in every aspect of the manufacturing and packing process. First, Yarra silica particles are synthesized using narrow tolerances for pore and particle size. Next, bonding of a proprietary hydrophilic ligand is tightly controlled and packing of each column requires validated recipes tested to high specifications to ensure excellent reproducibility.

Batch-to-Batch Reproducibility


Yarra 1.8µm SEC-X150



Column: Yarra 1.8µm SEC-X150
Dimensions: 150 x 4.6 mm
Part No.: 00F-4631-E0
Mobile Phase: 100 mM Sodium Phosphate in Water
pH 6.8 + 0.025 % NaN₃
Flow Rate: 0.3 mL/min
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
2. IgA (300 kDa)
3. IgG (150 kDa)
4. Ovalbumin (44 kDa)
5. Myoglobin (17 kDa)
6. Uridine

Batch-to-Batch Reproducibility

Yarra 3µm SEC-3000





Column: Yarra 3µm SEC-3000
Dimensions: 300 x 7.8 mm
Part No.: 00H-4513-K0
Mobile Phase: 100 mM Sodium Phosphate in Water
pH 6.8 + 0.025 % NaN₃
Flow Rate: 1 mL/min
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
2. IgA (300 kDa)
3. IgG (150 kDa)
4. Ovalbumin (44 kDa)
5. Myoglobin (17 kDa)
6. Uridine

Exceptional Separation Power - Insulin

The high resolution provided by Yarra 1.8 μ m SEC-X150 columns make them an excellent solution for the characterization of biotherapeutic insulin and any resulting aggregates. With two column lengths, you have the ability to further strengthen resolution with a 300 mm or speed up methods with a 150 mm.

Degraded Human Insulin




Conditions for all columns:

Column: Yarra 1.8 μ m SEC-X150
BEH125 SEC, 1.7 μ m
Dimensions: 150 x 4.6 mm

Mobile Phase: 1 g/L L-Arginine in Water/Acetic Acid/Acetonitrile (65:15:20)

Flow Rate: 0.2 mL/min
Detection: UV @ 280 nm
Sample: Human Insulin

Expanded Resolution of Recombinant Human Insulin



Column: Yarra 1.8 μ m SEC-X150
Dimensions: 300 x 4.6 mm
Part No.: 00H-4631-E0
Mobile Phase: Arginine (1 g/L) in Water/Acetic Acid/Acetonitrile (65:15:20)
Flow Rate: 0.4 mL/min

Temperature: Ambient
Detection: UV @ 276 nm
Sample: 1. Insulin Dimer
2. Insulin Monomer


Comparative separations may not be representative of all applications.

Exceptional Separation Power— Monoclonal Antibodies (mAb)




The large separation range of the Yarra SEC-X300 makes it ideal for high molecular weight aggregation studies of monoclonal antibodies. While the smaller MW focus of the SEC-X150 allows for excellent identification and separation of a mAb and its variable region fragments.

Mouse mAb (monoclonal antibody)




Monoclonal Antibody and Fragments



Expanded Separation Window


An expanded exclusion range was expertly developed and fine-tuned for both Yarra™ 1.8 µm products. Together they offer a wider applicability for popular biomolecule separations. The SEC-X150 focuses on smaller MW separations in the range of 1 kDa-450 kDa while the SEC-X300 targets those high MW compound mixes inside of 10 kDa-700 kDa.



Low MW Separation Mix

Column: Yarra 1.8 µm SEC-X150
Dimensions: 150 x 4.6 mm
Part No.: 00F-4631-E0
Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.03 % Na₃
Flow Rate: 0.35 mL/min
Temperature: Ambient
System: Agilent® 1260 UHPLC
Detection: UV @ 280 nm
Sample: 1. BSA Dimer (132 kDa)
2. BSA (66 kDa)
3. Ovalbumin (44 kDa)
4. Myoglobin Dimer (34 kDa)
5. Myoglobin (17 kDa)
6. Ribonuclease A (13.7 kDa)
7. Insulin (5800 Da)
8. Uridine (244 Da)

mAb with Dimer



Column: Yarra 1.8 µm SEC-X300
Dimensions: 150 x 4.6 mm
Part No.: 00F-4743-E0
Mobile Phase: 0.1 M Potassium phosphate in Water (pH 6.8)
Flow Rate: 0.35 mL/min


Temperature: Ambient
Detection: UV @ 280 nm
Sample: 1. Dimer
2. Monomer

Better Recovery with Greater Inertness




Phenomenex's proprietary surface chemistry provides an inertness that is hard to match by other GFC columns. The result is minimal adsorption of peptides, proteins, and other protein aggregates leading to more accurate quantitation.

Yarra® 1.8µm SEC-X150



Waters® BEH200 SEC 1.7µm



Conditions for all columns:

Columns: Yarra 1.8µm SEC-X150
BEH200 SEC 1.7µm

Dimensions: 150 x 4.6 mm

Mobile Phase: 50 mM Sodium Phosphate in Water pH 5.0/
0.1 M Sodium Sulfate

Flow Rate: 0.3 mL/min

Temperature: Ambient

Detection: UV @ 280 nm

Sample: 1. IgG Dimer
2. IgG Monomer
3. Albumin



New Bio-Inert Hardware

In addition to a highly inert particle chemistry, Yarra 1.8µm SEC-X150 columns utilize a Bio-Inert hardware to further ensure accurate recoveries.

Comparative separations may not be representative of all applications.

Fast GFC/SEC Methods on Your HPLC or UHPLC

Utilize the high performance of the new Yarra™ 1.8µm on the system(s) that you know and are comfortable with! The amazing separation power of the Yarra 1.8µm can easily be utilized on any HPLC and UHPLC system.



Waters® ACQUITY® UPLC®



Shimadzu® Nexera®



Agilent® 1200



Questions About LC System Compatibility?


Contact your local Phenomenex
representative for guidance!

Performance Upgrade




Upgrading to better performance from existing SEC columns is easy with the Yarra column line. Yarra 3 μm particle columns are a quick replacement for existing 5 μm columns, giving increased performance and better peak shape without method adjustment. However, large gains in analyses times, resolution and efficiency can be realized with the Yarra 1.8 μm SEC-X150 on the desired HPLC/UHPLC system.


Conventional 5 μm 2000 300 x 7.8 mm



Yarra 3 μm SEC-2000 300 x 7.8 mm



Yarra 1.8 μm SEC-X150 150 x 4.6 mm



Conditions for all columns:

Columns: Conventional 5 μm 2000

Yarra 3 μm SEC-2000

Yarra 1.8 μm SEC-X150

Dimensions: As noted in chromatograms

Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.025% NaN_3

Flow Rate: 1 mL/min

0.3 mL/min (Yarra 1.8 μm SEC-X150 150 x 4.6 mm)

Temperature: Ambient

Detection: UV @ 280 nm

Sample: 1. Thyroglobulin (669 kDa)

2. IgA (300 kDa)

3. IgG (150 kDa)

4. Ovalbumin (44 kDa)

5. Myoglobin (17 kDa)


6. Uridine

Comparative separations may not be representative of all applications.

Excellent Stability and Lifetime

The great care exhibited in both engineering and packing of Yarra™ columns can produce stable performance over longer lifetimes than other size exclusion columns on the market. By gaining extended lifetime with Yarra columns, you save significant time and money for your lab.


Increased Lifetime Under Extreme Conditions (10 % IPA and higher flow rate)



Column: Yarra 1.8 µm SEC-X150
Dimension: 150 x 4.6 mm
Part No.: 00F-4631-E0
Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.025 % NaN₃ + **10 % IPA**
Flow Rate: **0.4 mL/min**
Temperature: Ambient
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
 2. IgA (300 kDa)
 3. IgG (150 kDa)
 4. Ovalbumin (44 kDa)
 5. Myoglobin (17 kDa)
 6. Uridine

Minimal change
in retention
and
performance

Exceptional Stability Over 900 Injections






Column: Yarra 3 µm SEC-2000
Dimensions: 300 x 7.8 mm
Part No.: 00H-4512-K0
Guard Cartridge: AJO-4487
Guard Holder: KJO-4282, SecurityGuard Guard Cartridge Kit
Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.025 % NaN₃
Flow Rate: 1 mL/min
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
 2. IgA (300 kDa)
 3. IgG (150 kDa)
 4. Ovalbumin (44 kDa)
 5. Myoglobin (17 kDa)
 6. Uridine

High Performance Preparative SEC



Enjoy the same selectivity and ultra-high efficiency of Yarra 3 μm for your preparative gel filtration applications. Yarra SEC PREP features a 5 μm particle size version of each Yarra 3 μm chemistry in a 21.2 mm ID column for preparative purification, desalting, and characterization of biomolecules. Yarra 5 μm PREP is available at an affordable price while maintaining the high performance given with the analytical columns.



Conditions for both applications:



Columns: Yarra 5 μm SEC-2000 PREP
Dimensions: 300 x 21.2 mm
Part No.: 00H-4619-PO
Mobile Phase: 100 mM Sodium Phosphate in Water (pH 6.8)

Flow Rate: 5 mL/min
Temperature: Ambient
Detection: UV @ 280 nm (ambient)
Sample: 1. Bovine Serum Albumin (66 kDa)
 2. Myoglobin (17 kDa)

Select the Right Yarra Column for Your Application

With five Yarra™ phases covering a large exclusion window, the Yarra product line offers the perfect solution for HPLC/UHPLC size exclusion chromatography of biomolecules.

Molecular Weight (MW) Separation Ranges




Looking for a Yarra 5 µm Analytical Column?

We recommend using a 3 µm column to improve the resolution/sensitivity of your 5 µm method.

If method optimization is not an option, our BioSep™ 5 µm GFC/SEC columns are a guaranteed alternative to your current 5 µm column.

Visit www.phenomenex.com/biosep or contact your local representative for more details.



Ordering Information



Yarra Column Specifications

	SEC-X150	SEC-X300	SEC-2000	SEC-3000	SEC-4000	SEC-2000 PREP	SEC-3000 PREP	SEC-4000 PREP
Resin Type	Silica	Silica	Silica	Silica	Silica	Silica	Silica	Silica
Particle Size (µm)	1.8	1.8	3	3	3	5	5	5
Pore Size (Å)	150	300	145	290	500	145	290	500
MW Range (Da)	1K-450K	10K-700K	1K-300K	5K-700K	15K-1500K	1K-300K	5K-700K	15K-1500K
pH Range	1.5 - 8.5	1.5 - 8.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5
Typical Backpressure (psi)	3,250	3,250	1,300	1,300	1,000	800	800	800
Max Backpressure (psi)	7,000	7,000	3,000	3,000	1,700	1,500	1,500	1,500
Typical Efficiency	> 30,000 (150 x 4.6 mm)	> 30,000 (150 x 4.6 mm)	48,000 (300 x 7.8 mm)	48,000 (300 x 7.8 mm)	38,000 (300 x 7.8 mm)	30,000 (300 x 21.2 mm)	30,000 (300 x 21.2 mm)	30,000 (300 x 21.2 mm)
Max. Flow Rate (mL/min)	0.4	0.4	1.5	1.5	1.2	10	10	10

Yarra 1.8 µm SEC Bio-Inert Columns (mm)

	Analytical	Analytical
Phases	150 x 4.6 mm	300 x 4.6 mm

Yarra 1.8 µm SEC-X150 00F-4631-E0 00H-4631-E0

Yarra 1.8 µm SEC-X300 00F-4743-E0 00H-4743-E0



guarantee

If Yarra® analytical columns do not provide you with at least an equivalent separation as any other GFC column of similar porosity, type, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Yarra 1.8 µm SEC Stainless Steel Columns (mm)

	Analytical	Analytical	SecurityGuard ULTRA Cartridges***
Phases	150 x 4.6	300 x 4.6	3/pk
Yarra 1.8 µm SEC-X150	00F-4631-E0-SS	00H-4631-E0-SS	AJ0-9512
Yarra 1.8 µm SEC-X300	00F-4743-E0-SS	00H-4743-E0-SS	AJ0-9513
			For Stainless Steel Only

***SecurityGuard ULTRA cartridges require holder, Part No.: AJ0-9000

Yarra 3 µm SEC Columns (mm)

	Narrow Bore	Analytical	Analytical	SecurityGuard™ Cartridges (mm)
Phases	300 x 4.6	150 x 7.8	300 x 7.8	4 x 3.0*
				/10 pk
Yarra 3 µm SEC-2000	00H-4512-E0	00F-4512-K0	00H-4512-K0	AJ0-4487
Yarra 3 µm SEC-3000	00H-4513-E0	00F-4513-K0	00H-4513-K0	AJ0-4488
Yarra 3 µm SEC-4000	00H-4514-E0	—	00H-4514-K0	AJ0-4489
				for ID 4.6–7.8 mm

*SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282

Yarra 5 µm PREP SEC Columns (mm)

	Preparative	SecurityGuard™ Cartridges (mm)
Phases	300 x 21.2	15 x 21.2**
		/ea
Yarra 5 µm SEC-2000 PREP	00H-4619-P0	AJ0-8588
Yarra 5 µm SEC-3000 PREP	00H-4620-P0	AJ0-8589
Yarra 5 µm SEC-4000 PREP	00H-4621-P0	AJ0-8950
		for ID 18–29 mm

**PREP SecurityGuard™ Cartridges require holder, Part No.: AJ0-8223



Terms and Conditions

Offer on front cover of this brochure is valid only for addressee. Limit 6 columns at discounted price per company. Offer excludes columns with > 8mm ID. Offer must be specified at time of order and may not be combined with any other discount or promotional item. Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Yarra, BioSep, and SecurityGuard are trademarks of Phenomenex. TSKgel and Tosoh Bioscience are registered trademarks of Tosoh Corporation. ACQUITY, UPLC, and Waters are registered trademarks of Waters Corporation. Agilent is a registered trademark of Agilent Technologies, Inc. Shimadzu and Nexera are registered trademarks of Shimadzu Corporation.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Tosoh Corporation, or Waters Corporation. Comparative separations may not be representative of all applications.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPro, PREP or ULTRA holders, or to any cartridges.

© 2016 Phenomenex, Inc. All rights reserved.

NEW Yarra 1.8 µm SEC-X150 and X300

Replace Waters® BEH
1.7 µm SEC columns to:

- Save Money
- Minimize Adsorption
- Maximize Separation Range

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com