

# Fully Automated Sample Preparation for the Determination of Vitamin D on the Biomek® i7 Hybrid Using the Impact™ Protein Precipitation Plate Combined with Positive Pressure Solid Phase Extraction

Kerstin Thurow<sup>1</sup>, Heidi Fleischer<sup>1</sup>, Anna Bach<sup>1</sup>, Bhagya Wijayawardena<sup>2</sup>, Miranda Kheradmand<sup>2</sup>, Jenny Lewis<sup>3</sup>, and Dr. Bryan Tackett<sup>3</sup>

<sup>1</sup>Center for Life Science Automation, University Rostock, Friedrich-Barnewitz-Str. 8, 18119 Rostock, Germany

<sup>2</sup>Beckman Coulter Life Sciences, 5350 Lakeview Pkwy S Dr., Indianapolis, IN 46268 USA

<sup>3</sup>Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

## Overview

Quantitation of vitamin D<sub>2</sub> and D<sub>3</sub> (25(OH)D<sub>2</sub> and 25(OH)D<sub>3</sub>) is becoming increasingly important among researchers to understand the effects of varying amounts of vitamins. The manual sample preparation for LC-MS based vitamin D analysis is time consuming and prone to human errors. Consequently, automation of the complex sample preparation workflow is important to carry out vitamin D analysis in a high-throughput manner.

The protein precipitation and the separation of the metabolites from the vitamin D binding proteins is usually carried out by adding acetonitrile, methanol, as well as mixtures of acetonitrile and methanol or 2-propanol. The samples are then vortexed for 1 min, followed by centrifugation. Thereafter, supernatant is removed, and the solvent is evaporated in a gentle stream of nitrogen. Subsequent redissolution of the samples in a suitable solvent is required prior to analysis. In this application note, we aimed to develop an automated method for determining the vitamin D metabolites 25(OH)D<sub>2</sub> and 25(OH)D<sub>3</sub> without centrifugation. We used the Biomek i7 hybrid workstation with integrated Positive Pressure Extractor (Amplius GmbH, D). To avoid centrifugation, we used the Impact protein precipitation plate kit (Part No. [CE0-8201](#), Phenomenex, Torrance, US) in the workflow. Once the samples and the reagents are added to the Impact protein precipitation plate, the protein precipitation takes place in the plate. After mixing, the sample is moved through the filter by applying positive pressure or vacuum, separating vitamin D metabolites from the protein aggregates. The resulting filtrate can be further cleaned-up using solid phase extraction (SPE) Strata® C8 Cartridges with a sorbent mass of 30 mg/1 mL (Part No. [8B-S005-TAK](#), Phenomenex, Torrance, US).

We automated the complete workflow by integrating the Positive Pressure Extractor and the Impact protein precipitation plate. This automated workflow is less costly as it eliminates the centrifugation steps. We analyzed the processed samples using LC-MS and generated calibration curves to identify repeatability of the automated sample preparation (**Figure 2**). The low CV values (2 % -6 %) indicated high repeatability of the automated method. The results were comparable with manual sample preparation.

## Sample Preparation

<b>Pre-treatment:</b>	Transfer 500 µL Methanol/Water (60:40, v/v) followed by 200 µL Zinc Sulfate (0.2 M)
<b>Add:</b>	40 µL internal standard followed by 100 µL serum sample
<b>Shake:</b>	Plate for 3 min to enable suitable protein precipitation
<b>Incubate:</b>	Plate for 3 min at room temperature
<b>Apply:</b>	Positive pressure for 10 min (1500 bar) using the PP ALP
<b>Condition:</b>	Strata C8 cartridge with 250 µL of Methanol
<b>Equilibrate:</b>	Strata C8 cartridge with 500 µL Water
<b>Load:</b>	500 µL sample onto the Strata C8 Cartridge, apply pressure
<b>Wash:</b>	Strata C8 cartridge with 1 mL Methanol/Water (60:40, v/v)
<b>Dry:</b>	Strata C8 cartridge for 10 min
<b>Elute:</b>	Vitamin D components twice with 100 µL Methanol
<b>Measure:</b>	Samples using LC-MS or LC-MS/MS

**Figure 1.** Structures of 25(OH)D<sub>2</sub> and 25(OH)D<sub>3</sub>





Figure 2. Calibration Curves for (a) 25(OH)D<sub>2</sub> and (b) 25(OH)D<sub>3</sub>



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

**Australia**

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**🌐 All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Strata is a registered trademark, and Impact and BE-HAPPY are trademarks of Phenomenex.

Biomek and BECKMAN COULTER are registered trademarks of Beckman Coulter, Inc.

**Disclaimer**

Phenomenex is in no way affiliated with the Center for Life Science Automation, University of Rostock.

Biomek i-Series is for research use only and not intended or validated for use in the diagnosis of disease or other conditions.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat) to get in touch with one of our Technical Specialists