

APPLICATIONS

Utilizing a Rapid, Two-Step Method for the Clean-Up of Veterinary Drugs in Milk Using Strata[®]-X PRO Solid Phase Extraction (SPE)

Shahana Huq and Jenny Cybulski
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Analyzing trace veterinary drugs in milk is a growing area of concern meant to ensure the safety of food products. These drugs, which are used to treat animals, would ultimately be incorporated to humans by consumption of milk and pose health hazards. Any proposed solution needs to quickly and accurately quantify the residual drugs present in milk and related products before injection onto the LC column. When working with milk as a matrix, phospholipids from milk fat must be removed to reduce any ion suppression that could occur during LC-MS/MS analysis. To overcome these obstacles, Strata-X PRO, a solid phase extraction (SPE) product, offers a fast, two-step sample preparation method to remove phospholipids prior to MS analysis. This SPE product shows an improved solution over traditional protein precipitation methods and other types of SPE, due to clean up efficiency while maintaining a rapid and fast analysis time.

Materials and Methods

Reagents and Chemicals

Analytical reference standards were purchased from Cerilliant[®] Corporation (Round Rock, TX, and USA). All other chemicals were obtained from Sigma-Aldrich[®] (St. Louis, MO) and the analytes and internal standards were purchased from Cerilliant (Round Rock, TX). ALTA-DENA brand milk was used as a sample matrix for extraction. Ultrapure D.I. water was obtained from Sartorius[®] arium[®] comfort II, courtesy of Sartorius Corporation (Bohemia, NY).

SPE Protocol

Pre-treatment

To 1 mL of milk (spiked with analyte) add 3 mL of 0.2% Formic acid in Acetonitrile/Methanol (90:10) and mix or vortex for 15-20 seconds. Centrifuge for 5 minutes at 10,000 RPM and collect supernatant.

Cartridge: Strata-X PRO 60 mg/3 mL

Part No.: 8B-S536-UBJ

Load: Pass the pre-treated sample through the SPE cartridge and collect

Dry: Evaporate the extract to dryness under a gentle stream of nitrogen at room temperature

Reconstitute: The dried sample in 1 mL of initial mobile phase (0.1% Formic acid in Water/0.1% Formic acid in Methanol (95:5)) spiked with deuterated internal standard.

LC Conditions for Chromatogram

Column: Kinetex[®] 2.6 μ m Biphenyl 100 Å

Dimensions: 50 x 3.0 mm

Part No.: 00B-4622-Y0

SecurityGuard[™] ULTRA: AJ0-9208

Mobile Phase: A: 0.1% Formic acid in Water

B: 0.1% Formic acid in Methanol

Gradient: **Time (min)** **% B**

0 5

1.5 95

3 95

3.01 5

4.5 5

Flow Rate: 0.5 mL/min

Injection Volume: 5 μ L

Temperature: 45 °C

Instrument: Agilent[®] 1260

Detector: SCIEX Triple Quad[™] 4500 (ESI, +ve Ionization)

LC Conditions for Phospholipid Comparison

Column: Kinetex 2.6 μ m C18

Dimensions: 50 x 2.1 mm

Part No.: 00B-4462-AN

SecurityGuard ULTRA: AJ0-8782

Mobile Phase: A: 0.1% Formic acid in Water

B: 0.1% Formic acid in Methanol

Gradient: **Time (min)** **% B**

0 40

0.5 95

11.5 95

11.51 40

13 40

Flow Rate: 0.5 mL/min

Injection Volume: 2 μ L

Temperature: 45 °C

Instrument: Agilent 1260

Detector: SCIEX Triple Quad 4500 (ESI, +ve Ionization)

Sample: Phospholipid (Retention time in minute)

1. Lyso PC (2.25), MRM transition 496.4/184.2

2. PC-1 (4.14), MRM transition 760.7/184.2

3. PC-2 (4.6), MRM transition 786.8/184.2

Results and Discussion

Table 1. % Recovery and CV of Veterinary Drugs (50 ng/mL) in Milk using Strata-X PRO SPE

Peak No.	Analyte Name	Retention Time (min)	% Recovery	% CV	Q1	Q3
1	Sulfaguanidine	1.48	46	5	215	156.1
2	Lincomycin	2.07	92	5	407.1	126
3	Sulfadiazine	2.19	38	7	251	156
4	Cephapirin	2.22	76	7	424	292.1
5	Sulfamerazine	2.32	44	5	265.1	155.8
6	Sulfamethoxazole	2.36	53	13	254.1	156.1
7	Sulfamethizole	2.36	45	8	271.1	92
8	Cefalexin	2.39	66	4	348.2	174.2
9	Sulfamethazine	2.44	59	13	279.1	186.1
10	Cortisone	2.72	83	8	361.2	163.2
11	Cortisol	2.73	95	6	363.4	120.9
12	β -methasone	2.76	97	3	393.4	355.2
13	Prednisolone	2.81	92	10	361.2	147.2

APPLICATIONS

Figure 1. Chromatogram of Veterinary Drugs in Milk Following Clean-up by Strata[®]-X PRO using a Kinetex[®] Biphenyl LC Column





Figure 2. LC-MS/MS Total Ion Current (TIC) Comparison of Blank Milk Sample Processed with Standard Protein Precipitation (A) and Strata-X PRO SPE (B)



APPLICATIONS




In this technical note, Strata[®]-X PRO is compared to the protein precipitation method that is widely used for its simplicity, quick and minimal method development. In a parallel study, the TIC obtained from MS analysis, displays the high abundance of the phospholipids (**Fig. 2A**) in the sample extract from protein precipitation. The Strata-X PRO sample clean-up, on the other hand, depleted the majority of the phospholipids (**Fig. 2B**) resulting in absolute recovery of 13 analytes (**Fig. 1**) around 67% (**Table 1**) on the average. Failure to remove phospholipids in the sample results in matrix effect and eventually pre-mature LC column death due to continuous lipid build-up. Additionally, an increased instrument downtime and higher maintenance for the MS instrument is observed in the aftermath.

Conclusion

Though widely used, protein precipitation method co-extracts endogenous phospholipids that negatively affect chromatographic analysis resulting in skewed analytical data. A cleaner and yet a fast two-step method using Strata-X PRO, mitigates these effects, while ensuring the upkeep of the MS instrument. The analytical data elucidated in this study demonstrates effective removal of phospholipids and fats from milk matrix, maintaining the precision and accuracy of the assay, while Kinetex[®] Core-Shell LC columns were used in conjunction.


Ordering Information

Strata[®]-X PRO SPE

Format	Sorbent Mass	Part Number	Unit
Tube			
	10 mg	8B-S536-AAK	1 mL (100/box)
	30 mg	8B-S536-TAK	1 mL (100/box)
	30 mg	8B-S536-TBJ	3 mL (50/box)
	60 mg	8B-S536-UBJ	3 mL (50/box)
	200 mg	8B-S536-FBJ	3 mL (50/box)
	100 mg	8B-S536-ECH	6 mL (30/box)
	200 mg	8B-S536-FCH	6 mL (30/box)
	500 mg	8B-S536-HCH	6 mL (30/box)
96-Well Plate			
	10 mg/well	8E-S536-AGA	ea
	30 mg/well	8E-S536-TGA	ea
	60 mg/well	8E-S536-UGA	ea
96-Well Microelution Plate			
	2 mg/well	8M-S536-4GA	ea

Presston[™] 1000 Positive Pressure Manifold

Part No.	Description
AH1-7033	Presston 1000 Positive Pressure Manifold, 96-Well Plate


 Phenomenex warrants the Presston 1000 will be free of defects in materials and workmanship under normal installation, use, and maintenance for a period of 12 months following delivery. Please visit www.phenomenex.com/Presstonwarranty for complete warranty information.

BE-HAPPY[™]
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Kinetex[®] Core-Shell LC Columns

2.6 µm Minibore Columns (mm)						SecurityGuard [™] ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
Biphenyl	00A-4622-AN	00B-4622-AN	–	00D-4622-AN	00F-4622-AN	AJ0-9209
C18	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJ0-8782
						for 2.1 mm ID
2.6 µm MidBore [™] Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
Biphenyl	–	00B-4622-Y0	–	00D-4622-Y0	00F-4622-Y0	AJ0-9208
C18	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJ0-8775
						for 3.0 mm ID
2.6 µm Analytical Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
Biphenyl	–	00B-4622-E0	–	00D-4622-E0	00F-4622-E0	AJ0-9207
C18	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJ0-8768
						for 4.6 mm ID

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

Australia

t: +61 (0)2-9428-6444
 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
 t: +32 (0)2 511 8666 (Dutch)
 beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
 info@phenomenex.com

China

t: +86 400-606-8099
 cninfo@phenomenex.com

Denmark

t: +45 4824 8048
 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
 nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
 franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
 anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
 indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
 eireinfo@phenomenex.com

Italy

t: +39 051 6327511
 italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
 nzinfo@phenomenex.com

Norway

t: +47 810 02 005
 nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
 sginfo@phenomenex.com

Spain

t: +34 91-413-8613
 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
 twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
 ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
 info@phenomenex.com

**All other countries
 Corporate Office USA **

t: +1 (310) 212-0555
 info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Strata and Kinetex are registered trademarks and Preston, MidBore, Be-Happy, and SecurityGuard are trademarks of Phenomenex. Cerilliant is a registered trademark of Cerilliant Corporation. Sartorius and arium are registered trademarks of Sartorius Corporation. Sigma-Aldrich is a registered trademark of Sigma-Aldrich Co, LLC. Triple Quad is a trademark of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license. Agilent is a registered trademark of Agilent Technologies, Inc. Alta-Dena is a registered trademark of Dean Intellectual Property Services II, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.