

TN-1340

# Separation of Teriflunomide and its Organic Impurities per the Proposed USP Monograph in PF 48(5)

Zeshan Aqeel, Alyssa Nodland, and Bryan Tackett, PhD

Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

## Introduction

Teriflunomide is a pyrimidine synthesis inhibitor that inhibits the function of specific immune cells that have been implicated in Multiple Sclerosis (MS). It is related to Leflunomide, a drug used to treat rheumatoid arthritis. This study for Teriflunomide and its organic impurities is based on the newly proposed USP monograph where an end-capped octadecylsilyl silica gel stationary phase is used under gradient conditions. Because there is no existing USP monograph for this drug substance, a new monograph based on validated methods of analysis was proposed in PF 48(5). In this technical note, we report the separation of Teriflunomide and its related organic impurities using a Kinetex™ 2.6 μm C18 column compared to a Ascentis® Express 2.7 μm C18 column originally used in the proposed monograph.

System suitability per USP Monograph for the Teriflunomide Assay requires a tailing factor no more than (NMT) 2.0 and a percent relative standard deviation (%RSD) of NMT 0.85 % for six replicate injections. System suitability per USP Monograph for the Teriflunomide Related Organic Impurities requires resolution no less than (NLT) 1.5 between Teriflunomide related Compound B and Leflunomide Related Compound A, %RSD NMT 5.0 % each for Teriflunomide, Teriflunomide Related Compound B, and Leflunomide Related Compound A, and a Signal-to-Noise (S/N) ratio NLT 10.

According to USP General Chapter <621>, the configuration of the equipment employed may significantly alter the resolution, retention time, and relative retentions described. Differences in system dwell volume can have an impact on the results obtained for gradient methods. Monographs preferably include an isocratic step before the start of the gradient program so that an adaptation can be made to the gradient time points to take account of differences in dwell volume between the system used for analytical procedure development and that actually used for implementation. The system used in this study had a lower dwell volume (0.45 mL) than described in the monograph (1.1 mL)\*, so the time points stated in the gradient table must be replaced by the adapted time points calculated using the following equation:

$$t_c = t - \frac{(D - D_0)}{F}$$

$t_c$  = adapted time point (min)  
 $t$  = time point indicated in the monograph (min)  
 $D$  = Dwell volume (mL)  
 $D_0$  = Dwell volume used for development of the method (mL)  
 $F$  = flow rate (mL/min)

The adapted time points for the Ascentis Express 2.7 μm C18 column would be:

$$t_c = 2 - \frac{(0.45 - 1.1)}{1}$$

$$t_c = 2.65$$

Each time point in the gradient would be adjusted by +0.65 minutes and is shown in the gradient table.

Adjustments to column dimensions for gradient methods will be allowed provided that the L/dp ratio remains constant or within the range between -25 % to +50 % of the prescribed L/dp ratio indicated in the monograph. In this monograph, the indicated column length was 150 mm, and the particle size was 2.7 μm; therefore, the Kinetex 2.6 μm column used here would be an allowed adjustment. When the particle size is changed, the flow rate requires adjustment because smaller-particle columns will require higher linear velocities for the same performance. The flow rate is adjusted for particle size using the following equation:

$$F_2 = F_1 \times \frac{dc_2^2 \times dp_1}{dc_1^2 \times dp_2}$$

$F_1$  = flow rate indicated in the monograph (mL/min)  
 $F_2$  = adjusted flow rate (mL/min)  
 $dc_1$  = internal diameter of the column indicated in the monograph (mm)  
 $dc_2$  = internal diameter of the column used (mm)  
 $dp_1$  = particle size indicated in the monograph (μm)  
 $dp_2$  = particle size of the column used (μm)

The adjusted flow rate for the Kinetex 2.6 μm column would be:

$$F_2 = 1 \times \frac{4.6^2 \times 2.7}{4.6^2 \times 2.6}$$

$$F_2 = 1.04$$

A change in column dimensions, and thus in column volume, impacts the gradient volume which controls selectivity. Gradients are adjusted to the column volume by changing the gradient volume in proportion to the column volume. This applies to every gradient segment volume. The new gradient time for each gradient segment can be calculated using the following equation:

$$t_{g2} = t_{g1} \times \left(\frac{F_1}{F_2}\right) \times \left(\frac{L_2 \times dc_2^2}{L_1 \times dc_1^2}\right)$$

$t_{g1}$  = gradient time indicated in the monograph, or adjusted for dwell volume (min)  
 $t_{g2}$  = adjusted gradient time (min)  
 $F_1$  = flow rate indicated in the monograph (mL/min)  
 $F_2$  = adjusted flow rate (mL/min)  
 $L_1$  = column length indicated in the monograph (mm)  
 $L_2$  = new column length (mm)  
 $dc_1$  = internal diameter of the column indicated in the monograph (mm)  
 $dc_2$  = internal diameter of the column used (mm)

For the second gradient segment, the adjusted gradient time would be:

$$t_{g2} = 2.65 \times \left(\frac{1}{1.04}\right) \times \left(\frac{150 \times 4.6^2}{150 \times 4.6^2}\right)$$

$$t_{g2} = 2.55$$

A gradient adjustment factor can be calculated and used to determine the new gradient segment times using:

$$\text{Gradient adjustment factor} = \left(\frac{t_{g2}}{t_{g1}}\right)$$


$$\text{Gradient adjustment factor} = 0.96$$

The new gradient timetable for the Kinetex 2.6 μm C18 column is shown in the gradient table.

All requirements for System Suitability for Teriflunomide Assay and Organic Impurities were met by all columns.

All solutions were prepared as indicated in the USP Monograph for Teriflunomide. USP Teriflunomide RS and USP Teriflunomide Related Compound B RS were unavailable, so EDQM Teriflunomide CRS (Catalog No. Y0002247) and EDQM Teriflunomide Impurity B CRS (Catalog No. Y0002254) were substituted and purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026 F - 67081 Strasbourg (France). Leflunomide Related Compound A RS (Catalog No. 1357045) was purchased from USP.

Figure 1. Teriflunomide Structure



## LC Conditions

**Column:** Ascentis® Express 2.7 µm C18, 150 x 4.6 mm  
Kinetex™ 2.6 µm C18, 150 x 4.6 mm ([00F-4462-E0](#))

**Mobile Phase:** A: Acetonitrile / Buffer (10:90, v/v)  
B: Buffer / Acetonitrile (10:90, v/v)

Buffer: 3.85 g/L of Ammonium Acetate in Water. Adjust with glacial Acetic Acid to a pH of 5.5.

| Gradient: | Time (min) | Ascentis Express Adjusted Time (min) | Kinetex Adjusted Time (min) | %B |
|-----------|------------|--------------------------------------|-----------------------------|----|
|           | 0          | 0                                    | 0                           | 24 |
|           | 2          | 2.65                                 | 2.55                        | 24 |
|           | 12         | 12.65                                | 12.18                       | 77 |
|           | 13         | 13.65                                | 13.14                       | 24 |
|           | 17         | 17.65                                | 17                          | 24 |

**Flow Rate:** 1.0 mL/min (Ascentis Express)  
1.04 mL/min (Kinetex)

**Injection Volume:** 5 µL

**Temperature:** 40 °C


**Detector:** UV @ 249 nm

**System:** Waters® ACQUITY® H-Class UHPLC

**Table 1.** Preparation of Solutions


| Solution  | Composition   |
|---|---|
| <b>Diluent</b>  | Acetonitrile / Buffer (80:20, v/v)<br>Buffer: 3.85 g/L of Ammonium Acetate in Water. Adjust with glacial Acetic Acid to a pH of 5.5.  |
| <b>Standard Solution (Assay)</b>                      | 0.2 mg / mL of Teriflunomide CRS in Diluent.  |
| <b>Sample Solution (Assay and Organic Impurities)</b> | Same as Standard Solution.  |
| <b>Sensitivity Solution (Organic Impurities)</b>      | 0.1 µg/mL of Teriflunomide CRS in Diluent.  |
| <b>Standard Stock Solution A (Organic Impurities)</b> | 0.04 mg/mL of USP Leflunomide Related Compound A RS in Acetonitrile. Sonicate to dissolve.  |
| <b>Standard Stock Solution B (Organic Impurities)</b> | 0.04 mg/mL of Teriflunomide Related Compound B CRS in Acetonitrile. Sonicate to dissolve.<br>0.2 µg/mL each of Teriflunomide CRS, Teriflunomide Related Compound B CRS, and 0.02 µg/mL of USP Leflunomide Related Compound A RS, prepared as follows:<br>Transfer a suitable amount of Teriflunomide CRS to a suitable volumetric flask. Add Diluent equivalent to 20 % of the flask volume and sonicate to dissolve. Add a suitable volume of Standard Stock Solution A and Standard Stock Solution B and dilute with Diluent to volume. |
| <b>Standard Solution (Organic Impurities)</b>         | Transfer a suitable amount of Teriflunomide CRS to a suitable volumetric flask. Add Diluent equivalent to 20 % of the flask volume and sonicate to dissolve. Add a suitable volume of Standard Stock Solution A and Standard Stock Solution B and dilute with Diluent to volume.  |

**Figure 2.** Standard Solution - Assay



| Peak No. | Analyte       | Retention Time (min) | Area    | Area %RSD | Height    | Symmetry Factor |
|----------|---------------|----------------------|---------|-----------|-----------|-----------------|
| 1        | Teriflunomide | 5.09                 | 1659767 | 0.41      | 339289.67 | 1.59            |

N = 6 Injections




| Peak No. | Analyte       | Retention Time (min) | Area       | Area %RSD | Height    | Symmetry Factor |
|----------|---------------|----------------------|------------|-----------|-----------|-----------------|
| 1        | Teriflunomide | 5.36                 | 1638156.17 | 0.44      | 355644.67 | 1.51            |

N = 6 Injections




**Figure 3.** Standard Solution – Organic Impurities



| Peak No. | Analyte                          | Retention Time (min) | Area     | Area %RSD | Resolution | Symmetry Factor |
|----------|----------------------------------|----------------------|----------|-----------|------------|-----------------|
| 1        | Teriflunomide                    | 5.157                | 3723.5   | 0.352     | -          | 1.14            |
| 2        | Teriflunomide Related Compound B | 7.513                | 4246.167 | 1.041     | 4.498      | 1.098           |
| 3        | Leflunomide Related Compound A   | 7.866                | 6435.33  | 0.759     |            | 1.067           |

N = 3 Injections




| Peak No. | Analyte                          | Retention Time (min) | Area    | Area %RSD | Resolution | Symmetry Factor |
|----------|----------------------------------|----------------------|---------|-----------|------------|-----------------|
| 1        | Teriflunomide                    | 5.157                | 3723.5  | 0.352     | -          | 1.14            |
| 2        | Teriflunomide Related Compound B | 7.513                | 4241    | 1.067     | 4.502      | 1.097           |
| 3        | Leflunomide Related Compound A   | 7.866                | 6410.67 | 0.682     |            | 1.068           |


N = 3 Injections



Figure 4. Sensitivity Solution – Organic Impurities



| Peak No. | Analyte       | Retention Time (min) | Area | S/N Ratio | Symmetry Factor |
|----------|---------------|----------------------|------|-----------|-----------------|
| 1        | Teriflunomide | 5.421                | 608  | 16.28     | 1.24            |



| Peak No. | Analyte       | Retention Time (min) | Area | S/N Ratio | Symmetry Factor |
|----------|---------------|----------------------|------|-----------|-----------------|
| 1        | Teriflunomide | 5.162                | 593  | 15.42     | 1.08            |

### Conclusions

The USP General Chapter <621> now allows for adjustment of gradient methods. However, differences in system dwell volume and column dimensions (L and dp) must be properly accounted for to ensure performance of the method per the original monograph. Initial attempts to replicate the proposed USP monograph method for Teriflunomide using the original column (Ascentis Express 2.7  $\mu\text{m}$  C18) resulted in failure to meet system suitability requirements. After accounting for the lower dwell volume of the system used in this study by increasing the isocratic hold time and subsequent gradient time points, system suitability requirements were met with the original column.

Further adjustments to the flow rate and gradient timetable were required when using the Kinetex 2.6  $\mu\text{m}$  C18 column, with the L/dp ratio for the Kinetex column well within the allowable adjustment per USP <621>. Comparable results to the original column were obtained using the Kinetex 2.6  $\mu\text{m}$  C18 column, demonstrating that this column is a suitable alternative for the proposed USP Teriflunomide monograph for Assay and Organic Impurities. In this technical note we have illustrated the various adjustments that must be made to successfully adjust for system differences and column dimensions.



## Kinetex™ Ordering Information

| 2.6 µm Analytical Columns (mm) |                             |                             |                             |                             |                             |                             | SecurityGuard™<br>ULTRA Cartridges* |
|--------------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-------------------------------------|
| Phases                         | 30 x 4.6                    | 50 x 4.6                    | 75 x 4.6                    | 100 x 4.6                   | 150 x 4.6                   | 250 x 4.6                   | 3/pk                                |
| EVO C18                        | <a href="#">00A-4725-E0</a> | <a href="#">00B-4725-E0</a> | —                           | <a href="#">00D-4725-E0</a> | <a href="#">00F-4725-E0</a> | <a href="#">00G-4725-E0</a> | <a href="#">AJ0-9296</a>            |
| PS C18                         | <a href="#">00A-4780-E0</a> | <a href="#">00B-4780-E0</a> | —                           | <a href="#">00D-4780-E0</a> | <a href="#">00F-4780-E0</a> | <a href="#">00G-4780-E0</a> | <a href="#">AJ0-8949</a>            |
| Polar C18                      | <a href="#">00A-4759-E0</a> | <a href="#">00B-4759-E0</a> | —                           | <a href="#">00D-4759-E0</a> | <a href="#">00F-4759-E0</a> | —                           | <a href="#">AJ0-9530</a>            |
| Biphenyl                       | —                           | <a href="#">00B-4622-E0</a> | —                           | <a href="#">00D-4622-E0</a> | <a href="#">00F-4622-E0</a> | —                           | <a href="#">AJ0-9207</a>            |
| XB-C18                         | —                           | <a href="#">00B-4496-E0</a> | <a href="#">00C-4496-E0</a> | <a href="#">00D-4496-E0</a> | <a href="#">00F-4496-E0</a> | —                           | <a href="#">AJ0-8768</a>            |
| C18                            | <a href="#">00A-4462-E0</a> | <a href="#">00B-4462-E0</a> | <a href="#">00C-4462-E0</a> | <a href="#">00D-4462-E0</a> | <a href="#">00F-4462-E0</a> | —                           | <a href="#">AJ0-8768</a>            |
| C8                             | —                           | <a href="#">00B-4497-E0</a> | <a href="#">00C-4497-E0</a> | <a href="#">00D-4497-E0</a> | <a href="#">00F-4497-E0</a> | —                           | <a href="#">AJ0-8770</a>            |
| HILIC                          | —                           | <a href="#">00B-4461-E0</a> | <a href="#">00C-4461-E0</a> | <a href="#">00D-4461-E0</a> | <a href="#">00F-4461-E0</a> | —                           | <a href="#">AJ0-8772</a>            |
| Phenyl-Hexyl                   | —                           | <a href="#">00B-4495-E0</a> | <a href="#">00C-4495-E0</a> | <a href="#">00D-4495-E0</a> | <a href="#">00F-4495-E0</a> | —                           | <a href="#">AJ0-8774</a>            |
| F5                             | <a href="#">00A-4723-E0</a> | <a href="#">00B-4723-E0</a> | —                           | <a href="#">00D-4723-E0</a> | <a href="#">00F-4723-E0</a> | —                           | <a href="#">AJ0-9320</a>            |

for 4.6 mm ID

\*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**🌐 All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/phx-terms-and-conditions-of-sale](http://www.phenomenex.com/phx-terms-and-conditions-of-sale).

**Trademarks**

Kinetex, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation. Ascendis is a registered trademark of Merck KGaA, Darmstadt, Germany.

**Disclaimer**

Comparative separations may not be representative of all applications.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

*CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.*

Phenomenex is in no way affiliated with Waters Technologies Corporation or Merck KGaA, Darmstadt, Germany.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

