

Separation of Fluvoxamine Maleate and its Organic Impurities per USP Monograph

Lauren Nakasone, Zeshan Aqeel, and Bryan Tackett, PhD
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA


Overview

Fluvoxamine Maleate is a selective serotonin reuptake inhibitor (SSRI) in the brain. The development of a quick and efficient analysis of Fluvoxamine Maleate and its related organic impurities is of interest for generic drug manufacturers. In this application note, we report the separation of Fluvoxamine Maleate and its related organic impurities using a Luna 5 μm C8(2) column according to the USP monograph for Fluvoxamine Maleate.

System suitability per USP Monograph for the Fluvoxamine Maleate is a symmetry factor no more than (NMT) 2.0 and a percent relative standard deviation (%RSD) of NMT 2.0 % for Fluvoxamine Maleate. A resolution no less than (NLT) 3.0 between the Z-isomer and Fluvoxamine Maleate, and a resolution NLT 5.0 between Succinyl Fluvoxamine and the Z-isomer are also requirements of system suitability. All system suitability requirements for Fluvoxamine Maleate were met with the Luna C8(2) column.

All solutions were prepared as indicated in the USP Monograph for Fluvoxamine Maleate. USP Fluvoxamine Maleate RS (Catalog No. 1285909) was purchased from USP, and Maleic Acid was purchased from Sigma-Aldrich®.

Figure 1. Fluvoxamine Maleate



LC-UV Conditions

Column: Luna™ 5 μm C8(2) ([00G-4249-E0](#))

Dimension: 250 x 4.6 mm

Mobile Phase: Acetonitrile/Solution A (38:62, v/v)

Solution A: 8 g/L of Sodium 1-Pentanesulfonate and 1.1 g/L of Monobasic Potassium Phosphate in water. Adjust with Phosphoric Acid to a pH of 3.00 ± 0.05 .

Flow Rate: 1.7 mL/min (Isocratic)

Injection Volume: 20 μL

Temperature: 40 °C

Detector: UV @ 234 nm


System: Agilent® 1260 Binary UHPLC

Table 1. Preparation of Solutions

Solution	Composition
Standard Solution	0.05 mg/mL of USP Fluvoxamine Maleate RS in mobile phase
System Suitability Solution	Transfer 6 mg of USP Fluvoxamine Maleate RS to a 50 mL volumetric flask. Heat the sample at 120 °C for 10 min. Cool to room temperature and add 3.0 mL of 0.1 N Hydrochloric Acid. Heat the solution in a water bath for 10 min. Cool to room temperature, add 50 mg of USP Fluvoxamine Maleate RS, and dissolve in 25 mL of mobile phase. Dilute with mobile phase to volume.
Identification Solution	0.35 mg/mL of Maleic Acid in mobile phase.




Figure 2. Standard Solution




Peak No.	Analyte	Area %RSD	Symmetry Factor
1	Maleic Acid	-	-
2	Fluvoxamine Maleate	0.072	1.127
N = 6 Injections			

Figure 3. System Suitability Solution



Peak No.	Analyte	Resolution
1	Maleic Acid	-
2	Succinyl Fluvoxamine	-
3	Fluvoxamine Z-isomer	4.41
4	Fluvoxamine Maleate	-
5	Dealkyl benzyl Fluvoxamine	-
N = 6 Injections		

Figure 4. Identification Solution



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☉ All other countries/regions

Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. Sigma-Aldrich is a registered trademark of Merck KGaA, Darmstadt, Germany.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Inc. or Merck KGaA, Darmstadt, Germany.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

