

APPLICATIONS

Development of a Highly Selective Bioanalytical Quantitation Method for the LC-MS/MS Analysis of (R)-Amlodipine and (S)-Amlodipine Enantiomers in Human Plasma Using Lux[®] Cellulose-4 and Strata[®]-X

Rahul Baghla¹, Ian Moore¹, Sahana Mollah¹, Namrata Saxena², Phil Koerner², and Marc Jacob²

¹SCIEX, Redwood City, CA 94065;

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

In this technical note we report the development of an improved selective and rugged validated LC-MS/MS method for the quantitation of Amlodipine enantiomers in human plasma samples. Sample extraction utilized Strata-X polymeric solid phase extraction cartridges and chromatographic separation of (R)-Amlodipine and (S)-Amlodipine enantiomers was achieved on a polysaccharide-based chiral stationary phase, Lux Cellulose-4, under polar organic conditions. This method can be used in pre-clinical and bioequivalence studies for rapid analysis of both pharmacologically active and inactive enantiomers of Amlodipine in human plasma.


Introduction

Amlodipine is a calcium channel blocker and although (S)-Amlodipine is the pharmacologically active enantiomer, a racemic mixture of Amlodipine is used for therapeutic purposes. (S)-Amlodipine and (R)-Amlodipine exhibit different pharmacokinetics (PK) and pharmacodynamics (PD); and according to regulatory guidelines^{1,2} individual enantiomers should be measured in cases where enantiomers exhibit different pharmacokinetic or pharmacodynamic properties, or the exposure ratio of enantiomers is modified by a difference in the rate of absorption. If one enantiomer is pharmacologically active and the other is inactive or has low contribution to activity, it is sufficient to demonstrate bioequivalence (BE) for the active enantiomer; in this case (S)-Amlodipine. This validated LC-MS/MS method was developed to monitor concentrations of (S)-Amlodipine and (R)-Amlodipine in human plasma samples to evaluate the pharmaceutical equivalence for both racemic and (S)-Amlodipine formulations of Amlodipine.

Method Development

Several methods have been published for the separation of (S)-Amlodipine and (R)-Amlodipine in various biological matrices³⁻¹³. However, the major challenge in developing an adequate LC-MS/MS method was in achieving good peak shape for the Amlodipine enantiomers. Published methods using various chiral columns and mobile phase with various basic mobile phase additives such as ammonium hydroxide, trimethylamine and diethylamine, typically gave poor peak shapes and required long analysis times. During method development we tried various additives with 0.05% (v/v) ethanolamine providing the best peak shapes and signal intensity for amlodipine isomers. It is important to set up a rapid, selective and sensitive LC-MS/MS method for bioanalysis approach to ensure wide concentration range (50 pg/mL to 50 ng/mL) to cover PK and BE studies for all dosage forms ranging from 2.5 mg to 10 mg per day.

Structures of (S)-Amlodipine and (R)-Amlodipine



To minimize potential matrix interferences compromising quantitation due to ion suppression, Strata-X SPE tubes were used for sample preparation of the human plasma samples prior to LC-MS/MS analysis.

Reagents and Solutions

Analytical Standards of (S)-Amlodipine, (R)-Amlodipine, (S)-Amlodipine-d4 and (R)-Amlodipine-d4 were purchased from Toronto Research Chemicals, Canada with the purity of 99%. Acetonitrile, Water, Methanol, Ethanolamine, Formic acid and Isopropyl alcohol were purchased from Sigma Aldrich.

A stock solution of (S)-Amlodipine and (R)-Amlodipine enantiomers was prepared in methanol at 1 mg/mL. Working standard solutions at various concentrations for spiking in plasma were prepared with 50:50 (Methanol/Water). Working stock solution of internal standard was prepared at 100 ng/mL.

Sample Preparation

The extraction of Amlodipine enantiomers from human plasma was carried out using Strata-X cartridges per the SPE method described below:

Spiked Plasma Samples: 5 µL of internal standard working solution was added to 100 µL spiked human plasma (spiked at 0.050, 0.150, 20.0, 40.0, and 50.0 ng/mL) and mixed, 500 µL of 0.2% (v/v) Ethanolamine in water was added to sample and vortexed.

SPE Protocol:

Cartridge: Strata-X, 30 mg/1 mL tubes
Part Number: 8B-S100-TAK
Condition: 1 mL Methanol
Equilibrate: 1 mL Water
Load: Spiked plasma sample was loaded and allowed to pass through the cartridge at moderate speed
Wash 1: 1 mL Water (to remove polar interferences)
Wash 2: 1 mL of 20% Methanol in Water (to remove non-polar interferences)
Elute: 1 mL of 0.1% Formic Acid in Methanol
Dry: Take to dryness under a gentle stream of Nitrogen at 50 °C
Reconstitute: 100 µL of Mobile Phase

LC-MS/MS Conditions:

Column: Lux 3 µm Cellulose-4
Dimensions: 150 x 2.0 mm
Part No.: 00F-4490-B0
Mobile Phase: 0.05% Ethanolamine in Acetonitrile and Isopropyl Alcohol (96:4 v/v)
LC condition: Isocratic
Flow Rate: 0.3 mL/min
Run time: 4 min
Temperature: 25 °C
Injection Volume: 10 µL
HPLC System: SCIEX[®] ExionLC[™] AD HPLC (AB SCIEX PTE, Ltd.)
Mass Spectrometer: SCIEX QTRAP[®] 4500, ESI positive (AB SCIEX PTE, Ltd.)
Scan Type: MRM
IS Voltage: 5500V
Source Temperature: 300 °C
Collision Energy: 15V

QTRAP 4500 Compound Dependent Parameters

Compound	Precursor m/z	Fragment m/z
(R)-Amlodipine	409.3	237.9
(S)-Amlodipine	409.3	237.9
(R)-Amlodipine-d4	413.1	237.9
(S)-Amlodipine-d4	413.1	237.9

Results and Discussion

(S)-Amlodipine and (R)-Amlodipine enantiomers were analyzed on a Lux Cellulose-4 chiral column with the simple isocratic mobile phase conditions shown and detected using electrospray ionization (ESI) positive acquisition mode with multiple reaction monitoring (MRM) scan type using the SCIEX QTRAP® 4500. The optimized LC-MS/MS conditions and MS/MS compound dependent parameters are shown above. The chromatogram for the separation on the Lux Cellulose-4 column is shown in **Figure 1**; full baseline separation with excellent peak shape was achieved in less than 4 minutes.

The developed LC-MS/MS method was fully validated, and these results will be discussed below.

System suitability was assessed every day at the start of the experiment by injecting six replicates of analyte and internal standard mixture with concentration 1 ng/mL to ensure instrument performance at the start of each day. The criteria for system suitability was established at $\leq 5\%$ for the analyte-to-internal standard peak area, and retention time deviation for both the analytes and ISTD was less than 2 % for system suitability. Chromatograms of (R)-Amlodipine and (S)-Amlodipine and the respective deuterated internal standards is shown in **Figure 2**.

Calibration curves were generated by spiking analytes into blank human plasma at 8 different concentration levels ranging from 0.050 ng/mL to 50 ng/mL (**Figure 3**). The calibration curve was found to be linear over the specified range for both (R)-Amlodipine and (S)-Amlodipine enantiomers. Excellent linearity across the calibration range was exhibited, with the regression coefficient (r) greater than 0.99 for all calibration curves generated from the accuracy and precision sample batches.

For specificity and carryover studies eight different lots of blank plasma were extracted and analysed using the LC-MS/MS method and injected along with the corresponding lots of plasma spiked at the lower limit of quantitation (LLOQ) and working internal standard concentration. Carryover was evaluated by injecting the following sample sets: extracted plasma blank, LLOQ, ULOQ and extracted plasma blank. The percentage interference at the retention time of the analyte was monitored and calculated against the LLOQ area response and internal standard area response. Analytes and internal standard showed high degree of selectivity and specificity in the proposed method. Blank samples injected after ULOQ samples did not show quantifiable response ($<20\%$ of the LLOQ peak area response), which demonstrates that the proposed method is free of carryover. Chromatograms of extracted plasma blank and plasma spiked LLOQ QC sample are shown in **Figure 4**, and a summary of carryover data is presented in **Table 1**.

Table 1. Evaluation of Carry Over

Sample	Precursor m/z (R)-Amlodipine Peak Area	Fragment m/z (S)-Amlodipine Peak Area
Plasma Blank	22	66
LLOQ	3761	3706
ULOQ	6536117	6468086
Plasma Blank	82	119
Below 20% of LLOQ area response		

Figure 1.
Chromatogram for Amlodipine Enantiomers on Lux Cellulose-4




Figure 2.
Chromatogram for Amlodipine Enantiomers and Internal Standards




Figure 3.
Calibration Curves for (R)-Amlodipine (blue) and (S)-Amlodipine (purple)





Figure 4.
Chromatograms for Plasma Blank (left) and Plasma Spiked at LLOQ QC Sample (right); (m/z 409.3 \rightarrow 237.9)



Intraday and interday precision and accuracy was evaluated using 6 replicates of extracted quality control (QC) samples in human plasma at LLOQ, low QC (LQC), mid QC (MQC), and high QC (HQC) as follows: LLOQ QC (0.050 ng/mL), LQC (0.150 ng/mL), MQC (20.000 ng/mL) and HQC (40.000 ng/mL). According to the US FDA Regulatory guidelines, mean % nominal concentration at each QC level, other than LLOQ QC, must be within 85 % to 115 % and precision (% CV) should be $\leq 15\%$. Mean % of nominal concentration at LLOQ QC level must be within 80 % to 120 % and the precision (% CV) should be $\leq 20\%$. Intraday and interday precision and accuracy was studied and met both the precision and accuracy requirements (**Table 2**).

Table 2. Evaluation of Interday and Intraday Accuracy and Precision

Interday Precision and Accuracy, n = 6				Intraday Precision and Accuracy, n = 6		
(R)-Amlodipine						
Nominal Conc., ng/mL	Measured Conc., ng/mL	Precision (% CV)	Accuracy (%)	Measured Conc., ng/mL	Precision (% CV)	Accuracy (%)
0.050	0.048	6.00	96	0.047	8.66	94
0.150	0.151	4.05	101	0.149	3.47	99
20	20.399	3.11	102	20.576	2.46	103
40	41.922	2.79	105	42.410	2.74	106
(S)-Amlodipine						
0.050	0.048	9.67	96	0.047	8.51	94
0.150	0.152	3.10	101	0.151	4.64	101
20	20.812	2.89	104	20.693	2.38	103
40	42.599	3.19	107	42.279	2.98	106

A significant benefit associated with SPE for sample preparation is the removal of endogenous compounds present in plasma, such as proteins and phospholipids, prior to LC-MS/MS analysis. The presence of these matrix compounds in injected samples can impact analyte peak shapes, HPLC column lifetimes, and negatively affect accurate quantitation with MS detection due to potential ion suppression effects. The polymeric Strata[®]-X sorbent used in this method allows for both a simple water wash to remove polar endogenous compounds (sugars, salts, etc), as well as a stronger wash (20 % methanol in this case) to help remove some of the non-polar matrix compounds while retaining the analytes of interest (improved recoveries). To confirm that matrix effect was minimized for this method, the matrix effect was assessed by measuring the corresponding analyte peak area response in reconstituted matrix samples from 6 different plasma lots against the analyte mean peak area response in reference solution at 2 different concentration levels (low and high QC concentrations). Matrix effect was also studied for all 6 replicates and the developed method showed no matrix effect.

Extraction recoveries for both (R)-Amlodipine and (S)-Amlodipine at three concentration levels were determined by measuring the mean peak area response of 6 replicates of extracted QC samples and comparing to the mean peak area response of extracted blank matrix that was spiked after preparation with the nominal amount of both analytes at the low, medium and high QC levels. The extraction recoveries of analytes were uniform and reproducible, with 94.14 % extraction recovery for (R)-Amlodipine and 92.23 % for (S)-Amlodipine. This confirms that the stronger second solvent wash (20 % methanol) did not negatively impact recovery of (R)-Amlodipine and (S)-Amlodipine.

For this method, the upper level of quantitation (ULOQ) for (R)-Amlodipine and (S)-Amlodipine was found to be 50ng/mL. A dilution integrity test was performed by spiking 1.6x concentrations of ULOQ in plasma samples and diluted 4-fold and 2-fold with blank plasma. The average percentage recoveries for 4-fold and 2-fold dilutions for 6 replicate injections each were 94.90 % for (R)-Amlodipine and 96.40 % for (S)-Amlodipine.

Stability studies in human plasma were evaluated to cover a wide range of expected experimental and storage conditions. Analytes and internal standard were found to be stable on the benchtop for 6 hours, in the autosampler for 48 hours and for 3 freeze thaw cycles of stability studies. Sample stock solutions were also found to be stable when refrigerated at 2 °C to 8 °C for 4 days.

Conclusions

This technical note highlights the selective and sensitive LC-MS/MS method which was developed for the separation and quantitation of (R)-Amlodipine and (S)-Amlodipine enantiomers in human plasma using a Lux[®] Cellulose-4 chiral column and SCIEX QTRAP[®] 4500 mass spec system. A sample preparation method for extraction of Amlodipine from human plasma, with high analyte recoveries, was developed using polymeric Strata-X SPE tubes, reducing potential matrix interference from proteins and phospholipids that could impact accuracy of the quantitative method. The Lux Cellulose-4 column provided baseline chromatographic separation of the (R) and (S) enantiomers of Amlodipine with excellent peak shape, and the SCIEX QTRAP 4500 provides the sensitivity required to detect the Amlodipine enantiomers at therapeutic levels in human plasma. This method is unique in that it can quantify both pharmacologically active and inactive enantiomers of Amlodipine in human plasma samples with LLOQ of 0.050 ng/mL. We believe that this simple, sensitive and reproducible method for the quantitation of (R)-Amlodipine and (S)-Amlodipine in human plasma will have application in bioequivalence, bioavailability and preclinical studies. We anticipate that this method can be easily extended to other preclinical species with little or no modification.

Acknowledgement: the original SCIEX Technical Note can be found at <https://sciex.com>

References

- Guidance for Industry Bioanalytical Method Validation. US Department of Health and Human Services, US FDA, Rockville, MD, USA (2001).
- Guideline on bioanalytical method validation. European Medicines Agency, London, UK (2011).
- Jaivik V. Shah, Jignesh M. Parekh, Priyanka A. Shah, Priya V. Shah, Mallika Sanyal, Pranav S. Shrivastav; *J. Pharm. Anal.* 7(5), 309-316 (2017).
- Hye Hyun Yoo, Tae Kon Kim, Bong-Yong Lee, and Dong-Hyun Kim; *Mass Spectrom. Lett.* 2(4), 88-91 (2011).
- Katina S. S. Dossou, Patrick A. Etorh, Patrice Chiap, Bezhan Chankvetadze, Anne-Catherine Servais, Marianne Fillet, Jacques Crommen; *J. Sep. Sci.* 34 (15), 1772-1780 (2011).
- Massaroti, P., Moraes, L. A., Marchioretto, M. A., Cassiano, N. M., Bernasconi, G., Calafatti, S. A., Barros, F. A., Meurer, E. C., Pedrazzoli, J; *Anal. Bioanal. Chem.* 382(4), 1049-1054 (2005).
- Zou Q, Zhan Y, Ge Z, Wei P, Ouyang P; *Arzneimittelforschung* 59(8), 383-391 (2009).
- Alsarra IA; *J Chromatogr. Sci.* 47(10), 863-867 (2009).
- Hotha KK, Roychowdhury S, Mullangi R, Ravindranath LK; *Biomed Chromatogr.* 27(9), 1192-1199 (2013).
- Taguchi R, Naito T, Sato H, Kawakami J; *Ther Drug Monit.* 39(6), 625-631 (2017).
- Chan-Mei Lv, Chun-Min Wei, Fan-Long Bu, Rui Chen, Xiao-Lu Wang, Rong Li, Ben-Jie Wang, Rui-Chen Guo; *Pharmacology & Pharmacy* 4(2), 191-200 (2013).
- Somagoni J, Reddy S, Koorelli S, Manda S, Yamsani MR; *Pharm Anal Acta.* 2(5) (2011).
- Jiaqi Xie, Qi Tan, Lin Yang, Shenzhi Lai, Shaotan Tang, Changqun Cai, Xiaoming Chen; *Analytical Methods* 6(12), 4408-4413 (2014).

APPLICATIONS

Lux[®] Chiral LC & SFC Columns

Ordering Information

3µm Minibore, MidBore, and Analytical Columns (mm)									SecurityGuard [™] Cartridges (mm)	
Phases	50 x 2.0	150 x 2.0	100 x 3.0	150 x 3.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 2.0*	4 x 3.0*
									/10pk	/10pk
i-Cellulose-5	00B-4755-B0	00F-4755-B0	00D-4755-Y0	00F-4755-Y0	00B-4755-E0	00D-4755-E0	00F-4755-E0	00G-4755-E0	AJO-8631	AJO-8632
Cellulose-1	00B-4458-B0	00F-4458-B0	00D-4458-Y0	00F-4458-Y0	00B-4458-E0	00D-4458-E0	00F-4458-E0	00G-4458-E0	AJO-8402	AJO-8403
Cellulose-2	00B-4456-B0	00F-4456-B0	00D-4456-Y0	00F-4456-Y0	00B-4456-E0	00D-4456-E0	00F-4456-E0	00G-4456-E0	AJO-8398	AJO-8366
Cellulose-3	00B-4492-B0	00F-4492-B0	00D-4492-Y0	00F-4492-Y0	00B-4492-E0	00D-4492-E0	00F-4492-E0	00G-4492-E0	AJO-8621	AJO-8622
Cellulose-4	00B-4490-B0	00F-4490-B0	00D-4490-Y0	00F-4490-Y0	00B-4490-E0	00D-4490-E0	00F-4490-E0	00G-4490-E0	AJO-8626	AJO-8627
Amylose-1	00B-4729-B0	00F-4729-B0	00D-4729-Y0	00F-4729-Y0	00B-4729-E0	00D-4729-E0	00F-4729-E0	00G-4729-E0	AJO-9337	AJO-9336
Amylose-2	00B-4471-B0	00F-4471-B0	00D-4471-Y0	00F-4471-Y0	00B-4471-E0	00D-4471-E0	00F-4471-E0	00G-4471-E0	AJO-8471	AJO-8470

for ID: 2.0–3.0 mm 3.2–8.0 mm

Strata-X

Ordering Information

Format	Sorbent Mass	Part Number	Unit
Tube	30 mg	8B-S100-TAK**	1 mL (100/box)

*SecurityGuard Analytical Cartridges require holder, Part No.: KJO-4282

**Tab-less tubes available. Contact Phenomenex for details.



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Lux and Strata are registered trademarks of Phenomenex. SCIEX is a registered trademark of AB SCIEX Pte. Ltd. AB SCIEX[™] is being used under license.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Agilent Technologies, Inc. or Sigma-Aldrich.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
esinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

All other countries
Corporate Office USA 

t: +1 (310) 212-0555
info@phenomenex.com



phenomenex[®]
...breaking with traditionSM

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com