



## Irbesartan Related Substances Method per IP Monograph

Krishna Chaitanya Routhu<sup>1</sup>, Swetha Kotikalapudi<sup>1</sup>, and Bryan Tackett, PhD<sup>2</sup>

<sup>1</sup>India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India

<sup>2</sup>Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview

Irbesartan is a potent selective angiotensin II type 1 (AT1) receptor antagonist and is used in the treatment of hypertensive patients and diabetic nephropathy. The development of a quick and efficient analysis of Irbesartan and its related substances is of interest for India generic drug manufacturers. In this application note, we report the separation of Irbesartan and its related substances using a Luna™ 5 µm C18(2) column according to the Indian Pharmacopoeia (IP) Monograph for Irbesartan, which references use of a 25 cm x 4.0 mm stainless-steel column packed with octadecylsilane bonded to porous silica (5 µm).


System suitability per IP Monograph for the Irbesartan Related Substances is resolution of not less than (NLT) 3.0 between Irbesartan and Irbesartan Impurity A.

The results clearly show that the system suitability criteria (resolution) met and surpassed the minimum requirement. This illustrates that the Luna 5 µm C18(2) meet all requirements for system suitability as outlined in the Indian Pharmacopoeia monograph for Irbesartan.

All solutions were prepared as indicated in the IP Monograph for Irbesartan. The following certified reference standards (CRS) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: Allee Kastner CS 30026 F - 67081 Strasbourg (France):

- Y0001166, Irbesartan CRS
- Y0001156, Irbesartan Impurity A CRS

**Figure 1.** Irbesartan



### LC-UV Conditions

**Column:** Luna 5 µm C18(2)

**Dimensions:** 250 x 4.0 mm

**Part No.:** [00G-4252-D0](#)

**Mobile Phase:** A mixture of 67 volumes of buffer solution prepared by diluting 5.5 mL of Orthophosphoric Acid in 950 mL of Water, adjusted to pH 3.2 with Triethylamine and 33 volumes of Acetonitrile.

**Flow Rate:** 1 mL/min (Isocratic)

**Injection Volume:** 10 µL

**Temperature:** 25 °C

**Detector:** UV @ 220 nm


**System:** Waters® ACQUITY Arc® HPLC

**Table 1.** Preparation of Solutions


Solution	Composition
Test Solution	Dissolved 50 mg of Irbesartan RS in Methanol, dilute to 50.0 ml with Methanol.
Reference Solution (a)	Diluted 1.0 ml of the Test Solution to 20.0 ml with Methanol. Diluted 1.0 ml of this solution to 50.0 ml with Methanol.
Reference Solution (b)	Dissolved 5 mg of Irbesartan RS and 5 mg of Irbesartan impurity A CRS in Methanol and diluted to 10 mL with the same solvent. Diluted 1 mL of the solution to 10 mL with Methanol.




**Figure 2. Reference Solution (a)**



**Figure 3. Reference Solution (b)**



**Figure 3. Test Solution**



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™  
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/phx-terms-and-conditions-of-sales](http://www.phenomenex.com/phx-terms-and-conditions-of-sales).

**Trademarks**

Luna and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY Arc are registered trademarks of Waters Technologies Corporation.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

