


# Preparative Chromatography Meets UHPLC

**The Power of Kinetex Core-Shell Technology with  
the Performance of Piston Packed Prep Columns**

- **Ultra-High Efficiencies**
- **Ultra-Long Lifetimes**
- **NO Added Pressure**



# Two Award Winning Technologies



**NOW COMBINED** to give you unmatched purification performance for both HPLC and SFC applications.

## KINETEX<sup>®</sup> Core-Shell Technology

Kinetex Core-Shell Technology produces increased efficiencies over traditional, fully porous columns, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity, so labs can get even more out of their HPLC analyses!

The benefits of Kinetex Core-Shell Technology include:

- Increased efficiencies over traditional fully porous columns
- Seamless scalability from HPLC/UHPLC to Preparative LC
- Kinetex 5  $\mu\text{m}$  provides better performance than traditional fully porous 5 and 3  $\mu\text{m}$  materials



## AXIA<sup>™</sup> Preparative Column Packing Technology

An advanced preparative column packing and hardware design, Axia incorporates patented Hydraulic Piston Compression technology that offers increased sorbent bed density and eliminates media bed collapse as a source of premature column failure in preparative HPLC columns.

Unlike traditional column packing methods, the Axia packing method is completely automated and monitored by multiple sensors to allow for measurement and recording of all process parameters for every column. The result is a vastly improved packing process that offers the following benefits:


- Extended column lifetimes
- Improved reproducibility: Column-to-Column and Batch-to-Batch
- Efficiencies and peak symmetries on par with analytical separations
- Increased column stability under high flow rates



# Combining These Two Award Winning Technologies Gives You...


## Higher Efficiency!

Start with sharper peaks by taking advantage of the high efficiencies of Kinetex® 5 µm Axia™ preparative columns.



### Conditions for both columns:


**Columns:** Kinetex 5 µm XB-C18 Axia Packed  
XBridge 5 µm C18 Prep OBD  
**Dimensions:** 150 x 21.2 mm (Kinetex)  
150 x 19 mm (XBridge)  
**Mobile Phase:** Water/ Acetonitrile (50:50)  
**Injection Volume:** 10 µL



**Flow Rate:** 25 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**Sample:** 1. Uracil  
2. Acetophenone  
3. Toluene  
4. Naphthalene

## Excellent Loadability!


With narrower peak widths than fully porous columns across every sample load, Axia packed Kinetex 5 µm columns give you the capability of increased sample load and higher throughput for vastly improved purification performance and economics.



### Conditions for both columns:

**Columns:** Kinetex 5 µm C18 Axia Packed  
XBridge 5 µm C18 Prep OBD  
**Dimensions:** 50 x 21.2 mm (Kinetex)  
50 x 19 mm (XBridge)  
**Mobile Phase:** A: Water with 0.5 % Formic acid  
B: Acetonitrile with 0.5 % Formic acid  
**Gradient:**

Time (min)	% B
0	20
8	50
11	50



**Increased Sensitivity and Excellent Loadability!**

**Flow Rate:** 30 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**Sample:** 200 mg/mL in DMSO  
1. Doxepin (From 1 - 500 mg on-column)  
2. Amitriptyline (From 1 - 500 mg on-column)

Waters is a registered trademark of Waters Corp. OBD and XBridge are trademarks of Waters Corp.  
Phenomenex is in no way affiliated with Waters Corp. Comparative separations may not be representative of all applications.



# Preparative Technology Redefined

## Axia™ Technology Vs. Traditional Prep Column Packing

### Waters® “OBD” Patented Prep Column Packing Process:

In traditional slurry packing processes, like the Waters® OBD™ (Optimum Bed Density) Prep column packing approach, pressure on the packed bed is released when the column is removed from the column packing station to allow attachment of the endfitting.

This conventional packing process involves:

**Compression → Decompression → Recompression**

Several problems with this packing method are:

- Variability in column performance due to increased number of manual operations required for assembly
- Potential silica media damage during recompression
- Limited level of process control is based on traditional slurry packing technology





Diagram from Waters Corporation U.S. Patent No. 7,399,410

### Our Vastly Improved Patented Packing Process:

In contrast to the multi-step process required in conventional slurry packing, Axia packed preparative columns are packed using a single axial compression step. The ideal column bed density is custom calculated and automated for each specific media and column size. Computer control of the entire process ensures both proper bed density and column uniformity every time.

During the Axia packing process, the packing piston is locked in place, eliminating any decompression and recompression of the packed bed. This improves media integrity and column bed stability, and solves the primary lifetime and performance problems associated with conventional slurry packed preparative columns.



Axia patent: U.S. Patent No. 7,674,383

Waters is a registered trademark of Waters Corp. OBD is a trademark of Waters Corp. Phenomenex is in no way affiliated with Waters Corp.


# The Packing Difference You Can See!

## Traditional Prep Packing

Decompression and then recompression during packing can damage the media and lead to increased column-to-column variability, flow disturbances, and decreased column lifetimes.




\*SEM of Waters® OBD™ inlet frit



Crushed media or silica fines at frit surface after packing

## Axia Packing Technology

Highly tuned patented process and hardware eliminates potential decompression and ensures media stability and optimal bed density.



\*SEM of Axia inlet frit



Intact media at frit surface after packing

\*The images are believed to be representative, but individual columns may vary.

View an animated packing process comparison [www.AxiaPrep.com](http://www.AxiaPrep.com)


## Unmatched Column Reproducibility

The completely automated Axia packing system provides feedback control and infinite tuning of packing density for specific media characteristics such as mechanical strength and porosity. An optimum bed density can be consistently reproduced column-to-column. This directly translates into consistent efficiency and peak asymmetry measurements and decreases the column variability seen in traditionally packed preparative columns.

### Reproducible

#### Column-to-Column Efficiency


Average Efficiency (N) with Synergi™  
4 µm Hydro-RP 100 x 21.2 mm




### Reproducible

#### Column-to-Column Peak Asymmetry

Average Peak Asymmetry with Gemini®  
5 µm C18 50 x 21.2 mm



### Density Comparison of Packed Beds




Waters is a registered trademark of Waters Corp. OBD is a trademark of Waters Corp. Phenomenex is in no way affiliated with Waters Corp. Comparative separations may not be representative of all applications.


# Seamless Scalability from HPLC/UHPLC to PREP

The recent addition of the Kinetex® 5 µm in the Axia packed format (21.2 mm ID) makes it the first core-shell sorbent commercially available for small-scale preparative applications. Combine this with the added flexibility that the entire Kinetex core-shell line (1.3 µm, 1.7 µm, 2.6 µm and 5 µm) is fully scalable in retention and selectivity, makes transferring high performance HPLC/UHPLC methods to preparative and SFC applications, simple.


## Analytical method — Kinetex 2.6µm XB-C18



## Preparative scale-up and fraction collection — Kinetex 5µm XB-C18



## Analytical fraction analysis — Kinetex 2.6µm XB-C18



### Tip:



For more information on the power of Kinetex core-shell scalability, request technical note TN-1135 at:

[www.phenomenex.com/Kinetex/AxiaRequest](http://www.phenomenex.com/Kinetex/AxiaRequest)





# A Broad Spectrum of Column Selectivities



Kinetex® core-shell columns are available in a wide range of stationary phases, allowing you to optimize your separation for maximum resolution and loadability across HPLC, UHPLC, and Preparative HPLC and SFC applications.





**Endcapped C18 phase**  
Increased retention for polar basic compounds



**Protective isobutyl side chains**  
Increased retention of polar acidic compounds



**Phenyl-Hexyl**  
Greater retention and separation of aromatic hydrocarbons



**Pentafluorophenyl phase**  
Unique aromatic and polar selectivity

## Conditions for all columns:

**Column:** Kinetex 2.6 µm C18  
Kinetex 2.6 µm XB-C18  
Kinetex 2.6 µm Phenyl-Hexyl  
Kinetex 2.6 µm PFP

**Dimensions:** 50 x 2.1 mm

**Mobile Phase:** A: 0.1% Formic acid in Water  
B: 0.1% Formic acid in Acetonitrile

Gradient:	Time (min)	% B
	0.0	5
	0.2	5
	4.2	95
	4.21	5
	5.5	5

**Flow Rate:** 0.8 mL/min

**Temperature:** 30 °C

**Detection:** UV @ 254 nm (ambient)

**Sample:**

1. Pyridine	10. Triprolidine
2. Acetaminophen	11. Nortriptyline
3. Pindolol	12. Prednisolone
4. Quinidine	13. 3-Methyl-4-nitrobenzoic acid
5. Sulfathiazole	14. 2-Hydroxy-5-methylbenzaldehyde
6. Acebutolol	15. Diflunisal
7. Benzyl alcohol	16. Hexanophenone
8. Chlorpheniramine	
9. Phenol	

red-acids  
 blue-bases  
 white-neutrals


Comparative separations may not be representative of all applications. Columns are pH stable from 1.5-10 under isocratic conditions. Columns are pH stable 1.5-8.5 under gradient conditions.

# A New Era of Technical Support Services, Let Us Do the Work for You

PhenoLogix, our in-house application support lab, saves you time and money by screening multiple scout columns and solvent strategies for new purification methods or revalidating your current methods. We work together to make you successful by minimizing your process purification development time and optimizing your purification method.

## 1 Column Screening

- Normal Phase
- Reversed Phase
- Polar Organic
- SFC
- Chiral



## 2 Method Optimization Services

- Fast Turnaround
- Easy Method Transfer
- Continued Support

## 3 Preparative and Process Scale-Up

- Media Screening
- Small Scale Purification
- DAC Packing Assistance

Get  
started  
today

and let us exceed your expectations.

[www.phenomenex.com/phenologix](http://www.phenomenex.com/phenologix)



# Extend the Lifetime of Your Axia™ Preparative Column

## Use the SecurityGuard™ PREP column protection system!


- Extend preparative column lifetime by as much as 5x
- Protect column from samples that precipitate out of solution
- Protect column from contaminants
- Stable and leak-free up to 60 mL/min

The SecurityGuard PREP system was designed to effectively (and inexpensively) protect your valuable Prep columns from the damaging effects of mobile phase and sample chemical contaminants and particulates, without altering your chromatographic results.



## Forced Degradation Lifetime Study


### Axia packed column with SecurityGuard PREP cartridge after initial injection



#### Conditions

**Column:** Luna® 10 µm C18(2) Axia Packed  
**Dimension:** 50 x 21.2 mm  
**Part No.:** 00B-4253-P0-AX  
**Mobile Phase:** A: 0.1 % TFA in Water  
B: 0.1 % TFA in Water/Acetonitrile (25:75)  
**Gradient:** Linear 93:7 (A/B) to 100 % B over 5 minutes  
**Injection Volume:** 420 µL  
**Flow Rate:** 60 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 270 nm  
**Sample:** 1. Nadolol  
2. Metoprolol  
3. Propranolol


### Axia packed column with SecurityGuard PREP cartridge after multiple injections



Time to change  
the PREP cartridge



### Axia packed column after removing SecurityGuard column protection system



Original column  
performance maintained by  
using SecurityGuard PREP



# guarantee

If Axia™ packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase and dimensions, return the column with comparative data within 45 days for a **FULL REFUND**. Only applies to 21.2 mm ID columns.

## Kinetex® Analytical Columns

5 µm Columns (mm)		SecurityGuard ULTRA Cartridges*					SecurityGuard ULTRA Cartridges*
50 x 2.1		3/pk	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
<b>XB-C18</b>	00B-4605-AN	AJO-8782	00B-4605-E0	00D-4605-E0	00F-4605-E0	00G-4605-E0	AJO-8768
<b>C18</b>	00B-4601-AN	AJO-8782	00B-4601-E0	00D-4601-E0	00F-4601-E0	00G-4601-E0	AJO-8768
<b>PFP</b>	00B-4602-AN	AJO-8787	00B-4602-E0	00D-4602-E0	00F-4602-E0	00G-4602-E0	AJO-8773
<b>Phenyl-Hexyl</b>	00B-4603-AN	AJO-8788	00B-4603-E0	00D-4603-E0	00F-4603-E0	00G-4603-E0	AJO-8774
		for 2.1 mm ID					for 4.6 mm ID

2.6 µm Analytical Columns (mm)						SecurityGuard ULTRA Cartridges*
30 x 4.6 50 x 4.6 75 x 4.6 100 x 4.6 150 x 4.6						3/pk
<b>XB-C18</b>	—	00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	AJO-8768
<b>C18</b>	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJO-8768
<b>PFP</b>	00A-4477-E0	00B-4477-E0	00C-4477-E0	00D-4477-E0	00F-4477-E0	AJO-8773
<b>Phenyl-Hexyl</b>	—	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0	AJO-8774
						for 4.6 mm ID

2.6 µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges*
30 x 3.0 50 x 3.0 75 x 3.0 100 x 3.0 150 x 3.0						3/pk
<b>XB-C18</b>	00A-4496-Y0	00B-4496-Y0	00C-4496-Y0	00D-4496-Y0	00F-4496-Y0	AJO-8775
<b>C18</b>	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJO-8775
<b>PFP</b>	00A-4477-Y0	00B-4477-Y0	00C-4477-Y0	00D-4477-Y0	00F-4477-Y0	AJO-8780
<b>Phenyl-Hexyl</b>	—	—	—	00D-4495-Y0	00F-4495-Y0	AJO-8781
						for 3.0 mm ID

2.6 µm Minibore Columns (mm)						SecurityGuard ULTRA Cartridges*
30 x 2.1 50 x 2.1 75 x 2.1 100 x 2.1 150 x 2.1						3/pk
<b>XB-C18</b>	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
<b>C18</b>	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782
<b>PFP</b>	00A-4477-AN	00B-4477-AN	00C-4477-AN	00D-4477-AN	00F-4477-AN	AJO-8787
<b>Phenyl-Hexyl</b>	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJO-8788
						for 2.1 mm ID

1.7 µm MidBore Columns (mm)					SecurityGuard™ ULTRA Cartridges*
30 x 3.0 50 x 3.0 100 x 3.0					3/pk
<b>XB-C18</b>	00A-4498-Y0	00B-4498-Y0	00D-4498-Y0		AJO-8775
<b>C18</b>	—	00B-4475-Y0	00D-4475-Y0		AJO-8775
<b>PFP</b>	—	00B-4476-Y0	00D-4476-Y0		AJO-8780
					for 3.0 mm ID

1.7 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges*
30 x 2.1 50 x 2.1 100 x 2.1 150 x 2.1					3/pk
<b>XB-C18</b>	00A-4498-AN	00B-4498-AN	00D-4498-AN	00F-4498-AN	AJO-8782
<b>C18</b>	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJO-8782
<b>PFP</b>	00A-4476-AN	00B-4476-AN	00D-4476-AN	00F-4476-AN	AJO-8787
<b>Phenyl-Hexyl</b>	—	00B-4500-AN	00D-4500-AN	00F-4500-AN	AJO-8788
					for 2.1 mm ID



1.3 µm Columns (mm)	
50 x 2.1	
<b>C18</b>	00B-4515-AN

## Axia Packed Kinetex Preparative Columns

5 µm Axia™ Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges**
50 x 21.2 100 x 21.2 150 x 21.2 250 x 21.2					15 x 21.2 /ea
<b>XB-C18</b>	00B-4605-P0-AX	00D-4605-P0-AX	00F-4605-P0-AX	00G-4605-P0-AX	AJO-9145
<b>C18</b>	00B-4601-P0-AX	00D-4601-P0-AX	00F-4601-P0-AX	00G-4601-P0-AX	AJO-9145
<b>PFP</b>	00B-4602-P0-AX	00D-4602-P0-AX	00F-4602-P0-AX	00G-4602-P0-AX	AJO-9146
<b>Phenyl-Hexyl</b>	00B-4603-P0-AX	00D-4603-P0-AX	00F-4603-P0-AX	00G-4603-P0-AX	AJO-9147
					for 21.2 mm ID



## SecurityGuard™ Column Protection

Protect your Kinetex and Axia Packed preparative columns from the damaging effects of mobile phase and sample chemical contaminants with SecurityGuard ULTRA and SecurityGuard PREP

- Dramatically extends column lifetime and performance
- Virtually no change in chromatography
- Simple to use

\*SecurityGuard ULTRA cartridges require holder, Part No. AJO-9000

\*\*SecurityGuard PREP cartridges require holder, Part No. AJO-8223

SecurityGuard ULTRA			
Part No.	Description	Unit	Price
AJO-9000	SecurityGuard ULTRA Holder for UHPLC/ HPLC Columns 2.1 to 4.6mm ID	ea	



SecurityGuard PREP			
Part No.	Description	Unit	Price
AJO-8223	SecurityGuard PREP HPLC Guard Cartridge Holder Kit, 21.1mm ID, includes column coupler	ea	







KINETEX

# Preparative Chromatography Meets UHPLC

## Australia

t: 02-9428-6444  
f: 02-9428-6445  
auinfo@phenomenex.com

## Austria

t: 01-319-1301  
f: 01-319-1300  
anfrage@phenomenex.com

## Belgium

t: 02 503 4015 (French)  
t: 02 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: (800) 543-3681  
f: (310) 328-7768  
info@phenomenex.com

## Denmark

t: 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: 09 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: 01 30 09 21 10  
f: 01 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: 06021-58830-0  
f: 06021-58830-11  
anfrage@phenomenex.com

## India

t: 040-3012 2400  
f: 040-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: 01 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: 051 6327511  
f: 051 6327555  
italiainfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 001-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: 030-2418700  
f: 030-2383749  
nlinfo@phenomenex.com

## New Zealand

t: 09-4780951  
f: 09-4780952  
nzinfo@phenomenex.com

## Norway

t: 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: (800) 541-HPLC  
f: (310) 328-7768  
info@phenomenex.com

## Sweden

t: 08 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: 01625-501367  
f: 01625-501796  
ukinfo@phenomenex.com

## United States

t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com

## All other countries: Corporate Office USA



t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com



## www.phenomenex.com

Phenomenex products are available worldwide.  
For the distributor in your country, contact  
Phenomenex USA, International Department at  
international@phenomenex.com

## Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at  
[www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

## Trademarks

Phenomenex, Gemini, Kinetex, and Luna are registered trademarks, Axia, MidBore, Synergi,  
and SecurityGuard are trademarks of Phenomenex. Waters is a registered trademark, XBridge  
and OBD are trademarks of Waters Corporation.

## Disclaimer

Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be  
representative of all applications.

Axia is patented by Phenomenex. U.S. Patent No. 7,674,383

© 2013 Phenomenex. All rights reserved.

BR11940313\_W