

ULTIMATE HPLC/UHPLC COMBO

Kinetex + Luna Omega

**Gain Incredible Performance,
PEAK Reproducibility,
Versatile Selectivities and
Upgrade Throughput**



NEW



Luna Omega C18 (1.6 μ m, 3 μ m, 5 μ m)



**Kinetex and Luna Omega products with
the Next Level of Reproducibility**

 **phenomenex**
...breaking with traditionSM




Luna Omega

Cutting Edge Fully Porous Silica Particle

Luna® is one of the most recognized HPLC brands on the market, delivering high efficiency, ruggedness, reproducibility, and dependability for a wide range of analyses. The new Luna Omega builds upon this legacy with an innovative yet rugged silica particle architecture, designed and manufactured by Phenomenex based on more than 20 years of applied knowledge, invention, and customer experience.




Novel Design and Manufacturing Process

Within the novel manufacturing process of Luna Omega silica, we implement a proprietary processing technique to gain greater particle inertness, a stronger particle morphology, and more consistent porosity.



Thermal Modified Pore Structure

Most importantly, through our proprietary process, we eliminate micropores, further improving column efficiency, inertness, and reproducibility.




Kinetex

The Chosen Core-Shell Brand



Kinetex® Core-Shell Technology delivers dramatic improvements in efficiency over conventional fully porous media which can be leveraged to increase resolution, greatly improve productivity, reduce solvent consumption, and decrease costs. Whether you are running HPLC or UHPLC methods, the Kinetex core-shell family can deliver shockingly improved performance over the current column you are using.

Finely Tuned, Unique Core-Shell Manufacturing Process

Phenomenex designs, manufactures, and sells its very own silica and organo-silica core-shell particles. Using silica sol-gel processing techniques that incorporate nano-structuring technology, a durable, homogeneous porous shell is grown on a solid silica core to create a core-shell particle. The combination of a consistent, solid high density core along with proprietary column packing technologies ensures optimum bed structure and high column performance.




You can unleash the power of Kinetex Core-Shell Technology columns to achieve faster and better results just not possible on conventional fully porous materials. No matter what type of system you have or type of analysis you are performing, there is a Kinetex solution for you.



Gain Incredible Performance With Kinetex and Luna Omega


The undeniably high efficiency levels found in each Luna® Omega and Kinetex® column provide you with the potential of huge gains in method performance. While traditional silica and hybrid fully porous particles may claim high performance, when compared to Luna Omega or Kinetex, they fall short and prevent HPLC/UHPLC scientists from reaching their goals.

Efficiency Levels (plates/m)




UHPLC Performance – Cyclosporine Impurity Profile


Kinetex 2.6µm Polar C18



Luna Omega 1.6µm Polar C18



Conventional 1.7µm C18



Conditions for all columns same except where noted:

Columns: Kinetex 2.6 µm Polar C18
Luna Omega 1.6 µm Polar C18
Conventional Fully Porous 1.7 µm C18

Dimensions: 50 x 2.1 mm

Mobile Phase: Acetonitrile/Tert-butyl methyl ether/Water/Phosphoric acid (430:50:520:1)

Flow Rate: 0.30 mL/min


Temperature: 80 °C

Detection: UV @ 210 nm

Sample: Cyclosporine


Comparative separations may not be representative of all applications.

Efficiency Levels (plates/m)




HPLC Performance – Water Soluble Vitamins


Kinetex 2.6µm Polar C18



Luna Omega 3µm Polar C18



Core-Shell 2.7µm C18



Conditions for all columns same except where noted:

Columns: Kinetex 2.6 µm Polar C18
Luna Omega 3 µm Polar C18
Core-Shell 2.7 µm C18

Dimensions: 100 x 4.6 mm

Mobile Phase: A: 20 mM Potassium Phosphate
B: Methanol

Gradient:

Time (min)	% B
0	0
1	0
10	60

Gradient for conventional core-shell C18 column starts at 2 % B

Flow Rate: 1.2 mL/min

Temperature: Ambient


Detection: UV @ 210 nm

Sample: 1. Thiamine
2. Nicotinamide
3. Pyridoxal
4. Pyridoxine
5. Pantothenic Acid
6. Riboflavin

Gain Retention and Resolution With Enhanced Polar Phases

The polar modified functionalities of the Polar C18 and PS C18 stationary phases provide greater polar compound retention and resulting improvement in resolution values. Additionally, the advanced proprietary bonding technology used for both the Polar C18 and PS C18 ensures 100% aqueous stability as well as balanced retention for non-polar compounds.

Nicotine and Metabolites




Conditions for all columns:

Columns: Luna Omega 1.6 µm Polar C18
Core-Shell 1.6 µm C18
Dimension: 50 x 2.1 mm
Mobile Phase: A: 10 mM Ammonium Formate with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid
Gradient:

Time (min)	% B
0	2
3	90
3.1	2

Flow Rate: 0.4 mL/min
Temperature: 25 °C
Detection: MS/MS (SCIEX® API 4000™)
Sample: 1. Normicotine
2. 3-Hydroxycotinine
3. Nicotine
4. Cotinine
5. Anabasine

MMA and Succinic Acid



Conditions for all columns:

Columns: Luna Omega 1.6 µm PS C18
Conventional Fully Porous 2 µm C18
Dimension: 50 x 2.1 mm
Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid
Gradient:


Time (min)	% B
0	0
5	50
5.1	0
7	0

Flow Rate: 0.5 mL/min
Temperature: 22 °C
Detection: MS/MS (SCIEX API 4000™)
Sample: 1. Succinic acid
2. MMA

Optimal Selectivity for Complex Polar/Non-Polar Mixtures

Use the combined core-shell performance and dual polar/non-polar selectivity gains of the Kinetex 2.6 µm Polar C18 to easily expand the elution window and resolution values for your multi-compound methodologies that contain polar and non-polar compounds.

Columns: Kinetex 2.6 µm Polar C18
Dimension: 50 x 4.6 mm
Mobile Phase: A: Water
B: 0.1 % Formic Acid in Methanol
Gradient: 5-100 % B in 5 min, hold 1 min
Flow Rate: 0.7 mL/min
Temperature: Ambient
Detection: MS/MS (SCIEX API 4000™)
Sample: 206 Pesticides.
Find the full compound list online at www.phenomenex.com/Application/Detail/24017



Increase Retention with Alkaline Mobile Phase

The chemical stability and wide pH range (1-12) of Kinetex EVO C18 organo-silica particles allow you to use high pH mobile phases that are above a basic compound's pK_a , thus enhancing retention and potentially improving resolution and sensitivity values.

Conditions for all columns same except where noted:


Column: Kinetex 1.7 µm EVO C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4726-AN
Mobile Phase: A: 20 mM Potassium phosphate (pH 2.5)
B: Methanol
Mobile Phase: A: 20 mM Potassium phosphate (pH 10)
B: Methanol
Gradient:

Time (min)	% B
0	5
8	90


Flow Rate: 0.5 mL/min
Temperature: 30 °C
Detection: UV @ 210 nm
Sample: 1. Fluvoxamine

Comparative separations may not be representative of all applications.

Kinetex 1.7 µm EVO C18








Kinetex 1.7 µm EVO C18



Combined Method Development Flexibility

Expand your method development options by screening a combination of core-shell and fully porous particles with different selectivities which can increase retention, improve resolution, or even flip peak elution orders. The combination of available Luna® Omega stationary phases represent an outstanding tool set to start with for the separation of acids, bases, neutrals, or mixtures. At the same time, the Kinetex® core-shell EVO C18, Biphenyl, and F5 are outstanding orthogonal phases that can easily provide you with highly useful changes in elution windows and retention times.

Cold Medicine Ingredients



Conditions for all columns:

Columns: Kinetex 2.6µm Biphenyl
Luna Omega 3µm Polar C18
Kinetex 2.6µm EVO C18
Kinetex 2.6µm F5

Dimension: 100 x 4.6 mm

Mobile Phase: A: 20 mM Potassium Phosphate pH 3.5
B: Acetonitrile

Gradient:

Time (min)	% B
0	2
7	30
12	80

Flow Rate: 1 mL/min
Temperature: Ambient
Detection: UV @ 254 nm
Sample: 1. Maleic Acid
2. Phenylephrine
3. Acetaminophen
4. Doxylamine
5. Pyrilamine
6. Chlorpheniramine

7. Dextromethorphan
8. Bromopheniramine
9. 4-Nitrophenol
10. Acetylsalicylic Acid
11. Diphenhydramine

Improve Peak Shape for Bases

The unique surface of the Luna Omega PS C18 contains a positive charged ligand which aids in the retention of acidic compounds through ionic interactions, while also greatly improving the peak shape of basic compounds through ionic repulsion.

Conditions for all columns:

Column: Luna Omega 3µm PS C18
Conventional Fully Porous 3µm C18

Dimensions: 100 x 4.6 mm

Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid

Gradient: 10 - 95 % B in 10 min


Flow Rate: 1.5 mL/min

Temperature: Ambient


Detection: UV @ 254 nm

Sample: 1. Amitriptyline
2. Nortriptyline
3. Imipramine
4. Clomipramine

Luna Omega 3µm PS C18



Conventional Fully Porous 3µm C18




Comparative separations may not be representative of all applications.

Generating the Next Level

Of Reliability Through Advanced Process Optimization

Over the past three years, our scientists and engineers with the help of customers and Danaher colleagues, have optimized our processes to provide products that deliver very high levels of performance and newly achievable levels of reliability and reproducibility. This new advanced series of products and process optimization is called PEAK.



Luna® Omega C18

Reproducible and Scalable


By setting a new standard for reliability, the Luna Omega C18 spans UHPLC and HPLC with a scalable range of high-performance particle sizes that will ensure that your developed methods are easily transferred. From single compound identification to complex impurity profiles, the Luna Omega C18 will serve as a pillar for your lab to count on day in and day out.

Batch-to-Batch Reproducibility Study

In this example, we compared three batches of Luna Omega C18 using all three different particle sizes on a complex QC Pharmaceutical representative sample.

Conditions for all columns:
Mobile Phase: A: Water with 0.1 % Formic Acid
B: Acetonitrile with 0.1 % Formic Acid
Temperature: 30 °C
Detection: UV @ 254 nm
Injection Volume: 5 μL
Sample: 5 mg/mL of Chlorhexidine and Related Substances


Luna Omega 5 μm C18 Impurity Profile 3 Batch Comparison



Column: Luna Omega 5 μm C18
Dimension: 250 x 4.6 mm
Part No.: 00G-4785-E0
Gradient:

Time (min)	% B
0	2
2.5	2
52.5	35
55	35
57.5	2
62.5	2


Luna Omega 3 μm C18 Impurity Profile 3 Batch Comparison



Column: Luna Omega 3 μm C18
Dimension: 150 x 4.6 mm
Part No.: 00F-4784-E0
Gradient:

Time (min)	% B
0	2
1.5	2
31.5	35
34.5	35
36	2
42	2

Luna Omega 1.6 μm C18 Impurity Profile 3 Batch Comparison



Column: Luna Omega 1.6 μm C18
Dimension: 50 x 2.1 mm
Part No.: 00B-4742-AN
Gradient:

Time (min)	% B
0	2
0.5	2
10.5	35
11.5	35
12	2
14	2

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Axia, MidBore, SecurityGuard, Aeris, BE-HAPPY, and Synergi are trademarks, Luna, Gemini, Lux, Clarity, and Kinetex are registered trademarks of Phenomenex. SCIEX is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX™ is being used under license.

Disclaimer

Comparative separations may not be representative of all applications.

Axia column and packing technology is patented by Phenomenex.
U.S. Patent No. 7, 674, 383

Gemini and Kinetex EVO are patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures

© 2020 Phenomenex, Inc. All rights reserved.

