



Strata™  
XL  
from Phenomenex

Quick and Simple Cleanup for

# Medicinal Chemistry

strata™ XL

It's BIG.

 phenomenex®  
...breaking with tradition™



# Guaranteed Cleanup of Post Synthesis Reactions

**Strata™ -XL Giga™ Tubes are designed to make Medicinal Chemistry easy!**



*“With Strata Giga tubes I’ve been able to purify a compound containing a quaternary amine successfully. I would recommend Strata... for sure!”*

- Maxim Maheux, TransBIOTech

- **Quick and simple cleanup procedures**
- **Tubes available for a wide range of sample sizes; 1 mL up to 150 mL**
- **Pre-packed tubes eliminate the excess maintenance involved with cleaning and handling glass packed flash columns**
- **Large particle sizes for fast flow**
- **Robust polymeric material for pH and chemical resistance, including DMSO**

guarantee

If Strata SPE products do not perform as well or better than your current SPE product of similar phase, mass and size, send in your comparative data within 45 days and keep the Strata SPE product for FREE.

© 2010 Phenomenex, Inc. All rights reserved.

# Customize Your Cleanup in 2 Easy Steps

## step 1: Select a Sorbent

Choose from 5 different sorbent chemistries.

Cleanup Strategy	Sorbent Recommendation
Desalt neutral compounds	<b>Strata™-XL</b>
Cleanup basic compounds	<b>Strata-XL-C</b>
Cleanup strongly basic compounds	<b>Strata-XL-CW</b>
Cleanup acidic compounds	<b>Strata-XL-A</b>
Cleanup strongly acidic compounds	<b>Strata-XL-AW</b>

## step 2: Choose a Tube

Use the chart below to determine proper sorbent mass based on your sample size.


Strata-XL					
<b>Neutral Compound Loading Capacity</b>	185 mg	460 mg	935 mg	1.85 g	4.65 g
<b>Recommended Sorbent Mass</b>	2 g	5 g	10 g	20 g	50 g

Strata-XL-C, Strata-XL-CW, Strata-XL-A and Strata-XL-AW			
<b>Ionic Compound Loading Capacity</b>	200 mg	500 mg	1 g
<b>Recommended Sorbent Mass</b>	2 g	5 g	10 g

# Cleanup Basic Compounds

Achieve high recoveries of basic compounds with Strata™-XL-C and Strata™-XL-CW.

## Purify 1°, 2°, and 3° Amines with Strata-XL-C



Sulfonic Acid ligand offers selectivity for the retention of weaker bases

### Strata-XL-C for weakly basic compounds

#### Target Analytes

Basic products with a  $pK_a < 10.5$  (1°, 2°, and 3° amines)

#### Sorbent Chemistry

Benzenesulfonic acid conjugated to a large particle polymer which has been optimized for gravity flow

#### Particle Size


100  $\mu\text{m}$

#### pH Stability

1 - 14

#### Ionic Capacity

1 meq/g




# Cleanup Basic Compounds

Load reaction mixture



B Bases  
A Acids  
N Neutrals

Purify 1°, 2°, 3°, and 4° Amines with Strata<sup>™</sup>-XL-CW



Carboxylic Acid ligand offers selectivity for the retention of strong bases

## Strata-XL-CW for strongly basic compounds

### Target Analytes

Basic products with a  $pK_a > 8$   
(1°, 2°, 3°, and 4° amines)

### Sorbent Chemistry

Carboxylic acid group conjugated to a large particle polymer which has been optimized for gravity flow

### Particle Size

100  $\mu\text{m}$

### pH Stability

1 - 14

### Ionic Capacity

0.76 meq/g

Wash away acids and neutrals



Elute purified basic compounds




# Cleanup Neutral Compounds

## Ion-Exchange Method

Quickly and easily remove ionic species from Neutral compounds by performing an ion-exchange cleanup.

### Remove Basic AND Acidic Compounds with Strata®-ABW



Strong cation - exchange and weak anion - exchange groups allow neutral compounds to flow through the sorbent while ionic compounds are trapped in the sorbent


#### Strata ABW for neutral compounds

**Target Analytes** Ionically charged compounds, neutral compounds flow through sorbent

**Sorbent Chemistry** Sulfonic Acid and Amino groups for retention of basic and acidic compounds, respectively

**Particle Size** 55 µm

**pH Stability** 2.5 - 8.5



# Cleanup Neutral Compounds

## Reversed Phase/Desalting Method

Load reaction mixture



B	Bases
A	Acids
N	Neutrals


Wash away acids and bases



Elute purified neutral compounds



### Desalt Neutral Compounds with Strata<sup>TM</sup>-XL



Modified pyrrolidone ligand for retention of neutral compounds


#### Strata-XL for neutral compounds

<b>Target Analytes</b>	Neutral, hydrophobic compounds
<b>Sorbent Chemistry</b>	N-vinylpyrrolidone conjugated to a large particle polymer which has been optimized for gravity flow
<b>Particle Size</b>	100 $\mu$ m
<b>pH Stability</b>	1 - 14

# Cleanup Acidic Compounds

Achieve superior cleanup and high recoveries of acidic compounds with Strata™-XL-A and Strata™-XL-AW.

## Purify Weakly Acidic Compounds with Strata™-XL-A



Di-methylbutyl quaternary amine ligand offers selectivity for the retention of weak acids

### Strata-XL-A for weakly acidic compounds

**Target Analytes** Acidic compounds with a  $pK_a > 3$  (carboxylic acids)

**Sorbent Chemistry** Di-methylbutyl quaternary amine ligand conjugated to a large particle polymer which has been optimized for gravity flow

**Particle Size** 100  $\mu\text{m}$


**pH Stability** 1 - 14

**Ionic Capacity** 0.30 meq/g

1 Load reaction mixture

Strata™  
XL-A  
from Phenomenex

## Purify Strongly Acidic Compounds with Strata™-XL-AW



Di-amino ligand offers selectivity for the retention of strong acids

### Strata-XL-AW for strongly acidic compounds

**Target Analytes** Acidic compounds with a  $pK_a < 5$  (carboxylic and sulfonic acids)

**Sorbent Chemistry** Di-amino ligand conjugated to a large particle polymer which has been optimized for gravity flow

**Particle Size** 100  $\mu\text{m}$

**pH Stability** 1 - 14

**Ionic Capacity** 0.60 meq/g

2 Wash away bases and neutrals



3 Elute purified acidic compounds



B Bases  
A Acids  
N Neutrals



# Traditional Cleanup Solutions

Silica-based sorbents have been used to fractionate and cleanup synthetic organic compounds for many years. Strata silica-based products are made from the highest quality silica and each batch is quality tested for purity, recovery, and flow to ensure that your results are reproducible.

Compound to be Fractionated	Sorbent Recommendation
<b>Reversed Phase Solutions</b>	
Very hydrophobic compounds	<b>Strata® C8</b>
Hydrophobic compounds	<b>Strata C18-E</b>
Hydrophobic, aromatic compounds	<b>Strata Phenyl</b>
<b>Normal Phase Solutions</b>	
Polar compounds including those that contain a hydroxyl group	<b>Strata NH<sub>2</sub></b>
Polar compounds including those that contain an amino group	<b>Strata Silica</b>
<b>Additional Cleanup Solutions</b>	
Removal of aqueous residues	<b>Strata Sodium Sulfate</b>

<b>Strata ABW, Strata C8, Strata C18-E, Strata Phenyl, Strata NH<sub>2</sub>, and Strata Silica</b>					
<b>Maximum Loading Capacity of Neutral Compounds</b>	100 mg	250 mg	500 mg	1 g	3.5 g
<b>Sorbent Mass</b>	2 g	5 g	10 g	20 g	70 g

# Ordering Information

## Cleanup Solutions

### Ion-Exchange - Strata<sup>™</sup>-XL-C

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S044-KDG	20/box	
2 g	20 mL	8B-S044-KEG	20/box	
5 g	20 mL	8B-S044-LEG	20/box	
5 g	60 mL	8B-S044-LFF	16/box	
10 g	60 mL	8B-S044-MFF	16/box	

### Ion-Exchange - Strata-XL-AW

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S051-KDG	20/box	
2 g	20 mL	8B-S051-KEG	20/box	
5 g	20 mL	8B-S051-LEG	20/box	
5 g	60 mL	8B-S051-LFF	16/box	
10 g	60 mL	8B-S051-MFF	16/box	

### Ion-Exchange - Strata-XL-CW

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S052-KDG	20/box	
2 g	20 mL	8B-S052-KEG	20/box	
5 g	20 mL	8B-S052-LEG	20/box	
5 g	60 mL	8B-S052-LFF	16/box	
10 g	60 mL	8B-S052-MFF	16/box	

### Reversed Phase - Strata-XL

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S043-KDG	20/box	
2 g	20 mL	8B-S043-KEG	20/box	
5 g	20 mL	8B-S043-LEG	20/box	
5 g	60 mL	8B-S043-LFF	16/box	
10 g	60 mL	8B-S043-MFF	16/box	
20 g	60 mL	8B-S043-VFF	16/box	
50 g	150 mL	8B-S043-YSN	8/box	

### Ion-Exchange - Strata-XL-A

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	Available in August	20/box	*
2 g	20 mL	Available in August	20/box	*
5 g	20 mL	Available in August	20/box	*
5 g	60 mL	Available in August	16/box	*
10 g	60 mL	Available in August	16/box	*

\*Contact Phenomenex for more information.



# Ordering Information

## Traditional Cleanup Solutions

### Mixed Mode - Strata® ABW

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S030-KDG	20/box	
5 g	20 mL	8B-S030-LEG	20/box	
10 g	60 mL	8B-S030-MFF	16/box	
20 g	60 mL	8B-S030-VFF	16/box	
70 g	150 mL	8B-S030-ZSN	8/box	

### Reversed Phase - Strata C8

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S005-KDG	20/box	
5 g	20 mL	8B-S005-LEG	20/box	
10 g	60 mL	8B-S005-MFF	16/box	

### Reversed Phase - Strata C18-E

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
500 mg	12 mL	8B-S001-HDG	20/box	
1 g	20 mL	8B-S001-JEG	20/box	
2 g	12 mL	8B-S001-KDG	20/box	
5 g	20 mL	8B-S001-LEG	20/box	
10 g	60 mL	8B-S001-MFF	16/box	
20 g	60 mL	8B-S001-VFF	16/box	
20 g	150 mL	8B-S001-VSN	8/box	
50 g	150 mL	8B-S001-YSN	8/box	
70 g	150 mL	8B-S001-ZSN	8/box	



If Strata SPE products do not perform as well or better than your current SPE product of similar phase, mass and size, send in your comparative data within 45 days and keep the Strata SPE product for FREE.

Terms and Conditions  
Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

Trademarks  
Strata is a registered trademark of Phenomenex in the U.S. and other countries. Strata-X and Giga are trademarks of Phenomenex, Inc.

Strata-X is patented by Phenomenex

© 2010 Phenomenex, Inc. All rights reserved.

### Reversed Phase - Strata Phenyl

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
2 g	12 mL	8B-S006-KDG	20/box	

### Normal Phase - Strata NH<sub>2</sub>

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
500 mg	12 mL	8B-S009-HDG	20/box	
2 g	12 mL	8B-S009-KDG	20/box	
5 g	20 mL	8B-S009-LEG	20/box	
10 g	60 mL	8B-S009-MFF	16/box	
20 g	60 mL	8B-S009-VFF	16/box	
70 g	150 mL	8B-S009-ZSN	8/box	

### Normal Phase - Strata Silica

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
500 mg	12 mL	8B-S012-HDG	20/box	
1 g	12 mL	8B-S012-JDG	20/box	
2 g	12 mL	8B-S012-KDG	20/box	
5 g	20 mL	8B-S012-LEG	20/box	
10 g	60 mL	8B-S012-MFF	16/box	
20 g	60 mL	8B-S012-VFF	16/box	
50 g	150 mL	8B-S012-YSN	8/box	
70 g	150 mL	8B-S012-ZSN	8/box	

### Additional Cleanup - Strata Sodium Sulfate

Sorbent Mass	Tube Size	Part No.	Unit/Box	Price
1 g	6 mL	8B-S124-JCH	30/box	
1 g	20 mL	8B-S124-JEG	20/box	
5 g	20 mL	8B-S124-LEG	20/box	

\*Contact Phenomenex for more information.

**Australia**

t: 02-9428-6444  
f: 02-9428-6445  
auinfo@phenomenex.com

**Austria**

t: 01-319-1301  
f: 01-319-1300  
anfrage@phenomenex.com

**Belgium**

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

**Canada**

t: (800) 543-3681  
f: (310) 328-7768  
info@phenomenex.com

**Denmark**

t: 4824 8048  
f: 4810 6265  
nordicinfo@phenomenex.com

**Finland**

t: (09)4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

**France**

t: 01 30 09 21 10  
f: 01 30 09 21 11  
franceinfo@phenomenex.com

**Germany**

t: 06021-58830-0  
f: 06021-58830-11  
anfrage@phenomenex.com

**Ireland**

t: 01 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

**Italy**

t: 051 6327511  
f: 051 6327555  
italiainfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

**Mexico**

t: (55) 5018 3791  
f: (310) 328-7768  
tecnicomx@phenomenex.com

**Netherlands**

t: 030-2418700  
f: 030-2383749  
nlinfo@phenomenex.com

**New Zealand**

t: 09-4780951  
f: 09-4780952  
nzinfo@phenomenex.com

**Norway**

t: 81 00 20 05  
f: +45 4810 6265  
nordicinfo@phenomenex.com

**Puerto Rico**

t: (800) 541-HPLC  
f: (310) 328-7768  
info@phenomenex.com

**United Kingdom**

t: 01625-501367  
f: 01625-501796  
ukinfo@phenomenex.com

**All other countries:  
Corporate Office USA** 

t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com