

## Dramatically Improve Existing 5 $\mu\text{m}$ and 3 $\mu\text{m}$ Fully Porous Methods with Kinetex<sup>®</sup> 5 $\mu\text{m}$ Core-shell Technology

Jeff Layne, Phenomenex, Inc.  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

*With virtually no development effort or investment required chromatographers running traditional methods can instantly improve resolution, sensitivity, and productivity by simply replacing the 5  $\mu\text{m}$  and 3  $\mu\text{m}$  fully porous columns currently used with a Kinetex 5  $\mu\text{m}$  core-shell column. This newly introduced core-shell technology is able to achieve efficiency values that are equal or greater than fully porous 3  $\mu\text{m}$  columns, but operates at a pressure that is consistent with typical fully porous 5  $\mu\text{m}$  columns.*

### Introduction


Since the term “high performance liquid chromatography” was coined by Csaba Horvath in 1970, the dominant trend in HPLC media manufacturing has been the drive to create chromatographic media capable of delivering increasing efficiency values. Up until recently, this increase in column efficiency was achieved primarily through decreasing the mean diameter of the particles. Thus, from the 10  $\mu\text{m}$  media that was prevalent pre-1970s, we have seen a steady migration to 5  $\mu\text{m}$  which give efficiency values in the range of 90,000 -100,000 plates per meter, 3  $\mu\text{m}$  particles giving efficiencies of about 150,000 plates per meter, and more recently sub-2  $\mu\text{m}$  UHPLC (ultra high performance liquid chromatography) particles capable of achieving efficiency values in the range of almost 300,000 plates per meter. Of course, the drawback of these small particles was the accompanying increase in backpressure. In fact, the column backpressures generated by UHPLC particles was so great that specially designed UHPLC instrumentation, capable of operating at pressures of 600-1000 bar was required. When coupled with sub-2  $\mu\text{m}$  packed columns, these UHPLC systems were able to deliver amazing performance, but the high cost of these systems, typically costing 2-3 times as much as conventional HPLC systems, placed them out of reach of many potential end-users.

The development of core-shell media initially served as a counterpoint to the miniaturization of chromatographic media. In contrast to fully-porous silica particles which, as the name implies, consist of a porous silica microsphere, core-shell particles consist of an impermeable inner core surrounded by a layer of porous silica (**Figure 1**). We have recently released a core-shell product, Kinetex 5  $\mu\text{m}$ , which generates the same backpressure as traditional fully porous 5  $\mu\text{m}$  columns, but is able to generate efficiency values in the range of 170,000 -180,000 plates per meter. To put this into perspective, typical efficiencies for fully-porous 5  $\mu\text{m}$  packed columns are in the range of 90,000 -100,000 plates per meter and average 3  $\mu\text{m}$  fully-porous columns typically provide 150,000

-170,000 plates per meter. Thus, the new Kinetex 5  $\mu\text{m}$  core-shell media is able to deliver efficiencies that are as good as or better than 3  $\mu\text{m}$  media at the pressure of a 5  $\mu\text{m}$  particle.

Of course, efficiency values by themselves are only useful as bragging rights for HPLC column manufacturers. To the typical chromatographer, the true value of core-shell particles is the improved resolution, sensitivity, and productivity that they can provide when compared to fully porous particles of similar size. In this study, we present data generated using columns packed with Kinetex 5  $\mu\text{m}$  core-shell particles and show how they provide improved chromatographic performance when compared to fully-porous 5  $\mu\text{m}$  and 3  $\mu\text{m}$  packed columns.

**Figure 1.**  
Structure of Kinetex particle.



### Material and Methods

Unless stated otherwise, all chemicals and standards were obtained from Sigma Chemical Co. (St. Louis, Missouri). Standards for the analysis of paroxetine and its impurities were obtained from the European Directorate for the Quality of Medicines and Healthcare (EDQM). Solvents were purchased from EMD (San Diego, California). The Kinetex® 5 µm XB-C18 and C18 columns were obtained from Phenomenex (Torrance, California). Non-Phenomenex brand columns were purchased from the original equipment manufactures as noted in the figures. The column dimensions and running conditions for the various assays are detailed in the figure captions. Data from **Figures 2 and 3** were obtained using an Agilent 1200SL HPLC system equipped with a micro flow cell. All other chromatograms were obtained using an Agilent 1100 HPLC system equipped with a standard flow cell. Data was collected using ChemStation software (Agilent, Santa Clara, California).


### Results and Discussion

#### Kinetex 5 µm core-shell media provides improved resolution and sensitivity compared to traditional 5 µm fully porous media

As shown in **Figure 2**, Kinetex 5 µm core-shell particles can provide efficiency values that are as good or greater than conventional 5 µm and 3 µm columns, while operating at the same pressure as fully porous 5 µm columns (**Figure 3**). For demanding chromatographic methods such as stability-indicating assays and impurity profiles, the increased efficiency of Kinetex 5 µm can provide substantial improvements in both resolution and sensitivity when compared with columns packed with fully porous 5 µm

**Figure 2.**


Comparison of the efficiency of fully porous (F-P) 5 µm and 3 µm media with Kinetex 5 µm core-shell media. Running conditions: 65/35 Acetonitrile/Water; 1 ml/min; ambient temperature; efficiency determined from naphthalene peak.



Comparative separations may not be representative of all applications.

**Figure 3.**

Comparison of the backpressure of fully porous 5 µm and 3 µm media with Kinetex core-shell 5 µm media. Running conditions: 65/35 Acetonitrile/Water; 1 ml/min; ambient temperature; all columns 150 x 4.6 mm



**Figure 4** shows a stability-indicating assay developed for the antidepressant drug paroxetine. When the degraded drug sample is analyzed using a fully porous 5 µm packed column (Agilent®ZORBAX® 5 µm XDB-C18 250 x 4.6 mm) four principle degradants are identified and labeled as a-d (**Figure 4a**). There is also a small, unidentified degradant peak that is partially co-eluting with the main paroxetine peak ( $R_s = 0.98$ ; indicated in the zoomed-in view). If we run this same method and sample using the Kinetex 5 µm XB-C18 column (250 x 4.6 mm; **Figure 4b**), we first note that the peak height of all of the components is dramatically increased. For example, the peak height of the main paroxetine is 120 mAU on the fully-porous 5 µm ZORBAX column and 165 mAU on the Kinetex 5 µm XB-C18 column, a ~37 % increase in sensitivity. In addition, the unidentified impurity peak is completely resolved on the Kinetex 5 µm XB-C18 column ( $R_s = 2.23$ ; a 127 % increase). At least some of that increase in efficiency may also be attributable to differences in selectivity, not just the improved efficiency of the Kinetex 5 µm XB-C18 column, but a significant portion of that improved resolution can also be attributable to the peak widths being much narrower on the Kinetex 5 µm XB-C18 column.

Thus, without making any significant changes to the method itself, the Kinetex 5 µm XB-C18 column was able to provide a significant increase in resolution (127 %) and sensitivity (37 %) when compared to the fully porous 5 µm column.


In **Figure 5**, we present another comparison of the performance of a Kinetex core-shell 5 µm column versus a fully porous 5 µm

# TN-1132

## APPLICATIONS


**Figure 4.**  
Stability-indicating assay for paroxetine on a fully porous 5µm column and Kinetex core-shell 5µm column.

a. Agilent® ZORBAX® 5µm XDB-C18



App ID 20766

b. Phenomenex Kinetex® 5µm XB-C18



App ID 20767

**Conditions for both columns**

**Column:** Kinetex 5µm C18  
ZORBAX 5µm XDB-C18

**Dimensions:** 250 x 4.6 mm

**Mobile Phase:** A: TFA/Water/THF (5/900/100)  
B: TFA/Acetonitrile/THF (5/900/100)

Gradient: Time (min)	% B
0	10 %
40	80 %

**Flow Rate:** 1.5 mL/min

**Temperature:** 40 °C

**Detection:** UV @ 295 nm

<b>Sample:</b> 1. Impurity A	4. Unidentified
2. Impurity B	5. Paroxetine
3. Impurity C	6. Impurity D


Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc (Palo Alto, CA, USA). Phenomenex is not affiliated with Agilent Technologies. Comparative separations may not be representative of all applications.

# TN-1132


## APPLICATIONS

**Figure 5.**  
Impurity profile for omeprazole using a fully porous 5  $\mu\text{m}$  column and Kinetex 5  $\mu\text{m}$  core-shell column.

a. Waters® Symmetry® 5  $\mu\text{m}$  C18 250 x 4.6mm; Backpressure = 190 Bar



b. Phenomenex Kinetex® 5  $\mu\text{m}$  C18 250 x 4.6mm; Backpressure = 188 Bar



**Conditions for both columns**

**Column:** Kinetex 5  $\mu\text{m}$  C18  
Symmetry 5  $\mu\text{m}$  XDB-C18  
**Dimensions:** 250 x 4.6 mm  
**Mobile Phase:** A: mm Potassium Phosphate, pH 7  
B: Acetonitrile  
**Gradient:**

Time (min)	% B
0	5 %
40	70 %

**Flow Rate:** 1.2 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 295 nm  
**Sample:** Omeprazole degradants

packed column, in this case Waters Symmetry® 5  $\mu\text{m}$  C18. This is an impurity profile for the drug omeprazole. Overall, the selectivity between the two columns is fairly similar. However, the sensitivity, in terms of peak height response, is much greater on the Kinetex core-shell 5  $\mu\text{m}$  C18 column. This is immediately apparent if you look at the peak height response of the three impurity peaks in the red boxes – the peak heights are almost double on the Kinetex 5  $\mu\text{m}$  column as compared to the fully porous Symmetry 5  $\mu\text{m}$  C18 column.

Waters and Symmetry are registered trademarks of Waters Corporation (Milford, MA, USA). Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

# TN-1132

## APPLICATIONS

### Kinetex® 5 µm provides equivalent or better performance compared to a 3 µm column, but at a fraction of the pressure


In the previous section we showed that, with its improved efficiency, the Kinetex 5 µm core-shell column was able to provide significant improvements in resolution and sensitivity when compared to fully porous 5 µm columns. However, Kinetex 5 µm core-shell media is also an excellent alternative to fully porous 3 µm columns as it is able to provide equivalent or better efficiency compared to a 3 µm column (Figure 2), but operates at a much lower backpressure. This is well suited for pressure-limited operating systems

and also gives the chromatographer the opportunity to either (a) run at higher flow rates to decrease run time and increase productivity or (b) use a longer column to gain further resolution without introducing excessive backpressure to the system.

Figure 6 shows a stability-indicating assay developed for a proprietary API and run on a fully porous 3 µm packed column (CERI L-Column 3 µm ODS 150 x 4.6 mm) and on the Kinetex core-shell 5 µm C18 column (150 x 4.6 mm). The resolution between the two primary degradants of interest, indicated by the asterisk, are virtually identical on both columns ( $R_s = 1.59$  on the 3 µm fully po-

**Figure 6.**  
Stability-indicating assay for a proprietary API using a fully-porous 3µm column and Kinetex 5 µm C18.

A. CERI L-Column 3 µm ODS 150 x 4.6 mm




#### Conditions for both columns

**Column:** Kinetex 5 µm C18  
L-Column 3 µm ODS  
**Dimensions:** 150 x 4.6 mm  
**Mobile Phase:** A: Water with 0.1% TFA  
B: Acetonitrile with 0.1% TFA  
**Gradient:**

Time (min)	% B
0	10 %
20	70 %

  
**Flow Rate:** 1.0 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 210 nm  
**Sample:** Proprietary API degradants

B. Phenomenex Kinetex 5 µm C18 150 x 4.6 mm



L-Column is a trademark of CERI (Chemicals Evaluation and Research Institute, Japan). Phenomenex is not affiliated with CERI. Comparative separations may not be representative of all applications.

# TN-1132

## APPLICATIONS


rous column and 1.57 on the Kinetex 5  $\mu\text{m}$  core-shell column). However, the fully porous column is operating at almost double the backpressure that is generated by the Kinetex 5  $\mu\text{m}$  core-shell column (223 Bar versus 120 Bar). In addition, several of the smaller impurities, indicated by the arrows, are actually better resolved on the Kinetex 5  $\mu\text{m}$  core-shell column.

Thus far, we have focused primarily on LC/UV applications run using longer (150 mm and 250 mm) column formats. However, the benefits of Kinetex 5  $\mu\text{m}$  also extend to shorter lengths that are

typically associated with LC/MS analysis. In **Figure 7** we compare the performance of Kinetex 5  $\mu\text{m}$  core-shell media packed in a 50 x 2.1 mm format against a fully porous 3.5  $\mu\text{m}$  column of equivalent dimension (Waters® Sunfire™ 3.5  $\mu\text{m}$  C18). As is apparent from the chromatograms, the Kinetex core-shell 5  $\mu\text{m}$  C18 column provides much better performance than the fully porous 3.5  $\mu\text{m}$  column, with narrower peak widths and improved peak shape for these strongly basic analytes. For this method, the critical isobaric pair that needs to be resolved are 4,8- and 7,8-dimethyl IQx, which are extracted from the rest of the components in

**Figure 7.**  
LC/MS/MS analysis of heterocyclic amines extracted from salmon tissue using a fully porous 3.5  $\mu\text{m}$  column and the Kinetex 5  $\mu\text{m}$  core-shell column.


A. Waters® Sunfire™ 3.5  $\mu\text{m}$  C18 50 x 2.1 mm



**Conditions for both columns**

<b>Column:</b>	Kinetex 5 $\mu\text{m}$ C18	
	Sunfire 3.5 $\mu\text{m}$ C18	
<b>Dimensions:</b>	50 x 2.1 mm	
<b>Mobile Phase:</b>	A: 30 mM Ammonium formate	
	B: Acetonitrile	
<b>Gradient:</b>	<b>Time (min)</b>	<b>% B</b>
	0	5 %
	5	60 %
	5.1	100 %
	6	100 %
	7	5 %
<b>Flow Rate:</b>	1.2 mL/min	
<b>Temperature:</b>	Ambient	
<b>Detection:</b>	MS/MS (API 4000)	
<b>Sample:</b>	1. DMIP	10. PhIP
	2. Qx	11. AaC
	3. IQ	12. MeAaC
	4. 8-MeIQ	13. MeIQ
	5. 4,8-DiMe-Iqx	
	6. 7,8-DiMe-Iqx	
	7. Trp-P-2	
	8. Harman	
	9. Norharman	

B. Phenomenex Kinetex® 5  $\mu\text{m}$  C18 50 x 2.1 mm



Waters is registered trademark and Sunfire is a trademark of Waters Corporation (Milford, MA, USA). Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

# TN-1132

## APPLICATIONS

the bottom ion chromatogram. This critical isobaric pair is much better separated on the Kinetex® 5µm core-shell media than on the fully porous 3.5µm media. Also note that the backpressure of the Kinetex 5µm core-shell column is about one-half of the fully porous 3.5µm column.

### Conclusion

- Kinetex 5µm core-shell media is able to achieve efficiency values that are equal or greater than fully porous 3µm columns, but operates at a pressure that is consistent with typical fully porous 5µm columns
- For methods using fully porous 5µm columns, Kinetex 5µm core-shell media has the potential to provide improved chromatography, in terms of improved resolution and sensitivity, without the need for extensive method development.

- When compared to fully porous 3µm columns, Kinetex 5µm has the capacity to deliver equivalent or greater performance at about one-half of the pressure, allowing the chromatographer to take advantage of increasing flow rates to shorten analysis times or, alternatively, to use longer column lengths to further improve resolution. In addition, because of the low backpressure generated by the Kinetex 5µm columns, there is also the potential for increased column lifetime when compared to column packed with smaller particles.

### Ordering Information

5µm Columns (mm)		SecurityGuard™ ULTRA Cartridges*				SecurityGuard ULTRA Cartridges*	
		50 x 2.1	3/pk	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6
<b>XB-C18</b>	00B-4605-AN	AJO-8782	00B-4605-E0	00D-4605-E0	00F-4605-E0	00G-4605-E0	AJO-8768
<b>C18</b>	00B-4601-AN	AJO-8782	00B-4601-E0	00D-4601-E0	00F-4601-E0	00G-4601-E0	AJO-8768
<b>PFP</b>	00B-4602-AN	AJO-8787	00B-4602-E0	00D-4602-E0	00F-4602-E0	00G-4602-E0	AJO-8773
<b>Phenyl-Hexyl</b>	00B-4603-AN	AJO-8788	00B-4603-E0	00D-4603-E0	00F-4603-E0	00G-4603-E0	AJO-8774

for 2.1 mm ID

for 4.6 mm ID

\* SecurityGuard ULTRA cartridges require holder, Part No. AJO-9000.

### 5µm Axia Packed Preparative Columns (mm)

	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2
<b>XB-C18</b>	00B-4605-P0-AX	00D-4605-P0-AX	00F-4605-P0-AX	00G-4605-P0-AX
<b>C18</b>	00B-4601-P0-AX	00D-4601-P0-AX	00F-4601-P0-AX	00G-4601-P0-AX
<b>PFP</b>	00B-4602-P0-AX	00D-4602-P0-AX	00F-4602-P0-AX	00G-4602-P0-AX
<b>Phenyl-Hexyl</b>	00B-4603-P0-AX	00D-4603-P0-AX	00F-4603-P0-AX	00G-4603-P0-AX



If Kinetex core-shell columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

# TN-1132 APPLICATIONS

## Australia

t: 02-9428-6444  
f: 02-9428-6445  
auinfo@phenomenex.com

## Austria

t: 01-319-1301  
f: 01-319-1300  
anfrage@phenomenex.com

## Belgium

t: 02 503 4015 (French)  
t: 02 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: (800) 543-3681  
f: (310) 328-7768  
info@phenomenex.com

## Denmark

t: 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: 09 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: 01 30 09 21 10  
f: 01 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: 06021-58830-0  
f: 06021-58830-11  
anfrage@phenomenex.com

## India

t: 040-3012 2400  
f: 040-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: 01 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: 051 6327511  
f: 051 6327555  
italiainfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 001-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: 030-2418700  
f: 030-2383749  
nlinfo@phenomenex.com

## New Zealand

t: 09-4780951  
f: 09-4780952  
nzinfo@phenomenex.com

## Norway

t: 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: (800) 541-HPLC  
f: (310) 328-7768  
info@phenomenex.com

## Sweden

t: 08 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: 01625-501367  
f: 01625-501796  
ukinfo@phenomenex.com

## United States

t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com

## All other countries: Corporate Office USA

t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com

## Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

## Trademarks

Kinetex is a registered trademark of Phenomenex. SecurityGuard is a trademark of Phenomenex. Waters and Symmetry are registered trademarks of Waters Corporation. Sunfire is a trademark of Waters Corporation. Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc. L-Column is a trademark of CERI.

## Disclaimer

Phenomenex is in no way affiliated with Waters Corporation, Agilent Technologies, or CERI. Dimensions and chromatographic conditions are the same for all columns unless otherwise noted. Comparative separations may not be representative of all applications.

© 2012 Phenomenex, Inc. All rights reserved.