

APPLICATIONS

Tandem Digestion of Monoclonal Antibodies Using Novel Cysteine Proteases

Ivan Lebedev and Brian Rivera
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Monoclonal antibodies (mAbs) are well-established therapeutics, with a variety of different analytical methods associated with purity analysis. An analytical technique common for mAb purity is reversed phase LC (RPLC). This is a primary technique utilized for mAbs because of its relatively short analysis time, with high resolution, and ability to separation of hydrophobic variants such as oxidation, glycoforms, and lysine variants. Another utility of reversed phase chromatography can be used when analyzing fragments while using a reducing agent like dithiothreitol (DTT) to reduce interchain disulfides, yielding heavy chain and light chain.

Although the analysis of heavy and light chains is useful, the use of site-specific cysteine proteases is another approach that allows for better characterization of the antibody. Two common approaches to fragment generation are using the site-specific proteases IdeS and IgdE, which cleave below and above the hinge of IgG1 antibodies, respectively. This process allows for the isolation or characterization of specific regions of the antibody. Although analytical characterization by RPLC at the subunit level will provide more insight into sample heterogeneity, often further chemical digestion can be performed to obtain more detail. However, other variants - namely, disulfide variants - cannot be identified if this chemical reduction is performed.

As such, Faid and colleagues demonstrated the capabilities of digesting mAbs using both IgdE and IdeS as a method to identify free sulfhydryls.¹ Indeed, identification of free sulfhydryls has been limited mainly to colorimetric techniques; this tandem digestion allows for this identification of free sulfhydryls along with other variants commonly observed in subunit analysis. **Figure 1** shows NIST mAb that has been digested with IgdE and IdeS. Resulting fragments are observed, namely the Fc/2 and Fab, the former which gives insight on oxidation, glycoforms, and lysine variants, and the latter which gives insight on the hypervariable region. **Figure 1** inset shows the putative free sulfhydryl variants, which are later eluting to the Fc/2 fragment. It is important to note that there are some earlier eluting impurities, which are thought to be artifacts of the digestion itself.

Figure 2 shows a similar profile for trastuzumab, with free sulfhydryls also present in the sample. Later eluting variants may be deamidated variants, which have been observed with trastuzumab Fab. However, these might also be free sulfhydryl variants, which are also later eluting with the Fab fragment. Identification by high resolution mass spectrometry would be necessary to identify each variant. Once identified, the LC-UV method might then be appropriate for monitoring appropriately.

In summary, the use of cysteine proteases is an increasingly useful sample preparation technique to gain insight on mAb heterogeneity in conjunction with a wide pore core-shell LC column, bioZen™ 2.6 μm WidePore C4. The strategic partial proteolysis allows for characterization of specific regions of the mAb. By combining cysteine proteases IgdE and IdeS, one can gain insight on both conserved and hypervariable regions of the antibody, as well as identification of other variants, such as free sulfhydryls, all while gaining resolution and ruggedness associated with bioZen core-shell LC columns.

LC Conditions

- Column:** [bioZen 2.6 μm WidePore C4](#)
- Dimensions:** 100 x 2.1 mm
- Part No.:** [00D-4786-AN](#)
- Mobile Phase A:** 0.1 % TFA in Water
- Mobile Phase B:** 0.1 % TFA in Acetonitrile
- Gradient Program:** 30-40% B in 5 minutes (Fig 1-2)
- Flow Rate:** 0.8 mL/min
- Temperature:** 80°C
- Detection:** UV @ 214 nm
- Injection:** NIST, IdeS/IgdE Digested (0.5 mg/mL), Figure 1
Trastuzumab, IdeS/IgdE Digested (0.5 mg/mL), Figure 2





Figure 1. NIST mAb, IdeS and IgE Digested


Separation of mAb fragments generated by IdeS and IgE digested NIST mAb (RM 8671). Inset shows the early eluting hinge, as well as putative sulfhydryl variants. *Impurity peaks generated as a result of digestion.



App ID 26140

Figure 2. Trastuzumab, IdeS and IgE Digested

Separation of mAb fragments generated by IdeS and IgE digested trastuzumab. Inset shows separation of free sulfhydryl variants. Fab fragment also shows deamidated variants associated with trastuzumab Fab but may also contain potential sulfhydryl variants.



App ID 26141

APPLICATIONS

Reference

¹Faid, Valegh et al. "Middle-up analysis of monoclonal antibodies after combined IgD and IdeS hinge proteolysis: Investigation of free sulfhydryls." *Journal of pharmaceutical and biomedical analysis* vol. 149 (2018): 541-546. doi:10.1016/j.jpba.2017.11.046

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

<p>Australia t: +61 (0)2-9428-6444 auiinfo@phenomenex.com</p>	<p>India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com</p>	<p>Portugal t: +351 221 450 488 ptinfo@phenomenex.com</p>
<p>Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com</p>	<p>Ireland t: +353 (0)1 247 5405 eirinfo@phenomenex.com</p>	<p>Singapore t: +65 800-852-3944 sginfo@phenomenex.com</p>
<p>Belgium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com</p>	<p>Italy t: +39 051 6327511 italiainfo@phenomenex.com</p>	<p>Spain t: +34 91-413-8613 espinfo@phenomenex.com</p>
<p>Canada t: +1 (800) 543-3681 info@phenomenex.com</p>	<p>Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com</p>
<p>China t: +86 400-606-8099 cninfo@phenomenex.com</p>	<p>Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com</p>	<p>Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com</p>
<p>Denmark t: +45 4824 8048 nordicinfo@phenomenex.com</p>	<p>The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com</p>
<p>Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com</p>	<p>New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com</p>	<p>United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com</p>
<p>France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com</p>	<p>Norway t: +47 810 02 005 nordicinfo@phenomenex.com</p>	<p>USA t: +1 (310) 212-0555 info@phenomenex.com</p>
<p>Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com</p>	<p>Poland t: +48 22 104 21 72 pl-info@phenomenex.com</p>	<p>🌐 All other countries/regions Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com</p>

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen and BE-HAPPY are trademarks of Phenomenex.

Comparisons may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.