

Analysis of Alcohols in Hand Sanitizer and Impurities in Alcohol by GC-FID using Zebron™ ZB-624PLUS™ GC Columns

Dr. Ramkumar Dhandapani, Zandra Baja, Zara Jalai, and Dr. Bryan Tackett
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA


Overview

Alcohol-based hand sanitizers are used on a daily basis to kill microorganisms including bacteria. Alcohol-based sanitizers generally contains over 60% of alcohol. However, if the alcohols are not at the right concentrations, it becomes ineffective. Thus, quantitative estimation of alcohol % in hand sanitizer via GC-FID is a very important analysis for identifying and quantifying the alcohols and content in hand sanitizer for both quality and branding purposes. In addition, the quality of alcohol used as raw material is equally important to avoid adverse effects. In this application note, we have developed a GC-FID method for ascertaining the purity of raw material alcohol for impurity levels and quantitative method for % alcohol in hand sanitizer using a single method on a Zebron ZB-624PLUS GC Column.

GC-FID Conditions

Column: Zebron ZB-624PLUS
Dimension: 30 meter x 0.32 mm x 1.80 µm
Part No.: [7HM-G040-31](#)
Injection: Split 20:1 @ 200 °C, 1 µL
Recommended Liner: Zebron PLUS Z-Liner™ (Compatible with Agilent® & Thermo® GC instrument)
Liner Part No.: [AG2-0A03-05](#)
Carrier Gas: Helium @ 25 cm/sec (Constant Flow)
Oven Program: 36 °C for 12 min, 260 °C @ 10 °C/min for 15 min
Detector: FID
Detector Temperature: 280 °C

Figure 1. 50 ppm Acetaldehyde and 50 ppm methanol in 99.5% ethanol



Sample:

1. Acetaldehyde
2. Methanol
3. Ethanol

App ID 26182



Figure 2. 100 ppm Acetaldehyde and 100 ppm methanol in 99.5% ethanol




Figure 3. 200 ppm methanol in 99.5% ethanol




Figure 4. USP impurities in 99.5% ethanol




Figure 5. Analysis of % content of Alcohol in Hand Sanitizer




Table 1. Precision Data for Ethanol Content in Hand Sanitizer on a Zebron™ ZB-624_{PLUS} GC Column

Injection #	Ethanol	
	Retention Time (min)	Peak Area
1	4.625	4719
2	4.627	4702
3	4.633	4650
4	4.630	4714
5	4.627	4633
6	4.626	4675
% RSD	0.1%	0.8%

Results and Discussion

Due to increased demand for alcohol-based hand sanitizers, there is a need for a fast and efficient GC method to accurately quantify the % of alcohol in hand sanitizer. In addition, it is very important to make sure that the raw material alcohol used in the sanitizer manufacturing meets safety requirement for impurities. Presented in **Figure 1 to 4** are the various impurities that are possible in alcohol raw material. The method for analysis utilizes a modified USP method that is still within allowable adjustments. The method provides identification and quantification of impurities that can be present in alcohol. The same method parameters were extended to the analysis of alcohol-based sanitizer. The sample alcohol sanitizer was dissolved in Dimethylformamide to quantify % alcohol content as shown in **Figure 5**. The precision of the method with respect to peak retention and peak area are presented in **Table 1**. In addition to analysis of alcohols, this method provides a way to analyse % alcohol content in hand sanitizer and to quantify impurities on a single method without the necessity to change method parameters and GC column. This is possible due to the optimal selectivity of ZB-624_{PLUS} for volatile compounds like alcohols. In addition, the column has extensive cross-linkage through Engineered Self-Cross Linking (ESC™) and thermal stability of 300/320 °C max temperature to bake out contaminants.

Conclusion

Zebron ZB-624_{PLUS} provides optimal selectivity and reproducibility for the analysis of alcohol and alcohol-based sanitizers in a single method.



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Zebtron, Engineered Self Cross-linking, ESC, 624^{PLUS}, and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific Inc.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Inc, or Thermo Scientific.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

