


APPLICATION

Forced Degradation of Cefaclor Analysis using Kinetex® EVO C18 Core-Shell LC Column

Matthew Trass and Simon Lomas
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Cefaclor is a broad-spectrum antibiotic used to treat a wide-range of bacterial infections including lung, ear, throat, skin, and urinary tract infections. As part of a pharmaceutical stability indicating method, forced degradation studies are required to demonstrate specificity. The results are used to elucidate degradation pathways and the structures of degradation products, which then assist with formulation and packaging. To acquire as much information as possible from the degradation experiment, it is important to resolve as many degradation peaks as possible. Utilizing the unique selectivity and high efficiency of a Kinetex EVO C18 core-shell column we achieve excellent peak shape and resolution between the active and degradedants within this Cefaclor forced degradation. Additionally, the versatility of the 2.6 µm EVO C18 particle size, allows for excellent performance on both HPLC and UHPLC instrumentation.



LC/UV Conditions

Column: Kinetex 2.6 µm EVO C18
Dimensions: 100 x 4.6 mm
Part No.: 00D-4725-E0
Mobile Phase: A: 0.1% H₃PO₄ in Water
B: 0.1% H₃PO₄ in Acetonitrile
Gradient:

Time (min)	% B
0	5
6	40

Flow Rate: 1.75 mL/min
Temperature: Ambient
Detection: UV @ 254 nm
Sample: 1. Cefaclor

Sample Pretreatment

Reconstitute 1 mg of neat Cefaclor standard with 1 mL of 1 M HCl. Heat at 60 °C for 6 hours. Neutralize with 1 mL of 1 M NaOH. Refrigerate until ready to inject.



APPLICATION

Kinetex® EVO C18 Ordering Information

5 µm Minibore Columns (mm)				SecurityGuard™ ULTRA Cartridges [†]
30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJO-9298
for 2.1 mm ID				

5 µm MidBore™ Columns (mm)			SecurityGuard™ ULTRA Cartridges [†]
50 x 3.0	100 x 3.0	150 x 3.0	3/pk
00B-4633-YO	00D-4633-YO	00F-4633-YO	AJO-9297
for 3.0 mm ID			

5 µm Analytical Columns (mm)				SecurityGuard™ ULTRA Cartridges [†]
50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
00B-4633-EO	00D-4633-EO	00F-4633-EO	00G-4633-EO	AJO-9296
for 4.6 mm ID				

2.6 µm Minibore Columns (mm)				SecurityGuard™ ULTRA Cartridges [†]
30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00A-4725-AN	00B-4725-AN	00D-4725-AN	00F-4725-AN	AJO-9298
for 2.1 mm ID				

2.6 µm MidBore™ Columns (mm)			SecurityGuard™ ULTRA Cartridges [†]
50 x 3.0	100 x 3.0	150 x 3.0	3/pk
00B-4725-YO	00D-4725-YO	00F-4725-YO	AJO-9297
for 3.0 mm ID			

2.6 µm Analytical Columns (mm)			SecurityGuard™ ULTRA Cartridges [†]
50 x 4.6	100 x 4.6	150 x 4.6	3/pk
00B-4725-EO	00D-4725-EO	00F-4725-EO	AJO-9296
for 4.6 mm ID			

1.7 µm Minibore Columns (mm)			SecurityGuard™ ULTRA Cartridges [†]
50 x 2.1	100 x 2.1	150 x 2.1	3/pk
00B-4726-AN	00D-4726-AN	00F-4726-AN	AJO-9298
for 2.1 mm ID			

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

Australia
t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China
t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark
t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico
t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico
t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain
t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

**All other countries
Corporate Office USA** 
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide.
For the distributor in your country, contact Phenomenex USA,
International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at
<http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Kinetex is a registered trademark and SecurityGuard and MidBore are trademarks of Phenomenex.

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2015 Phenomenex, Inc. All rights reserved.



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.