

Improve Your European Pharmacopoeia (Ph. Eur.) and United States (USP) Monographs


- Reduce run times
- Achieve higher resolution
- Stay within Allowable Adjustments



High-throughput productivity is of critical importance for laboratories undertaking testing for generic drugs following the quality standards and test procedures of the United States Pharmacopeia (USP) and European Pharmacopoeia Monographs (Ph. Eur.). This guide provides analysts with solutions to standard USP and Ph. Eur. Monographs and also incorporates cutting edge Kinetex® core-shell LC columns to provide shorter separation times and improved resolution while meeting all the quality standards of the United States Pharmacopoeia and European Pharmacopoeia Monographs.

Table of Contents

European Pharmacopoeia (Ph. Eur.) Monographs for Generic Drugs.....	3-35
Allopurinol.....	4
Amlodipine Besylate.....	6
Atenolol.....	8
Carvedilol.....	10
Clarithromycin.....	12
Fluconazole.....	14
Fluoxetine Hydrochloride.....	16
Metoprolol Tartrate.....	18
Oxycodone Hydrochloride.....	20
Paroxetine Hydrochloride.....	22
Potassium Clavulanate.....	24
Pravastatin Sodium.....	26
Simvastatin.....	28
Tamsulosin Hydrochloride.....	30
Tramadol Hydrochloride.....	32
Trimethoprim.....	34
United States Pharmacopeia (USP) Monographs for Generic Drugs.....	36-53
Amlodipine Besylate.....	38
Clavulanate Potassium.....	40
Fluticasone Propionate.....	42
Ibuprofen.....	44
Lovastatin.....	46
Metformin Hydrochloride.....	48
Pravastatin Sodium.....	50
Trimethoprim.....	52
Ordering Information.....	54-58




**European
Pharmacopoeia
(Ph. Eur.)
Monographs
for Generic Drugs**

Allopurinol and Related Substances

Ph. Eur. monograph 0576

The related substances test of the Ph. Eur. Monograph 0576 outlines the separation of all relevant impurities from Allopurinol. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.




Ph. Eur. Monograph 0576 Details

Test Solution (a)	Dissolve 25.0 mg of Allopurinol CRS* in 2.5 mL of a 4 g/L solution of sodium hydroxide R and dilute immediately to 50.0 mL with the mobile phase
Reference Solution	(a) Dilute 2.0 mL of the test solution (a) to 100.0 mL with the mobile phase. Dilute 5 mL of this solution to 100.0 mL with the mobile phase. (b) Dissolve 5.0 mg of Allopurinol Impurity A CRS*, 5.0 mg of Allopurinol Impurity B CRS* and 5.0 mg of Allopurinol Impurity C CRS* in 5.0 mL of a 4 g/L solution of sodium hydroxide R and dilute immediately to 100.0 mL with the mobile phase. Dilute 1.0 mL of this solution to 100.0 mL with the mobile phase.
Column	
Size	250 x 4.6 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R (5 µm)
Mobile Phase	1.25 g/L solution of potassium dihydrogen phosphate R
Flow Rate	1.4 mL/min
Detection	Spectrophotometer @ 230 nm
Injection	20 µL (reference solution (a) and (b))
Run Time	Twice the retention time of Allopurinol
Elution Order	1. Impurity A 2. Impurity B 3. Impurity C 4. Allopurinol (about 10 min)

System Suitability

Reference Solution (b) Minimum resolution of 1.1 between peaks due to Impurities B and C

* Allopurinol CRS (A0350000), Allopurinol Impurity A CRS (A0350010), Allopurinol Impurity B CRS (A0350020) and Allopurinol Impurity C CRS (A0350030) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allee Kastner CS 30026F - 67081 STRASBOURG (France).



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1
Mobile Phase pH	± 0.2 units	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 0576 Details Table
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 0576 Details Table
Wavelength of Detector	No deviations permitted	230 nm (as specified)
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)
Column Temperature	± 10 °C	Ambient (as specified)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)
Column Length	± 70 %	250 mm (as specified)
Column Internal Diameter	± 25 %	4.6 mm (as specified)
Particle Size	-50 %	5 µm (as specified)
Flow Rate	± 50 %	1.4 mL/min (as specified)

Amlodipine Besylate and Related Substances

Ph. Eur. monograph 1491

The related substances test of the Ph. Eur. Monograph 1491 outlines the separation of all relevant impurities from Amlodipine Besylate. This method was studied and improvements were made to provide higher resolution (R_s) and a faster separation time within allowable adjustments.



Ph. Eur. Monograph 1491 Details

Reference Solution (b) Dissolve 2.5 mg of Amlodipine Impurity B CRS* and 2.5 mg of Amlodipine Impurity G CRS* in the mobile phase and dilute to 25 mL with the mobile phase. Dilute 1.0 mL of the solution to 10.0 mL with the mobile phase.

Column

Size	250 x 4.0 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R (5 μ m)
Temperature	30 °C
Mobile Phase	2.3 g/L solution of Ammonium acetate R, methanol R (30:70 V/V)
Flow Rate	1.5 mL/min
Detection	Spectrophotometer @ 237 nm
Injection	20 μ L
Run Time	Twice the retention time of Amlodipine

Relative Retention with Reference to Amlodipine (about 20 min)**

Impurity G	about 0.21
Impurity B	about 0.25

System Suitability

Reference Solution (b) Minimum resolution of 2.0 between peaks due to Impurities G and B

* Amlodipine Impurity B CRS (Y0001069) and Amlodipine Impurity G CRS (Y0001070) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).


** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1

Alternative Method Within Allowable Adjustments

Column: Luna® C18(2) 5 μ m Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: [00G-4252-E0](#)
Flow Rate: 1.5 mL/min
Sample: 1. Impurity G
 2. Impurity B
 3. Amlodipine

Elution Time of Last Peak: 8.9 min
 R_s Impurities G and B: 11.63





Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 1491 Details Table	As specified
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 1491 Details Table	As specified
Wavelength of Detector	No deviations permitted	237 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified
Column Temperature	± 10 °C	30 °C (as specified)	As specified
Stationary Phase	No change of the identity of the sub- stituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)	As specified
Column Length	± 70 %	250 mm (as specified)	As specified
Column Internal Diameter	± 25 %	4.6 mm (+15 %)	4.6 mm (+15 %)
Particle Size	-50 %	5 µm (as specified)	As specified
Flow Rate	± 50 %	1.5 mL/min (as specified)	As specified

Atenolol and Related Substances

Ph. Eur. monograph 0703

The related substances test of the Ph. Eur. Monograph 0703 outlines the separation of all relevant impurities from Atenolol. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Atenolol

Ph. Eur. Monograph 0703 Details

Reference Solution (a)	Dissolve 2 mg of Atenolol for system suitability CRS* (containing Impurities B, F, G, I and J) in 1 mL mobile phase
Column	
Size	125 x 4.0 mm
Stationary Phase	End-capped octadecylsilyl silica gel for chromatography R (5 µm)
Mobile Phase	Dissolve 1.0 g of sodium octanesulphonate R and 0.4 g of tetrabutylammonium hydrogen sulfate R in 1 L of a mixture of 20 volumes of tetrahydrofuran R, 180 volumes methanol R2 and 800 volumes of 3.4 g/L solution of potassium dihydrogen phosphate R; adjust the apparent pH to 3.0 with phosphoric acid R.
Flow Rate	0.6 mL/min
Detection	Spectrophotometer @ 226 nm
Injection	10 µL
Run Time	5 times the retention time of Atenolol
Relative Retention with Reference to Atenolol (about 8 min)**	
Impurity B	about 0.3
Impurity J	about 0.7
Impurity I	about 0.8
Impurity F	about 2.0 (pair of peaks)
Impurity G	about 3.5

System Suitability

Reference Solution (a)	Minimum resolution of 1.4 between peaks due to Impurities J and I
-------------------------------	---

* Atenolol for system suitability CRS (Y0001089) was purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allee Kastner CS 30026F - 67081 STRASBOURG (France).


** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1

Original Method as Described in the Monograph

Column: Luna® C18(2) 5 µm Fully Porous
Dimensions: 125 x 4.0 mm
Part No.: [00E-4252-DO](#)
Flow Rate: 0.6 mL/min
Sample: 1. Impurity B
 2. Impurity A
 3. Impurity J
 4. Impurity I
 5. Impurities D and E
 6. Impurity F
 7. Impurity G
 8. Impurity H

Elution Time of Last Peak: 33.3 min
Rs Impurities I and J: 2.0




App ID: 18825

Reduce run time by >10 min

Method 2

Faster and Higher Resolution Within Allowable Adjustments




Column: Kinetex® Core-Shell C18 2.6µm
Dimensions: 100 x 4.6mm
Part No.: 00D-4462-E0
Flow Rate: 0.6mL/min
Sample: 1. Impurity B
 2. Impurity A
 3. Impurity J
 4. Impurity I
 5. Impurities D and E
 6. Impurity F
 7. Impurity G
 8. Impurity H

Elution Time of Last Peak: 20.0 min
Rs Impurities I and J: 2.3

App ID: 18826

Method 3

Even Faster and Higher Resolution within Allowable Adjustments



Column: Kinetex Core-Shell C18 2.6µm
Dimensions: 100 x 4.6mm
Part No.: 00D-4462-E0
Flow Rate: 0.9mL/min
Sample: 1. Impurity B
 2. Impurity A
 3. Impurity J
 4. Impurity I
 5. Impurities D and E
 6. Impurity F
 7. Impurity G
 8. Impurity H

Elution Time of Last Peak: 11.9 min
Rs Impurities I and J: 3.72

App ID: 18827

Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2	Method 3
Mobile Phase pH	± 0.2 units	3 (as specified)	As specified	As specified
Concentration of Salts in Buffer	± 10%	As specified in Monograph 0703 Details Table	As specified	As specified
Composition of the Mobile Phase	± 30% of the minor solvent component relative or 2% absolute, whichever is the larger. No other component is altered by more than 10% absolute.	As specified in Monograph 0703 Details Table	As specified	As specified
Wavelength of Detector	No deviations permitted	226 nm (as specified)	As specified	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	10 µL (as specified)	As specified	As specified
Column Temperature	± 10°C	Ambient (as specified)	As specified	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	End-capped octadecylsilyl silica gel for chromatography (as specified)	As specified	As specified
Column Length	± 70%	125 mm (as specified)	100 mm (-20%)	100 mm (-20%)
Column Internal Diameter	± 25%	4.0 mm (as specified)	4.6 mm (+15%)	4.6 mm (+15%)
Particle Size	-50%	5 µm (as specified)	2.6 µm (-48%)	2.6 µm (-48%)
Flow Rate	± 50%	0.6 mL/min (as specified)	As specified	0.9 mL/min (+ 50%)

Carvedilol and Related Substances

Ph. Eur. monograph 1745

The Ph. Eur. Monograph 1745 outlines the separation of Carvedilol from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Ph. Eur. Monograph 1745 Details

Reference Solution

(b) Dissolve 5 mg of Carvedilol Impurity C CRS* in 5.0 mL of the mobile phase and dilute to 100.0 mL with the mobile phase. Dilute 4.0 mL of the solution to 100.0 mL with the mobile phase. Dilute 1.0 mL of this solution to 10.0 mL with the mobile phase.

(c) Dissolve 5 mg of Carvedilol for system suitability CRS* (containing Impurities A and D) in the mobile phase and dilute to 50.0 mL with the mobile phase.

Column

Size	150 x 4.6 mm
Stationary Phase	End-capped octylsilyl silica gel for chromatography R (5 µm)
Temperature	55 °C
Mobile Phase	Dissolve 1.77 g of potassium dihydrogen phosphate R in water and dilute to 650 mL with the same solvent; adjust to pH 2.0 with phosphoric acid R and add 350 mL of acetonitrile R
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 240 nm
Injection	20 µL
Run Time	6 times the retention time of Carvedilol

Relative Retention with Reference to Carvedilol (about 4 min)**


Impurity A	about 0.5
Impurity C	about 2.9
Impurity D	about 3.8

System Suitability

Reference Solution (b) Minimum resolution of 3.5 between peaks due to Impurity A and Carvedilol

*Carvedilol Impurity C CRS (Y0000103) and Carvedilol for system suitability CRS (Y0001426) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1
Mobile Phase pH	± 0.2 units	2.0 (as specified)
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 1745 Details Table
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 1745 Details Table
Wavelength of Detector	No deviations permitted	240 nm (as specified)
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)
Column Temperature	± 10 °C	55 °C (as specified)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	Octylsilyl silica gel for chromatography (as specified)
Column Length	± 70 %	150 mm (as specified)
Column Internal Diameter	± 25 %	4.6 mm (as specified)
Particle Size	-50 %	5 µm (as specified)
Flow Rate	± 50 %	1.0 mL/min (as specified)

Clarithromycin and Related Substances

Ph. Eur. monograph 1651

The Ph. Eur. Monograph 1651 outlines the separation of Clarithromycin from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Clarithromycin

Ph. Eur. Monograph 1651 Details

Reference Solution (d) Dissolve 15 mg of Clarithromycin for peak identification CRS* in 5 mL of acetonitrile and dilute to 10 mL with water

Column

Size 100 x 4.6 mm

Stationary Phase Octadecylsilyl silica gel for chromatography R (3.5 μm)

Temperature 40 °C

Mobile Phase
A: 4.76 g/L solution of potassium dihydrogen phosphate adjusted to pH 4.4 with dilute phosphoric acid
B: Acetonitrile

Gradient

Time (min)	%B
0 – 32 min	25 → 60
32– 34 min	60

Flow Rate 1.1 mL/min

Detection Spectrophotometer @ 205 nm

Injection 10 μL

Relative Retention with Reference to Clarithromycin (about 11 min)**

Impurity A	about 0.42
Impurity J	about 0.63
Impurity L	about 0.74
Impurity B	about 0.79
Impurity M	about 0.81
Impurity C	about 0.89
Impurity D	about 0.96
Impurity N	about 1.15
Impurity E	about 1.27
Impurity F	about 1.33
Impurity P	about 1.35
Impurity O	about 1.41
Impurity K	about 1.59
Impurity G	about 1.59
Impurity H	about 1.82

System Suitability

Peak-to-Valley Ratio Minimum 3.0, where H_p = height above the baseline of the peak due to Impurity D and H_v = height above the baseline of the lowest point of the curve separating this peak from the peak due to Clarithromycin in the chromatogram obtained with reference solution D

* Ph. Eur. Standard Clarithromycin for peak identification CRS Y0000321 was purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1

Original Method as Described in the Monograph

Achieve improved sensitivity and resolution using Kinetex Core-Shell Columns

Column: Kinetex® Core-Shell XB-C18 3.5µm
Dimensions: 100 x 4.6mm
Part No.: 00D-4744-E0
Flow Rate: 1.1 mL/min
Sample: 1. Impurity I
 2. Impurity A
 3. Impurity L
 4. Impurity B
 5. Impurity C
 6. Impurity D
 7. Clarithromycin
 8. Impurity E
 9. Impurity F
 10. Impurity P
 11. Impurity G
 12. Impurity H

Elution Time of Last Peak: 17.1 min
Peak-to-Valley Ratio: 10.8



App ID: 23503

Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1
Mobile Phase pH	± 0.2 units	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 1651 Details Table
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 1651 Details Table
Wavelength of Detector	No deviations permitted	205 nm (as specified)
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	1 µL (as specified)
Column Temperature	± 10 °C	40 °C (as specified)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)
Column Length	± 70 %	100 mm (as specified)
Column Internal Diameter	± 25 %	4.6 mm (as specified)
Particle Size	-50 %	3.5 µm (as specified)
Flow Rate	± 50 %	1.1 mL/min (as specified)

Fluconazole and Related Substances

Ph. Eur. monograph 2287

The Ph. Eur. Monograph 2287 outlines the separation of Fluconazole from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Fluconazole

Ph. Eur. Monograph 2287 Details

Reference Solution	(b) Dissolve 5.0 mg of Fluconazole for peak identification CRS* (containing Impurity A) in the mobile phase, sonicate if necessary, and dilute to 10 mL with the mobile phase
	(c) Dissolve 3.0 mg of Fluconazole Impurity B CRS* in the mobile phase, sonicate if necessary, and dilute to 100 mL with the mobile phase
	(d) Dissolve 3.0 mg of Fluconazole Impurity C CRS* in the mobile phase and dilute to 20 mL with the mobile phase

Column

Size	150 x 4.6 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R1 (5 µm)
Temperature	40 °C
Mobile Phase	Acetonitrile R, 0.63 g/L solution of ammonium formate R (14:86 V/V)
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 260 nm
Injection	20 µL
Run Time	3.5 times the retention time of Fluconazole

Relative Retention with Reference to Fluconazole (about 11 min)**

Impurity B	about 0.4
Impurity A	about 0.5
Impurity C	about 0.8

System Suitability

Reference Solution (a)	Minimum resolution of 3.0 between peaks due to Impurity C and Fluconazole
-------------------------------	---

* Fluconazole for peak identification CRS (Y0000558), Fluconazole Impurity B CRS (Y0000573) and Fluconazole Impurity C CRS (Y0000574) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1


Original Method as Described in the Monograph

Column: Luna® C18(2) 5 µm Fully Porous
Dimensions: 150 x 4.6 mm
Part No.: [00F-4252-E0](#)
Flow Rate: 1.0 mL/min

Sample: 1. Impurity B
 2. Impurity A
 3. Impurity C
 4. Fluconazole

Elution Time of Last Peak: 15.9 min

Rs Impurity C and Fluconazole: 4.39




Method 2

Faster Method Within Allowable Adjustments

Column: Luna C18(2) 3µm Fully Porous
Dimensions: 100 x 4.6 mm
Part No.: [00D-4251-E0](#)
Flow Rate: 1.5 mL/min
Sample: 1. Impurity B
 2. Impurity A
 3. Impurity C
 4. Fluconazole

Elution Time of Last Peak: 7.97 min

Rs Impurity C and Fluconazole: 3.55



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 2287 Details Table	As specified
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 2287 Details Table	As specified
Wavelength of Detector	No deviations permitted	260 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified
Column Temperature	± 10 °C	40 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)	As specified
Column Length	± 70 %	150 mm (as specified)	100 mm (-33.3 %)
Column Internal Diameter	± 25 %	4.6 mm (as specified)	As specified
Particle Size	-50 %	5 µm (as specified)	3 µm (-40 %)
Flow Rate	± 50 %	1.0 mL/min (as specified)	1.5 mL/min (+50 %)

Fluoxetine Hydrochloride and Related Substances

Ph. Eur. monograph 1104

The Ph. Eur. Monograph 1104 outlines the separation of Fluoxetine from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Fluoxetine Hydrochloride

Ph. Eur. Monograph 1104 Details

Reference Solution Dissolve 22 mg of Fluoxetine Hydrochloride CRS* in 10.0 mL of a 0.5 M sulfuric acid. Heat at about 85° C for 3 h. Allow to cool. The resulting solution contains considerable quantities of Impurity A and usually also contains 4-trifluoromethylphenol. To 0.4 mL of this solution add 28.0 mg of Fluoxetine hydrochloride CRS*, about 1 mg of Fluoxetine Impurity B CRS* and about 1 mg Fluoxetine Impurity C CRS*, then dilute to 25.0 mL with mobile phase.

Column

Size	150 x 4.6 mm
Stationary Phase	Octylsilyl silica gel for chromatography R (5 µm)
Mobile Phase	Mix 8 volumes of methanol R, 30 volumes of tetrahydrofuran R, and 62 volumes of a solution of trimethylamine R prepared as follows: to 10 mL of trimethylamine R, add 980 mL of water R, mix and adjust to pH 6.0 with phosphoric acid R (about 4.5 mL) and dilute to 1000 mL with water R.
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 215 nm
Injection	10 µL
Run Time	3 times the retention time of Fluoxetine

Relative Retention with Reference to Fluoxetine (about 10-18 min)**

Impurity A	about 0.24
Impurity B	about 0.27
Impurity C	about 0.90

System Suitability

Peak-to-Valley Ratio Minimum 11, where Hp = height above the baseline of the peak due to Impurity C and Hv = height above the baseline of the lowest point of the curve separating this peak from the peak due to Fluoxetine.

* Fluoxetine hydrochloride CRS (F0253000), Fluoxetine Impurity B CRS (F0253020) and Fluoxetine Impurity C CRS (F0253030) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).


** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1


Original Method as Described in the Monograph

Column: Luna® C8(2) 5 µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: 00G-4249-E0
Flow Rate: 1.0 mL/min
Sample: 1. Impurity A
 2. Impurity B
 3. Impurity C
 4. Fluoxetine

Elution Time of Last Peak: 10.4 min
Rs Impurity C and Fluoxetine: 2.66
Peak-to-Valley Ratio: >11



App ID 24199



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	6 (as specified)	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 1104 Details Table	As specified
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 1104 Details Table	As specified
Wavelength of Detector	No deviations permitted	215 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	10 µL (as specified)	As specified
Column Temperature	± 10 °C	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	Octylsilyl silica gel for chromatography (as specified)	As specified
Column Length	± 70 %	250 mm (as specified)	As specified
Column Internal Diameter	± 25 %	4.6 mm (as specified)	As specified
Particle Size	-50 %	5 µm (as specified)	As specified
Flow Rate	± 50 %	1.0 mL/min (as specified)	As specified

Metoprolol Tartrate and Related Substances

Ph. Eur. monograph 1028

The Ph. Eur. Monograph 1028 outlines the separation of Metoprolol Tartrate from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Metoprolol Tartrate

Ph. Eur. Monograph 1028 Details

Reference Solution (a) and Impurity G Dissolve 1.5 mg of Metoprolol Impurity A CRS*, 2.5 mg of Metoprolol Tartrate CRS* and 1.5 mg of Metoprolol impurity G in the mobile phase and dilute to 50.0 mL with the mobile phase.

Column

Size	150 x 3.9 mm
Stationary Phase	End-capped octadecylsilyl silica gel for chromatography R (5 µm)
Mobile Phase	Dissolve 3.9 g of ammonium acetate R in 810 mL of water R, add 2.0 mL of trimethylamine R, 10.0 mL of glacial acetic acid R, 3.0 mL of phosphoric acid R and 146 mL of acetonitrile R and mix
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 280 nm
Injection	20 µL
Run Time	3 times the retention time of Metoprolol
Elution Order	1. Impurity G 2. Impurity A 3. Metoprolol Tartrate

System Suitability

Reference Solution (a) Minimum resolution of 6.0 between peaks due to Impurity A and Metoprolol

* Metoprolol impurity A CRS (Y0000145) and Metoprolol Tartrate CRS (M1830000) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

Method 1

Original Method as Described in the Monograph within Allowable Adjustments

Column: Luna® C18(2) 5 µm Fully Porous

Dimensions: 150 x 4.6 mm


Part No.: 00F-4252-EO

Flow Rate: 1.0 mL/min

Sample: 1. Impurity A
2. Metoprolol Tartrate

Elution Time of Last Peak: 13.0 min


Rs Impurity A and Metoprolol: 10.10



Method 2 Faster Method Within Allowable Adjustments

Column: Kinetex® Core-Shell C18 2.6µm
Dimensions: 150 x 4.6mm
Part No.: [00F-4462-E0](#)
Flow Rate: 1.0mL/min
Sample: 1. Impurity A
 2. Metoprolol Tartrate

Elution Time of Last Peak: 7.30 min
Rs Impurity A and Metoprolol: 11.82




App ID 19014

Method 3 Even Faster Method within Allowable Adjustments

Column: Kinetex Core-Shell C18 2.6µm
Dimensions: 100 x 4.6mm
Part No.: [00D-4462-E0](#)
Flow Rate: 1.5mL/min
Sample: 1. Impurity A
 2. Metoprolol Tartrate

Elution Time of Last Peak: 3.30 min
Rs Impurity A and Metoprolol: 9.21



App ID 19013


Adjustments for Meeting System Suitability (European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2	Method 3
Mobile Phase pH	± 0.2 units	As specified	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 1028 Details Table	As specified	As specified
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 1028 Details Table	As specified	As specified
Wavelength of Detector	No deviations permitted	280 nm (as specified)	As specified	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified	As specified
Column Temperature	± 10°C	Ambient (as specified)	As specified	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	End-capped octadecylsilyl silica gel for chromatography (as specified)	As specified	As specified
Column Length	± 70 %	150 mm (as specified)	150 mm (as specified)	100 mm (-33 %)
Column Internal Diameter	± 25 %	4.6 mm (+18 %)	4.6 mm (+18 %)	4.6 mm (+15 %)
Particle Size	-50 %	5 µm (as specified)	2.6 µm (-48 %)	2.6 µm (-48 %)
Flow Rate	± 50 %	1.0 mL/min (as specified)	As specified	1.5 mL/min (+ 50 %)

Oxycodone Hydrochloride and Related Substances

Ph. Eur. monograph 1793

The Ph. Eur. Monograph 1793 outlines the separation of Oxycodone Hydrochloride from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Oxycodone Hydrochloride

Ph. Eur. Monograph 1793 Details

Test Solution	Dissolve 0.100 g of Oxycodone Hydrochloride CRS* in a 1 % V/V solution of dilute acetic acid R and dilute to 50.0 mL with the same solvent
Reference Solution	(a) Dissolve 20.0 mg of Oxycodone Impurity D CRS* in a 1.0 % V/V solution of dilute acetic acid R and dilute to 10.0 mL with the same solution (b) To 1.0 mL of the test solution, add 1 mL of reference solution (a) and dilute to 100.0 mL with a 1 % V/V solution of dilute acetic acid R. Dilute 1.0 mL of the solution to 10.0 mL with a 1.0 % V/V solution of dilute acetic acid R.
Column	
Size	150 x 4.6 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R (5 µm)
Temperature	40 °C
Mobile Phase	A: Mix 830 mL of a 1.1 g/L solution of sodium heptanesulfonate monohydrate R previously adjusted to pH 2.0 with a mixture of equal volumes of phosphoric acid R and water R, with 70 mL of acetonitrile R and 100 mL of methanol R B: Mix 600 mL of a 1.1 g/L solution of sodium heptanesulfonate monohydrate R previously adjusted to pH 2.0 with a mixture of equal volumes of phosphoric acid R and water R, with 150 mL of acetonitrile R and 250 mL of methanol R
Gradient	Time (min) %B 0 – 60 min 0 – 50
Flow Rate	1.5 mL/min
Detection	Spectrophotometer @ 230 nm
Injection	20 µL
Relative Retention with Reference to Oxycodone (about 24 min)**	
Impurity B	about 0.7
System Suitability	

Reference Solution (a) Minimum resolution of 3.0 between peaks due to Oxycodone and Impurity D

* Oxycodone Hydrochloride CRS (Y0000492) and Oxycodone Impurity D CRS (Y0000453) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1

Original Method as Described in the Monograph

Column: Luna® C18(2) 5 µm Fully Porous

Dimensions: 150 x 4.6 mm


Part No.: [00F-4252-E0](#)

Flow Rate: 1.5 mL/min

Sample: 1. Oxycodone
2. Impurity B

Elution Time of Last Peak: 34.51 min

Rs Oxycodone and Impurity D: 3.79




App. ID 21.485

Method 2**Faster Method Utilizing Core-Shell Technology**

Column: Kinetex® Core-Shell C18 5 µm
Dimensions: 150 x 4.6 mm
Part No.: [00F-4601-E0](#)
Flow Rate: 1.5 mL/min

Sample: 1. Oxycodone
 2. Impurity B

Elution Time of Last Peak: 23.30 min
Rs Oxycodone and Impurity D: 4.75

**Adjustments for Meeting System Suitability**


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (gradient elution)	Method 1	Method 2
Mobile Phase pH	No adjustment permitted	2 (as specified)	As specified
Concentration of Salts in Buffer	No adjustment permitted	As specified in Monograph 1793 Details Table	As specified
Composition of the Mobile Phase	Minor adjustments of the composition of the mobile phase and the gradient are acceptable, if the system suitability requirements are met, the principle peak(s) elute(s) within $\pm 15\%$ of the indicated retention time(s) and the final elution power of the mobile phase is not weaker in elution power than the prescribed composition	As specified in Monograph 1793 Details Table	As specified
Wavelength of Detector	No deviations permitted	230 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified
Column Temperature	$\pm 5^\circ\text{C}$	40 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C-18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)	As specified
Column Length	May be decreased, $\pm 70\%$	150 mm (as specified)	As specified
Column Internal Diameter	$\pm 25\%$	4.6 mm (as specified)	As specified
Particle Size	No adjustment permitted	5 µm (as specified)	As specified
Flow Rate	Acceptable when changing the column dimensions	1.5 mL/min (as specified)	As specified

Paroxetine Hydrochloride and Related Substances

Ph. Eur. monograph 2283

The Ph. Eur. Monograph 2283 outlines the separation of Paroxetine Hydrochloride from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Paroxetine Hydrochloride

Ph. Eur. Monograph 2283 Details

Solvent Mixture	Tetrahydrofuran R, water R (10:90 V/V)
Test Solution	Dissolve 50.0 mg of Paroxetine Hydrochloride (anhydrous) CRS* in the solvent mixture and dilute to 50.0 mL with the solvent mixture
Reference Solution	(a) Dilute 5.0 mL of the test solution to 50.0 mL with the solvent mixture (c) Dissolve 5.0 mg of anhydrous Paroxetine Hydrochloride Impurity C CRS* in 25 mL of tetrahydrofuran R and dilute to 50.0 mL with water R (f) Dissolve 2.5 mg of Paroxetine Impurity E CRS* in the solvent mixture, add 2.5 mL of the test solution and dilute to 100.0 mL with the solvent mixture (g) Dissolve 5 mg of Paroxetine Impurity A CRS* in the solvent mixture and dilute to 50 mL with the solvent mixture

Column

Size	250 x 4.6 mm
Stationary Phase	End-capped octadecylsilyl silica gel for chromatography R (5 µm)
Temperature	40 °C
Mobile Phase	A: Trifluoroacetic acid R, tetrahydrofuran R, water R (5:100:900 V/V/V) B: Trifluoroacetic acid R, tetrahydrofuran R, acetonitrile R (5:100:900 V/V/V)

Gradient	Time (min)	%B
	0 – 30 min	20
	30 – 50 min	20 → 80
	50 – 55 min	80
	55 – 60 min	80 → 20
60 – 65 min	20	

Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 295 nm
Injection	20 µL of the test solution and reference solutions

Relative Retention with Reference to Paroxetine (about 28 min)**

Impurity A	about 0.8
Impurity E	about 0.9
Impurity C	about 1.02

System Suitability

Reference Solution (b)	Minimum resolution of 3.5 between peaks due to Impurity E and Paroxetine
-------------------------------	--

* Paroxetine hydrochloride (anhydrous) CRS (Y0000578), Anhydrous Paroxetine Hydrochloride Impurity C CRS (Y0000579), Paroxetine Impurity E CRS (Y0000580) and Paroxetine Impurity A CRS (Y0000233) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).


** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1

Original Method as Described in the Monograph

Column: Luna® C8(2) 5 µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: [00G-4249-E0](#)
Flow Rate: 1.0 mL/min
Sample: 1. Impurity A
 2. Impurity E
 3. Paroxetine
 4. Impurity C

Elution Time of Last Peak: 42.5 min
Rs Impurity E and Paroxetine: 6.06
Impurity E Peak Height: 0.14 mAU




Method 2

Faster Method Utilizing Core-Shell Technology

Column: Kinetex® Core-Shell C8 5 µm
Dimensions: 250 x 4.6 mm
Part No.: [00G-4608-E0](#)
Flow Rate: 1.0 mL/min
Sample: 1. Impurity A
 2. Impurity E
 3. Paroxetine
 4. Impurity C

Elution Time of Last Peak: 40 min
Rs Impurity E and Paroxetine: 5.80
Impurity E Peak Height: 0.28 mAU



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (gradient elution)	Method 1	Method 2
Mobile Phase pH	No adjustment permitted	As specified	As specified
Concentration of Salts in Buffer	No adjustment permitted	As specified in Monograph 2283 Details Table	As specified
Composition of the Mobile Phase	Minor adjustments of the composition of the mobile phase and the gradient are acceptable, if the system suitability requirements are met, the principle peak(s) elute(s) within ± 15 % of the indicated retention time(s) and the final elution power of the mobile phase is not weaker in elution power than the prescribed composition	As specified in Monograph 2283 Details Table	As specified
Wavelength of Detector	No deviations permitted	295 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified
Column Temperature	± 5 °C	40° C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	End-capped octylsilyl silica gel for chromatography (as specified)	As specified
Column Length	± 70 %	250 mm (as specified)	As specified
Column Internal Diameter	± 25 %	4.6 mm (as specified)	As specified
Particle Size	No adjustment permitted	5 µm (as specified)	As specified
Flow Rate	Adjustment is acceptable when changing the column dimensions	1.0 mL/min (as specified)	As specified

Potassium Clavulanate and Related Substances

Ph. Eur. monograph 1140

The Ph. Eur. Monograph 1140 outlines the separation of Potassium Clavulanate from Amoxicillin. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Potassium Clavulanate

Ph. Eur. Monograph 1140 Details

Reference Solution (b) Dissolve 10 mg of Lithium Clavulanate CRS* and 10 mg of Amoxicillin Trihydrate CRS* in mobile phase A and dilute to 100 mL with mobile phase A

Column

Size 100 x 4.6 mm

Stationary Phase Octadecylsilyl silica gel for chromatography R (5 µm)

Temperature 40 °C

Mobile Phase
A: 7.8 g/L solution of sodium hydrogen phosphate R adjusted to pH 4.0 with phosphoric acid R
B: Mixture of equal volumes of methanol R and mobile phase A

Gradient

Time (min)	%B
0 – 4	0
4 – 15	0 → 50
15 – 18	50

Flow Rate 1.0 mL/min


Detection Spectrophotometer @ 230 nm

Injection 20 µL


System Suitability

Reference Solution (b) Minimum resolution of 13 between peaks due to Clavulanate (1st peak) and Amoxicillin (2nd peak)

* Amoxicillin Trihydrate CRS (A0800000) and Lithium Clavulanate CRS (L0720000) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).



App ID: 21.472



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (gradient elution)	Method 1	Method 2
Mobile Phase pH	No adjustment permitted	4 (as specified)	As specified
Concentration of Salts in Buffer	No adjustment permitted	As specified in Monograph 1140 Details Table	As specified
Composition of the Mobile Phase	Minor adjustments of the composition of the mobile phase and the gradient are acceptable, if the system suitability requirements are met, the principle peak(s) elute(s) within $\pm 15\%$ of the indicated retention time(s) and the final elution power of the mobile phase is not weaker in elution power than the prescribed composition	As specified in Monograph 1140 Details Table	As specified
Wavelength of Detector	No deviations permitted	230 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified
Column Temperature	$\pm 5^\circ\text{C}$	40 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)	As specified
Column Length	May be decreased, $\pm 70\%$	100 mm (as specified)	As specified
Column Internal Diameter	$\pm 25\%$	4.6 mm (as specified)	As specified
Particle Size	No adjustment permitted	5 µm (as specified)	As specified
Flow Rate	Acceptable when changing the column dimensions	1 mL/min (as specified)	As specified

Pravastatin Sodium and Related Substances

Ph. Eur. monograph 2059

The Ph. Eur. Monograph 2059 outlines the separation of Pravastatin from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Pravastatin Sodium

Ph. Eur. Monograph 2059 Details

Solvent Mixture	Methanol R, water R (9:11 V/V)
Test Solution	(a) Dissolve 0.1000 g of Pravastatin 1,1,3,3-tetramethylbutylamine CRS* in the solvent mixture and dilute to 100.0 mL with the solvent mixture (b) Dilute 10.0 mL of the test solution (a) to 100.0 mL with the solvent mixture
Reference Solution (a)	Dissolve the contents of a vial of Pravastatin Impurity A CRS* in 1.0 mL of test solution (b)
Column	
Size	150 x 4.6 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R (5 µm)
Temperature	25 °C
Mobile Phase	Glacial acetic acid R, trimethylamine R, methanol R, water R (1:1:450:550 V/V/V/V)
Flow Rate	1.3 mL/min
Detection	Spectrophotometer @ 238 nm
Injection	10 µL
Run Time	2.5 times the retention time of Pravastatin
Elution Order	1. Impurity A 2. Pravastatin

System Suitability

Minimum resolution of 7.0 between peaks due to Impurity A and Pravastatin

* Pravastatin 1,1,3,3-tetramethylbutylamine CRS (Y0000204) and Pravastatin Impurity A CRS (Y0000223) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

Method 1

Original Method as Described in the Monograph

Column: Luna® C18(2) 5 µm Fully Porous

Dimensions: 150 x 4.6 mm


Part No.: 00F-4252-E0

Flow Rate: 1.3 mL/min


Sample: 1. Impurity A
2. Pravastatin

Elution Time of Last Peak: 31.80 min

Rs Impurity A and Pravastatin: 11.91



APP ID 21471



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 2059 Details Table	As specified
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 2059 Details Table	As specified
Wavelength of Detector	No deviations permitted	238 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	10 µL (as specified)	As specified
Column Temperature	± 10 °C	25 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)	As specified
Column Length	± 70 %	150 mm (as specified)	As specified
Column Internal Diameter	± 25 %	4.6 mm (as specified)	As specified
Particle Size	-50 %	5 µm (as specified)	As specified
Flow Rate	± 50 %	1.3 mL/min (as specified)	As specified

Simvastatin and Related Substances

Ph. Eur. monograph 1563

The Ph. Eur. Monograph 1563 outlines the separation of Simvastatin from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Ph. Eur. Monograph 1563 Details

Solvent Mixture Mix 40 volumes of a 1.4 g/L solution of potassium dihydrogen phosphate R, adjusted to pH 4.0 with phosphoric acid R, and 60 volumes of acetonitrile R. Filter.

Reference Solution
(a) Dissolve 1.0 mg of Simvastatin CRS* and 1.0 mg of Lovastatin CRS* (Impurity E) in the solvent mixture and dilute to 50.0 mL with the solvent mixture
(d) Dissolve 5 mg of Simvastatin for peak identification CRS* (containing Impurities A, B, C, D, E, F, and G) in 5 mL of the solvent mixture

Column

Size 33 x 4.6 mm

Stationary Phase End-capped octadecylsilyl silica gel for chromatography R (3 µm)

Temperature 25°C

Mobile Phase
 A: Mix 50 volumes of acetonitrile R and 50 volumes of a 0.1 % V/V solution of phosphoric acid R
 B: 0.1 % V/V solution of phosphoric acid R in acetonitrile R

Gradient	Time (min)	%B
	0 – 4.5	0
	4.5 – 4.6	0 → 5
	4.6 – 8	5 → 95
	8.0 – 11.5	75

Flow Rate 3 mL/min

Detection Spectrophotometer @ 238 nm

Injection 5 µL

Relative Retention with Reference to Simvastatin (about 2.6 min)**

Impurity A	about 0.5
Impurities E + F	about 0.6
Impurity G	about 0.8
Impurities B + C	about 2.4
Impurity D	about 3.8

System Suitability


Minimum resolution of 4.0 between peaks due to Impurity E and Simvastatin

* Simvastatin CRS (S0650000), Lovastatin CRS (impurity E) (L0790000) and Simvastatin for peak identification CRS* (containing Impurities A, B, C, D, E, F, and G) (Y0001066) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1

Original Method as Described in the Monograph




Method 2 Alternative Method Within Allowable Adjustments

Column: Gemini® C18 3µm Fully Porous
Dimensions: 30 x 4.6 mm
Part No.: [00A-4439-E0](#)
Flow Rate: 3.0 mL/min
Sample: 1. Impurity A
 2. Impurity E/F
 3. Impurity G
 4. Simvastatin
 5. Impurity B/C
 6. Impurity D



Method 3 Faster Method Outside Allowable Adjustments

Column: Kinetex® Core-Shell C18 2.6µm
Dimensions: 50 x 4.6 mm
Part No.: [00B-4462-E0](#)
Flow Rate: 3.0 mL/min
Sample: 1. Impurity A
 2. Impurity E
 3. Impurity F
 4. Impurity G
 5. Simvastatin
 6. Impurity B/C
 7. Impurity D




Adjustments for Meeting System Suitability (European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (gradient elution)	Method 1	Method 2	Method 3
Mobile Phase pH	No adjustment permitted	As specified	As specified	As specified
Concentration of Salts in Buffer	No adjustment permitted	As specified in Monograph 1563 Details Table	As specified	As specified
Composition of the Mobile Phase	Minor adjustments of the composition of the mobile phase and the gradient are acceptable, if the system suitability requirements are met, the principle peak(s) elute(s) within $\pm 15\%$ of the indicated retention time(s) and the final elution power of the mobile phase is not weaker in elution power than the prescribed composition	As specified in Monograph 1563 Details Table	As specified	As specified
Wavelength of Detector	No deviations permitted	238 nm (as specified)	As specified	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	5 µL (as specified)	As specified	As specified
Column Temperature	$\pm 5^\circ\text{C}$	Ambient (as specified)	As specified	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	End-capped octadecylsilyl silica gel for chromatography (as specified)	As specified	As specified
Column Length	$\pm 70\%$	33 mm (as specified)	30 mm (-9%)	50 mm (+51%)
Column Internal Diameter	$\pm 25\%$	4.6 mm (as specified)	As specified	As specified
Particle Size	No adjustment permitted	3 µm (as specified)	As specified	2.6 µm (outside of allowed adjustments)
Flow Rate	Adjustment is acceptable when changing the column dimensions	3.0 mL/min (as specified)	As specified	As specified

Tamsulosin Hydrochloride and Related Substances

Ph. Eur. monograph 2131

The Ph. Eur. Monograph 2131 outlines the separation of Tamsulosin from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Ph. Eur. Monograph 2131 Details- Tamsulosin (A)

Reference Solution

(b) Dissolve 4 mg of Tamsulosin Impurity D CRS* and 4 mg Tamsulosin Hydrochloride CRS* in the mobile phase and dilute to 20.0 mL with the mobile phase. Dilute 2.0 mL of this solution to 20.0 mL with the mobile phase.

(c) Dissolve 4 mg of Tamsulosin Impurity H CRS* and 4 mg Tamsulosin Hydrochloride CRS* in the mobile phase and dilute to 20.0 mL with the mobile phase. Dilute 2.0 mL of this solution to 20.0 mL with the mobile phase.

Column

Size	150 x 4.6 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R (5 µm).
Temperature	40 °C
Mobile Phase	Dissolve 3.0 g of sodium hydroxide R in a mixture of 8.7 mL of perchloric acid R and 1.9 L of water R; adjust to pH 2.0 with 0.5 M sodium hydroxide and dilute to 2 L with water R; to 1.4 L of this solution, add 600 mL of acetonitrile R.
Flow Rate	1.3 mL/min
Detection	Spectrophotometer @ 225 nm
Injection	10 µL
Run Time	1.5 times the retention of Tamsulosin (about 6 min)

System Suitability


Reference Solution (b) Minimum resolution of 6.0 between peaks due to Impurity D and Tamsulosin

* Tamsulosin impurity D CRS* (Y0000651), Tamsulosin Impurity H CRS (Y0000652) and Tamsulosin Hydrochloride CRS (Y0000650) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

Method 1

Original Method as Described in the Monograph

Column: Kinetex® Core-Shell C18 5 µm
Dimensions: 150 x 4.6 mm
Part No.: 00F-4601-E0
Flow Rate: 1.3 mL/min
Sample: 1. Impurity B
 2. Tamsulosin
Elution Time of Last Peak: 5.47 min
Rs Impurity D and Tamsulosin: 11.78



Adjustments for Meeting System Suitability

(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1
Mobile Phase pH	± 0.2 units	2 (as specified)
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 2131 Details Table
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 2131 Details Table
Wavelength of Detector	No deviations permitted	225 nm (as specified)
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	10 µL (as specified)
Column Temperature	± 10 %	40 °C (as specified)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	Octadecylsilyl silica gel for chromatography (as specified)
Column Length	± 70 %	150 mm (as specified)
Column Internal Diameter	± 25 %	4.6 mm (as specified)
Particle Size	-50 %	5 µm (as specified)
Flow Rate	± 50 %	1.3 mL/min (as specified)

Ph. Eur. Monograph 2131 Details- Tamsulosin (B)

Reference Solution (c) Dissolve 4 mg of Tamsulosin Impurity H CRS* and 4 mg Tamsulosin Hydrochloride CRS* in the mobile phase and dilute to 20.0 mL with the mobile phase. Dilute 2.0 mL of this solution to 20.0 mL with the mobile phase.

Column

Size	150 x 4.6 mm
Stationary Phase	Octadecylsilyl silica gel for chromatography R (5 µm).
Temperature	40 °C
Mobile Phase	Dissolve 3.0 g of sodium hydroxide R in a mixture of 8.7 mL of perchloric acid R and 1.9 L of water R; adjust to pH 2.0 with 0.5 M sodium hydroxide and dilute to 2 L with water R; add 2 L of acetonitrile R.
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 225 nm
Injection	10 µL
Run Time	5 times the retention of Tamsulosin (about 2.5 min)

System Suitability


Reference Solution (c) Minimum resolution of 2.0 between peaks due to Tamsulosin and Impurity H

* Tamsulosin Impurity D CRS* (Y0000651), Tamsulosin Impurity H CRS (Y0000652) and Tamsulosin Hydrochloride CRS (Y0000650) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allee Kastner CS 30026F - 67081 STRASBOURG (France).

Method 1

Original Method as Described in the Monograph

Column: Kinetex® Core-Shell C18 5 µm
Dimensions: 150 x 4.6 mm
Part No.: [00F-4601-E0](#)
Flow Rate: 1.0 mL/min
Sample: 1. Impurity D
 2. Tamsulosin
 3. Impurity H
Elution Time of Last Peak: 3.01 min
Rs Tamsulosin and Impurity H: 15.37



Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1
Mobile Phase pH	± 0.2 units	2 (as specified)
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 2131 Details Table
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 2131 Details Table
Wavelength of Detector	No deviations permitted	225 nm (as specified)
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	10 µL (as specified)
Column Temperature	± 10 °C	40 °C (as specified)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography
Column Length	± 70 %	150 mm (as specified)
Column Internal Diameter	± 25 %	4.6 mm (as specified)
Particle Size	-50 %	5 µm (as specified)
Flow Rate	± 50 %	1.0 mL/min (as specified)

Tramadol Hydrochloride and Related Substances

Ph. Eur. monograph 1681

The Ph. Eur. Monograph 1681 outlines the separation of Tramadol from impurities. This method was studied and recommendations have been made to conform with the Ph. Eur. Monograph 1681 requirements.



Ph. Eur. Monograph 1681 Details

Test Solution	Dissolve 0.15 g of Tramadol Hydrochloride CRS* in the mobile phase and dilute to 100 mL with the mobile phase.
Reference Solution (b)	Dissolve 5 mg of Tramadol Impurity A CRS* in 4.0 mL of the test solution and dilute to 100 mL with the mobile phase.
Column	
Size	250 x 4.0 mm
Stationary Phase	End-capped base-deactivated octylsilyl silica gel for chromatography R (5 µm)
Temperature	25 °C
Mobile Phase	295 volumes of acetonitrile R and 705 volumes of a mixture of 0.2 mL of trifluoroacetic acid R and 100 mL of water R
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 270 nm
Injection	20 µL
Run Time	4 times the retention time of Tramadol
Relative Retention with Reference to Tramadol (about 5 min)**	
Impurity A	about 0.85 min

System Suitability


Reference Solution (b)	Minimum resolution of 2.0 between peaks due to Impurity A and Tramadol
-------------------------------	--

* Tramadol Hydrochloride CRS (Y0000155) and Tramadol Impurity A CRS (Y0000156) were purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1

Improved Resolution Within Allowable Adjustments



App ID 24098

Adjustments for Meeting System Suitability


(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1
Mobile Phase pH	± 0.2 units	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 1681 Details Table
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 1681 Details Table
Wavelength of Detector	No deviations permitted	270 nm (as specified)
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)
Column Temperature	± 10 °C	Ambient (as specified)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	Octylsilyl silica gel for chromatography (as specified)
Column Length	± 70 %	250 mm (as specified)
Column Internal Diameter	± 25 %	4.6 mm (+15 %)
Particle Size	-50 %	5 µm (as specified)
Flow Rate	± 50 %	1.0 ml/min (as specified)

Trimethoprim and Related Substances

Ph. Eur. monograph 0060

The Ph. Eur. Monograph 0060 outlines the separation of Trimethoprim from impurities. This method was studied and improvements were made to provide faster separations within allowable adjustments.



Trimethoprim

Ph. Eur. Monograph 0060 Details

Reference Solution (b) Dissolve the contents of a vial of Trimethoprim for system suitability CRS* (containing Impurity E) in 1 mL of mobile phase.

Column

Size	250 x 4.0 mm
Stationary Phase	Base-deactivated octadecylsilyl silica gel for chromatography R (5 µm)
Temperature	25 °C
Mobile Phase	Mix 30 volumes of methanol R and 70 volumes of a 1.4 g/L solution of sodium perchlorate R adjusted to pH 3.6 with phosphoric acid R.
Flow Rate	1.3 mL/min
Detection	Spectrophotometer @ 280 nm
Injection	20 µL loop injector
Run Time	11 times the retention time of Trimethoprim

System Suitability

Reference Solution (b) Minimum resolution of 2.5 between peaks due to Impurity E and Trimethoprim

* Trimethoprim for system suitability CRS (containing Impurity E) (Y0000684) was purchased from European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: 7 Allée Kastner CS 30026F - 67081 STRASBOURG (France).


Method 1

Original Method as Described in the Monograph

Column: Luna® C18(2) 5 µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: [00G-4252-E0](#)
Flow Rate: 1.3 mL/min
Sample: 1. Impurity E
 2. Trimethoprim

Elution Time of Last Peak: 6.28 min

Rs Impurity E and Trimethoprim: 2.92




Method 2

Faster Method Within Allowable Adjustments

Column: Kinetex® Core-Shell C18 5 µm
Dimensions: 250 x 4.6 mm
Part No.: [00G-4601-E0](#)
Flow Rate: 1.3 mL/min
Sample: 1. Impurity E
 2. Trimethoprim

Elution Time of Last Peak: 4.06 min

Rs Impurity E and Trimethoprim: 3.85



Adjustments for Meeting System Suitability

(European Pharmacopoeia 9.0, Chapter 2.2.46. Chromatographic separation techniques)


Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	3.6 (as specified)	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph 0060 Details Table	As specified
Composition of the Mobile Phase	± 30 % of the minor solvent component relative or 2 % absolute, whichever is the larger. No other component is altered by more than 10 % absolute.	As specified in Monograph 0060 Details Table	As specified
Wavelength of Detector	No deviations permitted	280 nm (as specified)	As specified
Injection Volume	May be decreased, provided detection and repeatability of the peak(s) to be determined are satisfactory.	20 µL (as specified)	As specified
Column Temperature	± 10 °C	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	Octadecylsilyl silica gel for chromatography (as specified)	As specified
Column Length	± 70 %	250 mm (as specified)	As specified
Column Internal Diameter	± 25 %	4.6 mm (+15 %)	4.6 mm (+15 %)
Particle Size	-50 %	5 µm (as specified)	As specified
Flow Rate	± 50 %	1.3 mL/min (as specified)	As specified



View our Webinar

**Understand How to Apply USP
Chapter <621> Allowable
Adjustments to Your USP
Pharmacopeia Method**

www.phenomenex.com/Webinars




**United States
Pharmacopeia
(USP)
Monographs
for Generic Drugs**

Amlodipine Besylate

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Amlodipine Besylate. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Amlodipine Besylate

USP Monograph: Amlodipine Besylate Details

pH 3.0 Buffer	Dissolve 7.0 of triethylamine in 800mL of water. Adjust with phosphoric acid to a pH of 3.0± 0.1, and dilute with water to 1 L.
System Suitability Solution	Dissolve about 5 mg of Amlodipine Besylate in 5 mL of hydrogen peroxide, and heat at 70° C for 45 minutes
Standard Preparation	Dissolve USP Amlodipine Besylate RS in mobile phase to obtain a concentration of 0.003 mg/mL
Test Solution	Dissolve 50 mg of Amlodipine Besylate in a 50 mL volumetric flask and dilute to volume with mobile phase

Column

Size	150 x 3.9 mm
Stationary Phase	L1: Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod
Mobile Phase	pH 3.0 Buffer, Methanol and acetonitrile (50:35:15)
Flow Rate	1.0 mL/min
Detection	Spectrophotometer @ 237 nm
Injection	10 µL

Relative Retention with Reference to Amlodipine*

Benzene Sulfonate	about 0.2
Impurity A	about 0.5

System Suitability

Minimum resolution of 4.5 between Amlodipine and Impurity A

* Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1

Standard Method Within Allowable Adjustments

Column: Luna® C18(2) 5µm Fully Porous


Dimensions: 150 x 4.6 mm

Part No.: [00F-4252-E0](#)

Flow Rate: 1.0 mL/min


Sample: 1. Impurity A
2. Amlodipine

Rs Impurity A and Amlodipine: 17.7



App ID: 24108

Reduce run times by >50% with
Kinetex Core-Shell Columns



Method 2

Faster and Higher Resolution Within Allowable Adjustments

Column: Kinetex® Core-Shell C18 5 μ m
Dimensions: 150 x 4.6 mm
Part No.: 00F-4601-E0
Flow Rate: 1.0 mL/min
Sample: 1. Impurity A
 2. Amlodipine

Rs Impurity A and Amlodipine: 18.4

App ID: 2410Z

Adjustments for Meeting System Suitability


Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	$\pm 10\%$	As specified	As specified
Composition of the Mobile Phase	$\pm 30\%$ Relative; cannot exceed $\pm 10\%$ Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	237 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	10 μ L (as specified)	As specified
Column Temperature	$\pm 10^\circ\text{C}$	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	150 mm (as specified)	As specified
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (+18%)	4.6 mm (+18%)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	5 μ m (as specified)	As specified
Flow Rate	$\pm 50\%$ (at given ID)	1.0 mL/min (as specified)	As specified

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25% to +50%

Clavulanate Potassium and Related Substances

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Clavulanate Potassium. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



USP Monograph: Clavulanate Potassium Details


Solution A	7.8 mg/mL of monobasic sodium phosphate in water. Adjust with phosphoric acid or 10 N sodium hydroxide to a pH of 4.4 ± 0.1 before final dilution.
Standard Solution	0.25 mg/mL of USP Clavulanate Lithium RS in water
System Suitability Solution	0.5 mg/mL of Amoxicillin dissolved in Standard Solution
Sample Solution	0.25 mg/mL of Clavulanate Potassium in water


Column

Size	30 x 4.0 mm
Stationary Phase	L1: Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod
Mobile Phase	Methanol and Solution A (1:19)
Flow Rate	2.0 mL/min
Detection	Spectrophotometer @ 220 nm
Injection	20 µL

System Suitability

Minimum resolution of 3.5 between Amoxicillin and Clavulanic Acid





Adjustments for Meeting System Suitability


Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	220 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	20 µL (as specified)	As specified
Column Temperature	± 10 °C	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	250 mm (-17 %)	250 mm (-17 %)
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (+15 %)	4.6 mm (+15 %)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	5 µm (as specified)	As specified
Flow Rate	± 50 % (at given ID)	2.0 mL/min (as specified)	As specified

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25 % to +50 %

Fluticasone Propionate and Related Substances

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Fluticasone Propionate. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Fluticasone Propionate

USP Monograph: Fluticasone Propionate Details

System Suitability Solution	Dissolve 2.0 mg of USP Fluticasone Propionate	
System Suitability Mixture	RS in 5 mL of Solution A using sonication. Add 5 mL of Solution C.	
Sample Solution	Dissolve 2.0 mg of Fluticasone Propionate in 5 mL of Solution A using sonication. Add 5 mL of Solution C.	
Column		
Size	250 x 4.6 mm	
Stationary Phase	5 µm, L1: Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod	
Mobile Phase	A: 0.5 mL of phosphoric acid in 1000 mL of acetonitrile B: 0.5 mL of phosphoric acid in 1000 mL of methanol C: 0.5 mL of phosphoric acid in 1000 mL of water	
Gradient	Time (min):	% (A/B/C)
	0	42/3/55
	40	53/3/44
	60	87/3/10
	70	87/3/10
	75	42/3/55
Flow Rate	1.0 mL/min	
Detection	Spectrophotometer @ 239 nm	
Injection	50 µL	

Relative Retention with Reference to Fluticasone Propionate*

Related Compound A	about 0.5
Related Compound B	about 0.75
Related Compound C	about 0.8
Related Compound D	about 0.95
Related Compound E	about 1.3

System Suitability

Minimum resolution of 0.6 between Related Compound B and Related Compound C.
 Minimum resolution of 1.5 between Related Compound D and Fluticasone Propionate.


* Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.

Method 1


Original Method as Described in the USP Monograph

- Column:** Luna® C18(2) 5 µm Fully Porous
- Dimensions:** 250 x 4.6 mm
- Part No.:** [00G-4252-E0](#)
- Flow Rate:** 1.0 mL/min
- Sample:**
 1. Related Compound A
 2. Related Compound B
 3. Related Compound C
 4. Related Compound D
 5. Fluticasone Propionate
 6. Related Compound E

- Elution Time of Last Peak:** 44.8 min
- Rs Related Compound B and Related Compound C:** 1.59
- Rs Related Compound D and Fluticasone Propionate:** 2.9



App ID: 24115




Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (gradient elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	Changes to gradient composition are not recommended	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	239 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	50 µL (as specified)	As specified
Column Temperature	± 10 °C	40 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	No deviations permitted	250 mm (as specified)	As specified
Column Internal Diameter	No deviations permitted	4.6 mm (as specified)	As specified
Particle Size	No deviations permitted	5 µm (as specified)	As specified
Flow Rate	No deviations permitted	1.0 mL/min (as specified)	As specified

Ibuprofen

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Ibuprofen. This method was studied and improvements were made to provide higher resolution (R_s) and a faster separation time within allowable adjustments.



USP Monograph: Ibuprofen Details

Resolution Solution	Prepare a solution in acetonitrile containing in each mL about 5 mg of Ibuprofen and 5 mg of Valerophenone
Test Preparation	Prepare a solution of Ibuprofen in acetonitrile containing about 5 mg per mL
Column	
Size	150 x 4.0 mm
Stationary Phase	5 μ m, L1: Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 μ m in diameter, or a monolithic rod
Temperature	30 °C \pm 0.5 °C
Mobile Phase	Prepare a suitable filtered mixture of water, previously adjusted with phosphoric acid to pH 2.5 and acetonitrile (1340:680).
Flow Rate	2.0 mL/min
Detection	Spectrophotometer @ 214 nm
Injection	5 μ L

Relative Retention with Reference to Ibuprofen*

Valerophenone	about 0.8
----------------------	-----------

System Suitability

Minimum resolution of 2.0 between Valerophenone and Ibuprofen

* Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1


Original Method Within Allowable Adjustments

Column: Luna[®] C18(2) 5 μ m Fully Porous
Dimensions: 150 x 4.6 mm
Part No.: [00F-4252-E0](#)
Flow Rate: 1.0 mL/min
Sample: 1. Valerophenone
 2. Ibuprofen

Elution Time of Last Peak: 3.2 min

R_s Valerophenone and Ibuprofen: 5.88





Adjustments for Meeting System Suitability


Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	214 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	5 µL (as specified)	As specified
Column Temperature	± 10 °C	30 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	150 mm (as specified)	As specified
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (+15 %)	4.6 mm (+15 %)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	5 µm (as specified)	As specified
Flow Rate	± 50 % (at given ID)	1.0 mL/min (-50 %)	1.0 mL/min (-50 %)

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25 % to +50 %

Lovastatin

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Lovastatin. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Lovastatin

USP Monograph: Lovastatin Details

System Suitability Solution	Dissolve USP Lovastatin RS and USP Lovastatin Related Compound A RS in acetonitrile to obtain a concentration of 2.0 µg/mL of each
Standard Solution	Dissolve USP Lovastatin RS in acetonitrile to obtain a concentration of about 2.0 µg/mL
Test Solution	Dissolve 25 mg of Lovastatin in a 25 mL volumetric flask and dilute to volume with acetonitrile, mix

Column

Size	250 x 4.6 mm
Stationary Phase	5 µm, L7: Octyl silane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod
Temperature	40 °C
Mobile Phase	Acetonitrile and 0.01 M Phosphoric acid (13:7)
Flow Rate	1.5 mL/min
Detection	Spectrophotometer @ 200 nm
Injection	10 µL

Relative Retention with Reference to Lovastatin*

Related Compound A	about 1.3
---------------------------	-----------

System Suitability

Minimum resolution of 6.0 between Lovastatin and Related Compound A

* Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1

Original Method within as Described in the USP Monograph

Column: Luna® C8(2) 5 µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: [00G-4249-E0](#)
Flow Rate: 1.5 mL/min
Sample: 1. Lovastatin
 2. Related Compound A

Elution Time of Last Peak: 10.9 min

Rs Lovastatin and Related Compound A: 12.33



App ID 24118

Method 2


Faster and Higher Resolution Within Allowable Adjustments

Column: Kinetex® Core-Shell C8 5µm
Dimensions: 250 x 4.6mm
Part No.: [00G-4608-E0](#)
Flow Rate: 1.5 mL/min
Sample: 1. Lovastatin
 2. Related Compound A

Elution Time of Last Peak: 5.6 min

Rs Lovastatin and Related Compound A: 9.92

Run time reduced by ~5 minutes with Kinetex Core-Shell Columns!



Adjustments for Meeting System Suitability


Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	200 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	10 µL (as specified)	As specified
Column Temperature	± 10 °C	40 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	L7 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	250 mm (as specified)	As specified
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (as specified)	As specified
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	5 µm (as specified)	As specified
Flow Rate	± 50 % (at given ID)	1.5 mL/min (as specified)	As specified

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25 % to +50 %

Metformin Hydrochloride

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Metformin Hydrochloride. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Metformin Hydrochloride

USP Monograph: Metformin Hydrochloride Details


System Suitability Stock Solution	0.25 mg/mL of Metformin Hydrochloride and 0.1 mg/mL of Melamine in water
System Suitability Solution	Transfer 1.0 mL of system suitability stock solution to a 50 mL volumetric flask, dilute with mobile phase to volume
Standard Stock Solution	0.2 mg/mL of USP Metformin Related Compound A RS in water
Standard Solution	0.001 mg/mL of USP Metformin Related Compound A RS in Mobile Phase from Standard stock solution
Sample Solution	5 mg/mL of Metformin Hydrochloride in mobile phase
Diluted Sample Solution	0.005 mg/mL of Metformin Hydrochloride in Mobile Phase from the Sample solution
Column	
Size	250 x 4.6 mm
Stationary Phase	L9: Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter
Mobile Phase	17 g/L of monobasic ammonium phosphate in water, adjusted with phosphoric acid to a pH 3.0
Flow Rate	1.0 - 1.7 mL/min
Detection	Spectrophotometer @ 218 nm
Injection	20 µL
Run Time	Not less than twice the retention time of Metformin
System Suitability	
Minimum resolution of 10 between Melamine and Metformin	

Method 1

Original Method as Described in the USP Monograph


Column: Luna® SCX 10µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: 00G-4401-E0
Flow Rate: 1.7 mL/min
Sample: 1. Melamine
 2. Metformin

Elution Time of Last Peak: 20.6 min
Rs Melamine and Metformin: 27.28



App ID 23671

Method 2 Faster Method Within Allowable Adjustments



Adjustments for Meeting System Suitability


Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	$\pm 10\%$	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	$\pm 30\%$ Relative; cannot exceed $\pm 10\%$ Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	218 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	20 μ L (as specified)	As specified
Column Temperature	$\pm 10^\circ\text{C}$	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L9 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	250 mm (as specified)	150 mm (-40%)
Column Internal Diameter	Can be adjusted so long as linear velocity if maintained	4.6 mm (as specified)	As specified
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	10 μ m (as specified)	5 μ m (as specified)
Flow Rate	$\pm 50\%$ (at given ID)	1.7 mL/min (as specified)	2.0 mL/min (+18)

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25% to +50%

Pravastatin Sodium

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Pravastatin Sodium. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



USP Monograph: Pravastatin Sodium Details

Diluent	Prepare a mixture of methanol and water (1:1)
Buffer pH 7.0	Prepare a 0.08 M phosphoric acid solution, adjust with triethylamine to pH 7.0, mix
Standard Solution*	Dissolve an accurately weighed quantity of USP Pravastatin 1,1,3,3-Tetramethylbutylamine RS in Diluent, and dilute quantitatively with Diluent to obtain a solution having a known concentration of about 1.25 µg of pravastatin 1,1,3,3-tetramethylbutylamine per mL
System Suitability Solution	Dissolve accurately weighed quantities of USP Pravastatin 1,1,3,3-Tetramethylbutylamine Rs and USP Pravastatin Related Compound A RS in Diluent to obtain a solution containing about 0.6 mg of USP Pravastatin 1,1,3,3 tetramethylbutylamine RS and 0.001 mg of USP Pravastatin Related Compound A RS per mL. (Note-USP Pravastatin Related Compound A RS is a sodium salt of 3α-hydroxisocompactin acid)
Test Solution*	Transfer about 50 mg of Pravastatin Sodium to a 100 mL volumetric flask, dissolve in and dilute with Diluent to volume, and mix

Column

Size	100 x 4.0 mm	
Stationary Phase	3 µm, L1: Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod	
Mobile Phase	Use variable mixtures of Solution A and Solution B as directed for: A. Prepare a filtered and degassed mixture of water, Buffer pH 7.0, and acetonitrile (52:30:10) B. Prepare a filtered and degassed mixture of acetonitrile, Buffer pH 7.0, and water (60:30:10)	
Gradient	Time	%B
	0 – 3.0 min	0
	3.0 – 26.5 min	0 → 100
	26.5 – 26.6 min	100 → 0
	26.6 – 30.0 min	0
Flow Rate	1.0 mL/min	
Detection	Spectrophotometer @ 238 nm	
Injection	10 µL	

Relative Retention with Reference to Pravastatin**

Related Compound A	about 1.1
---------------------------	-----------

System Suitability

Minimum resolution of 2.0 between Pravastatin and Pravastatin Related Compound A

*The Standard solution and the Test solution are maintained at 15° C until injected into the chromatograph

** Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1

Original Method Outside Allowable Adjustments

Column: Luna® C18(2) 3µm Fully Porous
Dimensions: 100 x 4.6mm
Part No.: [00D-4251-E0](#)
Flow Rate: 1.0mL/min
Sample: 1. Pravastatin
 2. Pravastatin Compound A

Elution Time of Last Peak: 11.5


Rs Pravastatin and Related Compound A: 13.74



App. ID: 241106

Method 2

Faster and Higher Resolution Outside Allowable Adjustments



Adjustments for Meeting System Suitability


Method Parameter	Allowed Adjustments (gradient elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified	As specified
Composition of the Mobile Phase	Changes to gradient composition are not recommended	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	238 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	10 μL (as specified)	As specified
Column Temperature	± 10 °C	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	No deviations permitted	100 mm (as specified)	As specified
Column Internal Diameter	No deviations permitted	4.6 mm (+15 %)	4.6 mm (+15)
Particle Size	No deviations permitted	3 μm (as specified)	2.6 μm (-13 %)
Flow Rate	No deviations permitted	1.0 mL/min (as specified)	As specified

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25 % to +50 %

Trimethoprim

USP

The related substances test of the USP monograph outlines the separation of all relevant impurities from Trimethoprim. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Trimethoprim

USP Monograph: Trimethoprim Details

Buffer Solution	Prepare a 10mM sodium perchlorate solution in water, adjust with phosphoric acid to pH 3.6, and mix
Resolution Solution	Dissolve accurately weighed quantities of USP Trimethoprim RS and Diaveridine; and dilute quantitatively with mobile phase to obtain a solution having known concentrations of about 10 µg per mL and 5 µg per mL, respectively
Test Solution	Transfer about 25.0mg of Trimethoprim to a 25 mL volumetric flask, dissolve in and dilute with mobile phase to volume, and mix

Column

Size	250 x 4.6 mm
Stationary Phase	L1: Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod
Mobile Phase	Prepare a filtered and degassed mixture of Buffer Solution and methanol (7:3)
Flow Rate	1.3 mL/min
Detection	Spectrophotometer @ 280 nm
Injection	20 µL

System Suitability


Minimum resolution of 2.5 between Trimethoprim and Diaveridine
 Relative standard deviation for replicate injections is not more than 2.0 %


Method 1

Original Method as Described in the USP Monograph

Column: Luna® C18(2) 5 µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: 00G-4252-E0
Flow Rate: 1.3 mL/min
Sample: 1. Trimethoprim
 2. Diaveridine

Elution Time of Last Peak: 5.5
Rs Trimethoprim and Diaveridine: 3.1





Method 2

Faster and Higher Resolution Utilizing Core-Shell Technology

Column: Kinetex® Core-Shell C18 5 μ m

Dimensions: 250 x 4.6 mm

Part No.: 00G-4601-E0

Flow Rate: 1.3 mL/min

Sample: 1. Trimethoprim

2. Diaveridine

Elution Time of Last Peak: 4.0 min

Rs Trimethoprim and Diaveridine: 3.2

Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	$\pm 10\%$	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	$\pm 30\%$ Relative; cannot exceed $\pm 10\%$ Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	280 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	20 μ L (as specified)	As specified
Column Temperature	$\pm 10^\circ\text{C}$	Ambient (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	250 mm (as specified)	As specified
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (as specified)	As specified
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	5 μ m (as specified)	As specified
Flow Rate	$\pm 50\%$ (at given ID)	1.3 mL/min (as specified)	As specified

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25% to +50%

Kinetex Ordering Information



5 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJ0-9298
F5	00A-4724-AN	00B-4724-AN	00D-4724-AN	00F-4724-AN	AJ0-9322
Biphenyl	00A-4627-AN	00B-4627-AN	00D-4627-AN	—	AJ0-9209
XB-C18	00A-4605-AN	00B-4605-AN	00D-4605-AN	—	AJ0-8782
C18	00A-4601-AN	00B-4601-AN	00D-4601-AN	00F-4601-AN	AJ0-8782
C8	—	00B-4608-AN	00D-4608-AN	—	AJ0-8784
Phenyl-Hexyl	—	00B-4603-AN	—	—	AJ0-8788

for 2.1 mm ID

5 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges [†]
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	00B-4633-YO	00D-4633-YO	00F-4633-YO	AJ0-9297
F5	00B-4724-YO	00D-4724-YO	00F-4724-YO	AJ0-9321
Biphenyl	00B-4627-YO	00D-4627-YO	00F-4627-YO	AJ0-9208
XB-C18	00B-4605-YO	00D-4605-YO	00F-4605-YO	AJ0-8775
C18	00B-4601-YO	00D-4601-YO	00F-4601-YO	AJ0-8775
C8	00B-4608-YO	00D-4608-YO	—	AJ0-8777
Phenyl-Hexyl	00B-4603-YO	00D-4603-YO	—	AJ0-8781

for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges [†]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJ0-9296
F5	00B-4724-E0	00D-4724-E0	00F-4724-E0	00G-4724-E0	AJ0-9320
Biphenyl	00B-4627-E0	00D-4627-E0	00F-4627-E0	00G-4627-E0	AJ0-9207
XB-C18	00B-4605-E0	00D-4605-E0	00F-4605-E0	00G-4605-E0	AJ0-8768
C18	00B-4601-E0	00D-4601-E0	00F-4601-E0	00G-4601-E0	AJ0-8768
C8	00B-4608-E0	00D-4608-E0	00F-4608-E0	00G-4608-E0	AJ0-8770
Phenyl-Hexyl	00B-4603-E0	00D-4603-E0	00F-4603-E0	00G-4603-E0	AJ0-8774

for 4.6 mm ID

5 µm Semi-Preparative Columns (mm)			SecurityGuard SemiPrep Cartridges ^{***}
Phases	150 x 10	250 x 10	3/pk
EVO C18	00F-4633-N0	00G-4633-N0	AJ0-9306
F5	—	00G-4724-N0	AJ0-9323
C18	00F-4601-N0	00G-4601-N0	AJ0-9278
Biphenyl	00F-4627-N0	00G-4627-N0	AJ0-9280

for 9-16 mm ID

3.5 µm Analytical Columns (mm)			SecurityGuard ULTRA Cartridges [†]
Phases	100 x 4.6	150 x 4.6	3/pk
XB-C18	00D-4744-E0	00F-4744-E0	AJ0-8768

for 4.6 mm ID

2.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
XB-C18	00B-4496-A0	00D-4496-A0	00F-4496-A0



Kinetex has earned a Gold Seal of Quality!

Learn more at www.phenomenex.com/Gold

[†]SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

^{***}SemiPrep SecurityGuard Cartridges require holder, Part No.: [AJ0-9281](#)



Kinetex Ordering Information (cont'd)



2.6 µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJ0-9298
Polar C18	00A-4759-AN	00B-4759-AN	—	00D-4759-AN	00F-4759-AN	AJ0-9530
F5	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJ0-9322
Biphenyl	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJ0-9209
XB-C18	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJ0-8782
C18	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJ0-8782
C8	00A-4497-AN	00B-4497-AN	00C-4497-AN	00D-4497-AN	00F-4497-AN	AJ0-8784
HILIC	00A-4461-AN	00B-4461-AN	00C-4461-AN	00D-4461-AN	00F-4461-AN	AJ0-8786
Phenyl-Hexyl	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJ0-8788

for 2.1 mm ID

2.6 µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	—	00B-4725-Y0	—	00D-4725-Y0	00F-4725-Y0	AJ0-9297
Polar C18	—	00B-4759-Y0	—	00D-4759-Y0	00F-4759-Y0	AJ0-9531
F5	—	00B-4723-Y0	—	00D-4723-Y0	00F-4723-Y0	AJ0-9321
Biphenyl	—	00B-4622-Y0	—	00D-4622-Y0	00F-4622-Y0	AJ0-9208
XB-C18	00A-4496-Y0	00B-4496-Y0	00C-4496-Y0	00D-4496-Y0	00F-4496-Y0	AJ0-8775
C18	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJ0-8775
C8	00A-4497-Y0	00B-4497-Y0	00C-4497-Y0	00D-4497-Y0	00F-4497-Y0	AJ0-8777
HILIC	00A-4461-Y0	—	—	—	00F-4461-Y0	AJ0-8779
Phenyl-Hexyl	—	00B-4495-Y0	—	00D-4495-Y0	00F-4495-Y0	AJ0-8781

for 3.0 mm ID

2.6 µm Analytical Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
EVO C18	—	00B-4725-E0	—	00D-4725-E0	00F-4725-E0	AJ0-9296
Polar C18	—	00B-4759-E0	—	00D-4759-E0	00F-4759-E0	AJ0-9532
F5	—	00B-4723-E0	—	00D-4723-E0	00F-4723-E0	AJ0-9320
Biphenyl	—	00B-4622-E0	—	00D-4622-E0	00F-4622-E0	AJ0-9207
XB-C18	—	00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	AJ0-8768
C18	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJ0-8768
C8	—	00B-4497-E0	00C-4497-E0	00D-4497-E0	00F-4497-E0	AJ0-8770
HILIC	—	00B-4461-E0	00C-4461-E0	00D-4461-E0	00F-4461-E0	AJ0-8772
Phenyl-Hexyl	—	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0	AJ0-8774

for 4.6 mm ID

1.7 µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk	
EVO C18	—	00B-4726-AN	00D-4726-AN	00F-4726-AN	AJ0-9298	
Biphenyl	—	00B-4628-AN	00D-4628-AN	00F-4628-AN	AJ0-9209	
XB-C18	00A-4498-AN	00B-4498-AN	00D-4498-AN	00F-4498-AN	AJ0-8782	
C18	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJ0-8782	
C8	00A-4499-AN	00B-4499-AN	00D-4499-AN	00F-4499-AN	AJ0-8784	
HILIC	00A-4474-AN	00B-4474-AN	00D-4474-AN	—	AJ0-8786	
Phenyl-Hexyl	—	00B-4500-AN	00D-4500-AN	00F-4500-AN	AJ0-8788	
F5	—	00B-4722-AN	00D-4722-AN	00F-4722-AN	AJ0-9322	

for 2.1 mm ID

1.7 µm MidBore™ Columns (mm)					SecurityGuard ULTRA Cartridges [†]
Phases	30 x 3.0	50 x 3.0	100 x 3.0	3/pk	
XB-C18	00A-4498-Y0	00B-4498-Y0	00D-4498-Y0	AJ0-8775	
C18	—	00B-4475-Y0	00D-4475-Y0	AJ0-8775	
C8	00A-4499-Y0	00B-4499-Y0	00D-4499-Y0	AJ0-8777	
HILIC	—	00B-4474-Y0	—	AJ0-8779	

for 3.0 mm ID

1.7 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
C18	00B-4726-AN	00D-4726-AN	00F-4726-AN

1.3 µm Minibore Columns (mm)		
Phases	30 x 2.1	50 x 2.1
C18	00A-4515-AN	00B-4515-AN

[†]SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

Luna Ordering Information (cont'd)



5 µm Axia™ Packed Preparative Columns (mm)								SecurityGuard™ Cartridges (mm)	
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	50 x 30	100 x 30	250 x 30	15 x 21.2**	15 x 30 *
								/ea	/ea
Silica(2)	—	00D-4274-PO-AX	00F-4274-PO-AX	00G-4274-PO-AX	—	—	00G-4274-UO-AX	AJ0-7229	AJ0-8312
C5	—	—	—	00G-4043-PO-AX	—	—	—	—	—
C8(2)	—	—	00F-4249-PO-AX	00G-4249-PO-AX	—	00D-4249-UO-AX	—	AJ0-7840	AJ0-8302
C18(2)	00B-4252-PO-AX	00D-4252-PO-AX	00F-4252-PO-AX	00G-4252-PO-AX	00B-4252-UO-AX	00D-4252-UO-AX	00G-4252-UO-AX	AJ0-7839	AJ0-8301
CN	—	—	—	00G-4255-PO-AX	—	—	00G-4255-UO-AX	AJ0-8220	AJ0-8311
Phenyl-Hexyl	—	—	00F-4257-PO-AX	00G-4257-PO-AX	—	—	00G-4257-UO-AX	AJ0-7841	AJ0-8303
NH ₂	—	—	00F-4378-PO-AX	00G-4378-PO-AX	—	—	—	AJ0-8162	AJ0-8309
PFP(2)	—	00D-4448-PO-AX	00F-4448-PO-AX	00G-4448-PO-AX	—	00D-4448-UO-AX	—	AJ0-8377	AJ0-8378
HILIC	—	00D-4450-PO-AX	00F-4450-PO-AX	00G-4450-PO-AX	—	—	00G-4450-UO-AX	AJ0-8829	AJ0-8830

for ID: 18-29 mm 30-49 mm

10 µm Axia™ Packed Preparative Columns (mm) (continued)						SecurityGuard Cartridges (mm)	
Phases	50 x 21.2	100 x 21.2	250 x 21.2	250 x 30	250 x 50	15 x 21.2**	15 x 30 *
						/ea	/ea
Silica(2)	—	—	00G-4091-PO-AX	00G-4091-UO-AX	00G-4091-V0-AX	AJ0-7229	AJ0-8312
C5	—	00D-4092-PO-AX	00G-4092-PO-AX	—	00G-4092-V0-AX	—	—
C8(2)	—	—	00G-4250-PO-AX	—	00G-4250-V0-AX	AJ0-7840	AJ0-8302
C18(2)	00B-4253-PO-AX	00D-4253-PO-AX	00G-4253-PO-AX	00G-4253-UO-AX	00G-4253-V0-AX	AJ0-7839	AJ0-8301
CN	—	—	00G-4300-PO-AX	—	—	AJ0-8220	AJ0-8311
Phenyl-Hexyl	—	—	00G-4285-PO-AX	00G-4285-UO-AX	—	AJ0-7841	AJ0-8303
NH ₂	—	—	00G-4379-PO-AX	—	—	AJ0-8162	AJ0-8309

for ID: 18-29 mm 30-49 mm

10 µm Analytical and Semi-Prep Columns (mm)			SecurityGuard Cartridges (mm)	
Phases	250 x 4.6	250 x 10	4 x 3.0*	10 x 10†
			/10 pk	/3 pk
Silica(2)	00G-4091-E0	00G-4091-N0	AJ0-4348	AJ0-7223
C8(2)	00G-4250-E0	00G-4250-N0	AJ0-4290	AJ0-7222
C18(2)	00G-4253-E0	00G-4253-N0	AJ0-4287	AJ0-7221
CN	00G-4300-E0	—	AJ0-4305	AJ0-7313
Phenyl-Hexyl	00G-4285-E0	00G-4285-N0	AJ0-4351	AJ0-7314
NH ₂	00G-4379-E0	00G-4379-N0	AJ0-4302	AJ0-7364
SCX	00G-4401-E0	00G-4401-N0	AJ0-4308	AJ0-7369

for ID: 3.2-8.0 mm 9-16 mm

15 µm Pilot Scale Columns (mm)	
Phases	250 x 4.6
C18(2)	00G-4273-E0
Phenyl-Hexyl	00G-4286-E0



*SecurityGuard™ Analytical Cartridges require holder, Part No.: [KJ0-4282](#)
†SemiPrep SecurityGuard Cartridges require holder, Part No.: [AJ0-9281](#)
**PREP SecurityGuard Cartridges require holder, Part No.: [AJ0-8223](#)
◆PREP SecurityGuard Cartridges require holder, Part No.: [AJ0-8277](#)

Gemini Ordering Information



3 µm Microbore, Minibore and MidBore™ Columns (mm)										SecurityGuard™ Cartridges (mm)	
Phases	50 x 1.0	20 x 2.0	30 x 2.0	50 x 2.0	100 x 2.0	150 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*	
C18	00B-4439-A0	00M-4439-B0	00A-4439-B0	00B-4439-B0	00D-4439-B0	00F-4439-B0	00B-4439-Y0	00D-4439-Y0	00F-4439-Y0	/10pk	AJ0-7596
C6-Phenyl	00B-4443-A0	—	00A-4443-B0	00B-4443-B0	00D-4443-B0	00F-4443-B0	00B-4443-Y0	00D-4443-Y0	00F-4443-Y0	/10pk	AJ0-7914
NX-C18	00B-4453-A0	00M-4453-B0	00A-4453-B0	00B-4453-B0	00D-4453-B0	00F-4453-B0	00B-4453-Y0	00D-4453-Y0	00F-4453-Y0	/10pk	AJ0-8367

for ID: 2.0-3.0 mm

3 µm Analytical Columns (mm)					SecurityGuard™ Cartridges (mm)	
Phases	30 x 4.6	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
C18	00A-4439-E0	00B-4439-E0	00D-4439-E0	00F-4439-E0	00G-4439-E0	AJ0-7597
C6-Phenyl	00A-4443-E0	00B-4443-E0	00D-4443-E0	00F-4443-E0	00G-4443-E0	AJ0-7915
NX-C18	—	00B-4453-E0	00D-4453-E0	00F-4453-E0	00G-4453-E0	AJ0-8368

for ID: 3.2-8.0 mm

5 µm Minibore and MidBore Columns (mm)										SecurityGuard™ Cartridges (mm)	
Phases	30 x 2.0	50 x 2.0	150 x 2.0	250 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	250 x 3.0	4 x 2.0*		
C18	00A-4435-B0	00B-4435-B0	00F-4435-B0	00G-4435-B0	00B-4435-Y0	00D-4435-Y0	00F-4435-Y0	00G-4435-Y0	/10pk	AJ0-7596	
C6-Phenyl	—	00B-4444-B0	00F-4444-B0	—	00B-4444-Y0	—	00F-4444-Y0	00G-4444-Y0	/10pk	AJ0-7914	
NX-C18	00A-4454-B0	00B-4454-B0	00F-4454-B0	—	00B-4454-Y0	00D-4454-Y0	00F-4454-Y0	00G-4454-Y0	/10pk	AJ0-8367	

for ID: 2.0-3.0 mm

5 µm Analytical Columns (mm)					SecurityGuard™ Cartridges (mm)	
Phases	30 x 4.6	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
C18	00A-4435-E0	00B-4435-E0	00D-4435-E0	00F-4435-E0	00G-4435-E0	AJ0-7597
C6-Phenyl	—	00B-4444-E0	00D-4444-E0	00F-4444-E0	00G-4444-E0	AJ0-7915
NX-C18	—	00B-4454-E0	00D-4454-E0	00F-4454-E0	00G-4454-E0	AJ0-8368

for ID: 3.2-8.0 mm

5 µm Semi-Prep Columns (mm)			SecurityGuard™ Cartridges (mm)	
Phases	150 x 10	250 x 10	10 x 10 [†]	
C18	00F-4435-N0	00G-4435-N0	AJ0-7598	
C6-Phenyl	—	00G-4444-N0	AJ0-9156	
NX-C18	00F-4454-N0	00G-4454-N0	AJ0-8369	

for ID: 9-16 mm

*SecurityGuard™ Analytical Cartridges require holder, Part No.: [KJ0-4282](#)
[†]SemiPrep SecurityGuard™ Cartridges require holder, Part No.: [AJ0-9281](#)
**PREP SecurityGuard™ Cartridges require holder, Part No.: [AJ0-8223](#)
◆PREP SecurityGuard™ Cartridges require holder, Part No.: [AJ0-8277](#)

Axia™ Packed Preparative Columns (mm)							SecurityGuard™ Cartridges (mm)	
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	50 x 30	75 x 30	15 x 21.2**	15 x 30.0*
5 µm							/ea	/ea
C18	00B-4435-P0-AX	00D-4435-P0-AX	00F-4435-P0-AX	00G-4435-P0-AX	00B-4435-U0-AX	—	AJ0-7846	AJ0-8308
C6-Phenyl	—	00D-4444-P0-AX	00F-4444-P0-AX	00G-4444-P0-AX	—	—	AJ0-9157	AJ0-9158
5 µm							/ea	/ea
NX-C18	00B-4454-P0-AX	00D-4454-P0-AX	00F-4454-P0-AX	00G-4454-P0-AX	00B-4454-U0-AX	00C-4454-U0-AX	AJ0-8370	AJ0-8371
10 µm							/ea	/ea
C18	—	00D-4436-P0-AX	00F-4436-P0-AX	00G-4436-P0-AX	—	—	AJ0-7846	AJ0-8308
10 µm							/ea	/ea
NX-C18	00B-4455-P0-AX	00D-4455-P0-AX	00F-4455-P0-AX	00G-4455-P0-AX	—	—	AJ0-8370	AJ0-8371

for ID: 18-29 mm 30-49 mm

Axia™ Packed Preparative Columns (mm) continued						SecurityGuard™ Cartridges (mm)	
Phases	100 x 30	150 x 30	250 x 30	100 x 50	150 x 50	250 x 50	15 x 30.0*
5 µm							/ea
C18	00D-4435-U0-AX	00F-4435-U0-AX	00G-4435-U0-AX	—	—	—	AJ0-8308
5 µm							/ea
NX-C18	00D-4454-U0-AX	00F-4454-U0-AX	00G-4454-U0-AX	—	—	—	AJ0-8371
10 µm							/ea
C18	00D-4436-U0-AX	00F-4436-U0-AX	00G-4436-U0-AX	—	00F-4436-V0-AX	00G-4436-V0-AX	AJ0-8308
10 µm							/ea
NX-C18	00D-4455-U0-AX	00F-4455-U0-AX	00G-4455-U0-AX	00D-4455-V0-AX	00F-4455-V0-AX	00G-4455-V0-AX	AJ0-8371

for ID: 30-49 mm



guarantee

If Phenomenex analytical columns do not provide at least an equivalent separation as compared to a competing column of same particle size, similar phase, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Improve Your European Pharmacopoeia (Ph. Eur.) and United States (USP) Monographs

- Reduce run times
- Achieve higher resolution
- Stay within Allowable Adjustments

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Trademarks

Kinetex, Luna, and Gemini are registered trademarks and Axia, MidBore, and SecurityGuard are trademarks of Phenomenex.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Gemini and Kinetex EVO are patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

© 2017 Phenomenex, Inc. All rights reserved.