

APPLICATIONS

Improving Selectivity for Intact Mass Spectrometry of Biotherapeutic mAbs

Dawn Chen, Ivan Lebedev, Brian Rivera, and Chad Eichman
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

With many regulatory characterization requirements in biotherapeutic development, one essential requirement is intact liquid chromatography-mass spectrometry (LC-MS) by reversed phase chromatography. This technique rapidly provides data for monoclonal antibody (mAb) primary sequence confirmation, impurity identification, and heterogeneity information. High resolution mass spectrometry (HRMS) is the most common detection for intact MS analyses and is the focus of many methods in the literature.

An often-overlooked consideration of intact MS analysis is chromatographic media selection. Subtle changes in selectivity are likely to occur, even with columns of the similar particle morphology and stationary phase (e.g. – C4) mainly due to differences in pore size distribution and surface chemistry. To demonstrate these differences, we tested two columns with a core-shell particle and butyl stationary phase: [bioZen™ 2.6 µm WidePore C4](#) and [BIOshell™ 2.7 µm IgG C4](#).

For the comparison application, we evaluated a commercial biotherapeutic, rituximab, under LC-UV conditions using trifluoroacetic acid (TFA) as the mobile phase. The performance of the two columns is shown in **Figure 1**. As evident, there are minimal differences in the LC-UV trace, which shows both columns are suitable for this specific mAb.

However, when looking at intact LC-MS analysis of NIST mAb, separation profile differences are observed (**Figure 2**). In this case, the peak of interest is a light chain impurity of NIST mAb that elutes prior to the intact, full mAb. While both columns provide separation of the impurity, the baseline separation on the bioZen WidePore C4 column avoids any spectral overlap in the MS detector. Therefore, in the case of NIST mAb analysis, the stationary phase of the bioZen WidePore C4 provides the necessary separation relative to the BIOshell IgG C4. While this difference does not occur in all mAb assessments (e.g. – rituximab), it clearly demonstrates the differences in similar column options.


LC-UV Conditions (Figure 1)

Columns: bioZen 2.6 µm WidePore C4
BIOshell 2.7 µm IgG C4
Dimensions: 100 x 2.1 mm
Part No.: [OOD-4786-AN](#) (WidePore C4)
Mobile Phase A: 0.1 % TFA in Water
Mobile Phase B: 0.1 % TFA in Acetonitrile
Gradient Program: 30-40 % B in 5 min
Temperature: 60 °C
Flow Rate: 0.8 mL/min
Detection: UV @ 280 nm
Injection: 2 µL
Sample: Rituximab, 1 mg/mL


LC-MS Conditions (Figure 2)

Columns: bioZen 2.6 µm WidePore C4
BIOshell 2.7 µm IgG C4
Dimensions: 100 x 2.1 mm
Part No.: [OOD-4786-AN](#) (WidePore C4)
Mobile Phase A: 0.1 % Formic acid in Water
Mobile Phase B: 0.1 % Formic acid in Acetonitrile
Gradient Program: 5-90 % B in 6.5 min
Temperature: 0.3 mL/min
Flow Rate: 80 °C
Detection: QTOF
Injection: 2 µL
Sample: NIST mAb, 1 mg/mL


Figure 1. Comparison of BIOshell™ IgG C4 with bioZen™ WidePore C4



**Figure 2. Intact LC-MS analysis of NIST mAb
bioZen WidePore C4**



BIOshell IgG C4



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 All other countries/regions Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen is a trademark of Phenomenex. BIOshell is a trademark of MilliporeSigma.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Millipore Sigma.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.