


# Biozen Nano LC Columns For Omics Analysis

- Optimize IDs
- Consistent Results
- Easy Connections






## 2 Advanced Particle Platforms

Both Biozen particle platforms were individually designed and built by Phenomenex to take advantage of integral levels of performance, ruggedness, and reproducibility for omics applications. Individually, each platform differs in the proprietary processing techniques used to control particle size and morphology.



- High Efficiency
- Excellent Inertness
- Increased Sensitivity
- Exceptional Quality and Robustness




Thermally Modified  
Fully Porous




Core-Shell Technology

## 3 Nano Chemistries and Growing!




### Biozen Peptide PS-C18

Excellent retention by combined positively charged surface ligand and C18 ligand.



### Biozen Peptide XB-C18

Overall retention of both acidic and basic peptides through C18 stationary phase with di-isobutyl side chains.



### Biozen Polar C18

Enhanced selectivity / retention for polar analytes without diminishing useful non-polar retention

### Material Characteristics

Biozen Nano Phases	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	pH Stability*	Temp* (°C)	Pressure (psi)
Peptide XB-C18	2.6, 5	100	200	1.5-9	90	10,000
Peptide PS-C18	3	100	260	1.5-8.5	60	10,000
Polar C18	3	100	260	1.5-8.5	60	10,000


\* Temperature and pH limits are dependent on method running parameters. Continuous use of Biozen columns at the maximum limit may compromise column longevity.

# Advanced Core-Shell Particle in Nano Format

The Biozen Nano Columns utilize core-shell particles with a highly consistent morphology that minimizes band broadening associated with diffusion and mass transfer, leading to higher efficiency and minimal peak widths, which enhances the separation for omics analysis.



## Core-Shell Particle Chemistry



## High Efficiency Core-Shell Particle

Using a rigorous core construction process, a uniform porous silica layer is grown around the spherical solid silica core. This unique combination of precise particle architecture and particle size provides dramatic leaps in performance.

Fully Porous		Biozen Nano Core-Shell	Average Efficiency Gain with Biozen Nano Core-Shell
5 µm	vs	5 µm	90 % Higher
3 µm	vs	2.6 µm	85 % Higher
1.7 µm	vs	2.6 µm	Equivalent Efficiency

## Better Performance than Fully Porous Particles

Core-Shell Technology provides extremely high efficiencies for omics analysis. Industry leading column packing technology in combination with high particle consistency produces highly reproducible columns that generate greater performance compared to fully porous particles. This ultra-high efficiency can be leveraged to achieve increased resolution, improved sensitivity, and higher productivity.






Illustration - not actual test data

# Zero Dead-Volume Nano LC Connections in a CLICK

Biozen Nano's fully integrated SecurityLINK™ fingertight fitting system simplifies your system connections while providing consistent performance through Torque Limiting Technology that prevents overtightening or undertightening making every connection leak-free.



## Fingertight

Easy fingertight installation, the "CLICK" feedback confirms connection is secure

## Click n' Done


Torque Limiting Technology prevents overtightening or undertightening

## Zero Dead-Volume


Sealing technology eliminates extra internal volume

## What is Torque Limiting Technology?

Once the perfect connection has been made through fingertightening, the SecurityLINK fitting offers a haptic "CLICK" to confirm that optimum torque has been reached. This ensures a consistent connection each and every time and prevents over or undertightening that may cause column performance issues.




## Biozen Nano Column with Integrated SecurityLINK Fingertight Fittings



VS.

## Nano Columns using Standard Nut Ferrule Fittings




# Robust Selectivity Portfolio for Improved IDs and Peak Shape

The total number of proteins and peptides that were identified using nano LC-MS analysis while using Biozen 2.6µm Peptide XB-C18, Biozen 3µm Polar C18, and Biozen 3µm Peptide PS-C18 columns.




## Column Selectivity Comparison

### Protein Identifications






### Peptide Identifications



## Column Peak Width Comparison

The core-shell based Peptide XB-C18 column had narrower peaks overall. The Polar C18 column had greater overall retention of the peptides. The Peptide PS-C18 column showed very short retention for the most hydrophilic IGN peptide despite retention of the peptide being short, the peak was extremely narrow.




# Improved IDs for Trap and Elute Analysis


Number of proteins and peptides that were identified on a nano LC-MS analysis of a digested HeLa sample using a Biozen 2.6µm Peptide XB-C18 column formatted in direct inject, trap and elute using a RP-1 trap and a RP-2 trap respectively.




**Peptide XB-C18 + Trap Selectivity Protein Identifications**




**Peptide XB-C18 + Trap Selectivity Peptide Identifications**




**Polar C18 + Trap Selectivity Protein Identifications**




**Polar C18 + Trap Selectivity Peptide Identifications**



**Peptide PS-C18 + Trap Selectivity Protein Identifications**



**Peptide PS-C18 + Trap Selectivity Peptide Identifications**




# Robust Performance for Low-Flow Analysis

Peak widths and tailing factors obtained from a mixture of 20 isotopically labeled peptides injected on columns packed with Biozen 2.6  $\mu\text{m}$  core-shell Peptide XB-C18, Thermo Fisher® Acclaim™ PepMap™ 100 nanoViper™ 3  $\mu\text{m}$  fully porous C18, and Waters® nanoEase® M/Z Peptide BEH 1.7  $\mu\text{m}$  fully porous C18 particles, respectively




## Minimum Peak Widths with Nano Core-Shell



## Reproducible Performance



Extracted ion chromatograms of injection 1 and 100 from a mixture of 20 isotopically labeled peptides run on a Biozen 2.6  $\mu\text{m}$  Peptide XB-C18 150 x 0.075 mm column.





# Robust Column Performance


Number of proteins and peptides that were identified on a nano LC-MS analysis of a digested HeLa sample using a Thermo Fisher® Acclaim™ PepMap™ 100 nanoViper™ 3µm C18, Waters® nanoEase® M/Z Peptide BEH 1.7µm C18, and Biozen 2.6µm Peptide XB-C18, in trap and elute mode with Thermo Fisher Acclaim PepMap nanoViper, Waters nanoEase M/Z Symmetry C18, and Nano Trap RP-1 (General RP) traps, respectively.




## Nano LC Column Manufacturer Comparison

Number of proteins and peptides that were identified on a nano LC-MS analysis of a digested HeLa sample using a Thermo Fisher Acclaim PepMap 100 nanoViper C18, Waters nanoEase M/Z Peptide BEH 1.7µm C18, and Biozen 2.6µm Peptide XB-C18, in direct inject mode.

### Protein ID Comparison



### Peptide ID Comparison





# Fingers Only Installation

## Leak-Free Connections in 2 Easy Steps

The Biozen Nano LC Column's integrated SecurityLINK™ fingertight fitting system simplifies your system connection process and provides consistent performance with Torque Limiting Technology that prevents column damaging overtightening. No tools required for installation ever again!



### STEP 1

#### Insert Biozen Nano LC Column...

SecurityLINK fitting into the column port.



### STEP 2

#### Slowly fingertighten the...


SecurityLINK fitting until the first "**CLICK**" is received.  
Your **Leak-Free** connection is now complete!



#### System Compatibility

Biozen Nano Columns with integrated SecurityLINK Fittings were designed to be compatible with the majority of systems that are available on the market today.

*Compatible with any System with  $1/16$  in. Ports.*



# Ordering Information



## Biozen Nano LC Columns with Integrated SecurityLINK™ Fitting



Phases	150 x 0.075 mm	250 x 0.075 mm	500 x 0.075 mm
Biozen 3 µm Peptide PS-C18	<a href="#">00F-4771-AW-21</a>	<a href="#">00G-4771-AW-21</a>	—
Biozen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-21</a>	<a href="#">00G-4768-AW-21</a>	<a href="#">00J-4768-AW-21</a>
Biozen 3 µm Polar-C18	<a href="#">00F-4782-AW-21</a>	<a href="#">00G-4782-AW-21</a>	—
Biozen 5 µm Peptide XB-C18	—	—	<a href="#">00J-4792-AW-21</a>


## Biozen Nano LC Columns with Open Fused-Silica Inlet Fitting



Phases	150 x 0.075 mm	250 x 0.075 mm	500 x 0.075 mm
Biozen 3 µm Peptide PS-C18	<a href="#">00F-4771-AW-11</a>	<a href="#">00G-4771-AW-11</a>	—
Biozen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-11</a>	<a href="#">00G-4768-AW-11</a>	—
Biozen 3 µm Polar-C18	<a href="#">00F-4782-AW-11</a>	<a href="#">00G-4782-AW-11</a>	—
Biozen 5 µm Peptide XB-C18	—	—	<a href="#">00J-4792-AW-11</a>


## High pH Fractionation Column

Fractionation Column		
Part No.	Description	Dimension
<a href="#">00F-4793-AN</a>	Biozen 3 µm High pH Fractionation Column	150 x 2.1 mm




## Nano Trap Columns

Trap Columns		
Phases	10 x 0.075 mm	Unit
RP-1 (General RP)	<a href="#">05N-4252-AW</a>	3/pk
RP-2 (Aqueous Stable RP)	<a href="#">05N-4754-AW</a>	3/pk




## Fittings

Trap Fittings		
Part No.	Description	Unit
<a href="#">AQO-7602</a>	PEEKlok™ fittings with 6-40 thread for 1/32" OD tubing (2 x fittings, 6 x ferrules and 1 x tightening tool)	ea
<a href="#">AQO-7603</a>	PEEKlok fittings with 6-32 thread for 1/32" OD tubing (2 x fittings, 6 x ferrules and 1 x tightening tool)	ea
<a href="#">AQO-7600</a>	PEEKlok fittings with 10-32 thread for 1/32" OD tubing with low profile hex head (2 x fittings, 6 x ferrules and 1 x wrench)	ea



## Trap Fitting Guide

Traps		
Threads per Inch	Pitch (inches)	Pitch (mm)
32	0.0313	0.794
40	0.025	0.635



### Caution

The installation of an improper nut could potentially cause cross-threading or damage to the port and fitting

Verify fit: Traps are available for 1/16" connections (10-32 thread) or with 1/32" connections (6-40 or 6-32 thread).

## BE-HAPPY™ GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

# Ordering Information



Nano Cartridge		
Part No.	Description	Dimension
<a href="#">00F-4768-AW-SX</a>	Biozen 2.6 µm Peptide XB-C18 Nano Cartridge for SCIEX	150 x 0.075 mm
<a href="#">00G-4768-AW-SX</a>	Biozen 2.6 µm Peptide XB-C18 Nano Cartridge for SCIEX	250 x 0.075 mm
<a href="#">00J-4768-AW-SX</a>	Biozen 2.6 µm Peptide XB-C18 Nano Cartridge for SCIEX	500 x 0.075 mm
<a href="#">00F-4782-AW-SX</a>	Biozen 3 µm Peptide Polar-C18 Nano Cartridge for SCIEX	150 x 0.075 mm
<a href="#">00G-4782-AW-SX</a>	Biozen 3 µm Peptide Polar-C18 Nano Cartridge for SCIEX	250 x 0.075 mm
<a href="#">00F-4771-AW-SX</a>	Biozen 3 µm Peptide PS-C18 Nano Cartridge for SCIEX	150 x 0.075 mm
<a href="#">00G-4771-AW-SX</a>	Biozen 3 µm Peptide PS-C18 Nano Cartridge for SCIEX	250 x 0.075 mm
<a href="#">00J-4792-AW-SX</a>	Biozen 5 µm Peptide XB-C18 Nano Cartridge for SCIEX	500 x 0.075 mm




## SCIEX® OptiFlow® Source for Nanoflow Chromatography


### Single Source for All Low-Flow Chromatography

- Select the right chromatography for the workflow and samples because switching is easy and fast
- Couple with Phenomenex micro and nano columns for that plug-and-play experience
- Phenomenex column cartridge for plug-and-play solution SecurityLINK™ tubing for connection to inlet short tubing at outlet to connect to electrode

Always be working in the sweet spot of sensitivity and robustness



# Introducing Biozen Nano LC



## Biozen Nano LC Columns

### For Omics Analysis

- Optimize IDs
- Consistent Results
- Easy Connections

#### Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

#### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

#### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

#### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

#### China

t: +86 400-606-8099  
cninfo@phenomenex.com

#### Czech Republic

t: +420 272 017 077  
cz-info@phenomenex.com

#### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

#### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

#### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

#### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

#### Hong Kong

t: +852 6012 8162  
hkinfo@phenomenex.com

#### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

#### Indonesia

t: +62 21 5010 9707  
indoinfo@phenomenex.com

#### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

#### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

#### Japan

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

#### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

#### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

#### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

#### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

#### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

#### Poland

t: +48 22 104 21 72  
pl-info@phenomenex.com

#### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

#### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

#### Slovakia

t: +420 272 017 077  
sk-info@phenomenex.com

#### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

#### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

#### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

#### Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

#### Thailand

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

#### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

#### USA

t: +1 (310) 212-0555  
info@phenomenex.com

#### 🌐 All other countries/regions Corporate Office USA

t: +1 (310) 212-0555  
info@phenomenex.com



#### www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

#### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Biozen, BE-HAPPY, and SecurityLINK are trademarks of Phenomenex. Thermo Fisher is a registered trademark and Acclaim, PepMap, and nanoViper are trademarks of Thermo Fisher Scientific Inc. Waters is a registered trademark and nanoEase is a trademark of Waters Technologies Corporation. PEEKLok and PEEKsil are trademarks of Trajan Scientific Australia Pty Ltd. SCIEX and OptiFlow are registered trademarks and SteadySpray is a trademark of AB SCIEX Pte. Ltd.

#### Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Thermo Fisher or Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.