

Assessment of Column Priming for sub-2 μm Size Exclusion Columns

Dr. Ivan Lebedev, Brian Rivera, and Dr. Bryan Tackett
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Size Exclusion Chromatography (SEC) is a technique for the separation of large molecules such as proteins and polymers based on their size in solution. SEC methods are commonly used for the quantitation of aggregates in a biotherapeutic sample. Consequently, any non-specific adsorption to the column media or hardware could compromise data, as peak areas may be underreported.

As a failsafe for aggregate and sample adsorption to the column, column “priming” is often included in an analytical method sequence. This includes several injections of a protein; commonly an inert protein like bovine serum albumin (BSA) or reference standard. This helps minimize the chance for sample adsorption to occur.

In this application, we demonstrate a priming protocol, wherein 10 subsequent injections of BSA are injected onto a freshly packed 150 x 4.6 mm ID column with sub-2 μm particles. **Table 1** shows a summary of total peak areas for monomer and dimer. **Figure 1** shows a chromatographic overlay.

Table 1. Retention Times for Protein Standards

Injection	Peak Area, Monomer	Peak Area, Dimer
1	10757707	462531
2	10774224	474567
3	10780383	475693
4	10775497	478940
5	10773430	480767
6	10779073	480131
7	10785349	481109
8	10785780	480578
9	10782926	481759
10	10771993	482644
average	10776636.2	477871.9
%CV	0.1%	1.3%


LC Conditions

Column: Biozen™ 1.8 μm dSEC-2, 200 Å
Part No.: [00F-4787-E0](#)
Dimensions: 150 x 4.6 mm
Mobile Phase: 200 mM Potassium Phosphate + 250 mM KCl, pH 6.2
Flow Rate: 0.35 mL/min
Injection Volume: 10 μL
Temperature: 25 °C
Detection: UV @ 280 nm
Sample: Bovine serum albumin (10 mg/mL)

Upon initial injection, peak areas for both monomer and dimer show almost complete recovery. Additionally, as shown in the inset for **Figure 1**, only slight shifts in retention time for monomer and dimer are observed. This indicates that the extent of priming necessary for the Biozen dSEC-2 chemistry is nominal.

In summary, although the practice of column priming is typically performed, column selection may minimize the need for extensive priming.

Figure 1. Chromatographic Overlays for BSA



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

- Australia**
t: +61 (0)2-9428-6444
auinfo@phenomenex.com
- Austria**
t: +43 (0)1-319-1301
anfrage@phenomenex.com
- Belgium**
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com
- Canada**
t: +1 (800) 543-3681
info@phenomenex.com
- China**
t: +86 400-606-8099
cninfo@phenomenex.com
- Czech Republic**
t: +420 272 017 077
cz-info@phenomenex.com
- Denmark**
t: +45 4824 8048
nordicinfo@phenomenex.com
- Finland**
t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com
- France**
t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com
- Germany**
t: +49 (0)6021-58830-0
anfrage@phenomenex.com
- Hong Kong**
t: +852 6012 8162
hkinfo@phenomenex.com
- India**
t: +91 (0)40-3012 2400
indiainfo@phenomenex.com
- Indonesia**
t: +62 21 5010 9707
indoinfo@phenomenex.com
- Ireland**
t: +353 (0)1 247 5405
eireinfo@phenomenex.com
- Italy**
t: +39 051 6327511
italiainfo@phenomenex.com
- Japan**
t: +81 (0) 120-149-262
jpinfo@phenomenex.com
- Luxembourg**
t: +31 (0)30-2418700
nlinfo@phenomenex.com
- Mexico**
t: 01-800-844-5226
tecnicomx@phenomenex.com
- The Netherlands**
t: +31 (0)30-2418700
nlinfo@phenomenex.com
- New Zealand**
t: +64 (0)9-4780951
nzinfo@phenomenex.com
- Norway**
t: +47 810 02 005
nordicinfo@phenomenex.com
- Poland**
t: +48 22 104 21 72
pl-info@phenomenex.com
- Portugal**
t: +351 221 450 488
ptinfo@phenomenex.com
- Singapore**
t: +65 800-852-3944
sginfo@phenomenex.com
- Slovakia**
t: +420 272 017 077
sk-info@phenomenex.com
- Spain**
t: +34 91-413-8613
espinfo@phenomenex.com
- Sweden**
t: +46 (0)8 611 6950
nordicinfo@phenomenex.com
- Switzerland**
t: +41 (0)61 692 20 20
swissinfo@phenomenex.com
- Taiwan**
t: +886 (0) 0801-49-1246
twinfo@phenomenex.com
- Thailand**
t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com
- United Kingdom**
t: +44 (0)1625-501367
ukinfo@phenomenex.com
- USA**
t: +1 (310) 212-0555
info@phenomenex.com
- 🌐 All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPYTM

GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions
Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.
Trademarks
Biozen and BE-HAPPY are trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2021 Phenomenex, Inc. All rights reserved.