

# APPLICATIONS



## LC/MS/MS Analysis of Synthetic Cathinones (Bath Salts) from Urine and Whole Blood using a Luna<sup>®</sup> Omega 1.6 μm C18 UHPLC Column and Strata<sup>®</sup>-X-Drug B Solid Phase Extraction (SPE)

Daniel Spurgin and Dr. Jeff Layne  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA


### Overview

Bath Salts also known as synthetic cathinones have recently emerged as popular drugs of abuse. Synthetic cathinones are phenylalkylamines derivatives and are often called bk-amphetamines due to a ketone attached at the beta position on the amino alkyl chain which is attached to the phenyl ring.<sup>1</sup> Similar to phenylethylamines like MDMA, they may possess both amphetamine-like properties and the ability to modulate serotonin, causing distinct psychoactive effects.<sup>2</sup> This application illustrates the separation of 5 common synthetic cathinones using a Luna Omega 1.6 μm C18 UHPLC column. Luna Omega UHPLC columns contain a unique silica which is modified using a proprietary, post-synthetic thermal treatment process to provide extraordinary mechanical strength and significantly greater inertness than traditional fully porous and hybrid materials. This greatly minimizes any secondary interactions that negatively affect peak shape, allowing for greater method accuracy.




**Figure 1.**  
Chemical Structures of Cathinones


#### Mephedrone




#### Butylone




#### Methylone



#### α-PVP



#### Ethylone



### Experimental Conditions

#### LC/MS/MS Conditions

**LC Column:** Luna Omega 1.6 μm C18  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4742-AN  
**Mobile Phase:** A: 0.1% Formic Acid in Water  
B: 0.1% Formic Acid in Acetonitrile  
**Gradient:**

Time (min)	% B
0	10
4	95
4.1	10
6	10

**Flow Rate:** 700 μL/min  
**Injection Volume:** 1 μL  
**Temperature:** 22 °C  
**Detection:** MS/MS  
**Detection System:** API 4000<sup>TM</sup> (SCIEX)  
**Analytes:** 1. Methylone  
2. Ethylone  
3. Butylone  
4. Mephedrone  
5. α-PVP

### Sample Preparation

#### Pretreatment

Urine	Add 2 mL of 100 mM sodium acetate buffer (pH 5.0) and 50 μL of internal standards (@10 ppm) to 2 mL of urine
Whole blood	1. Add 2 mL of ice cold methanol:acetonitrile (50:50) and 20 μL of internal standards (@10 ppm) and 2 mL of 100 mM sodium acetate to 1 mL of blood. 2. Centrifuge at 4700 rpm 10 °C for 5 min 3. Transfer supernatant for SPE

#### Solid Phase Extraction (SPE)


<b>Cartridge:</b>	Strata-X-Drug B
<b>Part No.:</b>	8B-S128-UCH
<b>Condition:</b>	1 mL each of methanol, DI water, and 100 mM sodium acetate
<b>Load:</b>	Load pretreated sample
<b>Weak Wash:</b>	2 mL of 100 mM sodium acetate (pH 5.0)
<b>Strong Wash:</b>	1 mL of methanol
<b>Dry Down:</b>	2 minutes at >10" hg
<b>Elute:</b>	3 mL of ethyl acetate: IPA: ammonium hydroxide (70:20:10)
<b>Evaporate:</b>	to 500 μL and add 100 μL of HCl:methanol (20:80), evaporate to dryness under nitrogen
<b>Reconstitute:</b>	100 μL of methanol

Q1	Q2	Analyte	Retention Time (min)
208.1	160.2	Methylone	1.56
222.3	174.1	Ethylone	1.69
222.3	174.1	Butylone	1.80
178.1	160.2	Mephedrone	1.84
232.4	91	α-PVP	2.15




# APPLICATIONS

## Extracted Ion Chromatogram of Bath Salts in Urine



APP ID 23946

## Extracted Ion Chromatogram of Bath Salts in Whole Blood



APP ID 23948



If the Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same particle size, phase, and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

# APPLICATIONS











## Luna® Omega Ordering Information

1.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges†
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>Polar C18</b>	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJO-9505
<b>PS C18</b>	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJO-9508
<b>C18</b>	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJO-9502

for 2.1 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

## Strata®-X-Drug B Ordering Information

Format	Sorbent Mass	Part Number	Unit
<b>Tube</b>			
	10 mg	8B-S128-AAK	1 mL (100/box)
	30 mg	8B-S128-TAK	1 mL (100/box)
	30 mg	8B-S128-TBJ	3 mL (50/box)
	60 mg	8B-S128-UBJ	3 mL (50/box)
	60 mg	8B-S128-UCH	6 mL (30/box)
	60 mg	8B-S128-UCL	6 mL (200/bag)
<b>Giga™ Tube</b>			
	100 mg	8B-S128-EDG	12 mL (20/box)
<b>96-Well Plate</b>			
	10 mg	8E-S128-AGB	2 Plates/box
	30 mg	8E-S128-TGB	2 Plates/box
	60 mg	8E-S128-UGB	2 Plates/box

### Reference:

- Zaitseva K, Katagi M, Tatsuno M et al. Recently abused b-keto derivatives of 3,4-methylenedioxyphenylalkylamines: a review of their metabolisms and toxicological analysis. *Forensic Toxicol* 2011.
- Prosser J, Nelson L. *The Toxicology of Bath Salts: A Review of Synthetic Cathinones*. *J. Med. Toxicol*. 2012



### www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

#### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Luna Omega and Strata are registered trademarks and SecurityGuard and Giga are trademarks of Phenomenex. API 4000 is a trademark of SCIEX Pte. Ltd. SCIEX is being used under license.

#### Disclaimer

Strata-X is patented by Phenomenex. U.S. patent No 7,119,145.  
For RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2017 Phenomenex, Inc. All rights reserved.

#### Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auiinfo@phenomenex.com

#### Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

#### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

#### Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

#### China

t: +86 400-606-8099  
f: +86 (0)22 2532-1033  
phen@agela.com

#### Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

#### Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

#### France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

#### Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

#### India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

#### Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eirinfo@phenomenex.com

#### Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

#### Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

#### Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

#### The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

#### New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

#### Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

#### Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

#### Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

#### Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

#### United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

#### USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

#### All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

