

APPLICATIONS

Assessment of Disulfide Variants of IgG2 Monoclonal Antibodies by Intact Reversed Phase

Laurence Brill, Ivan Lebedev, and Brian Rivera
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

A primary analytical technique for purity for recombinant proteins is reversed phase LC (RPLC). This simple, robust methodology has a relatively short analysis time, as well being high resolution, capable of separation of hydrophobic variants such as oxidation, glycoforms, and clipping. As such, it is a common analytical approach for monoclonal antibodies (mAbs), especially considering it can be implemented at both intact and subunit levels without extensive method optimization.

One application where intact reversed phase is commonly used is in the assessment of disulfide heterogeneity for IgG2 isotypes. Because their Fc-mediator effector function differ when compared to IgG1s, IgG2s may be a more appropriate therapeutic platform depending on the biotherapeutic mechanism of action. As it pertains to reversed phase chromatography, IgG2s uniquely have 4 interchain disulfides, which consequently yields unique arrangements with different structural properties. Each variant behaves differently chromatographically, and although orthogonal methods such as disulfide mapping are necessary for full characterization, intact reversed phase is still a useful tool for monitoring and assessing overall disulfide variant isoforms.

When developing a platform method for IgG1s, implementing a relatively shallow gradient and short run times is more feasible. Since physicochemical properties for most IgG1s are very similar, a gradient slope of 0.4 % B per column volume and a gradient length of 5 minutes might yield optimal results on a wide variety of different molecules. However, platform methods for IgG2 antibodies requires a wider window of organic to accommodate the differences in physicochemical properties between them. To ensure shallow gradient slope, the gradient length subsequently must be extended.

The resulting method uses a 0.6 % B per column volume gradient slope, with a 20-45 % B gradient program over 10 minutes. Of particular importance for this method is the maintenance of the high flow-rate; ensuring that the linear velocity is high allows for the shallowing of gradient slope without compromising the method run time. This is important for intact reversed phase methods, as temperature is a method that is optimized to ensure optimal separation, and aside from gradient slope, is the most critical method parameter to optimize for intact reversed phase methods.

Figure 1 shows the separation of panitumumab, a relatively hydrophobic IgG2. The inset shows different isoforms, putative disulfide variants, for panitumumab. Figure 2 shows a similar profile for denosumab, another IgG2 monoclonal antibody. Previous studies suggest that earlier eluting variants are IgG2-B variants, with IgG-A/B and IgG-A being more hydrophobic thus retaining longer¹. These would need to be confirmed via disulfide mapping.


In summary, purity analysis by reversed phase LC is a primary method for the analysis of recombinant proteins. For intact analysis of IgG2 isotypes, intact reversed phase is a relatively simple and high throughput method to understand the overall disulfide heterogeneity of disulfide isoforms.

HPLC Conditions

Column:	bioZen 2.6 µm WidePore C4
Dimension:	100 x 2.1 mm
Part No.:	00D-4786-AN
Mobile Phase A:	0.1 % TFA in Water
Mobile Phase B:	0.1 % TFA in Acetonitrile
Gradient Program:	20-45 % B in 10 minutes, 0.6 % B/CV
Flow-rate:	0.8 mL/min
Temperature:	80° C
Detection:	UV-Vis @ 280 nm
Injection:	2 µL
Samples:	As indicated (0.5 mg/mL)

Figure 1. Disulfide Variants of Panitumumab


Intact reversed phase analysis of panitumumab. The shallow 0.6% B per column volume gradient slope facilitates the separation. Putative variants include IgG2-B, IgG2-A/B, and IgG2-A, in order of elution. Shallowing of the gradient further might improve the separation for this relatively hydrophobic IgG2.



App ID 26266

Figure 2. Disulfide Variants of Denosumab

Intact reversed phase analysis of denosumab. Slightly earlier eluting than panitumumab, showing similar variants. These include IgG2-B, IgG2-A/B, and IgG2-A.



App ID 26267

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thainfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☎ **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen are registered trademarks of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.