



# 5 Frequently Asked Question About Lux<sup>®</sup> Chiral Columns






 **phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>

[phenomenex.com/Lux](http://phenomenex.com/Lux)

## Q1: How can I reverse enantiomeric elution order?

Because of the intrinsic nature of enantiomers, it is difficult to predict selectivity and therefore the relative analyte elution order. However, the Lux® polysaccharide chiral column portfolio contains a wide range of chiral stationary phases that can easily be screened for complementary or orthogonal selectivity. For instance, Lux i-Amylose-3 and i-Cellulose-5 have complementary but distinct selectivity in comparison to each other. In addition, the phase immobilization affords greater solvent flexibility, increasing the potential for enantioselectivity.



**Columns:** Lux 5  $\mu$ m i-Amylose-3  
Lux 5  $\mu$ m i-Cellulose-5  
Lux 5  $\mu$ m i-Amylose-1

**Dimensions:** 250 x 4.6 mm

**Part No.:** 00G-4779-E0  
00G-4756-E0  
00G-4762-E0

**Mobile Phase:** Water with 0.1 % Diethylamine/Acetonitrile (35:65)


**Flow Rate:** 1.0 mL/min

**Injection Volume:** 10  $\mu$ L (2 mg/mL)

**Detection:** UV @ 254 nm

**Sample:** 1. Diniconazole  
2. Diniconazole


### Enantiomers of Diniconazole




## Frequently Asked Chiral Questions

# Q2: Why does Phenomenex use DEA over TEA as the primary basic modifier in application notes?

Both DEA (diethylamine) and TEA (triethylamine) are widely published as good basic modifiers for improving peak shapes on polysaccharide-type chiral columns. We chose DEA for our initial screening data and have continued with it routinely to maintain consistency. TEA is also just as effective and commonly used by many customers successfully on our Lux<sup>®</sup> polysaccharide chiral columns.



### Enantiomers of Thalidomide



**Columns:** Lux 5  $\mu\text{m}$  i-Amylose-3  
**Dimensions:** 250 x 4.6 mm  
**Part No.:** 00G-4779-E0  
**Mobile Phase:** Acetonitrile with 0.1 % Diethylamine  
**Flow Rate:** 1.0 mL/min  
**Injection Volume:** 10  $\mu\text{L}$  (2 mg/mL)  
**Detection:** UV @ 254 nm  
**Sample:** 1. Thalidomide  
2. Thalidomide

**Need a chiral  
method right now?**


Start your chiral method search by an analyte  
name or structure search in seconds here:

[www.phenomenex.com/lux](http://www.phenomenex.com/lux)

## Frequently Asked Chiral Questions


### Q3: How does the aromaticity of chiral compounds affect selectivity?

In chiral compounds, the proximity of aromatic groups to the stereocenter is typically linked to the ease of enantiomeric resolution. For instance, the separations of enantiomers in which the aromatic functionality is 4 or more atoms away from the stereocenter can be challenging and chromatographically uncommon. Enantioselectivity is most effective when the distances between the aromatic group and stereocenter are equivalent in both chiral conformations. If the aromatic group of the compound has electron withdrawing groups like halogens or oxygen that will be more electron deficient and will interact more effectively with electron-rich aromatic groups of the chiral stationary phase.



App ID 26340

#### Enantiomers of Econazole






- Columns:** Lux® 5 µm i-Amylose-3
- Dimensions:** 250 x 4.6 mm
- Part No.:** 00G-4779-E0
- Mobile Phase:** Water with 5 mM Ammonium Acetate + 0.05 % Formic Acid/Acetonitrile (65:35)
- Flow Rate:** 1.0 mL/min
- Injection Volume:** 10 µL (2 mg/mL)
- Detection:** UV @ 254 nm
- Sample:** 1. Econazole  
2. Econazole


## Frequently Asked Chiral Questions

# Q4: What is the difference between Amylose and Cellulose polysaccharide backbones?

Both backbones form 3-dimensional helical structures that are well suited to provide grooves and cavities for potential steric interactions, Hydrogen bonding, dipole-dipole, and  $\pi$ - $\pi$  based interactions. However, the structure for amylose is considered more tightly coiled in comparison to the looser cellulose structure which may accommodate interactions differently. Practically, this difference can result in distinct selectivity for amylose and cellulose even if with an identical chiral selector, as the 3-dimensional orientation around the CSP will differ.



### Enantiomers of Napropamide



**Columns:** Lux<sup>®</sup> 5  $\mu$ m i-Amylose-3  
Lux 5  $\mu$ m i-Cellulose-5  
Lux 5  $\mu$ m i-Amylose-1

**Dimensions:** 250 x 4.6 mm

**Part No.:** 00G-4779-E0  
00G-4756-E0  
00G-4762-E0

**Mobile Phase:** Water with 5 mM Ammonium Acetate + 0.05% Formic Acid/Acetonitrile (53:47)

**Flow Rate:** 1.0 mL/min

**Injection Volume:** 10  $\mu$ L (2 mg/mL)


**Detection:** UV @ 254 nm

**Sample:** 1. Napropamide  
2. Napropamide

## Frequently Asked Chiral Questions

### Q5: What separation modes can Lux<sup>®</sup> polysaccharide chiral columns operate in?

Lux chiral columns can be used in reversed phase, normal phase, polar organic, polar ionic and SFC conditions. However, the immobilization and bonding technology used within the Lux i-Amylose-3, i-Cellulose-5 and i-Amylose-1 columns promotes stability in strong organic solvents, which affords you the ability to expand your chiral separation success with even more solvent systems.





**Terms and Conditions** - Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions) • **Trademarks** - Lux is a registered trademark, SecurityLINK, and SecurityGuard are trademarks of Phenomenex. • **Disclaimer** - Comparative separations may not be representative of all applications. SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362. • **CAUTION:** this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges. • **FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.** • © 2020 Phenomenex, Inc. All rights reserved.