


Customized SPE Solutions for Environmental Testing

Your samples aren't created the same, so don't settle for general solid phase extraction solutions. Finding the right SPE sorbent provides high recoveries, clean injections, and gains in sensitivity.




Polymeric SPE Options

Target Analyte	Strong Acids ($pK_a < 2$)	Weak Acids ($pK_a 2-4$)	Neutral Compounds	Weak Bases ($pK_a 8-10$)	Strong Bases ($pK_a > 10$)
	Reversed Phase & Weak Anion-Exchange	Reversed Phase & Strong Anion-Exchange	Reversed Phase	Reversed Phase & Strong Cation-Exchange	Reversed Phase & Weak Cation-Exchange
	Strata-X-AW	Strata-X-A	Strata-X	Strata-X-C	Strata-X-CW
Retention Mechanism	Di-amino	Quaternary Amine	N-Vinylpyrrolidone	Sulfonic Acid	Carboxylic Acid
Recommended Methods	<p>Condition: 1 mL Methanol</p> <p>Equilibrate: 1 mL Acidified Water</p> <p>Load: 2 mL Diluted Sample (pH 6-7)</p> <p>Wash: 1 mL 25 mM Ammonium Acetate* (pH 6-7)</p> <p>Elute Neutrals/Bases: 2 x 500 μL Methanol</p> <p>Elute Any Acid: 2x 500 μL 5% NH_4OH in Methanol</p> <p>Elute Weak Acids: 2x 500 μL 5% Formic Acid in Methanol</p>	<p>Condition: 1 mL Methanol</p> <p>Equilibrate: 1 mL Water</p> <p>Load: 2 mL Diluted Sample (pH 6-7)</p> <p>Wash: 1 mL 25 mM Ammonium Acetate* (pH 6-7)</p> <p>Elute Neutrals/Bases: 2x 500 μL Methanol*</p> <p>Wash: 1 mL Methanol*</p> <p>Elute Acids: 2x 500 μL 5% Formic Acid in Methanol</p>	<p>Condition: 1 mL Methanol</p> <p>Equilibrate: 1 mL Water</p> <p>Load: 2 mL Diluted Sample</p> <p>Wash: 1 mL 5-60% Methanol*</p> <p>Elute: 2x 500 μL 2% Formic Acid in Methanol or Acetonitrile</p>	<p>Condition: 1 mL Methanol</p> <p>Equilibrate: 1 mL Acidified Water</p> <p>Load: 2 mL Diluted Acidified Sample (pH 6-7)</p> <p>Wash: 1 mL 0.1 N HCl in Water</p> <p>Elute Neutrals/Acids: 2x 500 μL 0.1 N HCl in Methanol*</p> <p>Wash: 1 mL 0.1 N HCl in Methanol*</p> <p>Elute Bases: 2x 500 μL 5% NH_4OH in Methanol</p>	<p>Condition: 1 mL Methanol</p> <p>Equilibrate: 1 mL Water</p> <p>Load: 2 mL Diluted Sample (pH 6-7)</p> <p>Wash: 1 mL Water</p> <p>Elute Neutrals/Acids: 2 x 500 μL Methanol</p> <p>Elute Any Base: 2x 500 μL 5% Formic Acid in Methanol</p> <p>Elute Weak Bases: 2x 500 μL 5% NH_4OH in Methanol</p>

* Dry cartridge 5 - 10 mins prior to elution to remove residual solvents.

Matrix Options (Not limited to pictures below)



Recommended Loading Capacity

Strata-X Phase	Soil Extracts	Water (particulate-free)	Water (particulate-laden)	Mass (mg in tube)
Strata-X, X-C, X-CW, X-A, X-AW	N.A.	N.A.	N.A.	10 mg
	N.R.	N.R.	N.R.	30 mg
	N.R.	N.R.	N.R.	60 mg
	20 g	50 mL	25 mL	100 mg
	40 g	100 mL	50 mL	200 mg
	100 g	500 mL	100 mL	500 mg
Strata-XL, XL-C, XL-CW, XL-A, XL-AW	200 g	1 mL	200 mL	1 g
	N.A.	N.A.	N.A.	10 mg
	N.R.	N.R.	N.R.	30 mg
	N.R.	N.R.	N.R.	60 mg
	10 g	25 mL	13 mL	100 mg
	20 g	50 mL	25 mL	200 mg
50 g	250 mL	50 mL	500 mg	
100 g	500 mL	100 mL	1 g	

N.A. = Not Applicable (not commonly used)
N.R. = Not Recommended (may not provide expected results)

Wash and Elution Solvent Volumes

strata [®] Sorbent Mass	100mg	150mg	200mg	500mg	1g	2g	5g	10g
Practical Minimum Wash and Elution Volume 4 bed volumes	1 mL	1.5 mL	2 mL	5 mL	10 mL	20 mL	50 mL	100 mL
Recommended Wash and Elution Volume 8 bed volumes	2 mL	3 mL	4 mL	10 mL	20 mL	40 mL	100 mL	200 mL

Available Formats for Strata-X or Strata SPE

Giga [™] Tubes	Tubes	On-Line Columns and Cartridges
Large format tubes ideal for large volume samples	Excellent for most common extractions and available in a large variety of different volumes	Run samples in parallel with LC analysis for further time savings

Additional Strata Silica Based Phases for Your Analytes

- FL-PR (Florisil[®]) for the retention of polar and halogenated compounds, such as pesticides
- EPH (Extractable Petroleum Hydrocarbons) for the fractionation of aliphatic and aromatic hydrocarbons
- SDB-L (styrene divinylbenzene) hydrophobic and aromatic selectivity for reversed phase applications
- Si-1 (Silica) for strong polar selectivity
- C18-E for extraction of hydrophobic or polar organic analytes from aqueous matrices
- AI-N for the extraction of polar compounds from environmental samples
- NH₂ for strong polar selectivity and hydrogen bonding under normal phase conditions
- PAH (Polycyclic Aromatic Hydrocarbons) for polycyclic aromatic hydrocarbons and humic acids from water (as specified in EPA Method 550.1)



presston[™] 1000

Streamlined Sample Preparation with Presston 1000 Positive Pressure Manifold



www.phenomenex.com/Presston

Questions?
Email us at support@phenomenex.com



SPE Method Development Tool
Develop a customized solid phase extraction method in under 1 minute! www.phenomenex.com/MDTool



For more information and to order, go to www.phenomenex.com/StrataX