


Zebron PLUS GC Inlet Liners

Remarkably Inert

- Superior sample protection
- Easy to install
- A solution for all of your analyses

Safe and Simple Installation

No more broken or contaminated liners! Easy to open, easy to install Zebron PLUS liner packaging allows you to protect yourself from cuts or breakage while reducing any contamination that may occur during the installation process.



Choose Your Liner



We've made it easier than ever to find your ideal GC liner! Simply follow the selection chart below, use your current liner part number, or visit us online.

1. Find by GC System

| | |
|-------------------------|------|
| Agilent®..... | p. 6 |
| Thermo Scientific®..... | p. 6 |
| PerkinElmer®..... | p. 7 |
| Shimadzu® | p. 8 |

2. Find by Your Current Liner Part Number..... p. 9

3. Find Your Ideal GC Liner Online! Search by:



- Application
- Injection Type
- GC System
- Your Current Liner Part Number

www.phenomenex.com/FindLiner


Become A GC Liner Expert

- p. 10** Tips To Improve Your Analysis
- p. 12** Frequently Asked Questions
- p. 14** Meet The Zebtron Family

Deactivation Protects Your Samples

Zebron PLUS GC inlet liners undergo a unique deactivation process, resulting in a remarkably inert pathway that prevents sample adsorption and degradation for active compounds.

See The Difference: Zebron PLUS vs. Traditional Deactivation



Wool Enhances Your Analysis








Analysis of dirty samples, samples with wide boiling points, or samples with wide molecular weights can be enhanced by choosing a Zebron PLUS liner with quartz wool. Liners packed with wool can prevent non-volatile compounds from entering the column and also improve sample vaporization for a more robust, reproducible analysis.

Why Use Pre-Packed Liners?

Though the upfront cost of self-packing your liner may seem attractive, the time and headaches caused by resulting tailing or irreproducible peaks can be sizeable! Self-packed wool fibers commonly break during installation and any existing deactivation on the liner can also be scratched or damaged. Pre-packed Zebron PLUS liners undergo the deactivation process with the quartz wool already in place, which ensures that any active sites that form during packing are not exposed.

Liner Geometry Improves Performance

The geometry of a liner refers to the type and shape of a liner. Each geometry addresses a particular sample or injection type. Use the chart below to determine the best geometry for your sample.

| | Injection Type | | | Primary Challenge | | | | | | |
|--|----------------|-----------|--------|-------------------|------------------|----------------------|---------------|---|--------------------------|-----------------|
| | Split | Splitless | Direct | Trace Analysis | Active Compounds | Concentrated Samples | Dirty Samples | Gaseous Samples (e.g. headspace, purge and trap) | Wide Boiling Point Range | General Purpose |
| Straight  | ● | | | | | ● | ● | ● | | ● |
| Straight Z-Liner™  | ● | | | | | ● | ● | | | ● |
| Single Taper  | | ● | | ● | ● | | | | | ● |
| Single Taper with Wool  | ● | ● | | ● | | ● | ● | | | ● |
| Single Taper Z-Liner  | ● | ● | | ● | | | ● | | ● | ● |
| Direct Connect  | | | ● | ● | ● | | | | | ● |
| Single Baffle  | ● | ● | | ● | | | | | | |










It is important to ensure that your liner has the proper dimensions for your analysis. If the liner volume is too low, your sample may not fit in the liner after it is vaporized due to expansion. This can result in irreproducible peaks and carryover.

What's A Z-Liner?


Zebron PLUS liners with a Z-Liner geometry contain optimally-placed deactivated quartz wool, which is held in place by two tapered sections of glass inside the liner. This ensures that the wool remains in the correct position for injection, wipes the needle tip completely clean, and properly volatilizes the sample.


For Agilent® or Thermo Scientific® GC Systems

Ordering Information

| Description | Application | Dimensions ID x L (mm) | Part No. | Unit | Price |
|---|--------------------------------|------------------------|-------------|-------|-------|
| Direct Connect | Trace analysis | 4 x 78.5 | AG2-0A50-01 | Ea | |
|  | Splitless injections | | AG2-0A50-05 | 5/pk | |
| | | | AG2-0A50-25 | 25/pk | |
| Single Taper | Pesticides | 4 x 78.5 | AG2-0A10-01 | Ea | |
|  | | | AG2-0A10-05 | 5/pk | |
| | | | AG2-0A10-25 | 25/pk | |
| Single Taper Z-Liner™ | Semi-volatiles | 4 x 78.5 | AG2-0A13-01 | Ea | |
|  | Dirty samples | | AG2-0A13-05 | 5/pk | |
| | | | AG2-0A13-25 | 25/pk | |
| Single Taper with Wool | Semi-volatiles | 4 x 78.5 | AG2-0A11-01 | Ea | |
|  | | | AG2-0A11-05 | 5/pk | |
| | | | AG2-0A11-25 | 25/pk | |
| Straight | Volatiles | 4 x 78.5 | AG2-0A00-01 | Ea | |
|  | | | AG2-0A00-05 | 5/pk | |
| | | | AG2-0A00-25 | 25/pk | |
| Straight Z-Liner | Dirty samples | 4 x 78.5 | AG2-0A03-01 | Ea | |
|  | Volatiles | | AG2-0A03-05 | 5/pk | |
| | High initial oven temperatures | | AG2-0A03-25 | 25/pk | |
| Straight Single Baffle | Semi-volatiles | 1.8 x 71 | AG2-1F06-01 | Ea | |
|  | Pesticides | | AG2-1F06-05 | 5/pk | |
| | | | AG2-1F06-25 | 25/pk | |






Complete Your Analysis

| Easy Seals Inlet Base Seals | | | | 2/pk | | 10/pk | |
|---|----------------|--------------|--------------------------|----------|-------|----------|-------|
| Description | Injection Type | Groove Style | Inlet Hole Diameter (mm) | Part No. | Price | Part No. | Price |
| Easy Seals Gold Inlet Seal  | Splitless | Single | 0.8 | AG0-8619 | | AG0-8620 | |


| Description | | Diameter (mm) | (in.) | Includes GuideRight™ Hole | Part No. | Unit | Price |
|---|---|---------------|-------|---------------------------|----------|-------|-------|
| PhenoRed™ -400 Septa  | • Our most popular choice for low-bleed septa, rated to 400 °C | 9.5 | 3/8 | ✓ | AG0-7916 | 50/pk | |
| | • Designed and conditioned for high sensitivity | 11 | 7/16 | ✓ | AG0-7917 | 50/pk | |
| | • Durometer rating of 50 with typical injection life of 150 punctures | | | | | | |

For PerkinElmer® GC Systems

Ordering Information

| Description | Application | Dimensions ID x L (mm) | Part No. | Unit | Price |
|---|--------------------------------|---------------------------|-------------|-------|-------|
| Single Taper | Pesticides | 4 x 92 | AG2-2A10-01 | Ea | |
|  | | | AG2-2A10-05 | 5/pk | |
| | | | AG2-2A10-25 | 25/pk | |
| Single Taper Z-Liner™ | Semi-volatiles | 4 x 92 | AG2-2A13-01 | Ea | |
|  | Dirty samples | | AG2-2A13-05 | 5/pk | |
| | | | AG2-2A13-25 | 25/pk | |
| Straight | Volatiles | 4 x 92 | AG2-2A00-01 | Ea | |
|  | | | AG2-2A00-05 | 5/pk | |
| | | | AG2-2A00-25 | 25/pk | |
| Straight Z-Liner | Volatiles | 4 x 86.2 | AG2-2E03-01 | Ea | |
|  | Dirty samples | | AG2-2E03-05 | 5/pk | |
| | | | AG2-2E03-25 | 25/pk | |
| Straight Z-Liner | High initial oven temperatures | 4 x 92 | AG2-2A03-01 | Ea | |
|  | | | AG2-2A03-05 | 5/pk | |
| | | | AG2-2A03-25 | 25/pk | |

Complete Your Analysis

| Description | | Diameter (mm) | (in.) | Includes GuideRight™ Hole | Part No. | Unit | Price |
|---|--|---------------|-------|---------------------------|----------|-------|-------|
| PhenoRed™ -400 Septa  | <ul style="list-style-type: none"> Our most popular choice for low-bleed septa, rated to 400 °C Designed and conditioned for high sensitivity Durometer rating of 50 with typical injection life of 150 punctures | 9.5 | 3/8 | ✓ | AG0-7916 | 50/pk | |
| | | 11 | 7/16 | ✓ | AG0-7917 | 50/pk | |









More Inlet Consumables Are Available Online!

Need inlet seals, septa, or syringes? Explore hundreds of available parts online at


www.phenomenex.com/InletGC

For Shimadzu® GC Systems

Ordering Information

| Description | Application | Dimensions ID x L (mm) | Part No. | Unit | Price |
|---|--|---------------------------|---|---------------------|-------|
| Shimadzu 17A GC Systems | | | | | |
| Single Taper Z-Liner™  | Pesticides | 3.4 x 95 | AG2-3B13-01 AG2-3B13-05 AG2-3B13-25 | Ea 5/pk 25/pk | |
| Straight Z-Liner  | Volatiles Dirty samples High initial oven temperatures | 3.4 x 95 | AG2-3B03-01 AG2-3B03-05 AG2-3B03-25 | Ea 5/pk 25/pk | |
| Shimadzu 2010 GC Systems | | | | | |
| Single Taper Liner  | Pesticides | 3.4 x 95 | AG2-4B10-01 AG2-4B10-05 AG2-4B10-25 | Ea 5/pk 25/pk | |
| Single Taper Z-Liner  | Semi-volatiles Dirty samples | 3.4 x 95 | AG2-4B13-01 AG2-4B13-05 AG2-4B13-25 | Ea 5/pk 25/pk | |
| Straight Liner  | Volatiles | 3.4 x 95 | AG2-4B00-01 AG2-4B00-05 AG2-4B00-25 | Ea 5/pk 25/pk | |
| Straight Z-Liner  | Volatiles Dirty samples High initial oven temperatures | 3.4 x 95 | AG2-4B03-01 AG2-4B03-05 AG2-4B03-25 | Ea 5/pk 25/pk | |

Complete Your Analysis

| Description | Part No. | Unit | Price |
|---|----------|-------|-------|
| BTO® Silicone Septa Plug  | AG0-7517 | 50/pk | |

- Fits Shimadzu (9A, 14, 15A, 17A, 2010) and SRI injectors
- Rated to 400 °C

guarantee

If Zebron PLUS GC Inlet Liners do not provide you with equivalent or better separations as compared to any other GC inlet liners of the same geometry and comparable dimensions, return the liners with comparative data within 45 days for a FULL REFUND.

Liner Cross-Reference Guide

Liner Part Number






| Geometry | Instrument | ID (mm) | Length (mm) | Unit | Restek® Sky® | Restek | Agilent Ultra Inert | Agilent | Zebron™ PLUS |
|------------------------|-------------------------------|---------|-------------|-------|------------------|--------------|---------------------|------------|--------------|
| Straight | Agilent® / Thermo Scientific® | 4 | 78.5 | 5/pk | 23301.5 | 20782-214.5 | | 210-3003-5 | AG2-0A00-05 |
| Straight | Agilent / Thermo Scientific | 4 | 78.5 | 25/pk | 23301.25 | 20783-214.25 | | | AG2-0A00-25 |
| Straight Z-Liner™ | Agilent / Thermo Scientific | 4 | 78.5 | 5/pk | 23305.5 | | | 210-4004-5 | AG2-0A03-05 |
| Straight Z-Liner | Agilent / Thermo Scientific | 4 | 78.5 | 25/pk | 23305.25 | | | | AG2-0A03-25 |
| Single Taper | Agilent / Thermo Scientific | 4 | 78.5 | 5/pk | 23302.5 | | 5190-3162 | 5183-4695 | AG2-0A10-05 |
| Single Taper | Agilent / Thermo Scientific | 4 | 78.5 | 25/pk | 23302.25 | | 5190-3166 | 5183-4696 | AG2-0A10-25 |
| Single Taper with Wool | Agilent / Thermo Scientific | 4 | 78.5 | 5/pk | 23303.5 | | 5190-3163 | 5183-4693 | AG2-0A11-05 |
| Single Taper with Wool | Agilent / Thermo Scientific | 4 | 78.5 | 25/pk | 23303.25 | | 5190-3167 | 5183-4694 | AG2-0A11-25 |
| Single Taper Z-Liner | Agilent / Thermo Scientific | 4 | 78.5 | 5/pk | | 22984 | | 210-4022-5 | AG2-0A13-05 |
| Single Taper Z-Liner | Agilent / Thermo Scientific | 4 | 78.5 | 25/pk | | 22985 | | | AG2-0A13-25 |
| Direct Connect | Agilent / Thermo Scientific | 4 | 78.5 | 5/pk | | | | | AG2-0A50-05 |
| Direct Connect | Agilent / Thermo Scientific | 4 | 78.5 | 25/pk | | | | | AG2-0A50-25 |
| Straight | PerkinElmer® | 4 | 92 | 5/pk | | | | | AG2-2A00-05 |
| Straight | PerkinElmer | 4 | 92 | 25/pk | | | | | AG2-2A00-25 |
| Straight Z-Liner | PerkinElmer | 4 | 92 | 5/pk | 23450.5 | | | | AG2-2A03-05 |
| Straight Z-Liner | PerkinElmer | 4 | 92 | 25/pk | | | | | AG2-2A03-25 |
| Single Taper | PerkinElmer | 4 | 92 | 5/pk | | | | | AG2-2A10-05 |
| Single Taper | PerkinElmer | 4 | 92 | 25/pk | | | | | AG2-2A10-25 |
| Single Taper Z-Liner | PerkinElmer | 4 | 92 | 5/pk | | | | | AG2-2A13-05 |
| Single Taper Z-Liner | PerkinElmer | 4 | 92 | 25/pk | | | | | AG2-2A13-25 |
| Straight Z-Liner | PerkinElmer (PSS Injector) | 4 | 86.2 | 5/pk | | | | | AG2-2E03-05 |
| Straight Z-Liner | PerkinElmer (PSS Injector) | 4 | 86.2 | 25/pk | | | | | AG2-2E03-25 |
| Straight Z-Liner | Shimadzu® 17A | 3.4 | 95 | 5/pk | 23320.5 | | | | AG2-3B03-05 |
| Straight Z-Liner | Shimadzu 17A | 3.4 | 95 | 25/pk | | | | | AG2-3B03-25 |
| Single Taper Z-Liner | Shimadzu 17A | 3.4 | 95 | 5/pk | 23455.5 | | | | AG2-3B13-05 |
| Single Taper Z-Liner | Shimadzu 17A | 3.4 | 95 | 25/pk | | | | | AG2-3B13-25 |
| Straight | Shimadzu 2010 | 3.4 | 95 | 5/pk | 23318.5 | | | | AG2-4B00-05 |
| Straight | Shimadzu 2010 | 3.4 | 95 | 25/pk | 23318.25 | | | | AG2-4B00-25 |
| Straight Z-Liner | Shimadzu 2010 | 3.4 | 95 | 5/pk | 23320.5 | | | | AG2-4B03-05 |
| Straight Z-Liner | Shimadzu 2010 | 3.4 | 95 | 25/pk | | | | | AG2-4B03-25 |
| Single Taper | Shimadzu 17A / 2010 | 3.4 | 95 | 5/pk | 23321.5 | | | 8001-0104 | AG2-4B10-05 |
| Single Taper | Shimadzu 17A / 2010 | 3.4 | 95 | 25/pk | 23321.25 | | | | AG2-4B10-25 |
| Single Taper Z-Liner | Shimadzu 2010 | 3.4 | 95 | 5/pk | 23455.5 | | | | AG2-4B13-05 |
| Single Taper Z-Liner | Shimadzu 2010 | 3.4 | 95 | 25/pk | | | | | AG2-4B13-25 |
| Straight Single Baffle | Agilent PTV | 1.8 | 71 | 5/pk | 23432.10 (10/pk) | | | | AG2-1F06-05 |
| Straight Single Baffle | Agilent PTV | 1.8 | 71 | 25/pk | | | | | AG2-1F06-25 |

Find Your Liner Online!


Easily search by part numbers, application, injection mode, or system manufacturer for quick selection in under 1 minute!

www.phenomenex.com/FindLiner

Tips To Improve Your Analysis


| | |
|---|---|
|  | <p>Fronting Peaks</p> <p>Fronting peaks can occur when the expanded sample exceeds the liner's volume. Reducing your injection volume, increasing split, using a liner with a larger volume, or using a solvent with lower expansion volume can remedy this issue.</p> |
|  | <p>Ghost Peaks</p> <p>Septum bleed or cored septum fragments lodged in the liner can cause ghost peaks. Replacing the inlet liner or quartz wool can help to alleviate this problem. Be sure you are also using the correct septum that is temperature-rated for your analysis.</p> |
|  | <p>Tailing Peaks</p> <p>Contamination or active compounds can cause tailing. Replacing your liner can help to reduce inlet contamination and activity; ensure you are routinely maintaining the inlet for the best performance. If using a liner with wool, try switching to a liner without wool if your sample allows.</p> |
|  | <p>Split Peaks</p> <p>When samples do not vaporize in a single location, split peaks may occur. Try using a liner with wool positioned in the middle of the liner to minimize split peaks.</p> |
|  | <p>Loss of Early Eluting Compounds</p> <p>Early eluting compounds can be lost in cases of high inlet activity or thermal decomposition. Using a single taper liner can help to alleviate this problem. If using wool, try switching to a liner without wool if your sample allows.</p> |

Tips To Improve Your Analysis




Loss of Late Eluting Compounds

Late eluting compounds can be lost if they have a high molecular weight, which can reduce their volatility and result in inlet discrimination. Using a liner with wool can mitigate this problem.




Reproducibility Issues

If your liner contains wool, the wool may shift or become dislodged due to injection or pressure changes, resulting in irreproducible results. Using a Zebron™ PLUS liner with a Z-Liner™ geometry will prevent this shifting and improve reproducibility.



Low Sensitivity / Low Response

If the injector temperature is too low, inlet discrimination may occur resulting in low responses for late eluting and less volatile compounds. Increasing the injection temperature or using an on-column injection with a direct connect liner will help boost sensitivity.



Offset or Noisy Baseline

Contamination in the liner can result in a noisy baseline. Replace your liner regularly to prevent excess baseline noise.



Want More Tips?

Download a FREE copy of the complete GC Troubleshooting Guide online at

www.phenomenex.com/TroubleGC

Frequently Asked Questions

Q1

Why does the inlet require so much troubleshooting?

A

The inlet is the most common place where problems happen in GC. This is because:

- The inlet is subjected to very high temperatures.
- There are multiple consumables involved in the inlet that need routine maintenance.
- Some inlet modes, such as splitless injection, result in a slow process which means that target analytes and the inlet are able to interact with each other for some time before the sample enters the column.

When dealing with active analytes, a highly deactivated inlet is necessary to minimize problems that can occur in the inlet, such as asymmetric peak shape and adsorption of target analytes. Zebron™ PLUS liners are rigorously deactivated and tested for activity to ensure reliable results when working with highly active compounds.

Q2

How can a Zebron PLUS liner with wool improve my GC column lifetime?

A

Quartz wool in Zebron PLUS liners serves two main purposes:

1. Wool promotes homogenous vaporization of analytes and solvent in the inlet, and helps effectively transfer them to the head of the GC column.
 2. Wool acts as a filter, trapping non-volatile impurities and preventing them from proceeding further into the GC column.
-

Q3

Can I pack or re-pack my liner with wool?

A

Liners should not be packed or re-packed. When attempts are made to remove wool from the liner or re-pack them with new wool, scratches in the deactivation layer can be created. When wool is forcefully inserted or removed, these scratches can expose active sites, causing peak tailing and poor reproducibility.

Frequently Asked Questions

Q4

My sample is dirty and it produces a residue on my liner after 100 injections. Can I scrub/sonicate my liner and reuse it?

A

No. It is not recommended to scrub/sonicate your liners. Residue formation indicates contamination. More specifically, this indicates that the sample contains non-volatile impurities. In such cases, it is recommended to replace, rather than clean, your liner. Sonication in solvent can alter the liner's deactivation and scrubbing can create scratches and active sites. Both will result in chromatographic problems including and not limited to poor quantitation, analyte breakdown, poor reproducibility, and peak asymmetry.

Q5

Does the color of a glass liner represent deactivation?

A

The color of the liner does not determine its inertness. In fact, varying colors can be easily added to liners by using metal salts during manufacturing. The inertness of the liner stems from the deactivation process that the liner undergoes, rather than its color alone.

Q6

How often should I change my Zebron™ PLUS liner?

A

It is important to periodically check and maintain your liner. The frequency at which a liner must be changed can be hard to predict as it is entirely dependent on the sample matrix. In general, if you are performing a headspace injection, only the vapors enter the GC inlet so the liner will remain clean for months. However, with a neat matrix injection, it is necessary to check the liner at least twice a week to make sure that it is free of residues. Once a visible residue is noticed, it is time to change the liner.


Q7

I am using an inert column for my analysis and I still see peak tailing. What could be the reason for this?

A

For analysis of active compounds, it is extremely important to use both a highly inert column and a highly inert liner. The liner is the first place of potential analyte interaction during GC analysis, and it is important that your analytes are not adsorbed. Zebron PLUS liners undergo a rigorous deactivation process and are tested for inertness to ensure reliable results when working with highly active compounds such as underivatized acids and active bases.

Three Decades of GC Expertise



Try the NEW GC Column Finder!


Easily select a column by part number, manufacturer, industry, application or official method in under 1 minute.

www.phenomenex.com/FindGC



Septa

Durable, long lasting septa with GuideRight™ holes for less bent needles
See pp. 6-7



ZB-1XT SimDist

Robust simulated distillation



ZB-SemiVolatiles

Supreme inertness for environmental SVOCs, PBDEs, and PAHs



ZB-5MS^{PLUS}

The next generation of inertness for a wide variety of applications

2009

2010

2011

2012

2014

2015

2016

ZB-Drug-1

For drugs of abuse



ZB-Bioethanol

Fast, 5 minute results



Zebtron™ Easy Seals™

Easy, washerless, leak-tight performance
See pp. 6-7

ZB-CLPesticides Phases

Versatility for chlorinated pesticides by GC/ECD



Zebtron PLUS Inlet Liners

Inert, easy install liners for superior sample protection





www.phenomenex.com/Liners

- Superior sample protection
- Easy to install
- A solution for all of your analyses

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Z-Liner, WAXPLUS, Inferno, MultiResidue, Easy Seals, 5MSPLUS, GuideRight, PhenoRed, and Zebtron are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific, Inc. PerkinElmer is a registered trademark of PerkinElmer, Inc. Shimadzu is a registered trademark of Shimadzu Corporation. Restek and Sky are registered trademarks of Restek Corporation. Septa BTO is a registered trademark of Chromatography Research Supplies, Inc.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Thermo Fisher Scientific, PerkinElmer, or Restek Corporation.

© 2016 Phenomenex, Inc. All rights reserved.