


# APPLICATIONS

## Coupling Miniaturized SPE with bioZen™ XB-C18 Nano Core-Shell Technology for the Detection of SARS-CoV-2 Nucleocapsid Protein Viral Peptides in Nasopharyngeal Samples

Roxana Eggleston-Rangel<sup>1</sup>, Shahana Huq<sup>1</sup>, Dr. Jason Anspach<sup>1</sup>, Dr. Bryan Tackett<sup>1</sup>, Anthony Maus<sup>2</sup>, Ph. D., Jennifer Kemp<sup>2</sup>, Kari Gurtner<sup>2</sup>, Ravinder Singh<sup>2</sup>, Ph. D., and Stefan Grebe<sup>2</sup>, M.D., Ph. D.  
<sup>1</sup>Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA  
<sup>2</sup>Division of Clinical Biochemistry and Immunology, Department of Laboratory Medicine and Pathology, Mayo Clinic, 200 First St. SW Rochester, MN 55905

### Introduction

Amplification of genetic material by PCR (Polymerase Chain Reaction) has been the gold standard for viral detection in complex sample matrices. Although PCR methods for viral detection are well established, scientists still face many challenges such as dealing with sub-optimal sensitivity due to rapid viral mutation rate, low throughput due to long experimental steps and the accessibility of specialized reagents.<sup>1,2</sup> LC-MS methods have emerged as either complementary or as replacements for viral detection in complex matrices. LC-MS viral detection relies on the presence of peptides that are specific to the virus and not the virus host. For an LC-MS method to achieve the sensitivity needed for the detection of viral peptides, the sample cleanup method and the choice of LC column employed are key aspects of the task. Sample cleanup enhances the signal to noise ratio, which enables trace level detection of the viral peptides. SPE (Solid Phase Extraction) is the most effective sample extraction technique to help eliminate the interferences associated with sample matrix components, while selectively retaining the peptides of interest. Similarly, achieving the desired sensitivity is also dependent on the LC column format utilized. Nano LC has been long known to enhance sensitivity and allows for the user to scale down the chromatographic separation by decreasing the column inner diameter, resulting in a decrease in flow rate, an increase in analyte concentration and thus an increase in sensitivity.<sup>3-5</sup> While nano LC can increase sensitivity, the combination of narrow column inner diameters and Phenomenex bioZen core-shell technology results in minimal band broadening, a reduced diffusion path, an increase in peak height and thus an increase in sensitivity compared to fully porous column technology.<sup>6</sup> Microelution SPE can successfully be coupled to nano LC-MS workflows. Microelution is a miniaturized SPE format which is ideal when working with small sample volumes, low analyte concentration, and the need for an increase in analyte signal.<sup>7</sup> Thus, this technical note provides a miniaturized LC-MS and SPE method for the reliable and accurate detection of viral peptides at 1 fmol/μL using SARS-CoV-2 nucleocapsid protein (N-protein) as a model system, using Strata-X microelution SPE for sample cleanup and a bioZen 2.6 μm Peptide XB-C18 nano column coupled to a Sciex QTRAP 6500+ mass spectrometer for maximum sensitivity.



**Roxana Eggleston-Rangel**  
Application Scientist

Roxana likes to spend time with her dog and family. She has a German Shepherd mix named Cobi who was named after the official mascot of the 1992 Summer Olympics and not the basketball player. Besides the lab, you might find Roxana in old episodes of the Ghost Whisperer, MAD TV, The L word and others as she used to be a TV extra during her school years.

### Sample Preparation

#### Solid Phase Extraction

**Cartridge:** Strata®-X 96-Well Plate, 2 mg/well  
**Part No.:** [8M-S100-4GA](#)  
**Condition:** 200 μL Methanol  
**Equilibrate:** 200 μL Water  
**Load:** 400 μL diluted nasal matrix (200 μL sample with 200 μL 4% Phosphoric acid in water spiked with heavy peptide)  
**Wash 1:** 200 μL 2% Formic acid in water  
**Wash 2:** 200 μL Water  
**Elute:** 2x 40 μL 1% TFA in ACN:Water (75:25)  
**Dry:** Under N<sub>2</sub> at ≤ 40° C  
**Reconstitute:** Resuspend dried samples in 100 μL 0.1% Formic acid in water

### LC Conditions

**Column:** bioZen 2.6 μm Peptide XB-C18  
**Dimensions:** 150 x 0.075 mm  
**Part Number:** [00F-4768-AW-21](#)  
**Pressure (bar):** 200  
**Mobile Phase:** A: 0.1% Formic Acid in Water  
 B: 0.1% Formic Acid in Acetonitrile

Gradient	Time (min)	% B
	0	3
	30	40
	35	70
	40	80
	41	3
	50	3

**Flow Rate:** 250 nL/min  
**Temperature:** 25 °C  
**LC System:** NanoLC™ 425 (SCIEX)  
**Detection:** MS/MS  
**Detector:** 6500 QTRAP® (SCIEX)  
**Injection Volume:** 1 μL

### MS Parameters

**CUR:** 30  
**IS:** 2500  
**GS1:** 35  
**GS2:** A: 0  
**IHT:** 150  
**CAD:** Low  
**EP:** 10

### Materials and Methods

Nasopharyngeal swab samples were collected in PBS buffer followed by tryptic digest. Digested samples were spiked with synthetic peptides (**Table 1**) followed by serial dilution. The synthetic isotopically labeled peptides were labeled at the C-terminal of arginine (R) or lysine (K) (**Table 1**) and were spiked for quantification purposes. Microelution SPE purification was performed as shown.

# APPLICATIONS

**Table 1.**  
Multiple Reaction Monitoring (MRM) optimized transitions parameters.

ID	Q1 Mass (Da)	Q3 Mass (Da)	RT (min)	sMRM Dwell (ms)	DP	CE	CXP
GQGVPIINTSSPDDQIGYYR+3y16+2.heavy	731.015	925.433	17.1	208.871	70	31	50
GQGVPIINTSSPDDQIGYYR+3y16+2.light	727.678	920.429	17.1	208.871	70	31	50
GQGVPIINTSSPDDQIGYYR+3y9+2.light	727.678	563.762	17.1	208.871	70	31	50
NPANNAIIVLQLPQGTTLPK+3y8.heavy	690.06	849.492	20.9	224.491	80	40	40
NPANNAIIVLQLPQGTTLPK+3y8.light	687.388	841.478	20.9	224.491	80	40	40
NPANNAIIVLQLPQGTTLPK+3b9.light	687.388	865.453	20.9	224.491	80	40	40
DGIWVATEGALNTPK+3y7.heavy	564.973	708.413	22.3	155.65	80	36	10
DGIWVATEGALNTPK+3y7.light	687.388	865.453	20.9	224.491	80	40	40
DGIWVATEGALNTPK. +3y5.light	562.302	572.34	22.3	155.65	80	35	18
AYNVTQAFGR +2y6. heavy	569	689.4	16.3	209.429	50	24	12
AYNVTQAFGR +2y6. light	564	679.4	16	250	50	24	12
AYNVTQAFGR +2y6. light	564	892.6	16.3	209.429	50	26	15
QQTVTLLPAADLDDFSK+2y10.heavy	935.6	1086.8	22	129.08	80	40	20
QQTVTLLPAADLDDFSK+2y10.light	931.6	1078.5	22	129.08	80	40	20
QQTVTLLPAADLDDFSK+2y11.light	931.6	1191.6	22	129.08	70	37	20

## Data Processing


Sciex OS version (2.0.0) and Skyline (version 20.1.0.31) software were used for quantitative and qualitative purposes, respectively.

## Results and Discussion

Trypsin digested Nasopharyngeal swab samples (containing five spiked synthetic SARS-CoV-2 Nucleocapsid Protein specific peptides as in **Table 1**) were spiked with their peptide counterpart labeled at the C-terminal of arginine (R) or lysine (K) (**Table 1**) to allow for quantification using the peak area ratio of the heavy labeled peptide to the unlabeled (light) peptide. The chemistry of the bioZen™ 2.6 µm Peptide XB-C18 allowed for adequate retention of the viral peptides, providing great resolution among the five analytes. Furthermore, the nano dimensionality of the column together with miniaturized SPE resulted in an ultrapure clean extract with minimum background noise to yield a significantly lower level of quantitation such as 1 fmol/µL.

In order to determine the dynamic range of the assay, a calibration curve constructed for the extracted samples over a concentration range of 1 to 200 fmol/µL was analyzed. All five peptides showed linearity across the entire concentration range (**Figure 3**). This demonstrates that each peptide can be reliably quantitated down to a concentration of at least 1 fmol/µL. The signal to noise (S/N) ratios were calculated at the lowest concentration tested at 1 fmol/µL and were based upon the relative response against a blank matrix extract. S/N ratios ranging 112 to 1293 (**Table 2**) were observed. QC samples were tested for replicate extraction at 3 different levels of concentration (**Table 3**). Accuracy and CV were calculated and ranged between 93-118% and <7%, respectively, which are within acceptable industry standards.

**Figure 1.**  
Chromatographic separation of the 5 SARS-CoV-2 synthetic peptides.




App ID 26169


# APPLICATIONS

**Figure 2.**  
Extracted ion chromatograms of the individual MRM transitions.


**GGQVPINTNSSPDDQIGYYR**




**NPANNAIVLQLPQGTTLPK**




**DGIWVATEGALNTPK**



**AYNVTQAFGR**




**QQTVTLLPAADLDDFSK**




**Figure 3.**  
Calibration curves for each peptide.


**QGVPIINTNSSPDDQIGYYR**




**NPANNAIVLQLPQGTTLPK**




**DGIWVATEGALNTPK**



**AYNVTQAFGR**



**QQTVTLLPAADLDDFSK**



**Table 2.**  
Signal to noise ratio for individual peptides at 1 fmol/ $\mu$ L.

Peptide ID	S/N
AYNVTQAFGR+2y6.light	408
DGIWVATEGALNTPK+3y7.light	460
QQTVTLLPAADLDDFSK+2y10.light	1293
NPANNAIVLQLPQGTTLPK+3y8.light	323
QGVPIINTNSSPDDQIGYYR+3y16+2.light	112

# APPLICATIONS

**Table 3.**  
Precision and accuracy for QC samples of each peptide.

## AYNVTQAFGR

Expected Conc. (fmol/ $\mu$ L)	Sample	Replicates (N)	% CV	Accuracy
5	QC1-1	4	3.7	96
40	QC1-2	4	1.2	102
100	QC1-3	4	2.1	96

## DGIIVVATEGALNTPK

Expected Conc. (fmol/ $\mu$ L)	Sample	Replicates (N)	% CV	Accuracy
5	QC1-1	4	4	118
40	QC1-2	4	3.9	115
100	QC1-3	4	1.2	113

## GQGVPINTNSSPDDQIGYYR

Expected Conc. (fmol/ $\mu$ L)	Sample	Replicates (N)	% CV	Accuracy
5	QC1-1	4	0.6	110
40	QC1-2	4	6	107
100	QC1-3	4	4	93

## QQTVTLLPAADLDDFSK

Expected Conc. (fmol/ $\mu$ L)	Sample	Replicates (N)	% CV	Accuracy
5	QC1-1	4	7	104
40	QC1-2	4	5	102
100	QC1-3	4	5	110

## Conclusion

The bioZen™ Peptide XB-C18 nano core-shell column chemistry retains the SARS-CoV-2 viral peptides adequately, providing great resolution between these analytes. By employing Strata®-X microelution SPE, samples were extracted more efficiently and in a much shorter amount of time than conventional SPE. This selective extraction and sensitive chromatography yielded signal/noise ratios  $\geq 112$  at 1 fmol/ $\mu$ L. The selectivity at the lowest concentration tested of 1 fmol/ $\mu$ L, suggests that the detection range for these peptides can easily be expanded to a much lower concentration range to obtain the lowest level of quantitation. The dynamic range of microelution SPE tested demonstrates an acceptable linearity with  $r^2$  values  $\geq 0.995$  for all peptides. Also, industry standards call for an accuracy of 80-120% and a precision of <20% at LLOQ and <15% across the calibration curve, which were well met by the presented methodology. The data presented here demonstrates that using microelution SPE in conjunction with bioZen 2.6  $\mu$ m Peptide XB-C18 nano core-shell columns, peptides from the COVID-19 causing SARS-CoV-2 virus can be detected with great resolution and sensitivity.

## References

- Orooji, Yasin et al. "An Overview on SARS-CoV-2 (COVID-19) and Other Human Coronaviruses and Their Detection Capability via Amplification Assay, Chemical Sensing, Biosensing, Immunosensing, and Clinical Assays." *Nano-micro letters* vol. 13,1 (2021): 18. doi:10.1007/s40820-020-00533-y .
- Islam, Khursheed UI, and Jawed Iqbal. "An Update on Molecular Diagnostics for COVID-19." *Frontiers in cellular and infection microbiology* vol. 10 560616. 10 Nov. 2020, doi:10.3389/fcimb.2020.560616
- Wilm, M. & Mann, M. Analytical properties of the nano-electrospray ion source. *Anal. Chem.* 68, 1–8 (1996).
- Gatlin, C. L., Kleemann, G. R., Hays, L. G., Link, A. J. & Yates, J. R. III Protein identification at the low femtomole level from silver-stained gels using a new fritless electrospray interface for liquid chromatography–microspray and nanospray mass spectrometry. *Anal. Biochem.* 263, 93–101 (1998).
- Bian, Y., Zheng, R., Bayer, F.P. et al. Robust, reproducible and quantitative analysis of thousands of proteomes by micro-flow LC–MS/MS. *Nat Commun* 11, 157 (2020)
- <https://phenomenex.blog/2017/07/25/porous-core-shell-tech/>
- Justyna Plotka-Wasyłka, Natalia Szczepańska, Miguel de la Guardia, Jacek Namieśnik, Miniaturized solid-phase extraction techniques, *TrAC Trends in Analytical Chemistry*, Volume 73, 2015, Pages 19–38 (<http://www.sciencedirect.com/science/article/pii/S0165993615002010>)

## Ordering Information

### bioZen™ Nano LC Columns with Integrated SecurityLINK™ Fitting

Phases	150 x 0.075	250 x 0.075	500 x 0.075
bioZen 3 µm Peptide PS-C18	<a href="#">00F-4771-AW-21</a>	<a href="#">00G-4771-AW-21</a>	–
bioZen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-21</a>	<a href="#">00G-4768-AW-21</a>	–
bioZen 3 µm Peptide Polar-C18	<a href="#">00F-4782-AW-21</a>	<a href="#">00G-4782-AW-21</a>	–
bioZen 5 µm Peptide XB-C18	–	–	<a href="#">00J-4792-AW-11</a>

### bioZen Nano LC Columns with Open Fused-Silica Inlet/Outlet Fitting

Phases	150 x 0.075	250 x 0.075	500 x 0.075
bioZen 3 µm Peptide PS-C18	<a href="#">00F-4771-AW-11</a>	<a href="#">00G-4771-AW-11</a>	–
bioZen 2.6 µm Peptide XB-C18	<a href="#">00F-4768-AW-11</a>	<a href="#">00G-4768-AW-11</a>	–
bioZen 3 µm Peptide Polar-C18	<a href="#">00F-4782-AW-11</a>	<a href="#">00G-4782-AW-11</a>	–
bioZen 5 µm Peptide XB-C18	–	–	<a href="#">00J-4792-AW-21</a>

### Trap Columns

Trap Columns	10 x 0.075 mm	Unit
RP-1	05N-4252-AW	3/pk
RP-2	05N-4754-AW	3/pk

### Trap Fittings

Trap Fittings	Part No.	Description	Unit
	<a href="#">A00-7602</a>	PEEKlok™ fittings with 6-40 thread for 1/32" OD tubing (2 x fittings, 6 x ferrules and 1 x tightening tool)	ea
	<a href="#">A00-7603</a>	PEEKlok fittings with 6-32 thread for 1/32" OD tubing (2 x fittings, 6 x ferrules and 1 x tightening tool)	ea
	<a href="#">A00-7600</a>	PEEKlok fittings with 10-32 thread for 1/32" OD tubing with low profile hex head (2 x fittings, 6 x ferrules and 1 x wrench)	ea

### Strata®-X Microelution 96-Well Plates

96-Well Plates (ea)	
Phase	2 mg
Strata X-AW	<a href="#">8M-S038-4GA</a>
Strata X-A	<a href="#">8M-S123-4GA</a>
Strata-X	<a href="#">8M-S100-4GA</a>
Strata-X-C	<a href="#">8M-S029-4GA</a>
Strata-X-CW	<a href="#">8M-S035-4GA</a>

### Strata-X Microelution Method Development 96-Well Plates

Part No.	Description	Unit
<a href="#">KS0-9528</a>	<b>Strata-X Peptide Screening</b> Strata-X-CW 2 mg/well (6 rows) Strata-X-A 2 mg/well (6 rows)	ea
<a href="#">KS0-9529</a>	<b>Strata-X Method Development</b> Strata-X-C 2 mg/well (3 rows) Strata-X-AW 2 mg/well (3 rows) Strata-X-CW 2 mg/well (3 rows) Strata-X-A 2 mg/well (3 rows)	ea

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Australia**

t: +61 (0)2-9428-6444  
 auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
 anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
 t: +32 (0)2 511 8666 (Dutch)  
 beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
 info@phenomenex.com

**China**

t: +86 400-606-8099  
 cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
 cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
 nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
 nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
 franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
 anfrage@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
 indiainfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
 eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
 italiainfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
 nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
 tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
 nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
 nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
 nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
 pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
 ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
 sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
 sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
 espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
 nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
 swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
 twinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
 ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
 info@phenomenex.com

☎ **All other countries/regions**  
**Corporate Office USA**  
 t: +1 (310) 212-0555  
 info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**Terms and Conditions**

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Strata is a registered trademark and bioZen, and SecurityLINK are trademarks of Phenomenex. PeekLok is a trademark of Trajan Scientific Australia Pty Ltd. QTRAP is a registered trademark and nanoLC is a trademark of AB SCIEX Pte Ltd. AB SCIEX is being used under license.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.