



TN-1253

## Subunit Analysis of Monoclonal Antibodies Using Biozen™ 3.6 μm Intact XB-C8

Helen Whitby, PhD<sup>2</sup>, M Christina Malinao<sup>1</sup>, Brian Rivera<sup>1</sup>, Chad Eichman, PhD<sup>1</sup>, and Bryan Tackett, PhD<sup>1</sup>

<sup>1</sup>Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

<sup>2</sup>Phenomenex Ltd., Queens Avenue, Hurdsfield Ind. Est., Macclesfield, Cheshire SK10 2BN UK

### Introduction

Monoclonal antibodies (mAbs) are the most common protein therapeutic with an estimated 180 or more mAb drugs on the global market by 2025. Fragmentation of monoclonal antibodies is a Critical Quality Attribute (CQA) which must be measured for all biotherapeutics to assess the purity and integrity of the protein. Fragmentation can occur during expression but also upon storage or metabolism *in vivo*. The fragmentation pattern of a mAb can offer information on its stability making this an essential part of the characterization process. Fragments are formed at sites where the peptide bonds are prone to cleavage and although the factors which influence which sites are most susceptible can be complex, some general fragmentation patterns are typically observed. In this technical note, we employ a forced reduction approach to identify and report scFc and F(ab')<sub>2</sub> profiles for two common mAbs through an IdeS digestion. We then look at the fragmentation patterns of a Dithiothreitol (DTT) reduction, which provides separated heavy chain (HC) and light chains (LC).

Impurity profiling and characterization of intact biological subunits present many challenges because of a need to identify small differences between variants. Biozen 3.6 μm Intact XB-C8 and C4 columns are packed with specialized wide pore core-shell particles which facilitate faster mass transfer kinetics for large biomolecules, providing narrower, taller peaks and greater resolving power when analyzing HC / LC, Fc / Fab regions or isoforms, compared to other commercially available fully porous silica-based column solutions designed for mAb analysis. The introduction of high efficiency UHPLC instruments offering low system dispersion coupled with highly efficient core-shell particles has offered significant improvements in the resolution and characterization of antibody fragments.

### Sample Preparation

Cetuximab was treated with IdeS enzyme (1 μL per 50 μg mAb) and incubated for 60 min at 37 °C. Rituximab was treated with 10 mM DTT at pH 5.5 for 30 min.

### LC Conditions

<b>Columns:</b>	Biozen 3.6 μm Intact XB-C8
	AdvanceBio RP-mAb 3.5 μm SB-C8
<b>Dimensions:</b>	150 x 2.1 mm
<b>Part No.:</b>	<a href="#">00F-4766-AN</a>
<b>Mobile Phase:</b>	A: 0.1 % Trifluoroacetic Acid in Water
	B: 0.1 % Trifluoroacetic Acid in Acetonitrile
<b>Gradient:</b>	<b>Time (min)</b> <b>%B</b>
	0                      40
	1                      40
	13                     60
<b>Flow Rate:</b>	0.5 mL/min
<b>Injection Volume:</b>	1 μL
<b>Temperature:</b>	80 °C
<b>Detection:</b>	UV @ 280 nm

### Results and Discussion

For the IdeS digestion of Cetuximab, the F(ab')<sub>2</sub> was identified along with minor amounts of scFc using a Biozen Intact XB-C8 column (**Figure 1**). Under standard reversed phase gradient conditions, we observed excellent peak shape of the F(ab')<sub>2</sub>. Analysis of this sample under the same running conditions with the AdvanceBio® RP-mAb SB-C8 column was unable to detect any scFc fragment. Significantly less tailing was also observed in the F(ab')<sub>2</sub> fragments with the Biozen Intact XB-C8 column, which allows for improved integration and more accurate results. Notably, both applications were run at an elevated temperature to improve the peak shape of the fragments.

Analysis of both HC and LC variants was also assessed using a DTT reduction of Rituximab (**Figure 2**). After employing the same method conditions, superior separation was observed for the HC variants on the Biozen Intact XB-C8 column. Moreover, this column was able to detect minor amounts of post-translational modifications (HC variants). When using the Advance Bio RP-mAb SB-C8 column the HC variants were not fully resolved.

**Figure 1.** Subunit Analysis of an IdeS Digestion of Cetuximab.





Figure 2. HC / LC Analysis of DTT-reduced Rituximab.



**Conclusions**

Both chemistries evaluated (Biozen Intact XB-C8 and AdvanceBio RP-mAb SB-C8) are bonded with a C8 ligand. However, the combination of the core-shell particle and the sterically hindered XB-C8 ligand offer better peak shape compared with a fully porous C8 phase. The superior peak shape and resolution observed with the Biozen Intact XB-C8 column was present in both the IdeS digestion of Cetuximab and the DTT reduction of Rituximab. Optimizing the performance of this subunit analysis is pivotal as it serves as a critical quality attribute that needs to be monitored to assess the integrity of the mAb. The wide pore core-shell particle morphology of the Biozen Intact XB-C8 column allows for faster mass transfer kinetics of antibodies in and out of the stationary phase, which reduces the diffusion pathway and improves chromatographic efficiency. The benefits of producing more efficient, narrower peaks is seen with improved resolution for subunit analysis.

**Biozen Ordering Information**

Biozen Columns (mm)								Biocompatible Guard Cartridges		
	50 x 2.1	100 x 2.1	150 x 2.1	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	for 2.1 mm	for 4.6 mm	Holder
								/3pk	/3pk	ea
Biozen 2.6 µm WidePore C4	<a href="#">00B-4786-AN</a>	<a href="#">00D-4786-AN</a>	<a href="#">00F-4786-AN</a>	<a href="#">00B-4786-E0</a>	<a href="#">00D-4786-E0</a>	<a href="#">00F-4786-E0</a>	<a href="#">00G-4786-E0</a>	<a href="#">AJ0-9816</a>	<a href="#">AJ0-9818</a>	<a href="#">AJ0-9000</a>
Biozen 3.6 µm Intact XB-C8	<a href="#">00B-4766-AN</a>	<a href="#">00D-4766-AN</a>	<a href="#">00F-4766-AN</a>	<a href="#">00B-4766-E0</a>	—	<a href="#">00F-4766-E0</a>	—	<a href="#">AJ0-9812</a>	<a href="#">AJ0-9814</a>	<a href="#">AJ0-9000</a>



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/phx-terms-and-conditions-of-sale](http://www.phenomenex.com/phx-terms-and-conditions-of-sale).

**Trademarks**

Biozen and BE-HAPPY are trademarks of Phenomenex. AdvanceBio is a registered trademark of Agilent Technologies, Inc.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

