

The Effect of Column Hardware on the Analysis of Synthetic Oligonucleotides by LC-MS

Ivan Lebedev, PhD¹, Brian Rivera¹, Helen Whitby, PhD², and Bryan Tackett, PhD¹

¹Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

²Phenomenex Ltd., Queens Ave., Hurdsfield Ind. Est., Macclesfield, SK10 2BN, UK

Overview


Therapeutic oligonucleotides represent a recent breakthrough in the pharmaceutical industry. Improved reliability to synthetic routes coupled with new drug delivery technologies has given a heightened level of interest in these clinically interesting targets.

Chromatographically, they present a challenge due to their polar nature and phosphodiester backbone. As such, ion-pair reversed phase chromatography (IP-RPLC) is necessary, typically using a mobile phase consisting of an alkylamine ion pair and a perfluorinated alcohol such as Hexafluoro-2-propanol.

Although optimizing the IP-RPLC conditions can improve chromatographic separation, this does not address oligonucleotide propensity to chelate with trace metal in the stainless steel of traditional HPLC column hardware. These non-specific interactions are often associated with analyte loss, poor peak shape, and irreproducibility. In this application note, we demonstrate the positive effect BioTi, a multi-layered titanium tube which provides a bio-inert flow path, has on chromatography and electrospray ionization of two synthetic oligos compared to traditional stainless steel LC columns.

When analyzing a double-stranded RNA by traditional column hardware, presumed nonspecific interactions are observed (**Figure 1**). This results in a bimodal peak for the later eluting anti-sense strand, with full scan and MS/MS data confirming this to be the same analyte. Conversely, the use of bio-inert hardware drastically improves resolution of the sense and anti-sense strand peaks. MS data in **Figure 2** confirms that BioTi hardware mitigates the effects of non-specific binding and subsequently addresses the bimodal chromatographic artifact. We also observe greater sensitivity when run on the BioTi hardware, showing its utility for oligonucleotide characterization using relatively low concentrations of Hexafluoro-2-propanol and ion-pairing reagents.

This improvement in ionization is also observed with single stranded nucleic acids. The analysis of Nusinersen shows a marked improvement with a full order of magnitude in sensitivity when the method was run on a Biozen Oligo column (packed BioTi) compared with a stainless steel column packed to a similar packing efficiency, and with the same batch of media (**Figure 3**). In summary, by minimizing non-specific interactions, the use of BioTi hardware not only serves to improve chromatographic performance and reproducibility of oligonucleotide separations, but also maximizes analyte recovery and sensitivity enabling both consistent quantitation and characterization.



LC Conditions – Figure 1 and 2

Columns: Biozen 2.6 µm Oligo (BioTi™)
Clarity™ 2.6 µm Oligo-XT (stainless steel)

Dimensions: 100 x 2.1 mm

Part No.: [00D-4790-AN](#) (Biozen)
[00D-4746-AN](#) (Clarity)

Mobile Phase: A: 4 mM Triethylamine in Water + 12.5 mM Hexafluoro-2-propanol
B: 4 mM Triethylamine in Methanol + 12.5 mM Hexafluoro-2-propanol

Gradient:	Time (min)	%B
	0	5
	2	5
	16	30
	16.1	95
	20	95
	20.1	5

Flow Rate: 0.3 mL/min

Injection Volume: 2 µL (12.5 ng)

Temperature: 55 °C

Instrument: Shimadzu® LC-20A Prominence®

Detection: TOF-MS

Detector: SCIEX® TripleTOF® 6600

LC Conditions – Figure 3

Columns: Biozen 2.6 µm Oligo (BioTi)
Clarity 2.6 µm Oligo-XT (stainless steel)

Dimensions: 100 x 2.1 mm

Part No.: [00D-4790-AN](#) (Biozen)
[00D-4746-AN](#) (Clarity)

Mobile Phase: A: 10 mM Hexylamine in Water + 12.5 mM Hexafluoro-2-propanol
B: 10 mM Hexylamine in Methanol + 12.5 mM Hexafluoro-2-propanol

Gradient:	Time (min)	%B
	0	25
	2	25
	16	75
	16.1	95
	20	95
	20.1	25

Flow Rate: 0.3 mL/min

Injection Volume: 2 µL (12.5 ng)

Temperature: 55 °C

Instrument: Shimadzu LC-20A Prominence

Detection: TOF-MS

Detector: SCIEX TripleTOF 6600



Figure 1. siRNA on a Clarity™ 2.6 µm Oligo-XT Column (Standard Stainless Steel Hardware).




Figure 2. siRNA on a Biozen™ 2.6 µm Oligo Column (BioTi™ Bio-inert Hardware).





Figure 3. Nusinersen on Biozen 2.6 µm Oligo Column (Bio-inert) Compared to Clarity 2.6 µm Oligo-XT Column (Stainless Steel).



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

**🌐 All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Biozen, BioTi, Clarity, and BE-HAPPY are trademarks of Phenomenex. Shimadzu and Prominence are registered trademarks of Shimadzu Corporation. SCIEX and TripleTOF are registered trademarks of AB SCIEX Pte. Ltd.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Shimadzu Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

