



Separation of Quercetin and its Organic Impurities per USP Monograph

Alyssa Nodland, Zandra Baja, and Bryan Tackett, PhD
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview


Quercetin is a potent antioxidant flavonoid found in fruits and vegetables. The development of a quick and efficient analysis of Quercetin and its related organic impurities is of interest for generic drug manufacturers. In this application note, we report the separation of Quercetin and its related organic impurities using a Luna Omega 5 μ m C18 compared with a Inertsil 5 μ m ODS-3 column according to the USP monograph for Quercetin, which specifies a column containing L1 (C18) packing and dimensions of 250 x 4.6 mm.

System suitability per USP Monograph for Quercetin is the same for both Assay and Organic Impurities: Resolution no less than (NLT) 1.5 between Kaempferol and Isorhamnetin, a column efficiency NLT 2000 theoretical plates, and a percent relative standard deviation (%RSD) of no more than (NMT) 2.0% for Quercetin. All system suitability requirements for Quercetin Assay and Organic Impurities were met by both the Luna Omega and Inertsil columns. Other columns, including a Luna 5 μ m C18(2) column, a InertClone™ 5 μ m ODS-3 column, and a Prodigy™ 5 μ m ODS-3 column were also tested but did not meet acceptance criteria for impurity content according to the USP monograph for Quercetin.

It is important to note that Retention Time %RSD is included because there was variability between injections across all columns tested. The column temperature was not specified in the USP monograph, so it was assumed to be ambient. To maintain a reasonably constant column temperature during sample analysis, the HPLC column was placed inside the column oven with the oven door closed and the oven turned off.

All solutions were prepared as indicated in the USP Monograph for Quercetin. USP Quercetin RS (Catalog No. 1592409), USP Kaempferol RS (Catalog No. 1594900), and USP Isorhamnetin RS (Catalog No. 1591800) were purchased from USP.

Figure 1. Quercetin



LC-UV Conditions

Column: Luna™ Omega 5 μ m C18 (00G-4785-E0)
Inertsil® 5 μ m ODS-3

Dimensions: 250 x 4.6 mm

Mobile Phase: Methanol / Water / Phosphoric Acid
(100:100:1, v/v/v)

Flow Rate: 1.0 mL/min (Isocratic)

Injection Volume: 20 μ L

Temperature: Ambient

Detector: UV @ 370 nm (Assay) or 270 (Organic Impurities)

LC System: Agilent® 1260

Table 1. Preparation of Solutions

Solution	Composition
System Suitability Solution – Assay and Organic Impurities	0.02 mg/mL each of USP Quercetin RS, USP Kaempferol RS, and USP Isorhamnetin RS in Methanol.
Standard Solution – Assay and Organic Impurities	Transfer 5 mg of USP Quercetin RS to a 25-mL volumetric flask and add 10 mL of Methanol to dissolve. Add 10 mL of water, mix, and dilute with Methanol to volume.



Figure 2. Standard Solution – Assay and Organic Impurities







Figure 3. System Suitability Solution – Assay and Organic Impurities



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sales.

Trademarks

Luna, InertClone, Prodigy, and BE-HAPPY are trademarks of Phenomenex. Inertsil is a registered trademark of GL Sciences, Inc. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with GL Sciences, Inc. or Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

