


## Chiral Separation of the Drug Product Lacosamide on a Lux® 5 µm Amylose-1 Column According to Ph. Eur. Monograph 2992

Zeshan Aqeel, Dr. Bryan Tackett, and Anania Assefa Amenu  
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501-1430, USA



### Overview


Lacosamide is a chiral drug which is commercialized as pure enantiomer. It is an antiepileptic drug approved in the USA and Europe. The current trend is to commercialize drugs formulated in their enantiopure form, which means containing only the enantiomer responsible for the pharmaceutical activity. The S-enantiomer of Lacosamide, designated as Impurity A, is pharmaceutically inactive and therefore undesirable to have in the presence of the active R-Lacosamide.

In this application note, we report the enantiomeric separation between Lacosamide and its enantiomer (Impurity A) using the Lux Amylose-1 chiral column according to the Ph. Eur. Monograph 2992.

System suitability per Ph. Eur. Monograph 2992 for enantiomeric purity of Lacosamide is for resolution of not less than (NLT) 3.0. The separation of the enantiomers was achieved and peaks corresponding to the Impurity A and Lacosamide had a resolution of 4.7, exceeding the system suitability requirement (**Figure 2**). The higher resolution indicates the powerful chiral recognition ability of Lux 5 µm Amylose-1 chiral stationary phase column.

All reference solutions were prepared as indicated in Ph. Eur. Monograph 2992 for Lacosamide. Lacosamide Certified Reference Standard (CRS), (Catalog No Y0001982) and Lacosamide Impurity A CRS (catalog no. Y0001972) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) Council of Europe.

**Figure 1.**  
Lacosamide



### LC-UV Conditions


**Column:** Lux 5 µm Amylose-1  
**Dimensions:** 150 x 4.6 mm  
**Part No.:** [00F-4732-E0](#)  
**Mobile Phase:** Water, 2-Propanol, Heptane (3:100:900, v/v/v)  
**Pressure (bar):** 81  
**Flow Rate:** 1 mL/min  
**Injection:** 20 µL  
**Temperature:** 25 °C  
**Detector:** UV @ 215 nm  
**System:** Agilent® 1260 Binary UHPLC

**Table 1.** Preparation of Test and Reference Solutions

Solution	Composition
Test Solution	Dissolve 50.0 mg of the substance to be examined in the mobile phase (MP) and dilute to 50.0 mL with MP
Reference Solution (a)	Dissolve 1 mg/mL of Lacosamide Impurity A CRS in MP and dilute to 10.0 mL with MP.
Reference Solution (b)	Dissolve 20 mg of the substance to be examined in MP, add 1.0 mL of reference solution (a) and dilute to 20.0 mL with MP.
Reference Solution (c)	Dilute 1.0 mL of the test solution to 100.0 mL with MP. Dilute 1.0 mL of this solution to 20.0 mL with MP.




**Figure 2.** System Suitability Solution (Reference Solution (b))



App ID 26580

Peak	Time	Area	Height	Width	Resolution
Impurity A	16.5	225	5.5	0.6	4.7
Lacosamide	23.1	37231	556.6	99.4	

**Figure 3.** Reference Solution (c): 0.5 µg/mL of Lacosamide CRS



App ID 26581



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

**Australia**

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**  
t: +1 (310) 212-0555  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Lux is a registered trademark and BE-HAPPY is a trademark of Phenomenex.

Agilent is a registered trademark of Agilent Technologies, Inc.

**Disclaimer**

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

