

# CASE STUDY

## USP Miconazole Nitrate Assay and Organic Impurities

Phil Koerner, PhD and Bryan Tackett, PhD  
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

### Introduction

Miconazole Nitrate belongs to the imidazole group of antifungal agents. The development of a quick and efficient analysis of Miconazole Nitrate and its related organic impurities is of interest for generic drug manufacturers. The USP monograph for Miconazole Nitrate Assay and Organic Impurities, which became official on May 1<sup>st</sup>, 2020, uses the Kinetex™ 2.6 µm, 100 x 4.6 mm Phenyl-Hexyl column. This method was revised by USP as part of the USP monograph modernization efforts and subject to the typical USP review and comment process before becoming official. The monograph revision for Miconazole Nitrate replaced both the original titration method (Assay) and the original TLC procedure (Organic Impurities) with HPLC. It is worth noting that, while the official USP monograph does not specify a column by name (simply indicated as “L11, 2.6 µm, 100 x 4.6 mm”), the draft which was published in USP-PF 41(4) in 2015 does indicate that the Kinetex 2.6 µm, 100 x 4.6 mm Phenyl-Hexyl column was the column that was used in the elaboration of the monograph.

Since the USP does not publish chromatograms in the official monograph methods, we wanted to demonstrate in this case study that the Kinetex Phenyl-Hexyl column gave the desired results and to share the chromatographic results with those interested in this USP monograph method for Miconazole Nitrate.

System suitability per USP Monograph for the Miconazole Nitrate Assay is a resolution no less than (NLT) 1.5 between Miconazole Nitrate and Miconazole Related Compound F, a symmetry factor no more than (NMT) 2.0 for Miconazole Nitrate, and a percent relative standard deviation (%RSD) of NMT 0.73 % for Miconazole Nitrate. All system suitability requirements for Miconazole Nitrate Assay were met using the Kinetex Phenyl-Hexyl column (Figure 1 and 2).

### Key Concepts:


- Applicability of Core-shell Phenyl-Hexyl for USP Miconazole Nitrate Assay and Organic Impurities.
- Weak analyte absorbance at low UV and concentration differences between sample/solvent can lead to baseline shifts.

System suitability per USP Monograph for the Miconazole Nitrate Organic Impurities is a resolution NLT 1.5 between Miconazole Related Compound C and Miconazole Related Compound I, a resolution NLT 1.5 between Miconazole Related Compound I and Econazole Nitrate, a resolution NLT 1.5 between Miconazole Related Compound F and Miconazole Nitrate, as well as a %RSD NMT 3.0 % for Miconazole Nitrate. The separation of Miconazole Nitrate and its organic impurities was achieved with all requirements for System Suitability for Organic Impurities met using the Kinetex Phenyl-Hexyl column (Figure 3).

All solutions were prepared as indicated in the USP Monograph for Miconazole Nitrate. USP Miconazole Nitrate RS (Catalog No. 1443500), USP Econazole Nitrate RS (Catalog No. 1231808), USP Miconazole Related Compound C RS (Catalog No. 1443431), USP Miconazole Related Compound F RS (Catalog No. 1443464), and USP Miconazole Related Compound I RS (Catalog No. 1443497) were purchased from USP.




Figure 1. Miconazole Nitrate



Results and Discussion

Figure 2. Standard Solution - Assay



Application Methods

LC-UV Conditions

Column: Kinetex™ 2.6 µm Phenyl-Hexyl

Dimensions: 100 x 4.6 mm

Part No.: [00D-4495-E0](#)

Mobile Phase: A: Methanol/Water/1 M Triethylammonium

Acetate (30:70:1, v/v/v)

B: Acetonitrile/Methanol/1 M Triethylammonium

Acetate (25:75:0.2, v/v/v)

| Gradient: Time (min) | %B |
|----------------------|----|
| 0                    | 30 |
| 5                    | 30 |
| 10                   | 56 |
| 27                   | 56 |
| 30                   | 75 |
| 35                   | 75 |
| 36                   | 30 |
| 40                   | 30 |

Flow Rate: 0.8 mL/min


Injection Volume: 10 µL

Temperature: 40 °C

Detector: UV @ 215

System: Waters® ACQUITY® I-Class UPLC

Figure 3. System Suitability Solution – Assay



In both Figure 2 and 3, there was a rise in the baseline after 10 minutes that continued and then increased in magnitude after 27 minutes. The baseline returned after 36 minutes when the mobile phase composition is returned to the initial conditions. This rise was also seen in the solvent blank (not shown), so it was not due to the composition of the prepared solutions for this application. The rising baseline is due to the increasing % of Methanol in the mobile phase during the gradient. Methanol has some absorbance (is not as UV transparent as Acetonitrile) at 215 nm, so as the % Methanol increases the UV absorbance also increases, resulting in the rising baseline.

Table 1. Preparation of Solutions

| Solution                               | Composition  |
|--|--|
| Diluent – Assay                        | Methanol/Water (70:30, v/v)  |
| Standard Solution – Assay              | 0.1 mg/mL USP Miconazole Nitrate RS in Diluent   |
| System Suitability Solution – Assay    | 0.1 mg/mL of USP Miconazole Nitrate RS and 6 µg/mL of USP Miconazole Related Compound F RS in Diluent  |
| Standard Solution – Organic Impurities | 1.2 µg/mL each of USP Miconazole Nitrate RS, USP Econazole Nitrate RS, USP Miconazole Related Compound C RS, USP Miconazole Related Compound F RS, and USP Miconazole Related Compound I RS in Diluent |



Have questions or want more details on implementing this method? We would love to help! Visit [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat) to get in touch with one of our Technical Specialists

The baseline rise due to the increase in Methanol is much more noticeable when the analyte concentrations are much lower than in the Assay.





Figure 4. Standard Solution – Organic Impurities



| Peak No.       | Analyte                       | Resolution | Area %RSD |
|----------------|-------------------------------|------------|-----------|
| 1              | Miconazole Related Compound C | -          | -         |
| 2              | Miconazole Related Compound I | 2.417      | -         |
| 3              | Econazole Nitrate             | 1.867      | -         |
| 4              | Miconazole Related Compound F | -          | -         |
| 5              | Miconazole Nitrate            | 2.2        | 2.592     |
| N=6 Injections |                               |            |           |

Figure 4 also shows the rise in the baseline due to the increasing % of Methanol in the mobile phase during the gradient, but to what seems a greater extent. This, however, is the same rise that was observed in Figure 2 and 3 but is much more noticeable because the analyte concentrations are much lower compared to the Assay and there is a difference in scaling.

### Key Learnings:

- This method was revised by USP as part of the USP monograph modernization efforts and subject to the typical USP review and comment process before becoming official.
- The results obtained using the Kinetex™ 2.6 μm Phenyl-Hexyl 100 x 4.6 mm column demonstrate that the Kinetex 2.6 μm Phenyl-Hexyl column met all system suitability requirements as published in the USP monograph.

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/phx-terms-and-conditions-of-sale](http://www.phenomenex.com/phx-terms-and-conditions-of-sale).

**Trademarks**

Kinetex is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

**Disclaimer**

Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.  
© 2023 Phenomenex, Inc. All rights reserved.

