

APPLICATIONS

Comparison of Two High-Performance Particle Morphologies in the Separation of Hydrochlorothiazide and Chlorothiazide

Zeshan Aqeel, J Preston, Phil Koerner, and Simon Lomas
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

In this application, we build upon previous work and explore the high-performance separation of hydrochlorothiazide and chlorothiazide with two different LC column particle morphologies. Under the same system and method conditions we compared a traditional Fully Porous 5 µm C18 to both a thermally modified Fully Porous Luna Omega 5 µm C18 and a Core-Shell Kinetex 5 µm C18 of identical dimensions. The Luna Omega and Kinetex columns satisfied the analysis requirement for a Hydrochlorothiazide Sulfide tablet under USP guidelines and were compared to the traditional Fully Porous column originally referenced in the same guidelines.

Compound: Hydrochlorothiazide
pK_{a1}: 7.9, pK_{a2}: 9.2
LogP: -0.07

USP Monograph: Hydrochlorothiazide Tablet Assay

Standard Solution 0.15 mg/mL of USP Hydrochlorothiazide RS in *Mobile Phase*

System Suitability Solution 0.015 mg/mL of Chlorothiazide RS and 0.015 mg/mL Hydrochlorothiazide RS in *Mobile Phase**

* Note: A volume of Acetonitrile not exceeding 10 % of the total volume of solution may be used to dissolve the USP Reference Standard (RS)

Column

Size Method 1: 250 x 4.6 mm, Method 2: 250 x 4.6 mm, Method 3: 250 x 4.6

Stationary Phase Method 1: Symmetry® 5 µm C18, Method 2: Kinetex 5 µm C18, Method 3: Luna Omega 5 µm C18

Temperature 30 °C

Mobile Phase Acetonitrile and 0.1 M Monobasic Sodium Phosphate (1:9). Adjust with Phosphoric Acid to a pH of 3.0 ± 0.1

Isocratic Isocratic: (1:9, A:B)
Total Run Time: 30 min

Flow Rate 2.0 mL/min

Detector UV @ 254 nm

Injection Volume 20 µL of System Suitability solution and Standard solution

System Suitability – System Suitability solution and Standard solution

Sample: Standard solution and System Suitability solution:

- Resolution (Rs): NLT 2.0 between Chlorothiazide and Hydrochlorothiazide for System Suitability solution
- Relative Standard Deviation: NMT 1.5 % for Standard solution (5 replicate injections)

Method: Standard solution Overlay


Kinetex 5 µm C18, Luna Omega 5 µm C18, and Symmetry 5 µm C18

Column: Symmetry 5 µm C18 (grey)
Kinetex 5 µm C18 (orange)
Luna Omega 5 µm C18 (blue)

Dimensions: 250 x 4.6 mm

Flow Rate: 2.0 mL/min

Sample: 1. Hydrochlorothiazide



App ID: 25877


APPLICATIONS

Method 1

System Suitability solution: Hydrochlorothiazide Symmetry[®] 5 µm C18

Column: Symmetry 5 µm C18 Traditional Fully Porous
Dimensions: 250 x 4.6 mm
Flow Rate: 2.0 mL/min
Sample: 1. Chlorothiazide
2. Hydrochlorothiazide

Standard Solution	Symmetry C18	Percent RSD
Resolution	6.34	0.190
Chlorothiazide RT	5.75	0.038
Hydrochlorothiazide RT	7.26	0.029




App ID: 25878

Method 2

System Suitability solution: Hydrochlorothiazide Kinetex[®] 5 µm C18

Column: Kinetex 5 µm C18 Core-Shell
Dimensions: 250 x 4.6 mm
Part No.: [00G-4601-E0](#)
Flow Rate: 2.0 mL/min
Sample: 1. Chlorothiazide
2. Hydrochlorothiazide

Standard Solution	Kinetex C18	Percent RSD
Resolution	8.11	0.055
Chlorothiazide RT	4.15	0.020
Hydrochlorothiazide RT	5.18	0.026




App ID: 25879

Method 3

System Suitability solution: Hydrochlorothiazide Luna[®] Omega 5 µm C18

Column: Luna Omega 5 µm C18 Thermally Modified Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: [00G-4785-E0](#)
Flow Rate: 2.0 mL/min
Sample: 1. Chlorothiazide
2. Hydrochlorothiazide

Standard Solution	Luna Omega C18	Percent RSD
Resolution	7.56	0.172
Chlorothiazide RT	6.77	0.043
Hydrochlorothiazide RT	8.58	0.038



App ID: 26147

APPLICATIONS

Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2 & 3
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10%	As specified	As specified
Composition of the Mobile Phase	± 30% Relative; cannot exceed ± 10% Absolute adjustment; cannot be reduced to zero	As specified	As specified
Wavelength of Detector	No deviations permitted	254 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	20 µL (as specified)	10 µL (Allowed)
Column Temperature	± 10 °C	30 °C (Allowed)	30 °C (Allowed)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	250 mm (as specified)	250 mm (As specified)
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (as specified)	4.6 mm (As specified)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25% and +50%*	5 µm (as specified)	5 µm (As specified)
Flow Rate	± 5.0 % (at given ID)	2.0 mL/min (as specified)	2.0 mL/min (As specified)

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25% to +50%.

Allowable Column Adjustments: L/dp Ratio -25 % to 50 %

Column	Length (mm)	ID (mm)	dp (µm)	L/dp	Allowable Range 37,500-75,000
Original	250	4.6	5	50,000	37,500 – 75,000
Alternative	250	4.6	5	50,000	ALLOWED

Method Summary and Comparison

	Method 1	Method 2	Method 3
Column	Symmetry® 5 µm C18	Kinetex® 5 µm C18	Luna® Omega 5 µm C18
System Suitability Chlorothiazide Average Rt	5.75 min	4.15 min	6.77 min
System Suitability Hydrochlorothiazide Average Rt	7.26 min	5.18 min	8.58 min
System Suitability Hydrochlorothiazide Average Rs	6.34	8.112	7.562
System Suitability Hydrochlorothiazide Average Peak Height	71.1	134.53	71.32
System Suitability Hydrochlorothiazide Average Peak Area	715.80	722.15	722.98
System Suitability Hydrochlorothiazide Peak Height RSD (n=5)	0.183 %	0.081 %	0.154 %
System Suitability Hydrochlorothiazide Peak Area RSD (n=5)	0.316 %	0.102 %	0.142 %
Backpressure (Bar)	350	346	326

Conclusion

Both the Kinetex and Luna Omega 5 µm C18 met the system suitability requirement for this monograph method, and provided improved resolution and method performance in comparison to a traditional fully porous C18. This application demonstrated the use of allowable adjustments to improve method performance without requiring a re-validation. Also, this application illustrates selectivity difference based on particle morphology and the benefit of a reproducible and robust particle starting point in method development and lifecycle management.

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

● **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

BE-HAPPYTM
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex and Luna are registered trademarks and BE-HAPPY is a trademark of Phenomenex. Symmetry is a registered trademark of Waters Technologies Corporation.

Disclaimer

Phenomenex is not affiliated with Waters Corporation. Comparative chromatograms may not be representative of all separations.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.