

REDUCE PROTOCOL TIME BY AT LEAST **40%** WITH A REVOLUTIONARY 3- AND 2-STEP SPE

NO EQUILIBRATION NEEDED. NO CONDITIONING NEEDED.
HIGH RECOVERIES



APPLICATION BOOK

Advantages



Isolate your analytes of interest while removing interferences with the highly targeted **Strata-X PRO SPE**

- Ultra-clean extracts
- Helps with the concentration of samples for better chromatographic results
- Solvent switching for GC or LC compatibility
- Extends column lifetime and improves chromatographic results

Strata-X PRO SPE	Increase Column Lifetime	Remove Particulates	Remove Proteins	Remove Phospholipids	De-salt	Solvent Switching	Specifically Extract Target Analyte	Concentrate
•	•	•	•	•	•	•	•	•

A separation process that is used to remove compounds from a mixture, based on their physical and chemical properties. Analytical laboratories use solid phase extraction to concentrate and purify samples for analysis from a wide variety of matrices.

Sample Preparation Tools and Resources



Syringe Filter Finder Tool



Search Hundreds of Applications



Sample Preparation Support at Your Fingertips



SPE Method Development Tool



Sample Preparation Basics Overview

www.phenomenex.com/sampleprep

Getting Started

1

Prepare Your Sample

Plasma/Serum

If the analyte of interest is an acid, 2% phosphoric acid can be used (20 μ L 85% H_3PO_4 to 1 mL of plasma (or serum) to disrupt the drug-protein interaction. If the analyte of interest is basic, 0.1M sodium hydroxide can be used to disrupt the drug-protein interaction. After addition of acid or base, the sample should be vortexed for 20-30 seconds followed by centrifugation. The supernatant is now ready for further analysis. Alternatively, a zinc sulfate or strong organic solvent can be used to disrupt the protein interaction.

Whole Blood

There are several pre-treatment strategies that can be followed for whole blood. If the target analyte is present in red blood cells, a hemolysis step is necessary.

a. Hemolysis: To 0.2 mL whole blood (spiked with analytes and internal standard) in a 1.2 mL centrifuge tube, add 400 μ L of 2% zinc sulfate/80% methanol. Vortex for 10-20 seconds followed by centrifugation at 14,000 rpm for 10 minutes. Collect the supernatant for further analysis.

Preparation of zinc sulfate/methanol: Into a 100 mL volumetric flask add 20 mL water and 3.6 g $ZnSO_4 \cdot 7H_2O$. After the solution is clear and the salt crystals have dissolved, add 100% methanol. Refrigerate the solution at 2-8 $^{\circ}C$ for 7 days.

b. Osmotic breakdown: To 1 mL of whole blood add internal standard and 4 mL of distilled water. Mix/vortex and let stand for 5 minutes. Centrifuge at 670g for 10 minutes and discard the pellet. Adjust the pH of the supernatant accordingly with the addition of a buffer solution.

c. Sonication: Add 3-6 mL of appropriate pH buffer (such as potassium phosphate) to 1 mL of whole blood and sonicate for 15 minutes at room temperature.

Urine

Enzymatic hydrolysis is necessary in case of conjugated forms (sulfated or glucuronide) of the analyte present. Enzymatic hydrolysis requires specific pH (pH 4-5) and temperature ranges. An acid or base hydrolysis can be performed as well, depending on the stability of the compound.

a. Enzymatic hydrolysis: To 500 μ L sample (spiked with analyte and internal standard) add 100 μ L acidic buffer (see below) and 100 μ L β -glucuronidase. Vortex 5-6 seconds. Incubate in a water bath at 63 $^{\circ}C$ for 30 minutes. Transfer sample to a 96-well collection plate or autosampler vial. Seal and centrifuge for 10 minutes at 2,000 rpm.

Preparation of acidic buffer (1.0 M acetate buffer, pH 4.0): Dissolve 3.0g of glacial acetic acid and 4.1 g of sodium acetate in a 1L volumetric flask.

b. Base hydrolysis: To 1 mL urine (spiked with analyte and internal standard) add 100 μ L 10 N KOH. Mix, vortex, and hydrolyze for 20 minutes at 60 $^{\circ}C$. Cool and adjust pH to 3.5- 4.0 (by adding 200 μ L glacial acetic acid).

c. Acid hydrolysis: To 1 mL urine add 0.25 mL HCl in a screw capped test tube. Screw the tube top on loosely and heat in a boiling water bath for 60 minutes. Adjust to pH 7 (or as needed) with 1.0 N NaOH.

Saliva

No hydrolysis is required for oral fluids and the generic protocol used for plasma/serum pretreatment may be followed.

Tissue

Homogenize with organic or aqueous solvent depending upon analyte solubility. Settle, decant, centrifuge or filter supernatant. Perform direct Matrix Solid Phase Dispersion (MSPD) extraction on tissue.

Getting Started *(cont'd)*

2 Find the Correct Volumes

Sample Matrix	Strata-X PRO Sorbent Mass	Sample Size	Recommended Sorbent Wash and Elution Volumes		
			Strata-X PRO Sorbent Mass	Practical Minimum Wash and Elution Volume	Recommended Wash and Elution Volumes
			4 bed volumes		8 bed volumes
Blood, serum, plasma	30 mg	250 µL	10 mg	100 µL	200 µL
Urine	30 mg	1 mL	30 mg	300 µL	600 µL
Filtered tissue homogenates	60 mg	100 mg	60 mg	600 µL	1.2 mL
Oral Fluid	30 mg	500 µL	100 mg	1 mL	2 mL
Environmental Samples	Strata-X PRO Sorbent Mass	Sample Size	200 mg	2 mL	4 mL
Water (particulate-free) drinking	500 mg	100 - 400 mL	500 mg	5 mL	10 mL
Water (particulate-laden) rivers, runoff, etc.	500 mg	100 - 400 mL			
Soil extracts	500 mg	100 g			

3 Determine Your Method



2-Step Protocol

Non-retentive SPE method to help achieve the fastest extraction.



Load

1 mL Pre-treated sample/0.1 % Formic acid in Acetonitrile (1:4)
Apply 5" Hg vacuum until all tubes or wells have cleared

Elute

75 µL Water/0.1 % Formic acid in Acetonitrile (1:4)
Apply 5" Hg vacuum until all tubes or wells have cleared

Protocols are written for 30 mg/1 mL tubes, adjust based on sorbent size.



3-Step Protocol

Rapid protocol to reduce matrix effects and increase recovery of polar analytes.



Load

500 µL Pre-treated sample/buffer* (1:1)
Apply 2-5" Hg vacuum until liquid is no longer visible above top frit

Wash

600 µL 5 % Methanol in Water

Elute

600 µL 0.1 % Formic acid in Acetonitrile/Methanol (90:10)
Apply 2-5" Hg vacuum for 1 minute

Protocols are written for 30 mg/1 mL tubes, adjust based on sorbent size.

Select a buffer that maximizes the hydrophobicity of the analytes. For example, if an analyte is basic, dilute with a base.


Cleaner Extractions

Strata-X PRO uses cutting-edge technology to reduce traditional SPE protocol time and still results in high recoveries of the analytes. Its broad range of target applications makes it a great selection for labs that need a clean, fast, and efficient sample preparation method.

SPE Protocol

- 96-Well Plate:** Strata-X PRO, 30 mg/well
- Part No.:** [8E-S536-TGA](#)
- Load:** 400 μ L Plasma/0.1 % Formic acid in Water (1:1)
- Wash:** 1 mL 5 % Methanol in Water
- Dry:** 1 minute at 5" Hg
- Elute:** 1 mL 0.1 % Formic acid in Acetonitrile/Methanol (90:10)
- Dry Down:** 1 minute at 5" Hg
- Reconstitute:** 200 μ L 5 % Methanol in Water

Recovery from Human Plasma



For extremely polar analytes, Strata-X PRO provides higher recoveries!

Ultra-Quick Clean-Up

Strata-X PRO accomplishes what was unimaginable. Its 3- and 2- rapid clean-up system reduces protocol time by 40%. When working with milk as a matrix, phospholipids from milk fat must be removed to reduce any ion suppression that could occur during LC-MS/MS analysis for veterinary drugs. This ultra-quick SPE works for complex matrices and applications from a variety of industries: Forensics/Toxicology, Environmental, Clinical, Bio/Pharma, and Food Safety and Quality.

SPE Protocol

Pre-treatment

To 1 mL of milk add 3 mL of 0.2% Formic acid in Acetonitrile/Methanol (90:10) and mix or vortex for 15-20 seconds. Centrifuge for 5 minutes at 10,000 RPM and collect supernatant.

Cartridge: Strata-X PRO, 60 mg/ 3 mL

Part No.: [8B-S536-UBJ](#)

Load: Pass the pre-treated sample through the SPE cartridge and collect

Dry: Evaporate the extract to dryness under a gentle stream of nitrogen at room temperature

Reconstitute: The dried sample in 1 mL of initial mobile phase (0.1% Formic acid in Water/0.1% Formic acid in Methanol (95:5)) spiked with deuterated internal standard.


% Recovery and % CVs for Veterinary Drugs from Milk Using Strata-X PRO

Peak No.	Analyte Name	Retention Time (min)	% Recovery	% CV	Q1	Q3
1	Sulfaguanidine	1.48	46	5	215	156.1
2	Lincomycin	2.07	92	5	407.1	126
3	Sulfadiazine	2.19	38	7	251	156
4	Cephapirin	2.22	76	7	424	292.1
5	Sulfamerazine	2.32	44	5	265.1	155.8
6	Sulfamethoxazole	2.36	53	13	254.1	156.1
7	Sulfamethizole	2.36	45	8	271.1	92
8	Cefalexin	2.39	66	4	348.2	174.2
9	Sulfamethazine	2.44	59	13	279.1	186.1
10	Cortisone	2.72	83	8	361.2	163.2
11	Cortisol	2.73	95	6	363.4	120.9
12	β -methasone	2.76	97	3	393.4	355.2
13	Prednisolone	2.81	92	10	361.2	147.2

So good you'll think it's made up.

Ultra-Quick Clean-Up (cont'd)


Veterinary Drugs from Milk




Comparative Phospholipid Trace of Clean-Up Methods

Phospholipid Trace After Protein Precipitation

LC Conditions for Phospholipid Comparison



Phospholipid Trace After Strata[®]-X PRO Extraction



Additional Applications

Strata-X Pro



[Strata-X PRO Solid Phase Extraction User Guide](#)

[Utilizing a Rapid, Two-Step Method for the Clean-Up of Veterinary Drugs in Milk Using Strata[®]-X PRO Solid Phase Extraction \(SPE\)](#)

[Comparison of Tradition Solid Phase Extraction \(SPE\) and Strata[®]-X PRO SPE for the Extraction of Acids, Neutrals, and Bases from Plasma](#)

[A Faster SPE Method for Pain Relievers Using Strata-X PRO from Plasma](#)

[High Recovery of Peptides Using Strata-X PRO Solid Phase Extraction](#)

[Fast Extraction of Barbiturates from Serum Using Strata-X PRO Solid Phase Extraction](#)

SPE Guides and Resources



[Solid Phase Extraction for Clinical Research](#)

[Solid Phase Extraction for Food Samples](#)

[Clinical Resources for Toxicology Analysis](#)

SPE Applications






[A Fast and Effective Quantitation Method for Uracil, 5,6-Dihydrouracil, and 5-Fluorouracil from Human Serum by LC-MS/MS](#)

[A Rapid and Robust Sample Preparation Method for Quantitation of Nicotine from Oral Fluid](#)

Download Applications at phenomenex.com/StrataXPRO

Ordering Information

Strata[®]-X PRO SPE

Format	Sorbent Mass	Part Number	Unit
Tube			
	10 mg	8B-S536-AAK	1 mL (100/box)
	30 mg	8B-S536-TAK	1 mL (100/box)
	30 mg	8L-S536-TAK	1 mL (100/box)
	30 mg	8B-S536-TBJ	3 mL (50/box)
	60 mg	8B-S536-UBJ	3 mL (50/box)
	200 mg	8B-S536-FBJ	3 mL (50/box)
	100 mg	8B-S536-ECH	6 mL (30/box)
	200 mg	8B-S536-FCH	6 mL (30/box)
	500 mg	8B-S536-HCH	6 mL (30/box)
96-Well Plate			
	10 mg/well	8E-S536-AGA	ea
	30 mg/well	8E-S536-TGA	ea
	60 mg/well	8E-S536-UGA	ea
96-Well Microelution Plate			
	2 mg/well	8M-S536-4GA	ea

BE-HAPPY[™]

guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Round Well Collection Plates (polypropylene)

Part No.	Well Bottom	Well Volume	Unit	Suggested Sealing Mats
AH1-7025	Round	1 mL	50/pk	AHO-8631 AHO-8632
AHO-8636	Round	2 mL	50/pk	AHO-8633 AHO-8634

Square Well Collection Plates (polypropylene)

Part No.	Well Bottom	Well Volume	Unit	Suggested Sealing Mats
AHO-7192	Conical	350 µL	50/pk	AHO-8597 AHO-8598 AHO-8199 AHO-7195
AHO-7193	Conical	1 mL	50/pk	AHO-8597 AHO-8598 AHO-8199 AHO-7195
AHO-7194	Conical	2 mL	50/pk	AHO-8597 AHO-8598 AHO-8199 AHO-7195

Round Well Collection Plate (Low Bind)

Part No.	Well Bottom	Well Volume	Unit	Suggested Sealing Mats
AH1-7036	Conical	2 mL	120/pk	AHO-8633 AHO-8634

Round Well Sealing Mats

Part No.	Description	Material	Unit
AHO-8631	Pierceable, 7 mm diameter	Silicone	50/pk
AHO-8632	Pre-Slit, 7 mm diameter	Silicone	50/pk
AHO-8633	Pierceable, 8 mm diameter	Silicone	50/pk
AHO-8634	Pre-Slit, 8 mm diameter	Silicone	50/pk
AHO-7362	Sealing Tap Pad	—	10/pk

Square Well Sealing Mats

Part No.	Description	Material	Unit
AHO-8597	Pierceable	Silicone	50/pk
AHO-8598	Pre-Slit	Silicone	50/pk
AHO-8199	Pierceable	Santoprene [™]	100/pk
AHO-7195	Pierceable	Ethylene Vinyl Acetate (EVA)	50/pk
AHO-7362	Sealing Tap Pad	—	10/pk

Redefine Your Solid Phase Extraction with Strata-X PRO!

Ordering Information *(cont'd)*

Sample Processing

Streamline your 96-well plate processing for easier sample preparation with a pneumatic positive pressure manifold.

Pneumatic Handling

Consistent Flow Rates

Safe and Easy-to-Use



Presston 1000
Part No.: [AH1-7033](#)

Do More with Presston™ 1000

Safer Lab Environment

Use both hands to move manifold shield, ensuring enhanced safety.

Determine Your Operating Pressure

Conveniently monitor and maintain a consistent pressure.

3 Easily Load Samples

Moveable locator plate makes sample loading and cleaning easy.

Never Lose Pressure

Always maintain a tight seal between the manifold and 96-well plate.

Sleek, Low Profile Design

Width: 11", Depth: 15", Height: 14.8"

Know Your Step

Simply move the "SPE Procedure Indicator" to the correct step to stay on track of your extraction.

Presston™ 1000 Positive Pressure Manifold

Part No.	Description
AH1-7033	Presston 1000 Positive Pressure Manifold, 96-Well Plate



Phenomenex warrants the Presston 1000 Positive Pressure Manifold against defects in materials and workmanship under normal installation, use, and maintenance for a period of 12 months following delivery.

Please visit www.phenomenex.com/presstonwarranty for complete warranty information.




**Your solution to
increased productivity.**

Learn more at www.phenomenex.com/Presston


Expand Your Solid Phase Extraction Portfolio



3 Unique Sorbent Platforms

 <p>Underline reversed phase polymer with matrix removal technology offers a faster, cleaner way to perform SP.</p> <div style="display: flex; justify-content: space-around;"> <div data-bbox="272 864 400 987"> </div> <div data-bbox="445 864 572 987"> </div> </div>	 <p>Polymeric sorbent available in reversed phase and ion-exchange capabilities for wide range of applications.</p> <div data-bbox="751 857 879 981"> </div>	 <p>Silica-based SPE sorbent provides a reliable and clean extracts with high recoveries for target analytes across all sample matrices.</p> <div data-bbox="1149 857 1276 981"> </div>
---	---	--

Identify Your SPE Retention Mechanism



STRATA[®]X[®] PRO

A Rapid Solid Phase Extraction Solution

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thainfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Strata and Kinetex are registered trademarks and SecurityGuard, BE-HAPPY, and Presston are trademarks of Phenomenex. SCIEX is a registered trademark, Triple Quad and API 4000 are trademarks of AB SCIEX Pte. Ltd. Santoprene is a trademark of Exxon Mobil Corporation.

Disclaimer

Comparative separations may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.