

Fully Automated Sample Preparation for the Determination of Vitamin D on the Biomek® i7 Hybrid Using Advanced Sample Cleanup with Phree™ Phospholipid Removal Cartridges

Kerstin Thurow¹, Heidi Fleischer¹, Anna Bach¹, Bhagya Wijayawardena², Miranda Kheradmand², Jenny Lewis³, and Dr. Bryan Tackett³

¹Center for Life Science Automation, University Rostock, Friedrich-Barnewitz-Str. 8, 18119 Rostock, Germany

²Beckman Coulter Life Sciences, 5350 Lakeview Pkwy S Dr., Indianapolis, IN 46268 USA

³Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

Researchers often quantify the level of vitamin D using LC-MS based methods. The sample preparation for these analytical methods, however, are labor intensive and prone to human error. In this application note we evaluate the feasibility of automating the sample prep for vitamin D determination on the Biomek i7 hybrid workstation. We used Phree phospholipid removal cartridges (Part No. [8B-S133-TAK](#), Phenomenex, Torrance, US) for the protein precipitation prior to the sample cleanup by solid phase extraction.

The quantitative determination of 25(OH)D₂ and 25(OH)D₃ is becoming increasingly important to understand the effects of overdosing and underdosing of vitamins. Due to the complicated matrix, the sample preparation for vitamin D analysis is extensive, including protein precipitation and extraction steps. The protein precipitation and the separation of the metabolites from the vitamin D-binding proteins is usually carried out by adding acetonitrile, methanol, mixtures of acetonitrile and methanol, or 2-propanol followed by a centrifugation step. In addition to the proteins, there are other substances in biological matrices, in particular phospholipids (components of the cell membrane), whose highly ionic nature interferes with the ionization in the mass spectrometer. In addition, phospholipids can cause reduced column lifetime as they build up on the stationary phase, which leads to changes in the separation properties and an increase in the baseline noise. Therefore, the removal of the phospholipids is an important step in the sample preparation. The main methods of removing phospholipids are solid phase extraction (SPE) with strong cation exchangers, and liquid-liquid extraction (LLE). However, performing these methods increases the time it takes to complete the analytical measurement process. To reduce sample preparation time, we automated a sample preparation method for vitamin D determination that uses Phree phospholipid removal cartridges and solid phase extraction.

The use of Phree phospholipid removal cartridges for combined protein precipitation and sample clean-up reduced the number of processing steps and the processing time. The high recovery rate is indicative of efficient removal of phospholipids (**Figure 2**). The low CV values (2-9 %) of the automated workflow indicates high repeatability.

Sample Preparation

Pre-treatment: Transfer 850 µL Acetonitrile

Add: Transfer 50 µL internal standard followed by 100 µL serum sample

Mix: Solution with multiple aspiration/dispensing

Incubate: Samples at room temperature for 4 min

Apply: Positive pressure for 100 sec (500 mbar)

Measure: Samples using LC-MS or LC-MS/MS

Figure 1. Structures of 25(OH)D₂ and 25(OH)D₃





Figure 2. Calibration Curves for (a) 25(OH)D₂ and (b) 25(OH)D₃



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Phree and BE-HAPPY are trademarks of Phenomenex. Biomek and BECKMAN COULTER are registered trademarks of Beckman Coulter, Inc.

Disclaimer

Phenomenex is in no way affiliated with the Center for Life Science Automation, University of Rostock.

Biomek i-Series is for research use only and not intended or validated for use in the diagnosis of disease or other conditions.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists