

Reproducible Extraction of 1,4-Dioxane by EPA Method 522 using Strata™ Activated Carbon in Comparison to Competitor Brand

Eduardo Morales¹, Agustin Pierri, PhD¹, Ramkumar Dhandapani, PhD², Philip Leong², and Bryan Tackett, PhD²

¹Weck Laboratories, Inc., 14859 Clark Ave, Hacienda Heights, CA 91745 USA

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

EPA method 522 is a popular testing method for 1,4-Dioxane by GC-MS. This analyte is challenging to get good recovery with traditional sample preparation. Here, EPA method 522 is performed using activated carbon for the extraction of 1,4-Dioxane. While dealing with polar small molecules like 1,4-Dioxane, traditional sample preparation using Solid Phase Extraction (SPE) might not be efficient. Since the analyte is polar, a specialized extraction technique is needed. Strata Activated Carbon, which contains porous carbon with defined pore volume and high surface area, is an appropriate choice for polar analyte extraction and recovery. In this application note we have extracted 1,4-Dioxane from a 100 mL water sample using Strata Activated Carbon 400 mg pass through cartridges, which are convenient for small sample volume and demonstrate extraction-to-extraction and batch-to-batch reproducibility.

GC-MS Conditions

Column:	Zebtron™ ZB-624 ^{PLUS} ™		
Dimension:	30 meter x 0.25 mm x 1.40 μm		
Part No.:	7HG-G040-27		
Injection:	Splitless, 0.5min @ 175 °C, 1 μL		
Recommended Liner:	Zebtron PLUS Z-Liner™ (Compatible with Agilent® & Thermo Scientific® GC instrument)		
Liner Part No.:	AG2-0A11-05		
Carrier Gas:	Helium @ 1 mL/min (Constant Flow)		
Oven Program:	Ramp(°C/min)	Temp (°C)	Time(min)
	-	35	1.0
	12	100	0.0
	25	230	0.0
Detector:	MSD		
SIM Ions (m/z):	58, 64, 88, 96		
Transfer Line Temperature:	300 °C		

Sample Preparation

Condition:	Strata Activated Carbon, 400 mg pass through cartridge (Part No. CS0-9210) or Waters® Sep-Pak® AC2 Plus 400 mg short cartridge with 4 mL Methylene Chloride, followed by 4 washes of 5 mL Methanol
Equilibrate:	Cartridges with 2 washes of 5 mL Water
Load:	100 mL water sample spiked with internal standard onto cartridges
Dry:	Cartridges for 10 min
Elute:	With 2 mL Methylene Chloride and repeat elution 1 more time (4 mL total volume)
Water Removal:	Pass eluent through Sodium Sulfate tubes, 1 g/6 mL (Part No. 8B-S124-JCH) and wash 2 times with 2 mL Methylene Chloride
Evaporate:	Solvent to approximately 0.1 mL at 45 °C under a gentle stream of Nitrogen
Adjust:	Final volume to 0.5 mL with Methylene Chloride

Figure 1. Strata Activated Carbon Extraction Cartridges




Figure 2. Batch-to-Batch Reproducibility of 1,4-Dioxane Extraction from Water Sample using the Strata™ Activated Carbon Cartridge




Figure 3. Comparison of 1,4-Dioxane-d8 % Recovery from a Water Sample




Figure 4. Comparison of 1,4-Dioxane-d8 Recovery Variability from a Water Sample




Figure 5. Comparison of 1,4-Dioxane % Recovery from a Water Sample





Figure 4. Comparison of 1,4-Dioxane Recovery Variability from a Water Sample



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☉ All other countries/regions

Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Strata, Zebtron, 624PLUS, Z-Liner, and BE-HAPPY are trademarks of Phenomenex. Waters and Sep-Pak are registered trademarks of Waters Technologies Corporation. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific, Inc.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Waters Technologies Corporation, Agilent Technologies, Inc. or Thermo Fisher Scientific, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

