

Fast Separation of Chlorinated Pesticides Using Zebron™ GC Columns

Kory Kelly
Phenomenex, Inc., 411 Madrid Avenue, Torrance, CA 90501, USA

- Fast analysis of chlorinated pesticides in under 10 minutes results in shorter cycle times and improved productivity
- Baseline separation of all analytes using Zebron ZB-MultiResidue™-1 and -2 columns provides greater confidence in identification

Introduction


Pesticides are classes of chemicals that are used to treat or control outbreaks of pests, especially insects. Food production as well as safety and health organizations heavily rely on pesticides to increase food yields and / or control disease outbreak spread by insects. Some pesticides have proven hazardous to other animals or the environment and have either been restricted or banned. The need to monitor products, especially foods, for pesticides is essential as more pesticides are being discovered to have adverse effects.

One class of commonly used pesticides is chlorinated pesticides. These compounds are commonly detected using EPA Method 8081. This method determines the concentrations of various organochlorine pesticides in extracts from solid and liquid matrices, including food products, using gas chromatography with electron capture detectors (ECD). The following applications provide two options for analyzing chlorinated pesticides by EPA Method 8081. Two different sets of columns are shown that provide separation of all analytes while still keeping retention times less than ten minutes to allow for faster cycle times and higher lab throughput.

Methods

If a compound is identified using a single column, the compounds presence should be confirmed on a second column. Simultaneous identification and verification can be performed using a dual-column configuration as shown in **Figure 1**. Further experimental conditions are shown with each application.



Figure 1.
Example of a Dual-Column Configuration.



Discussion

Figure 2 shows a dual-column separation using Zebron MultiResidue-1 and -2 columns. Each chromatogram shows complete resolution of all 20 analytes that are typically analyzed within EPA Method 8081 with run times less than 10 minutes. The chromatogram also includes surrogates to ensure that they don't co-elute with analytes. The two main advantages of using a Zebron MultiResidue-1 and -2 for these analyses are that all compounds have baseline resolution as well as the very short retention times. The improved resolution allows for more robust methods. The faster retention times result in short cycle times and increased lab productivity.

Figure 2.
Dual-Column Analysis of Chlorinated Pesticides using Zebron MultiResidue-1 and -2 Columns.



Conditions the same except where noted:

Injection: Splitless @ 250 °C, 1 µL
Carrier Gas: Helium @ 2.7 mL/min (constant flow)
Oven Program: 90 °C for 0.5 min to 200 °C at 35 °C/min to 340 °C at 25 °C/min for 1.5 min
Detector: Electron Capture (ECD) @ 350 °C



Sample: 1 TCMX (Surr)	10 trans-Chlordane (gamma)	17 Endosulfan II
2 1-Bromo-2-nitrobenzene (IS)	11 cis-Chlordane (alpha)	18 Endrin aldehyde
3 a-BHC	12 Endosulfan I	19 4,4'-DDT
4 g-BHC (Lindane)	13 4,4'-DDE	20 Endosulfan sulfate
5 b-BHC	14 Dieldrin	21 Methoxychlor
6 d-BHC	15 Endrin	22 Endrin ketone
7 Heptachlor	16 4,4'-DDD	23 Decachlorobiphenyl (Surr)
8 Aldrin		
9. Heptachlor epoxide		

TN-2039

APPLICATIONS

Figure 3 shows separation of the same chlorinated pesticides using an alternative set of Zebron™ ZB-35 and ZB-XLB columns. These two chromatograms also show separation of all analytes and surrogates in under ten minutes. The short run time again shows advantages of shorter cycle times as well as increased lab productivity

Figure 3.
Dual-Column Analysis of Chlorinated Pesticides using Zebron ZB-35 and ZB-XLB Columns.



Conditions the same except where noted:

Injection: Splitless @ 250 °C, 1 µL
Carrier Gas: Helium @ 2.6 mL/min (constant flow)
Oven Program: 110 °C for 0.5 min to 190 °C @ 18 °C/min to 330 °C @ 35 °C/min for 2 min
Detector: Electron Capture (ECD) @ 335 °C

Sample:	1 TCMX (Surr)	10 a-Chlordane	17 4,4'-DDT
	2 a-BHC	11 Endosulfan I	18 Endrin aldehyde
	3 g-BHC (Lindane)	12 4,4'-DDE	19 Endosulfan sulfate
	4 b-BHC	13 Dieldrin	20 Methoxychlor
	5 Heptachlor	14 Endrin	21 Endrin ketone
	6 d-BHC	15 4,4'-DDD	22 Decachlorobiphenyl (IS)
	7 Aldrin	16 Endosulfan II	
	8 Heptachlor epoxide		
	9 g-Chlordane (trans)		

Conclusions

In this application note, two column pairs were exhibited for organochlorine pesticide analysis. The Zebron ZB-XLB and ZB-35 are traditionally used for this analysis but the Zebron ZB-MultiResidue-1 and -2 pair show additional benefits of improved resolution. Both pairs were tested using a dual-column analysis technique that provides separation of all compounds in less than ten minutes. Very short run times for all columns tested allow for fast cycle times resulting in improved productivity for environmental labs under production deadlines. The Zebron ZB-MultiResidue columns also offer improved baseline resolution of all compounds. This could result in more consistent quantitation and greater method robustness by minimizing the possibility of misidentified peaks and removing quantitation complications due to closely eluting peaks.

TN-2039

APPLICATIONS

Ordering Information

Zebtron ZB-35 GC Columns

ID(mm)	df(μ m)	Temp. Limits °C	Part No.
15-Meter			
0.25	0.25	50 to 340/360	7EG-G003-11
0.25	0.50	50 to 340/360	7EG-G003-17
0.53	1.00	50 to 340/360	7EK-G003-22
30-Meter			
0.25	0.25	50 to 340/360	7HG-G003-11
0.25	0.50	50 to 340/360	7HG-G003-17
0.32	0.25	50 to 340/360	7HM-G003-11
0.53	0.50	50 to 340/360	7HK-G003-17
0.53	1.00	50 to 340/360	7HK-G003-22
60-Meter			
0.25	0.25	50 to 340/360	7KG-G003-11
0.32	0.25	50 to 340/360	7KM-G003-11

ZB-35 Test Mix Part No.: AG0-5156

Zebtron ZB-XLB GC Columns

ID(mm)	df(μ m)	Temp. Limits °C	Part No.
15-Meter			
0.25	0.25	30 to 340/360	7EG-G019-11
20-Meter			
0.18	0.18	30 to 340/360	7FD-G019-08
30-Meter			
0.25	0.25	30 to 340/360	7HG-G019-11
0.25	0.50	30 to 340/360	7HG-G019-17
0.32	0.25	30 to 340/360	7HM-G019-11
0.32	0.50	30 to 340/360	7HM-G019-17
0.53	1.50	30 to 320/340	7HK-G019-28
60-Meter			
0.25	0.25	30 to 340/360	7KG-G019-11

ZB-XLB Test Mix Part No.: AG0-7578

Zebtron ZB-MultiResidue™ GC Columns (MR-1)

ID(mm)	df(μ m)	Temp. Limits °C	Part No.
30-Meter			
0.25	0.25	-60 to 320/340	7HG-G016-11
0.32	0.25	-60 to 320/340	7HM-G016-17
0.53	0.50	-60 to 320/340	7HK-G016-17

Zebtron ZB-MultiResidue™ GC Columns (MR-2)

ID(mm)	df(μ m)	Temp. Limits °C	Part No.
30-Meter			
0.25	0.20	-60 to 320/340	7HG-G017-10
0.32	0.25	-60 to 320/340	7HM-G017-11
0.53	0.50	-60 to 320/340	7HK-G017-17



If Zebtron GC columns do not provide you with equivalent separations as compared to any other GC column of the same phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Trademarks

Zebtron and MultiResidue are trademarks of Phenomenex, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

© 2009 Phenomenex, Inc. All rights reserved.

TN-2039

APPLICATIONS

Australia

t: 02-9428-6444
f: 02-9428-6445
auinfo@phenomenex.com

Austria

t: 01-319-1301
f: 01-319-1300
anfrage@phenomenex.com

Belgium

t: +31 (0)30-2418700
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: (800) 543-3681
f: (310) 328-7768
info@phenomenex.com

Denmark

t: 4824 8048
f: 4810 6265
dkinfo@phenomenex.com

France

t: 01 30 09 21 10
f: 01 30 09 21 11
franceinfo@phenomenex.com

Germany

t: 06021-58830-0
f: 06021-58830-11
anfrage@phenomenex.com

Ireland

t: 01 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: 051 6327511
f: 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Netherlands

t: 030-2418700
f: 030-2383749
nlinfo@phenomenex.com

New Zealand

t: 09-4780951
f: 09-4780952
nzinfo@phenomenex.com

Puerto Rico

t: (800) 541-HPLC
f: (310) 328-7768
info@phenomenex.com

United Kingdom

t: 01625-501367
f: 01625-501796
ukinfo@phenomenex.com

All other countries: Corporate Office USA

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com.