

APPLICATION

Analyzing β -GoneTM β -Glucuronidase Removal Centrifuge Tube Recovery and Clean Up

Matthew Brusius, Jessica Detsch, and Eric Chapa
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Matt Brusius
Product Manager,
Sample Preparation

Matt Brusius is an avid ice hockey player. He likes skating backwards and taking slapshots from the point.

Introduction

This technical note describes how β -Gone β -Glucuronidase Removal 2mL Centrifuge Tubes can be used for the analysis of most comprehensive drugs research panel suites prepared from a urine hydrolysate. By selectively removing the β -glucuronidase, column lifetime is improved versus a standard "dilute-and-shoot" protocol while achieving exceptional recovery for all classes of compounds.

Materials and Methods

All reagents and solvents were HPLC or analytical grade. Analyses were performed using an API 4000TM LC-MS/MS (SCIEX[®], Framingham, MA).

Sample Preparation

Enzymatic Hydrolysis

Campbell Enzyme Hydrolysate was prepared as follows:

- 1) Add 10 μ L of analyte spike (1 μ g/mg) to 200 μ L urine
- 2) Dilute with 100 μ L of 0.1 M ammonium acetate buffer and add 40 μ L of Campbell Science β -Glucuronidase Enzyme Solution (Campbell Part No.: DR2102)
- 3) Add 400 μ L 0.1 % formic acid in water. Proceed to β -Gone Centrifuge Tube Protocol.

IMCSzyme[®] Hydrolysate Mix was prepared as follows:

- 1) Add 10 μ L of analyte spike (1 μ g/mg) to 140 μ L urine
- 2) Dilute with 80 μ L of IMCS buffer (IMCS Part No.: 04-EZ-RHB-20) and add 30 μ L of IMCS Enzyme (IMCS Part No.:04-E1F-010). Proceed to β -Gone Centrifuge Tube Protocol.

β -Gone Centrifuge Tube Protocol

- 1) Tap tube prior to use to consolidate sorbent
- 2) Combine 200 μ L of urine hydrolysate with 133 μ L of 0.1 % formic acid in methanol
- 3) Load the diluted sample into the centrifuge tube
- 4) Cap and invert the tube 10 times. Vortex for 30 seconds
- 5) Centrifuge tube at 14000 rpm for 10 min (16000 xg RCF)
- 6) Collect supernatant for analysis

Table 1. Recovery Values for IMCS and Campbell Enzyme Prepared Samples

Analyte	IMCSzyme Prepared Samples		Campbell Enzyme Prepared Samples	
	Average % Recovery	%CV (n=4)	Average % Recovery	%CV (n=4)
6-MAM	108	10	105	5
7-Aminoclonazepam	108	5	109	6
α -Hydroxyalprazolam	105	5	100	3
Alprazolam	106	5	103	7
Amitriptyline	104	7	102	3
Amphetamine	105	2	102	4
Benzoyllecgonine	110	4	98	4
Buprenorphine	100	2	103	5
Carisoprodol	110	5	98	5
Citalopram	104	5	104	2
Codeine	103	9	100	6
Cotinine	112	9	99	2
Diazepam	103	6	99	2
EDDP	104	4	105	4
Fentanyl	106	6	101	3
Fluoxetine	101	5	104	3
Gabapentin	112	6	98	3
Hydrocodone	110	8	103	5
Hydromorphone	110	10	106	3
Imipramine	103	5	101	2
Lorazepam	102	2	97	4
MDA	111	7	93	3
MDEA	109	6	98	3
MDMA	110	5	96	2
Meperidine	102	5	103	0
Meprobamate	106	4	104	4
Methadone	105	7	103	3
Methamphetamine	106	6	102	1
Methylphenidate	95	6	101	2
Morphine	115	8	101	6
Naloxone	112	13	100	7
Norbuprenorphine	106	4	114	6
Nordiazepam	106	6	100	4
Norfentanyl	97	4	109	5
Norhydrocodone	108	8	102	4
Noroxycodone	99	10	99	7
Nortriptyline	105	6	102	2
O-Desmethyltramadol	105	4	100	3
Oxazepam	103	7	97	7
Oxycodone	111	5	103	7
Oxymorphone	108	8	104	1
Paroxetine	102	7	107	3
PCP	102	6	102	2
Pregabalin	133	15	92	1
Ritalinic Acid	106	7	97	1
Tapentadol	102	6	94	2
Temazepam	105	6	98	3
THC-COOH	100	7	96	4
Tramadol	100	5	104	4
Zolpidem	109	5	108	4

Results and Discussion


Table 1 shows the average % recoveries and % CVs for the comprehensive drugs research panel suite for both the IMCS (recombinant enzyme) and the Campbell (non-recombinant/natural enzyme) solutions. In general, all compounds show an absolute recovery above 90%, and only pregabalin prepared with IMCS shows a % CV of 15 or greater, all other compounds are less than or equal to a % CV of 10.

Figure 1 monitors the change in backpressure over the course of 200 20 μ L injections between a urine sample that is prepared by dilution and one that is prepared with the β -Gone Centrifuge Tubes. The samples prepared with the β -Gone Centrifuge Tubes maintain a constant backpressure over the span of these injections while the urine samples that were prepared by "dilute-and-shoot" trend toward a higher backpressure, and in the case of the natural enzyme, pressure out before the completion of 200 injections. This shows that the sample is cleaner after being subjected to the β -Gone Centrifuge Tubes and could help improve column lifetime versus the standard "dilute-and-shoot" protocol.

HPLC Conditions

Column:	Kinetex® 2.6 μ m Phenyl Hexyl 100Å												
Dimensions:	50 x 3.0 mm												
Part No.:	00B-4495-Y0												
Mobile Phase:	A: 0.1 % Formic acid in Water B: 0.1 % Formic acid in Methanol												
Flow Rate:	0.6 mL/min												
Gradient:	<table border="0"> <thead> <tr> <th>Time (min)</th> <th>% B</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>5</td> </tr> <tr> <td>4</td> <td>95</td> </tr> <tr> <td>5.5</td> <td>95</td> </tr> <tr> <td>5.51</td> <td>5</td> </tr> <tr> <td>7</td> <td>5</td> </tr> </tbody> </table>	Time (min)	% B	0	5	4	95	5.5	95	5.51	5	7	5
Time (min)	% B												
0	5												
4	95												
5.5	95												
5.51	5												
7	5												
Detector:	MS/MS API 4000™												

Figure 1.
Column Lifetime: β -Gone Centrifuge Tubes vs. Dilute-and-Shoot



Conclusion:

By using the β -Gone β -Glucuronidase 2 mL Centrifuge Tubes, one can obtain acceptable recoveries for a large drugs of abuse suite, while improving column lifetime in comparison to a dilute-and-shoot procedure.

Ordering Information

β-Gone™ β-Glucuronidase Removal Products

Part No.	Description	Unit
8B-S139-TAK	1 mL Tubes, Recombinant Enzyme	100/Box
8B-S322-DAK	1 mL Tubes, Non-Recombinant Enzyme	100/Box
8E-S139-TGA	96-Well Plate, Recombinant Enzyme	1/Box
8E-S322-DGA	96-Well Plate, Non-Recombinant Enzyme	1/Box
8N-S323-TUK	2 mL Centrifuge Tubes, Recombinant and Non-Recombinant Enzyme	100/Box
8E-S323-TGA	96-Well Plate Plus 30 mg/well, Recombinant/Non-Recombinant Enzyme	1/Box
8B-S323-UGA	96-Well Plate Plus 60 mg/well, Recombinant/Non-Recombinant Enzyme	1/Box



Kinetex® Core-Shell LC Columns

2.6 μm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges†
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
Phenyl-Hexyl	00B-4495-Y0	00D-4495-Y0	00F-4495-Y0	AJ0-8781 for 3.0 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJ0-9000

Vacuum Manifolds

Part No.	Description	Unit
12-Position Vacuum Manifold for Tubes*		
AH0-6023	12-Position Vacuum Manifold Set, complete assembly	ea
24-Position Vacuum Manifold for Tubes*		
AH0-6024	24-Position Vacuum Manifold Set, complete assembly	ea
96-Well Plate Manifold		
AH0-8950	96-Well Plate Manifold, Universal with vacuum gauge	ea



* Manifolds include: Vacuum-tight glass chamber, vacuum gauge assembly, polypropylene lid with gasket, male and female luers and yellow end plugs, stopcock valves, collection rack assemblies, polypropylene needles, lid support legs. Waste container included with 12-position manifold.

Presston™ 1000 Positive Pressure Manifold

Part No.	Description
AH1-7033	Presston 1000 Positive Pressure Manifold, 96-Well Plate



Phenomenex warrants the Presston 1000 will be free of defects in materials and workmanship under normal installation, use, and maintenance for a period of 12 months following delivery. Please visit www.phenomenex.com/Presstonwarranty for complete warranty information.



Australia

t: +61 (0)2-9428-6444
 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
 t: +32 (0)2 511 8666 (Dutch)
 beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
 info@phenomenex.com

China

t: +86 400-606-8099
 cninfo@phenomenex.com

Denmark

t: +45 4824 8048
 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
 nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
 franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
 anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
 indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
 eireinfo@phenomenex.com

Italy

t: +39 051 6327511
 italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
 tecnicomx@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
 nzinfo@phenomenex.com

Norway

t: +47 810 02 005
 nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
 pl-info@phenomenex.com

Portugal

t: +351 221 450 488
 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
 sginfo@phenomenex.com

Spain

t: +34 91-413-8613
 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
 twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
 ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
 info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
 info@phenomenex.com

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark of Phenomenex. β -Gone, SecurityGuard, MidBore, BE-HAPPY, and Presston are trademarks of Phenomenex. IMCSzyme is a registered trademark of Integrated Micro-Chromatography Systems, LLC. SCIEX is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEXTM is being used under license.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.