

# Increased Sensitivity of THC and Metabolites using $\beta$ -Gone™ and Kinetex™ 2.6 $\mu$ m C18 Column by LC-MS/MS

Dawn Chen, PhD, Shahana Wahab Huq, Zara Jalali, and Bryan Tackett, PhD  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

## Overview

THC, or tetrahydrocannabinol, is the chemical responsible for the psychological effects of cannabis. THC and its metabolites are tested in urine samples but need to undergo hydrolysis in the sample pre-treatment prior to extraction. The Dilute-and-Shoot method simply dilutes the sample after pre-treatment hydrolysis and the sample is then injected onto an LC-MS system. This method does not remove any matrix interferences that might exist which could clog an LC column due to exposure to a higher percentage of organic material during the gradient of the mobile phases in a chromatographic run.  $\beta$ -Gone targets and removes  $\beta$ -glucuronidase from hydrolyzed urine samples without requiring additional method development.

In this application note, we show improved sensitivity by an increase in intensity and signal to noise (S/N) ratio of THC and its metabolites using the  $\beta$ -Gone centrifuge tubes compared to the Dilute-and-Shoot method. The high efficiency Kinetex 2.6  $\mu$ m C18 LC column was used for chromatographic separation.

## Sample Preparation – Dilute-and-Shoot

**Sample Pre-treatment:** Combine 200  $\mu$ L of urine sample spiked with 20  $\mu$ L analyte mixture, 60  $\mu$ L hydrolysis buffer, and 20  $\mu$ L of IMCSzyme® RT enzyme. Incubate at room temperature for 15 minutes.

**Dilute:** 20  $\mu$ L from above 10 times with 0.1 % Formic acid in water.

**Inject:** 10  $\mu$ L

## Sample Preparation – $\beta$ -Gone

**Sample Pre-treatment:** Combine 200  $\mu$ L of urine sample spiked with 20  $\mu$ L analyte mixture, 60  $\mu$ L hydrolysis buffer, and 20  $\mu$ L of IMCSzyme RT enzyme. Incubate at room temperature for 15 minutes.

**Combine:** 200  $\mu$ L from above with 133  $\mu$ L 0.1 % Formic acid in methanol.

**Load:** Sample in  $\beta$ -Gone 2 mL centrifuge tube (Part No.: [8N-S323-TUK](#)). Mix by inverting tube 10 times.

**Centrifuge:** For 10 min at 14800 rpm.

**Inject:** 10  $\mu$ L supernatant.

## LC Conditions

**Column:** Kinetex 2.6  $\mu$ m C18

**Dimensions:** 50 x 2.1 mm

**Part No.:** [00B-4462-AN](#)

**Mobile Phase:** A: 0.1 % Formic acid in water

B: 0.1 % Formic acid in methanol

Gradient:	Time (min)	%B
	0.25	68
	5.25	70
	7.75	80
	7.95	100
	8.95	100
	9.15	68
	11.0	68

**Flow Rate:** 3 mL/min

**Injection Volume:** 10  $\mu$ L

**Temperature:** 40 °C

**LC System:** Agilent® 1260 Infinity

**Detection:** MS/MS

**Detector:** SCIEX® API 4000 QTRAP®

## MS Conditions

**Ion Source:** Positive or Negative

**Source Temperature:** 600 °C

**GS1:** 50

**GS2:** 50

**CUR:** 10


**IS:** +5500 or -4500

**Table 1.** MRM Transitions


Peak No.	Analyte	Q1 (m/z)	Q2 (m/z)
1	D9 THC-OH 1	331	193.1
1	D9 THC-OH 2	331	201.1
2	D9 THC-COOH 1	343	245
2	D9 THC-COOH 2	343	191
2	D9 THC-COOH 3	343	299.4
3	D9 THC	315.2	193.1
3	D9 THC ql 1	315.2	245
3	D9 THC ql 3	315.2	259.1
4	D8 THC	315.2	192.95
4	D8 THC ql 1	315.2	135.1



**Figure 1.** Separation of THC and Metabolites using Dilute-and-Shoot Sample Preparation and the Kinetex™ 2.6 µm C18 Column



**Figure 2.** Separation of THC and Metabolites using β-Gone™ Sample Preparation and the Kinetex 2.6 µm C18 Column



**Table 2.** Intensities and Signal to Noise Ratios of THC and Metabolites

		Dilute-N-Shoot	β-Gone	β-Gone/Dilute-and-Shoot
<b>D9 THC-OH</b>	Intensity	82525.8	528999.7	6.4
	% CV (n=3)	-	4.8	-
	S/N	218.7	484.3	2.2
<b>D9 THC-COOH</b>	Intensity	33604.2	196110.5	5.8
	% CV (n=3)	-	4.5	-
	S/N	33.9	135.5	4.0
<b>D9 THC</b>	Intensity	32885.1	140385.0	4.3
	% CV (n=3)	-	3.1	-
	S/N	20.0	76.9	3.8
<b>D8 THC</b>	Intensity	3587.9	7132.4	2.0
	% CV (n=3)	-	6.5	-
	S/N	3.5	12.0	3.4



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Kinetex,  $\beta$ -Gone, and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX and QTRAP are registered trademarks of AB SCIEX Pte. Ltd. IMCSyme is a registered trademark of Integrated Micro-chromatography Systems.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc. or Integrated Micro-chromatography Systems.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.  
© 2022 Phenomenex, Inc. All rights reserved.

