

TN-1326

Defining Robustness of Pesticide and Mycotoxin Analysis in Cannabis Matrices

Diana Tran¹, Karl A. Oetjen, PhD¹, Scott Krepich¹, Matthew Noestheden, PhD², and Bryan Tackett, PhD³

¹AB Sciex LLC, 500 Old Connecticut Path, Framingham, MA 01701, USA

²AB Sciex LP, 71 Four Valley Drive, Concord, Ontario L4K 4V8, Canada

³Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Introduction

Sample clean-up is an important step for high-throughput LC-MS/MS analyses. The more contaminants that are removed during the clean-up step, the longer the LC-MS/MS will be able to maintain the required sensitivity. Unfortunately, cannabis matrices contain high concentrations of cannabinoids, waxes, terpenes, and other secondary metabolites which present a significant analytical challenge. These compounds have the potential to interfere with the analysis of pesticides, making it difficult to meet the ng/g sensitivity levels required by most recreational United States regulations and Canadian regulations.

In this technical note, the unique selectivity and increased retention for both polar and non-polar analytes of the Luna Omega Polar C18 column and the robustness of the SCIEX® Triple Quad™ 6500+ system was evaluated by injecting a cannabis flower extract 830 times with no system maintenance. The cannabis flower was spiked with a mixture of commonly monitored cannabis pesticides and the peak area of these pesticides was monitored over time, with and without internal standard correction.

Sample Preparation

A 1:10 dilution was performed using 5 g of homogenized cannabis flower extracted in 50 mL of 0.1 % Formic Acid in Acetonitrile. Extracts were winterized at -20 °C for 2 hours before filtration with 0.2 µm PTFE syringe filters. The extract was fortified with an analytical pesticide mixture and vortexed before being dispensed into equal 1 mL aliquots to be stored at 4 °C prior to LC-MS/MS analysis.

LC Conditions

Column: Luna™ Omega 3 µm Polar C18

Dimensions: 150 x 3.0 mm

Part No.: [00F-4760-YO](#)

Mobile Phase: A: 5 mM Ammonium Formate + 0.1 % Formic Acid in Water

B: 5 mM Ammonium Formate + 0.1 % Formic Acid in Methanol

Gradient: Time (min)	%B
0	55
0.5	55
2.5	80
8.5	90
12.5	100
16.5	100
17	0
20	0

Flow Rate: 0.8 mL/min

Injection Volume: 2 µL

Temperature: 40 °C

LC System: SCIEX ExionLC™

Detection: MS/MS

Detector: SCIEX Triple Quad 6500+

Results and Discussion

Cannabis flower extracts are a particularly challenging matrix. Very few LC-MS/MS robustness studies have been conducted with this matrix, without MS system maintenance over a prolonged duration. When determining instrument stability using this type of robustness test, normalizing the analyte peak area to an internal standard (IS) area can be misleading, as the response from the internal standard and the native pesticide(s) are likely to change proportionately. Therefore, the IS ratio will stay consistent across many injections, as shown for Carbofuran (**Figure 1**, left), inaccurately suggesting ideal system performance despite the harsh conditions employed in this study.

The true measure of instrument robustness must be an evaluation of the uncorrected peak area as a function of time. Without MS system maintenance and given the conditions of this study, a decrease in peak area may be expected, as observed when the raw Carbofuran area is plotted (**Figure 1**, right).

However, these data show that the SCIEX Triple Quad 6500+ system, coupled with the Luna Omega Polar C18 column, achieves sensitivity that meets regulatory limits and reliably detects pesticides of interest in a complex matrix. These features persist over the analysis of 830 cannabis samples without cleaning the MS system.

An example of this robust sensitivity can be seen with Acequinocyl, which is hydrolytically unstable. It has poor ionization efficiency and coelutes with numerous cannabinoids late in the gradient. For all 830 injections, Acequinocyl was detected at a concentration 40x lower than Oregon Regulatory Limits (**Figure 2**, top). Additionally, Avermectin B1a, which is known for its thermal lability, was detected at a concentration 10x lower than Oregon Regulatory Limits after the 830 matrix injections (**Figure 2**, bottom).

It should be noted that although this study using this instrument shows that system maintenance isn't required during the study, regular maintenance of the column, LC system, and MS system should be done.



Figure 1. 830 Replicate Injections of Cannabis Flower Matrix Using a Luna™ Omega Polar C18 Column.





Figure 2. Stable Pesticide Peak Areas Across 830 Injections Using a Luna Omega Polar C18 Column.



Conclusion

The reality of analyzing a highly contaminating matrix is an inevitable decrease in sensitivity. In this technical note, it is shown that the way instrument robustness data is organized and presented can fail to capture changes in sensitivity over time. It is therefore important to assess both the ion ratio reproducibility (Figure 2) and the raw peak area reproducibility (Figure 1), as this will inform practical considerations in a testing lab such as how often an MS system must be cleaned to maintain sensitivity requirements.

Luna™ Omega Ordering Information

3 μm MidBore™ Columns (mm)	SecurityGuard™ Cartridges (mm)			
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0* /10pk
Polar C18	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJ0-7600
PS C18	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJ0-7605
C18	00B-4784-Y0	00D-4784-Y0	00F-4784-Y0	AJ0-7611
SUGAR	—	—	00F-4775-Y0	AJ0-4496

for ID: 2.0 – 3.0 mm

*SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)

Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna, SecurityGuard, Midbore, and BE-HAPPY are trademarks of Phenomenex. SCIEX is a registered trademark and ExionLC and Triple Quad are trademarks of AB SCIEX Pte. Ltd.

Disclaimer

Comparative separations may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

