

# APPLICATIONS

## THC Remediation from CBD Extract to Demonstrate Running Flash Chromatography on a Prep HPLC Instrument

J Preston, Sean Orłowicz  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

### Overview

Flash chromatography is often used for quick purifications during intermediate steps of pharmaceutical product synthesis. These intermediate purifications allow for higher yields of subsequent steps and ultimately more end product. Unlike LC columns, flash cartridges are considered disposable and as a result require less cleaning and maintenance since they are often disposed of after 1-3 purification runs.

The 2018 Farm Bill<sup>1</sup> effectively legalized hemp, with many restrictions that challenge hemp processors. Passed into law in December of 2018, this bill removed hemp from the definition of marijuana in the Controlled Substances Act<sup>2</sup> and provided an important and clearer definition of “hemp” in federal law. In Section 297A<sup>1</sup>, “The term ‘Hemp’ means the plant *Cannabis sativa* L, and any part of that plant, including the seeds thereof and all derivatives, extracts, cannabinoids, isomers, acids, salts, and salts of isomers, whether growing or not, with a delta-9 tetrahydrocannabinol concentration of not more than 0.3 percent on a dry weight basis”. This new definition places added responsibilities onto hemp cultivators and processors to meet the requirements to remain less than 0.3% Δ9-THC in their products. In addition to cultivation aspects, some common extraction techniques can concentrate cannabinoids, yielding extracts with greater Δ9-THC content. A popular technique for selectively removing Δ9-THC from these complex mixtures is large scale chromatography.

When compared to Prep HPLC instruments, Flash instruments typically have basic functions, including low pressure pumps, a fraction collector and a UV-VIS detector. This simple setup allows for easier maintenance and a lower purchase cost. Likewise, flash cartridges are available in many common stationary phase chemistries suitable for reversed phase, normal phase and even HILIC separations. The latest high-end Flash instruments have been equipped with pumps that reach higher pressures and include sophisticated detectors. These systems are capable of running Prep HPLC columns. A system that can run both Flash columns and HPLC columns is very versatile tool for a Prep Chromatography Lab. In most cases, a lab that already has a Prep HPLC (not designed as a Flash Instrument) can use their system to run Flash Columns. The work presented here will describe the adjustments necessary to run Flash columns on any typical Prep HPLC system in order to gain this extra level of versatility.

### Materials and Methods

#### Reagents and Chemicals

All chemicals and reagents were obtained from the Sigma-Aldrich Company (St. Louis, MO, USA).

Water purification via Sartorius® arium® Comfort II (Goettinger, Germany). CBD Sample is a hemp extract provided by local commercial supplier. This material contained the maximum amount of THC allowed before it is considered a schedule 1 drug.

#### Instrument Descriptions

CHEETAH™ MP Series Flash Chromatography System, Bonna-Agela Technologies OCTOPUS Purification System, Bonna-Agela Technologies, Shimadzu® LC-20AP pump, Agilent® 1100 HPLC

### Column Descriptions for Analytical, Prep and Flash

**Flash column:** Claricep™ Screw-on Spherical C18

**Dimensions:** I-Series, 20g/column

**Part No.:** SN-SO230020-0

**HPLC Prep column:** Luna® 5 µm C18(2)

**Dimensions:** 100 x 21.2 mm

**Part No.:** 00D-4252-P0-AX

**HPLC Analytical column:** Luna® 5 µm C18(2)

**Dimensions:** 100 x 4.6 mm

**Part No.:** 00D-4252-E0

### Conditions for Flash, Analytical, and Prep

**Columns:** As listed above

**Mobile Phase:** Mobile Phase: Isocratic: 10% Water with 0.1% Acetic Acid / 90% Acetonitrile

**Detection:** UV @ 220 nm

Hemp Extract, 500 mg/mL in 1-1-1

**Sample:** Tetrahydrofuran-Acetonitrile-Water

**Injection:** 100 µL (unless specified)

**Flow Rate:** 20 mL/min (unless specified)

**Temperature:** Ambient

**Injection Volume:** 100 µL (unless specified)

### Conditions for Fraction Analysis

**Column:** Kinetex® 2.6 µm C18

**Dimensions:** 150 x 4.6 mm

**Part No.:** 00F-4462-E0

**Mobile Phase:** Isocratic: 24% Water with 0.1% Formic Acid / 76% Methanol

**Flow Rate:** 0.75mL/min

**Injection Volume:** 5 µL

**Temperature:** Ambient

**Detector:** UV @ 228 nm

**Sample:** Fractions directly from the fraction collector

## Instrumentation for Flash Chromatography



Traditional flash chromatography instruments contain a low pressure pump, column holder, basic UV detector and fraction collector.  
(Ex. CHEETAH™ MP Series Flash Chromatography system)



High pressure flash instruments have been commercialized in the last 10 to 15 years and are capable of running both prep HPLC columns and Flash Chromatography columns. They utilize a pump that is capable of higher pressures, valves to switch between flash and prep columns, improved UV and other types of detectors, along with a fraction collector.  
(Ex. OCTOPUS Purification system)



Prep HPLC systems have been around for many years. They have always had the capabilities of the recently available high pressure flash systems. The main difference is the fittings used for prep HPLC columns and flash columns. If we could use suitable adaptors and fittings, we can use a prep HPLC system to run flash columns.



Comparing the low pressure flash, high pressure flash/prep and prep instruments above, all of the systems have a pump, detector and fraction collector. If a fraction collector is added to an analytical HPLC system, we could use it to run flash or prep chromatography. The maximum flow rate is the most significant limitation.

## Caution About Pressure

Another significant difference between prep HPLC and flash columns is the pressure limitations. Prep HPLC columns are designed for high pressure applications. The typical pressure limits for steel prep HPLC columns can range between 100 bar and 400 bar. Flash chromatography columns are a low pressure technique and these columns are not designed for high pressures. The typical pressure limit for plastic flash columns is usually below 20 bar.

The initial instruments designed for flash chromatography had pumps that could only function in the pressure range suitable for flash columns. The recent instruments (capable of using prep HPLC columns) are capable of over pressuring the flash column. If an HPLC system (prep or even analytical) was used with a flash column, over pressuring the column must also be a concern.

The Luer fittings on flash columns are typically rated to 35 bar. A Luer fitting can fail due to pressure as a partial leak or a total separation of the male and female Luer components. HPLC systems are designed with safety features such as leak sensors and low pressure limit settings that can turn off the eluent flow when a Luer fitting fails. Some individual Luer fittings, when assembled, can hold substantially more than 35 bar. If this happens and a high pressure instrument is being used, the plastic flash column body could be subjected to pressure much higher than their rated pressure limit. When the flash column body fails due to pressure, there could be a leak at the seal on the inlet side where the column is closed or the body could rupture. The leaking body is easily handled but the ruptured body could be a significant safety issue and should be addressed.

## Fittings




The main differences between prep HPLC and Flash chromatography columns is the type of fittings used for the column inlet and outlet. Prep columns use high pressure fittings in combination with either 1/16" or 1/8" stainless steel or PEEK tubing. Flash columns use Luer fittings, with the inlet typically being a Luer lock fitting and the outlet being either a Luer/Threaded or Snap lock fittings.




There are fittings and adaptors available that make Luer fittings work with either 1/16" or 1/8" high pressure tubing. There are other fittings and adaptors available that makes low pressure tubing work with high pressure fittings.


**Scalability: Analytical, Prep and Flash**



**Figure 1**  
Chromatogram of CBD extract on analytical scale  
Luna<sup>®</sup> 5  $\mu$ m C18(2) 100 Å, LC Column, 100 x 4.6mm  
Flow Rate = 1 mL/min  
Injection Volume = 5  $\mu$ L




**Figure 2**  
Chromatogram of CBD extract on prep scale  
Luna<sup>®</sup> 5  $\mu$ m C18(2) 100 Å, LC Column, 100 x 21.2mm  
Flow Rate = 20 mL/min  
Injection Volume = 100  $\mu$ L



**Figure 3**  
Chromatogram of CBD extract on flash  
Claricep<sup>™</sup> Screw-on Spherical C18, 20-35  $\mu$ m, 100Å,  
Flash column, 20g  
Flow Rate = 20 mL/min  
Injection Volume = 100  $\mu$ L

### Fractions Collected from Axia<sup>®</sup>

From a single injection of the crude sample run on an Axia prep column, 13 different fractions were collected. See Figure 4. Each fraction was evaluated with analytical scale HPLC methodology. See Figure 5.




**Figure 4**  
Chromatogram showing fraction collection  
From CBD extract on 100 x 21.2 mm, Axia prep column  
Flow Rate = 20 mL/min  
Injection Volume = 100  $\mu$ L




**Figure 5**  
Chromatograms of collected fractions  
From CBD extract on 100 x 21.2mm, Axia prep column  
Flow Rate = 20 mL/min  
Injection Volume = 100  $\mu$ L

### Fractions Collected from Flash

From a single injection of the crude sample run on a flash column, 4 different fractions were collected. See Figure 6. Each fraction was evaluated with analytical scale HPLC methodology. See Figure 7.



**Figure 6**  
Chromatogram showing fraction collection  
From CBD extract on 20g Flash column  
Flow Rate = 20 mL/min  
Injection Volume = 100  $\mu$ L



**Figure 7**  
Chromatograms of collected fractions  
From CBD extract on 20g Flash column  
Flow Rate = 20 mL/min  
Injection Volume = 100  $\mu$ L

## Results and Discussion

One of the best features of chromatography is scalability. Provided the instrumentation and column characteristics can be maintained, chromatographic performance can be scaled from small diameter columns to larger and larger columns. This principle is demonstrated here by scaling a CBD extract sample from a 4.6mm ID diameter analytical scale to a 21.2mm ID diameter prep scale. See Figure 1 and Figure 2. The principle of scaling can be applied to Flash chromatography as well, but the differences in column characteristics will cause changes in the performance. See Figure 1 and Figure 3.

In this case of THC removal from a CBD sample, a flash column easily provides enough resolution to perform this goal. See Figure 5 and Figure 6. The Axia packed prep column provides chromatography that is essentially equivalent to the analytical column. With the prep column, 13 fractions were collected and the components were isolated as significantly pure fractions. See Figure 4 and Figure 5.

## Conclusion

Any chromatography system that can pump eluent can be used with Flash Chromatography columns. There are fittings and adapters that can connect the different types of tubing and components. This allows for greater versatility. One word of caution is to be aware of the operating pressure when using plastic flash cartridges.

Flash Chromatography can be very useful when the resolution of peaks is not very demanding. Prep Chromatography is very powerful. If a separation can be done at the analytical scale it can scaled to prep.

Chromatography (HPLC or Flash) can be used for THC remediation from hemp extracts. This work was done with acetonitrile but the chromatography can easily be converted to another strong eluent such as ethanol.



# APPLICATIONS

## Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/ChatNow](http://www.phenomenex.com/ChatNow).

### Australia

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

### Poland

t: +48 (12) 881 0121  
pl-info@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

### Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**  
t: +1 (310) 212-0555  
info@phenomenex.com

## [www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

## References

1. 2018 Farm Bill, Part of the Agricultural Act of 2018, Pub. L. 115-334
2. Controlled Substances Act, 21 U.S.C. § 802(16).