



## TN-1311

# Micro-flow Bottom-up Proteomics or Detection of Peptides at Trace Levels Using a Kinetex™ XB-C18 Column

Roxana Eggleston-Rangel, Jason Anspach, PhD, Namrata Saxena, and Bryan Tackett, PhD  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Introduction

Proteomics research is often challenged by the complex nature of its samples. The identification and quantification of peptides of interest can prove particularly difficult when dealing with trace amounts of samples, small volumes, and a variety of complex biological sample matrices. While mass spectrometry (MS) coupled to liquid chromatography (LC) has proven to be a powerful tool in the area of discovery proteomics, miniaturization of liquid chromatography columns can further overcome the sensitivity challenges associated with smaller sample volumes and complex sample matrices. Reducing the inner-diameter of the LC column results in reduction in the chromatographic diffusion. This, in turn, can increase the method's sensitivity up to 60X in comparison to analytical dimensions, allowing for more efficient MS/MS sampling and thus a higher number of unique peak identifications. Another way to improve ion sensitivity is to select a core-shell particle. These particles consist of a compact non-porous core surrounded by a porous silica layer and bonded with the stationary phase primary and secondary ligands. By using core-shell technology, the peptide diffusion path is significantly reduced. This results in sharper chromatographic peaks leading to increased ion intensities and higher peak capacities. Most of the proteomics studies done on core-shell particles are focused on analytical flow chromatography which unfortunately does not truly display the benefits of core-shell particles. Therefore, in this technical note we have compared and investigated sub-2 μm fully porous particles with core-shell particles in a micro-LC format to demonstrate the benefits of core-shell particles in peptide separation performance as well as relative peptide and protein identifications in a commercially available HeLa cell lysate tryptic digest.

### LC Conditions

<b>Columns:</b>	Kinetex 2.6 μm XB-C18 (core-shell) Waters® nanoEase M/Z Peptide BEH 1.7 μm C18 (sub-2 μm fully porous) Thermo Scientific® Acclaim™ PepMap™ 100 2.0 μm C18 (fully porous)	
<b>Dimensions:</b>	50 x 0.3 mm (Kinetex, Waters, Thermo Scientific) 150 x 0.3 mm (Kinetex)	
<b>Part Nos.:</b>	<a href="#">00B-4496-AC</a> (50 mm) <a href="#">00F-4496-AC</a> (150 mm)	
<b>Mobile Phase:</b>	A: 0.1 % Formic Acid in Water B: 0.1 % Formic Acid in Acetonitrile	
<b>Gradient:</b>	<b>Time (min)</b>	<b>%B</b>
	0	2
	20	40
	20.1	80
	25	80
	25.1	2
	30	2
<b>Flow Rate:</b>	10 μL/min	
<b>Injection:</b>	1 μL	
<b>Temperature:</b>	40 °C	
<b>Detector:</b>	Q Exactive™ Plus Orbitrap™	
<b>System:</b>	NanoLC™ 425 (SCIEX®)	
<b>Detection:</b>	MS	

### MS/MS Conditions

<b>Source Type:</b>	Ion Max Source with HESI-II probe and 50 μm ID needle
<b>Scan Type:</b>	Full MS
<b>Polarity:</b>	Positive
<b>Resolution:</b>	70,000
<b>Scan Range:</b>	400 to 1650 m/z
<b>AGC Target:</b>	3e6
<b>Maximum IT:</b>	15 ms
<b>Scan Type:</b>	dd-MS <sup>2</sup>
<b>Resolution:</b>	35,000
<b>AGC Target:</b>	1e5
<b>Maximum IT:</b>	45 ms
<b>CE:</b>	29
<b>Isolation Window:</b>	1.2 m/z




## Results and Discussion

Bottom-up proteomics LC-MS/MS workflows rely on the fragmentation by parent ions to provide MS2 spectra containing enough peptide sequence information that can be assigned to a particular protein. Through increasing the parent ion intensity in the MS1 spectrum, higher fragment ion intensities (MS2s) will result (**Figure 2**). This provides high quality spectra which, in turn, generate larger numbers of unique peptide and protein identifications of any proteomics sample. It has been successfully established that core-shell particles reduce analytes' peak widths due to reduced diffusion in them. Narrower, sharper peaks and high ion intensities are the result. In contrast, sub-2  $\mu\text{m}$  fully porous particle columns have also been known to provide the same effect in ion intensity and peak shapes as core-shell particles, but with the disadvantage of generating much higher pressures due to the fully porous smaller particle size (**Table 1**).

As was expected, the core-shell column and sub-2  $\mu\text{m}$  columns produced Total Ion Chromatograms (TICs) with a higher maximum ion intensity than the fully porous particle column (**Figure 1**). Also, the average peak widths were narrower in the core-shell column and the sub-2  $\mu\text{m}$  column with almost identical increased peak capacities than in the fully porous column (**Table 1**). However, the sub-2  $\mu\text{m}$  column's pressure was twice as much in comparison to the core-shell particle column's pressure, restricting the modification of flow rates if required. Moreover, the relative number of peaks that were repeatedly sequenced was significantly higher for the fully porous column due to its higher average peak widths, which leads to increased redundant MS2 fragment data generation (**Table 1**).

**Figure 3** shows the relative number of unique peptides and proteins identified when using the core-shell, sub-2  $\mu\text{m}$ , and the fully porous column. Again, as was expected, the core-shell and sub-2  $\mu\text{m}$  columns resulted in a higher number of protein and peptide identifications. It is worth pointing out that the absolute number of unique groups found using the core-shell column, loading only 200 ng of digested protein, was 411. To improve protein IDs, a longer column can be used. Here, we also compared the number of protein IDs obtained with a 5 and a 15 cm long core-shell Kinetex™ 2.6  $\mu\text{m}$  XB-C18 column. The number of IDs increased by almost 50% by using a column three times longer (**Table 2**).


**Figure 1.** Total Ion Chromatograms of 200 ng of HeLa Cell Tryptic Digest. Maximum Ion Intensity is Shown.




**Table 1.** Average Peak Widths, Peak Capacities and Corresponding Column Pressures (bar) for a 20 Minute Gradient Run of 200 ng of HeLa Lysate Tryptic Digest.

Column	Average Peak Width	Peak Capacity	Column Pressure
Kinetex™ 2.6 μm XB-C18 (core-shell)	0.55	29	172
nanoEase M/Z Peptide BEH 1.7 μm C18 (sub-2 μm)	0.64	29	345
Acclaim™ PepMap™ 100 2.0 μm C18 (fully porous)	0.75	21	170

**Figure 2.** MS/MS Fragmentation of a 19 Amino Acid Long HeLa Peptide with m/z 934.9887+2, when Using a Kinetex XB-C18 Core-shell Column (Top), a nanoEase M/Z Peptide BEH C18 sub-2 μm Column (Middle) and an Acclaim PepMap 100 C18 Fully Porous Column (Bottom).



**Figure 3.** Relative Percent of Unique Peptide Identifications, Repeatedly Sequenced Peaks, and Unique Protein Identifications Found in a 20 min Gradient Run of 200 ng of HeLa Lysate Tryptic Digest.



**Table 2.** Number of Unique Protein Identifications Found Using a 5 cm Long Kinetex 2.6 µm XB-C18 Column and a 5 cm Long Kinetex 2.6 µm XB-C18 Column with 200 ng of HeLa Lysate Tryptic Digest and the % Protein ID Increase Found by Using a Longer Column.

Column	Protein Groups
5 cm Kinetex 2.6 µm XB-C18	411
15 cm Kinetex 2.6 µm XB-C18	901
% Protein ID Increase	45.62%

## Conclusions

Column miniaturization (reduction of column inner diameter) can be a powerful tool to increase ion intensities and thus the number of peptides and proteins identified in a particular sample. Apart from the smaller column inner diameters, core-shell technology can also further improve the peak shapes. This leads to improved ion intensities which enhances the number of unique IDs that can be detected in MS, as compared to what can be achieved through a conventional fully porous particle column. Another major advantage of core-shell technology is that the generated back pressure is almost half of what is generated with a sub-2 µm particle column. This is highly desirable in a micro-LC setting where complex sample components tend to accumulate in the column, further increasing the back pressures and leading to reduced column lifetimes.

In summary, it can be implied that miniaturized core-shell technology provides a useful medium to improve the identification of peptides and proteins in bottom-up proteomics workflows, while still maintaining the practical back pressure in the column.



## Kinetex™ Micro LC Columns Ordering Information

2.6 µm Micro LC Columns (mm)						
Phases	30 x 0.3	50 x 0.3	100 x 0.3	150 x 0.3	50 x 0.5	150 x 0.5
XB-C18	<a href="#">00A-4496-AC</a>	<a href="#">00B-4496-AC</a>	<a href="#">00D-4496-AC</a>	<a href="#">00F-4496-AC</a>	<a href="#">00B-4496-AF</a>	<a href="#">00F-4496-AF</a>
Biphenyl	—	<a href="#">00B-4622-AC</a>	—	<a href="#">00F-4622-AC</a>	<a href="#">00B-4622-AF</a>	—
C18	<a href="#">00A-4462-AC</a>	<a href="#">00B-4462-AC</a>	—	<a href="#">00F-4462-AC</a>	<a href="#">00B-4462-AF</a>	—
EVO C18	—	<a href="#">00B-4725-AC</a>	—	<a href="#">00F-4725-AC</a>	<a href="#">00B-4725-AF</a>	—
F5	—	<a href="#">00B-4723-AC</a>	<a href="#">00D-4723-AC</a>	<a href="#">00F-4723-AC</a>	<a href="#">00B-4723-AF</a>	—



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

🌐 **All other countries/regions**  
**Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Kinetex and BE-HAPPY are trademarks of Phenomenex. Thermo Scientific is a registered trademark and Acclaim, PepMap, Orbitrap, and Q Exactive are trademarks of Thermo Fisher Scientific. SCIEX is a registered trademark and nanoLC is a trademark of AB SCIEX Pte. Ltd. Waters is a registered trademark of Waters Technologies Corporation.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Thermo Fisher Scientific, Inc. or Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

