

# THANK YOU

for choosing Strata-X PRO  
Solid Phase Extraction (SPE) Products



Phenomenex products are available worldwide.

[www.phenomenex.com/mysupport](http://www.phenomenex.com/mysupport)

QUALITY  
MANAGEMENT SYSTEM  
CERTIFIED BY DIN EN  
- ISO 9001:2015 -

IL66990319\_W

# STRATA<sup>®</sup>X PRO

A Rapid Solid Phase Extraction Solution

## User Guide



 phenomenex<sup>®</sup>  
*breaking with tradition™*

# STRATA<sup>®</sup>X PRO

A Rapid Solid Phase Extraction Solution

1

## Prepare Your Sample

### Plasma/Serum

If the analyte of interest is an acid, 2% phosphoric acid can be used (20 µL 85% H<sub>3</sub>PO<sub>4</sub> to 1 mL of plasma (or serum) to disrupt the drug-protein interaction. If the analyte of interest is basic, 0.1 M sodium hydroxide can be used to disrupt the drug-protein interaction. After addition of acid or base, the sample should be vortexed for 20-30 seconds followed by centrifugation. The supernatant is now ready for further analysis. Alternatively, a zinc sulfate or strong organic solvent can be used to disrupt the protein interaction.

### Whole Blood

There are several pre-treatment strategies that can be followed for whole blood. If the target analyte is present in red blood cells, a hemolysis step is necessary.

- Hemolysis:** To 0.2 mL whole blood (spiked with analytes and internal standard) in a 1.2 mL centrifuge tube, add 400 µL of 2% zinc sulfate/80% methanol. Vortex for 10-20 seconds followed by centrifugation at 14,000 rpm for 10 minutes. Collect the supernatant for further analysis.

**Preparation of zinc sulfate/methanol:** Into a 100 mL volumetric flask add 20 mL water and 3.6 g ZnSO<sub>4</sub> · 7 · H<sub>2</sub>O. After the solution is clear and the salt crystals have dissolved, add 100% methanol. Refrigerate the solution at 2-8 °C for 7 days.

- Osmotic breakdown:** To 1 mL of whole blood add internal standard and 4 mL of distilled water. Mix/vortex and let stand for 5 minutes. Centrifuge at 670g for 10 minutes and discard the pellet. Adjust the pH of the supernatant accordingly with the addition of a buffer solution.
- Sonication:** Add 3-6 mL of appropriate pH buffer (such as potassium phosphate) to 1 mL of whole blood and sonicate for 15 minutes at room temperature.

### Urine

Enzymatic hydrolysis is necessary in case of conjugated forms (sulfated or glucuronide) of the analyte present. Enzymatic hydrolysis requires specific pH (pH 4-5) and temperature ranges. An acid or base hydrolysis can be performed as well, depending on the stability of the compound.

- Enzymatic hydrolysis:** To 500 µL sample (spiked with analyte and internal standard) add 100 µL acidic buffer (see below) and 100 µL β-glucuronidase. Vortex 5-6 seconds. Incubate in a water bath at 63 °C for 30 minutes. Transfer sample to a 96-well collection plate or autosampler vial. Seal and centrifuge for 10 minutes at 2,000 rpm.

**Preparation of acidic buffer (1.0 M acetate buffer, pH4.0):** Dissolve 3.0 g of glacial acetic acid and 4.1 g of sodium acetate in a 1 L volumetric flask.

- Base hydrolysis:** To 1 mL urine (spiked with analyte and internal standard) add 100 µL 10 N KOH. Mix, vortex, and hydrolyze for 20 minutes at 60 °C. Cool and adjust pH to 3.5- 4.0 (by adding 200 µL glacial acetic acid).
- Acid hydrolysis:** To 1 mL urine add 0.25 mL HCl in a screw capped test tube. Screw the tube top on loosely and heat in a boiling water bath for 60 minutes. Adjust to pH 7 (or as needed) with 1.0 N NaOH.

### Saliva

No hydrolysis is required for oral fluids and the generic protocol used for plasma/serum pretreatment may be followed.

### Tissue

Homogenize with organic or aqueous solvent depending upon analyte solubility. Settle, decant, centrifuge or filter supernatant. Perform direct Matrix Solid Phase Dispersion (MSPD) extraction on tissue.

2

## Find the Correct Volumes

Sample Matrix	Sorbent Mass	Strata-X PRO
Blood, serum, plasma	30mg	250µL
Urine	30mg	1 mL
Filtered tissue homogenates	60mg	100mg
Oral Fluid	30 mg	500µL
Environmental Samples	Sorbent Mass	Strata-X PRO
Water (particulate-free) drinking	500mg	100 - 400mL
Water (particulate-laden) rivers, runoff, etc.	500mg	100 - 400mL
Soil extracts	500mg	100g

**STRATA<sup>®</sup>X PRO**  
A Rapid Solid Phase Extraction Solution

## Recommended Sorbent Wash and Elution Volumes

Strata-X PRO Sorbent Mass	Practical Minimum Wash and Elution Volume	Recommended Wash and Elution Volumes
	4 bed volumes	8 bed volumes
10 mg	100µL	200µL
30 mg	300µL	600µL
60 mg	600µL	1.2 mL
100 mg	1 mL	2 mL
200 mg	2 mL	4 mL
500 mg	5 mL	10 mL

3

## Determine Your Method

## 2-Step Protocol

Non-retentive SPE method to help achieve the fastest extraction.

**Load**

1 mL Pre-treated sample/0.1 % Formic acid in Acetonitrile (1:4)  
Apply 5" Hg vacuum until all tubes or wells have cleared.

**Elute**

75 µL Water/0.1% Formic acid in Acetonitrile (1:4)  
Apply 5" Hg vacuum until all tubes or wells have cleared.

Protocols are written for 30 mg/1 mL tubes, adjust based on sorbent size.

## 3-Step Protocol

Rapid protocol to reduce matrix effects and increase recovery of polar analytes.

**Load**

500 µL Pre-treated sample/buffer\* (1:1)  
Apply 2-5" Hg vacuum until liquid is no longer visible above top frit.

**Wash**

600 µL 5% Methanol in Water.

**Elute**

600 µL 0.1% Formic acid in Acetonitrile/Methanol (90:10)  
Apply 2-5" Hg vacuum for 1 minute.

Protocols are written for 30 mg/1 mL tubes, adjust based on sorbent size

\*Select a buffer that maximizes the hydrophobicity of the analytes. For example, if the analyte is basic, dilute with a base.

## Presston™ 1000 Positive Pressure Manifold

Part No.	Description
AH1-7033	Presston 1000 Positive Pressure Manifold, 96-Well Plate

Phenomenex warrants the Presston 1000 will be free of defects in materials and workmanship under normal installation, use, and maintenance for a period of 12 months following delivery. Please visit [www.phenomenex.com/Presstonwarranty](http://www.phenomenex.com/Presstonwarranty) for complete warranty information.



## Accessories

Part No.	Description	Unit
<b>Collection Plates (deep well, polypropylene)</b>		
AH0-7192	96-Well Collection Plate, 350µL/well	50/pk
AH0-7193	96-Well Collection Plate, 1 mL/well	50/pk
AH0-7279	96-Well Collection Plate, 1 mL/well Round, 7 mm	50/pk
AH0-7194	96-Well Collection Plate, 2 mL/well	50/pk
AH0-8635	96-Well Collection Plate, 2 mL Square/Round-Conical	50/pk
AH0-8636	96-Well Collection Plate, 2 mL Round/Round, 8 mm	50/pk
AH1-7036	96-Well Collection Plate, 2 mL Low Protein Binding	120/pk

## Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which can be viewed at [www.phenomenex.com/TermAndConditions](http://www.phenomenex.com/TermAndConditions)

## Trademarks

Strata is a registered trademark and Presston is a trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.