

Application Guide

# CLINICAL SAMPLE PREPARATION

SIMPLE | FAST | CLEAN

Revision: 0

PHEN-RUO-00054

 **phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>

 [www.phenomenex.com/ClinicalSP](http://www.phenomenex.com/ClinicalSP)

# Table of Contents

Sample Preparation Selection Chart .....	p. 3
Sample Pre-treatment Overview .....	p. 4
<b><math>\beta</math>-Glucuronidase Removal</b> .....	p. 5
<b>Phospholipid Removal Solutions</b> .....	p. 6
<b>Simplified Liquid Extraction (SLE)</b> .....	pp. 7-11
Corticosteroids from Plasma .....	p. 8
Vitamin A and E from Human Serum .....	p. 9
NSAIDs from Plasma.....	p. 10
Acids, Neutrals, and Bases from Urine .....	p. 11
<b>Solid Phase Extraction (SPE)</b> .....	pp. 12-17
Urinary Catecholamines.....	p. 13
Underivatized Methylmalonic Acid (MMA) from Plasma.....	p. 14
Comprehensive Drug Research Panel from Oral Fluid Collection Devices.....	p. 15
Expanded Comprehensive Drug Research Panel from Oral Fluid Collection Devices .....	p. 16-17
Ordering Information .....	pp. 18-19

# Select the Appropriate Sample Prep Technique for Your Key Requirements



If Phenomenex products in this brochure do not provide at least equivalent separation as compared to other products of the same phase and dimensions, return the product with your comparative data within 45 days for a FULL REFUND.

# Sample Pre-treatment

Due to their nature, bioanalytical samples often require a pre-treatment step prior to further cleanup.



## Plasma/Serum

If the analyte is an acid, 2 % phosphoric acid can be used (20  $\mu$ L 85 % Phosphoric acid to 1 mL of plasma or serum) to disrupt the drug-protein interaction. If the analyte is basic, 0.1 M sodium hydroxide can be used to disrupt the drug-protein interaction. After addition of acid or base, the sample should be vortexed for 20-30 seconds followed by centrifugation. The supernatant is now ready for further analysis.



## Whole Blood

There are several pre-treatment strategies that can be followed for whole blood. If the target analyte is present in red blood cells, a hemolysis step is necessary.

### Hemolysis

To 0.2 mL whole blood (spiked with analytes and internal standard (I.S.) in a 1.2 mL centrifuge tube, add 400  $\mu$ L of 2 % Zinc sulfate/80 % Methanol. Vortex for 10-20 seconds followed by centrifugation at 14,000 rpm for 10 minutes. Collect the supernatant for further analysis.

Preparation of zinc sulfate/methanol: Into a 100 mL volumetric flask add 20 mL water and 3.6 g Zinc sulfate  $\cdot$  7 H<sub>2</sub>O. After the solution is clear and the salt crystals have dissolved, add 100 % Methanol. Refrigerate the solution at 2-8 °C for 7 days.

### Osmotic Breakdown

To 1 mL of whole blood add internal standard and 4 mL of distilled water. Mix/vortex and let stand for 5 minutes. Centrifuge at 670 g for 10 minutes and discard the pellet. Adjust the pH of the supernatant accordingly with the addition of a buffer solution.

### Sonication

Sonicate 1 mL whole blood for 15 minutes at room temperature. Add 3-6 mL of an appropriate pH buffer (such as potassium phosphate buffer). Mix/vortex. Let stand for 5 minutes. Centrifuge at 670 g for 15 minutes. Analyze supernatant.



## Urine

Enzymatic hydrolysis is necessary in case of conjugated forms (sulfated or glucuronide) of the analyte present and requires specific pH (~ 4-5) and temperature ranges. Depending on the compound's stability, an acid or base hydrolysis can be performed as well.

### Enzymatic Hydrolysis

To 500  $\mu$ L sample (spiked with analyte and I.S.) add 100  $\mu$ L acidic buffer (see below) and 20  $\mu$ L  $\beta$ -glucuronidase. Vortex 5-6 seconds. Incubate in a water bath at 63 °C for 30 minutes. Transfer sample to a 96-well collection plate or autosampler vial. Centrifuge for 10 minutes at 2,000 rpm. Preparation of acetic buffer (1.0 M acetate buffer, pH4.0): Dissolve 3.0 g of glacial acetic acid and 4.1 g of sodium acetate in a 1L volumetric flask.

### Base Hydrolysis

To 1 mL urine (spiked with analyte and I.S.) add 100  $\mu$ L 10 N Potassium hydroxide. Mix, vortex, and hydrolyze for 20 minutes at 60 °C. Cool and adjust pH to 3.5- 4.0 (by adding 200  $\mu$ L glacial acetic acid).

### Acid Hydrolysis

To 1 mL urine add 0.25 mL Hydrogen chloride in a screw capped test tube. Screw the tube top on loosely and heat in a boiling water bath for 60 minutes. Adjust to pH 7 (or as needed) with 1.0 N Sodium hydroxide.



## Saliva

No hydrolysis is required for oral fluids and the generic protocol used for plasma/serum pre-treatment may be followed.




## Tissue

Homogenize with organic or aqueous solvent depending upon analyte solubility. Settle, decant, centrifuge or filter supernatant. Perform direct Matrix Solid Phase Dispersion (MSPD) extraction on tissue.

β-Glucuronidase is a large enzyme that is present in urine samples after a urine hydrolysis is performed. In order to protect columns, β-Gone β-Glucuronidase Removal Products selectively remove the β-Glucuronidase enzyme and increase the sensitivity of analysis.

- Increase HPLC/UHPLC column lifetime
- Increase sensitivity of analysis
- Save time with a simple two-step method



## Increase Analyte Sensitivity


### β-Gone Protocol:

- 1) Dilute 200 μL of Urine Hydrolysate with 133 μL of 0.1 % Formic acid in Methanol
- 2) Load diluted sample onto β-Gone 96-Well Plate (Part No. 8E-S322-DGA) and apply 2-5 psi using a positive pressure manifold or a vacuum manifold
- 3) Collect eluent and inject 10 μL for analysis

## Recovery

Analyte	% Average Recovery	% CV
Benzoylceognine	109	3
Buprenorphine	93	6
Codeine	109	4
Lorazepam	79	5
Methamphetamine	106	3
Norbuprenorphine	109	5
PCP	102	3

## Norbuprenorphine: β-Gone vs. Dilute-and-Shoot



## LC-MS/MS Conditions

**Column:** Kinetex® 2.6 μm Biphenyl  
**Dimensions:** 50 x 3.0mm  
**Part No.:** 00B-4622-Y0  
**Guard:** SecurityGuard™ ULTRA Biphenyl Cartridge: AJO-9208  
**Mobile Phase:** A: 0.1% Formic acid in Water  
 B: 0.1% Formic acid in Methanol  
**Gradient:**

Time (min)	% B
0.01	10
1	10
4	100
5	100
5.01	10
6	10

**Flow Rate:** 0.7 mL/min  
**Detection:** MS/MS (SCIEX, API 4000™)


# Phospholipid Removal



Traditional protein precipitation does not remove phospholipids that are present in biological samples, such as plasma and whole blood, and phospholipids are a primary source of ion suppression in LC-MS analysis. Easily remove phospholipids and proteins in under 5 minutes using Phree Phospholipid Removal Solutions.

- Remove phospholipids
- Remove proteins
- No method development required

## Phospholipid Profiles




### HPLC Conditions

**Plasma Cleanup:** 100  $\mu$ L plasma plus 300  $\mu$ L Acetonitrile with 1% Formic acid  
**Column:** Kinetex<sup>®</sup> 2.6  $\mu$ m C18 100  $\text{\AA}$   
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4462-AN  
**Guard:** SecurityGuard<sup>™</sup> ULTRA C18 Cartridge: AJO-8782  
**Mobile Phase:** A: 0.1% Formic acid in Water  
 B: 0.1% Formic acid in Methanol  
**Gradient:**

Time (min)	% B
0	60
0.5	95
15.5	95
15.51	60
19.5	60

**Flow Rate:** 0.4 mL/min  
**Detection:** Mass Spectrometer (MS) @ 425°C; 184 amu  
**Temperature:** 22°C

~ 25x Zoom



Phospholipid profile monitored using  $m/z$  184-184

## Removal of Phospholipids from 100 $\mu$ L Plasma

	Lyso 1	Lyso 2	PC 1	PC 2	PC 4
<b>Phenomenex Phree</b>	>99.9 %	>99.9 %	>99.9 %	>99.9 %	>99.9 %
<b>Agilent Captiva ND Lipids</b>	28.9 %	36.4 %	>99.9 %	99.9 %	>99.9 %
<b>Supelco HybridSPE</b>	97.8 %	98.4 %	>96.3 %	>99.7 %	>80.9 %

**Lyso 1:** 1-Palmitoyl-2-OH-sn-glycero-phosphocholine, (16:0)  $m/z$  496-184

**Lyso 2:** 1-Oleoyl-2-OH-sn-glycero-phosphocholine, (18:1)  $m/z$  522-184

**PC 1:** 1-Palmitoyl-2-Oleoyl-sn-glycero-phosphocholine, (16:0, 18:1)  $m/z$  760-184

**PC 2:** 1-Stearoyl-2-Lindoleoyl-sn-glycero-phosphocholine, (18:0, 18:2)  $m/z$  786-184

**PC 4:** 1-Oleoyl-2-Linoleoyl-sn-glycero-phosphocholine, (18:1, 18:2)  $m/z$  784-184


Phenomenex is not affiliated with Agilent Technologies or Sigma-Aldrich Co. Comparative chromatograms may not be representative of all separations.

# Simplified Liquid Extraction

Simplified Liquid Extraction (SLE) is a faster, easier, and more reliable way to perform liquid-liquid extractions (LLE).

- Eliminate interferences from your analysis without extensive method development
- Achieve consistent, reliable results from lot-to-lot
- Available as a high quality synthetic sorbent (Novum) or as a traditional diatomaceous earth sorbent (Strata® DE)

## An Easy, Automatable Procedure



## Determine Which SLE Sorbent is Right for Your Extraction



Synthetic	<b>Sorbent Type</b>	Diatomaceous Earth
Lot-to-lot consistency and reproducibility	<b>Advantages</b>	Cost effective and large volume capabilities
Ethyl Acetate, MTBE	<b>Extraction Solvents</b>	Hexane, DCM, MTBE Ethyl Acetate
MINI 96-Well Plates, MAX 96-Well Plates	<b>Plate Formats</b>	200 µL 96-Well Plates, 400 µL 96-Well Plates
1 cc, 3 cc, 6 cc, 12 cc	<b>Tube Formats</b>	12 cc and 60 cc

SLE sorbent selections are dependent on extraction solvents and sample volume.

# Extraction of Corticosteroids

## from Plasma Using Strata DE SLE

We developed a method using Strata DE SLE for a wide range of corticosteroid compounds from plasma, which are then analyzed by LC-MS/MS. All compounds in the suite provided recovery greater than 90% with the exception of Triamcinolone. Triamcinolone is the most polar compound in the suite and is simply too hydrophilic to be extracted by DCM. Acceptable recoveries can be obtained by changing to ethyl acetate as an elution solvent. All compounds show a % CV of less than 12%.

### Pre-treatment

Dilute 100 µL of spiked plasma (125 ng/mL) with 200 µL of Water

### SLE Protocol

<b>96-Well Plate:</b>	Strata DE SLE 400 µL 96-Well Plate
<b>Part No.:</b>	8E-S325-5GB
<b>Load:</b>	Pre-treated sample onto plate (apply vacuum or positive pressure to pull/push sample into sorbent)
<b>Wait:</b>	5 minutes
<b>Elution:</b>	3x 600 µL Dichloromethane (DCM) or 3x 600 µL Ethyl Acetate
<b>Apply:</b>	Vacuum or apply positive pressure at 5-10" Hg for 10 seconds
<b>Dry Down:</b>	Sample under slow stream of Nitrogen at 30 °C
<b>Reconstitute:</b>	200 µL Acetonitrile/Water (20:80)

### LC-MS/MS Conditions

**Column:** Kinetex® 2.6 µm Biphenyl  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4622-AN  
**Guard Cartridge:** SecurityGuard™ ULTRA Biphenyl  
**Guard Part No.:** AJ0-9209  
**Mobile Phase:** A: 0.1% Formic acid in Water  
 B: 0.1% Formic acid in Acetonitrile

Gradient: Time (min)	% B
0	20
3	95
3.5	95
3.51	20
6	20

**Injection Volume:** 5 µL  
**Flow Rate:** 0.5 mL/min  
**Detection:** MS/MS (SCIEX API 4000™), ESI+

### Recovery Values and % CVs

Elution Solvent:	Dichloromethane		Ethyl Acetate	
	% Recovery	% RSD (n=4)	% Recovery	% RSD (n=10)
<b>β-Methasone</b>	92	4	98	6
<b>Cortisone</b>	96	10	96	8
<b>Coritcosterone</b>	92	3	74	10
<b>Cortisone Acetate</b>	90	12	112	12
<b>Triamcinolone</b>	13	8	92	9
<b>Prednisone</b>	94	7	93	10
<b>Testosterone</b>	95	5		

\*Testosterone was not extracted using Ethyl Acetate

**High recoveries and low RSD values using Strata DE!**

Learn more about Strata DE at [www.phenomenex.com/StrataDE](http://www.phenomenex.com/StrataDE)





# Extraction of NSAIDs

## from Plasma Using Novum SLE

In this application, we use the Novum SLE MAX 96-Well Plates to extract an eight compound NSAID suite from plasma.

### Sample Pre-treatment

Dilute 200 µL of human plasma with 200 µL of 1 % Formic acid in water. Vortex 3-5 seconds.

### SLE Protocol

**96-Well Plate:** Novum SLE MAX 96-Well Plate

**Part No.:** 8E-S138-5GA

**Load:** 400 µL of pre-treated sample onto plate and apply a short pulse of vacuum (~5" Hg) for 5-10 seconds or until sample has completely entered the sorbent.

**Elute:** 2x 900 µL of 10 % Ethyl acetate in dichloromethane (DCM) and allow the solvent to elute by gravity (~5 minutes) and collect eluant. At the completion of the second aliquot, apply vacuum at 5" Hg for 30 seconds to complete the extraction.


**Dry Down:** Evaporate the final extract to complete dryness under a slow stream of nitrogen at 40 °C for about 60 minutes.

**Reconstitute:** 200 µL of Methanol/Water (10:90) by vortexing the plate at 1400 rpm for 2 minutes.

### Recovery

Analyte	% Average Recovery	% CV
Ibuprofen	82	6.70
Diclofenac	79	3.20
Naproxen	96	2.80
Ketoprofen	96	3.10
Mefenamic Acid	77	12.6
Flurbiprofen	82	12.0
Sulindac	87	10.1
Salicylic Acid	93	7.30

High recoveries across all NSAIDs!



# Acids, Neutrals, and Bases

## from Urine Using Novum SLE

In this application, we will show how a specific pH manipulation can lead to extraction conditions of a relatively hydrophobic acid (THC-COOH) along with more polar bases (buprenorphine and norbuprenorphine) and neutrals (barbiturates). We developed a SLE application for acids, neutrals, and bases from a urine matrix containing  $\beta$ -glucuronidase followed by two LC-MS/MS methods.


### Pre-treatment

To 200  $\mu$ L of urine, add 25  $\mu$ L of  $\beta$ -Glucuronidase Enzyme, 25  $\mu$ L Ammonium Acetate Buffer (100 mM, pH 4), 180  $\mu$ L Ammonium Bicarbonate Buffer (100 mM, pH 9) and 20  $\mu$ L Internal Standard (1  $\mu$ g/mL). Final total volume is 450  $\mu$ L.

### SLE Protocol

<b>96-Well Plate:</b>	Novum Max SLE 96-Well Plate
<b>Part No.:</b>	8E-S138-5GA
<b>Load:</b>	Pre-treated sample and pulse vacuum at 5" Hg for 2-3 seconds, or until the sample completely enters the sorbent bed. Wait for 6 minutes.
<b>Elute:</b>	2x 900 $\mu$ L Ethyl acetate and elute by gravity. Apply 5" vacuum at end of elution to collect residual solvent from tips in collection plate.
<b>Dry Down:</b>	Under a gentle stream of nitrogen at 30°C.
<b>Reconstitute</b>	For ESI+ samples (THC-COOH, Buprenorphine and Norbuprenorphine) reconstitute in Methanol/0.1% Formic acid in water (1:4). For ESI- samples (Barbiturates) reconstitute in Methanol/1% NH <sub>4</sub> OH in water (1:4).


### ESI+ Chromatogram (Buprenorphine/Norbuprenorphine/THC)



### Positive Mode LC-MS/MS Conditions

<b>Column:</b>	Kinetex® 2.6 $\mu$ m Biphenyl										
<b>Dimensions:</b>	50 x 2.1 mm										
<b>Part No.:</b>	00B-4622-AN										
<b>Guard:</b>	SecurityGuard™ ULTRA Biphenyl Cartridge: AJO-9209										
<b>Mobile Phase:</b>	A: 0.1% Formic acid in Water B: 0.1% Formic acid in Acetonitrile										
<b>Gradient:</b>	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>B (%)</th> </tr> </thead> <tbody> <tr><td>0</td><td>5</td></tr> <tr><td>5</td><td>100</td></tr> <tr><td>5.1</td><td>5</td></tr> <tr><td>7</td><td>5</td></tr> </tbody> </table>	Time (min)	B (%)	0	5	5	100	5.1	5	7	5
Time (min)	B (%)										
0	5										
5	100										
5.1	5										
7	5										
<b>Injection:</b>	4 $\mu$ L										
<b>Flow Rate:</b>	0.5 mL/min										
<b>Temperature:</b>	Ambient										
<b>Detection:</b>	MS/MS (SCIEX API 4000™)										
<b>Sample:</b>	1. Norbuprenorphine 2. Buprenorphine 3. THC-COOH										

### ESI- Chromatogram (Barbiturates Mix)




### Negative Mode LC-MS/MS Conditions

<b>Column:</b>	Kinetex 2.6 $\mu$ m EVO C18 100Å																
<b>Dimensions:</b>	50 x 2.1 mm																
<b>Part No.:</b>	00B-4725-AN																
<b>Guard:</b>	SecurityGuard ULTRA EVO-C18 Cartridge: AJO-9298																
<b>Mobile Phase:</b>	A: 10 mM Ammonium bicarbonate, pH 9 B: Acetonitrile																
<b>Gradient:</b>	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>B (%)</th> </tr> </thead> <tbody> <tr><td>0</td><td>5</td></tr> <tr><td>2</td><td>15</td></tr> <tr><td>5</td><td>20</td></tr> <tr><td>5.01</td><td>60</td></tr> <tr><td>6</td><td>60</td></tr> <tr><td>6.1</td><td>5</td></tr> <tr><td>7.5</td><td>5</td></tr> </tbody> </table>	Time (min)	B (%)	0	5	2	15	5	20	5.01	60	6	60	6.1	5	7.5	5
Time (min)	B (%)																
0	5																
2	15																
5	20																
5.01	60																
6	60																
6.1	5																
7.5	5																
<b>Injection:</b>	4 $\mu$ L																
<b>Flow Rate:</b>	0.5 mL/min																
<b>Temperature:</b>	Ambient																
<b>Detection:</b>	MS/MS (SCIEX API 4000)																
<b>Sample:</b>	1. Phenobarbital 2. Butalbital 3. Pentobarbital 4. Amobarbital 5. Amobarbital-DS 6. Secobarbital																

Solid Phase Extraction (SPE) is a very targeted form of sample preparation that allows you to isolate your analyte of interest, while removing any interfering compounds that may be in your sample.

- Targeted analyte extraction for cleaner analysis
- Concentration of samples for better chromatographic results
- Solvent switching for GC or LC compatibility

### Solid Phase Extraction General Protocol



### A Choice for Every Analyte

OPTION  
**01**

#### Acidic Compounds

Strong Acids  
( $pK_a < 2$ )

Weak Acids  
( $pK_a 2-4$ )



OPTION  
**02**

#### Neutral Compounds

Neutral Compounds



OPTION  
**03**

#### Basic Compounds

Weak Bases  
( $pK_b 8-10$ )

Strong Bases  
( $pK_b > 10$ )



Learn more, visit  
[www.phenomenex.com/SPE](http://www.phenomenex.com/SPE)

# Urinary Catecholamines

## Using Strata-X Microelution SPE

Metanephrine and normetanephrine are both metabolites of epinephrine and norepinephrine. In this application, Strata-X-CW Microelution SPE 96-Well Plates were used in conjunction with a Kinetex® Biphenyl HPLC column in order to resolve an interference that coelutes with 3-Methoxytyramine on a standard C18 HPLC column, while reaching low limits of quantification for specific urinary catecholamines, metanephrine, and normetanephrine.

### Urine Pre-treatment

500 µL of urine was diluted with 500 µL of 50 mM Ammonium acetate buffer, (pH 7). Urine was pre-spiked from 10 ng/mL to 63 pg/mL with metanephrine, normetanephrine, and 3-methoxytyramine (standards provided by Cerilliant®).


### SPE Method

<b>Microelution 96-Well Plate:</b>	Strata-X-CW Microelution 96-Well Plate, 2 mg/well
<b>Part No.:</b>	8M-S035-4GA
<b>Condition:</b>	200 µL Methanol
<b>Equilibrate:</b>	200 µL 50 mM Ammonium acetate buffer, pH 7
<b>Load:</b>	1 mL of pre-treated sample
<b>Wash 1:</b>	200 µL of 50 mM Ammonium acetate buffer, pH 7
<b>Wash 2:</b>	200 µL Acetonitrile/IPA (1:1)
<b>Elute:</b>	2x 25 µL of Water/Acetonitrile/Formic acid (85:10:5)
<b>Injection:</b>	Dilute eluent with 100 µL of 0.1 % Formic acid in water (Metanephrine-D3 internal standard was included at 1ng/mL)


### Recovery Values from 10 ng/mL to 63 pg/mL

Analyte Concentration (ng/mL)	Average % Recovery	% CV (n=6)
<b>Metanephrine</b>		
10	102	5
1	102	3
0.5	99	2
0.25	99	3
0.125	97	3
0.063	94	6
<b>Normetanephrine</b>		
10	100	10
1	87	12
0.5	110	10
0.25	89	9
0.125	110	13
0.063	108	15
<b>Amobarbital</b>		
10	91	3
1	89	6
0.5	95	4
0.25	86	5
0.125	87	6
0.063	92	7

### Chromatogram of unresolved interference for 3-methoxytyramine using a HPLC C18 column



### Representative Chromatogram of Urinary Catecholamines



### LC-MS/MS Conditions

<b>Column:</b>	Kinetex 5 µm Biphenyl										
<b>Dimensions:</b>	50 x 4.6 mm										
<b>Part No.:</b>	00B-4627-E0										
<b>Guard:</b>	SecurityGuard™ ULTRA Biphenyl Cartridges: AJ0-9207										
<b>Mobile Phase:</b>	A: 0.1 % Formic acid in Water B: 0.1 % Formic acid in Methanol										
<b>Gradient:</b>	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>B (%)</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>5</td> </tr> <tr> <td>3</td> <td>90</td> </tr> <tr> <td>3.1</td> <td>5</td> </tr> <tr> <td>5</td> <td>5</td> </tr> </tbody> </table>	Time (min)	B (%)	0	5	3	90	3.1	5	5	5
Time (min)	B (%)										
0	5										
3	90										
3.1	5										
5	5										
<b>Injection:</b>	30 µL										
<b>Flow Rate:</b>	0.7 mL/min										
<b>Temperature:</b>	Ambient										
<b>Detection:</b>	MS/MS (SCIEX API 4000™)										

# Underivatized Methylmalonic Acid (MMA)

## from Plasma Using Strata-X SPE

Methylmalonic acid (MMA) is a small dicarboxylic acid. This hydrophilic molecule can present chromatographic challenges both in achieving adequate retention under reversed phase conditions as well as resolution from the isomeric/isobaric species such as succinic acid, especially at low analyte concentrations. To combat these challenges, many published LC-MS/MS methods require a sample derivatization step, however, this step can add time to the overall analysis. Presented is a fast, reproducible LC-MS/MS method to analyze underivatized MMA by utilizing a unique Luna® Omega 1.6 µm PS C18 UHPLC column. The method runtime is 5 minutes including column re-equilibration. For the sample preparation procedure we used Strata-X-AW Solid Phase Extraction (SPE) to produce a clean sample from plasma. Analyte detection was performed using negative mode electrospray ionization of a triple quadrupole MS.

### Sample Pre-treatment

Combine 0.5 mL of 1 % aqueous acetic acid and 50 µL of internal standard with 100 µL blank, standard, or sample.


### SPE Protocol

<b>Cartridge:</b>	Strata-X-AW 30 mg/1 mL
<b>Part No.:</b>	8B-S038-TAK
<b>Condition:</b>	1 mL of Methanol
<b>Equilibrate:</b>	1 mL of 1 % Acetic acid in Water
<b>Load:</b>	Pre-treated sample (see above)
<b>Wash:</b>	0.5 mL of Methanol/Water (50:50)
<b>Dry:</b>	5 to 10 minutes at max vacuum (or apply positive pressure using Preston™ 100 Positive Pressure Manifold)
<b>Elute:</b>	2x 0.6 mL 2 % Ammonium hydroxide in Methanol
<b>Dry Down:</b>	Evaporate solvent to dryness @ 45-50 °C under a gentle stream of nitrogen
<b>Reconstitute:</b>	200 µL of mobile phase A (0.1 % Formic acid in Water)

### Analyte Recovery

Sample Name	Spike (nmol/L)	Average Concentration N=3 (nmol/L)	% CV	% Recovery
Prespiked 250 nmol/L in plasma	250	696	9.44	114
Prespiked 750 nmol/L in plasma	750	1157	2.00	102
Extracted unspiked plasma	0	385	3.01	N/A

### MMA Representative Chromatogram



Peaks in order of elution: plasma interference (0.81 min), succinic acid (1.00 min), methyl-D3-malonic acid (1.20 min), and methylmalonic acid (1.23 min).

### LC Conditions

<b>Analytical Column:</b>	Luna Omega 1.6 µm PS C18												
<b>Dimensions:</b>	50 x 2.1 mm												
<b>Part No.:</b>	00B-4752-AN												
<b>Guard:</b>	SecurityGuard™ ULTRA PS C18 Cartridges: AJO-9508												
<b>Mobile Phase:</b>	A: 0.1% Formic acid in Water B: 0.1% Formic acid in Acetonitrile												
<b>Gradient:</b>	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>B (%)</th> </tr> </thead> <tbody> <tr><td>0.01</td><td>2</td></tr> <tr><td>2</td><td>90</td></tr> <tr><td>3</td><td>90</td></tr> <tr><td>3.01</td><td>2</td></tr> <tr><td>5</td><td>2</td></tr> </tbody> </table>	Time (min)	B (%)	0.01	2	2	90	3	90	3.01	2	5	2
Time (min)	B (%)												
0.01	2												
2	90												
3	90												
3.01	2												
5	2												
<b>Injection:</b>	5 µL												
<b>Flow Rate:</b>	0.4 mL/min												
<b>Temperature:</b>	40 °C												

### MS/MS Conditions

<b>Detector:</b>	SCIEX 4000 QTRAP®
<b>Mode:</b>	Negative Ionization Mode
<b>Scan Type:</b>	MRM
<b>Curtain Gas (CUR):</b>	10.0 psi
<b>Collision Gas (CAD):</b>	Medium
<b>IonSpray Voltage (IS):</b>	-4500 V
<b>Temperature (TEM):</b>	600 °C
<b>Ion Source Gas 1 (Gas1):</b>	50 psi
<b>Ion Source Gas 2 (Gas2):</b>	50 psi
<b>Interface Heater (Ihe):</b>	On

# Comprehensive Drug Research Panel

## from Oral Fluid Collection Devices Using Strata-X SPE

Drug testing in oral fluid has steadily gained popularity due to the easy and non-intrusive sample collection procedure. The oral fluid collection device provides a buffer solution that contains a number of antibacterial agents and surfactants that act to prevent bacterial growth and increase the analytes stability during the sample transit to testing laboratories. The buffer solution poses many chromatography challenges, such as ion suppression. Here, we present a fast sample preparation procedure to reduce the effects of the device's buffer solution while maintaining good recovery of analytes using a Strata-X-C SPE method for different classes of analytes from a comprehensive drug research panel. For a more selective SPE solution, our dual cartridge SPE method is recommended (page 16).

### Sample Pre-treatment

Transfer 1 mL of oral fluid collected on an applicator tip in its preservative buffer. Leave it for 2 hours followed by centrifugation for 15 minutes at 600 g. Remove 0.5 mL of supernatant and combine with 1 mL of 1% Formic acid. Vortex for 5-10 seconds.

### SPE Protocol


<b>96-Well Plate:</b>	Strata-X-CW 30 mg/well
<b>Part No.:</b>	8E-S035-TGB
<b>Condition:</b>	1 mL Methanol
<b>Equilibrate:</b>	1 mL DI Water
<b>Load:</b>	Pre-treated Sample
<b>Wash 1:</b>	1 mL 1% Formic acid in DI Water
<b>Wash 2:</b>	1 mL DI Water
<b>Dry:</b>	5-6 minutes at maximum vacuum (20" Hg or higher)
<b>Elution:</b>	2x 500 µL (2 aliquots of 500 µL) Methylene chloride/Isopropanol/30% Ammonium hydroxide in Water (80:18:2)
<b>Dry Down:</b>	Evaporate to dryness under nitrogen at 45-50°C
<b>Reconstitute:</b>	200 µL of 0.1% Formic acid in Water/0.1% Formic acid in Methanol (90:10)

**Analyzing additional analytes?**  
Use our Strata-X double cartridge solution for increased clean-up and increased sensitivity across all classes of drug compounds.  
pg. 16

Analyte Concentration (ng/mL)	Average % Recovery
6-MAM	79
α-Hydroxyalprazolam	88
Alprazolam	79
Amphetamine	94
Benzoyllecgonine	89
Carisoprodol	95
Citalopram	81
Cocaine	84
Codeine	89
Diazepam	74
EDDP	72
Fentanyl	74
Hydrocodone	102
Hydromorphone	98
Meperidine	84
Mephedrone	84
Meprobamate	92
Methamphetamines	81
Naloxone	83
Norbuprenorphine	87
Nordiazepam	75
Norfentanyl	85
Norhydrocodone	89
Noroxycodone	85
Oxycodone	84
PCP	84
Phenobarbital	83
Secobarbital	82
Tapentadol	83
Temazepam	86
Tramadol	87

# Expanded Comprehensive Drug Research Panel

## from Oral Fluid Collection Devices Using Double Cartridge Strata-X SPE



### SPE Method

Step	Basic analyte extraction	Acidic analyte extraction
<b>Cartridge:</b>	Strata-X-C, 30 mg/3 mL	Strata-X-A, 30 mg/3 mL
<b>Part No.:</b>	8B-S029-TBJ	8B-S123-TBJ
<b>Condition:</b>	1 mL 100% Methanol	1 mL 100% Methanol
<b>Equilibrate:</b>	1 mL DI Water	1 mL DI Water
<b>Load:</b>	Combine 0.5 mL of pre-treated sample with 1 mL 1% Formic acid, mix/vortex 5-10 sec and load on Strata-X-C	Combine 0.5 mL of pre-treated sample with 1 mL 1% Ammonium hydroxide, mix/vortex 5-10 sec and load on Strata-X-A
<b>Weak Wash:</b>	1 mL DI Water	1 mL DI Water
<b>Strong Wash:</b>	1 mL Acetone/Water (50:50)	1 mL Acetone/Water (50:50)
<b>Dry down:</b>	3-4 minutes at maximum vacuum (15" Hg or higher)	3-4 minutes at maximum vacuum (15" Hg or higher)
<b>Elute:</b>	2x 500 µL Methanol/Acetonitrile/Ammonium hydroxide (5:5:2)	2x 500 µL Methanol/Acetonitrile/Formic acid (50:50:5)
<b>Dry down:</b>	Evaporate to dryness under a gentle steam of Nitrogen at 45-50 °C	Evaporate to dryness under a gentle steam of Nitrogen at 45-50 °C
<b>Reconstitute:</b>	With 125 µL initial Mobile Phase	With 125 µL initial Mobile Phase
<b>Combine into a single sample vial</b>		

Note: If not testing for THC-COOH, lorazepam and other select barbiturates, use Strata-X-C only.




# Expanded Comprehensive Drug Research Panel (cont'd)

## from Oral Fluid Collection Devices Using Double Cartridge Strata-X SPE

### Comprehensive Drug Research Panel Chromatogram

### Positive ESI Panel




**Column:** Kinetex 2.6  $\mu$ m Biphenyl  
**Dimensions:** 50 x 3.0 mm  
**Part No.:** 00B-4622-Y0  
**Guard:** SecurityGuard™ ULTRA Biphenyl Cartridge: AJO-9208  
**Mobile Phase:** A: 0.1% Formic acid in Water  
                   B: 0.1% Formic acid in Methanol  
**Gradient:**

Time (min)	% B
0	10
4	95
5.5	95
5.51	10
7.5	10

  
**Injection:** 10  $\mu$ L  
**Flow Rate:** 0.5 mL/min  
**Temperature:** Ambient  
**Detector:** MS/MS (SCIEX API 5000™) ESI+

### Representative LC-MS Chromatogram of Buffer Solution

See the reduction in buffer interferences!



# Ordering Information

Ordering Information



## β-Gone β-Glucuronidase Removal

Part No.	Description	Unit
8B-S139-TAK	1 mL Tubes, Recombinant Enzyme	100/Box
8B-S322-DAK	1 mL Tubes, Non-Recombinant Enzyme	100/Box
8E-S139-TGA	96-Well Plate, Recombinant Enzyme	1/Box
8E-S322-DGA	96-Well Plate, Non-Recombinant Enzyme	1/Box
8N-S323-TUK	2 mL Centrifuge Tubes, Recombinant and Non-Recombinant Enzyme	100/Box



## Strata DE SLE

Strata DE Diatomaceous Earth (SLE)		
Part No.	Description	Unit
8E-S325-FGB	Strata DE SLE 200 µL 96-Well Plate	2/pk
8E-S325-5GB	Strata DE SLE 400 µL 96-Well Plate	2/pk
8B-S325-KDG	Strata DE SLE 12 cc Tubes	20/pk
8B-S325-VFF	Strata DE SLE 60 cc Tubes	16/pk



## Phree Phospholipid Removal

Part No.	Description	Unit
8B-S133-TAK	Phree Phospholipid Removal 1 mL Tube	100/box
8E-S133-TGB	Phree Phospholipid Removal 96-Well Plates	2/box



If Phenomenex products in this brochure do not provide at least equivalent separation as compared to other products of the same phase and dimensions, return the product with your comparative data within 45 days for a FULL REFUND.



## Novum SLE

Part No.	Description	Unit
8E-S138-FGA	Novum SLE MINI 96-Well Plate	1/Box
8E-S138-5GA	Novum SLE MAX 96-Well Plate	1/Box
8B-S138-FAK	Novum SLE 1 cc tubes	100/Box
8B-S138-5BJ	Novum SLE 3 cc tubes	50/Box
8B-S138-JCH	Novum SLE 6 cc tubes	30/Box
8B-S138-KDG	Novum SLE 12 cc tubes	20/Box

## Strata-X SPE Tubes



Tubes	1 mL (100/box)		3 mL (50/box)			6 mL (30/box)		
	30 mg	60 mg	60 mg	200 mg	500 mg	100 mg	200 mg	500 mg
Strata-X	8B-S100-TAK	8B-S100-UAK	8B-S100-UBJ	8B-S100-FBJ	8B-S100-HBJ	8B-S100-ECH	8B-S100-FCH	8B-S100-HCH
Strata-X-C	8B-S029-TAK	—	8B-S029-UBJ	8B-S029-FBJ	8B-S029-HBJ	8B-S029-ECH	8B-S029-FCH	8B-S029-HCH
Strata-X-CW	8B-S035-TAK	—	8B-S035-UBJ	8B-S035-FBJ	8B-S035-HBJ	8B-S035-ECH	8B-S035-FCH	8B-S035-HCH
Strata-X-A	8B-S123-TAK	—	8B-S123-UBJ	8B-S123-FBJ	8B-S123-HBJ	8B-S123-ECH	8B-S123-FCH	8B-S123-HCH
Strata-X-AW	8B-S038-TAK	—	8B-S038-UBJ	8B-S038-FBJ	8B-S038-HBJ	8B-S038-ECH	8B-S038-FCH	8B-S038-HCH
Strata-XL	8B-S043-TAK	—	8B-S043-UBJ	8B-S043-FBJ	8B-S043-HBJ	8B-S043-ECH	8B-S043-FCH	8B-S043-HCH
Strata-XL-C	8B-S044-TAK	—	8B-S044-UBJ	8B-S044-FBJ	8B-S044-HBJ	8B-S044-ECH	8B-S044-FCH	8B-S044-HCH
Strata-XL-CW	8B-S052-TAK	—	8B-S052-UBJ	8B-S052-FBJ	8B-S052-HBJ	8B-S052-ECH	8B-S052-FCH	8B-S052-HCH
Strata-XL-A	8B-S053-TAK	—	8B-S053-UBJ	8B-S053-FBJ	8B-S053-HBJ	8B-S053-ECH	8B-S053-FCH	8B-S053-HCH
Strata-XL-AW	8B-S051-TAK	—	8B-S051-UBJ	8B-S051-FBJ	8B-S051-HBJ	8B-S051-ECH	8B-S051-FCH	8B-S051-HCH

## Strata-X SPE 96-Well Plates

96-Well Plates (2/Box)			
Phase	10 mg	30 mg	60 mg
Strata-X-AW	8E-S038-AGB	8E-S038-TGB	8E-S038-UGB
Strata-X-A	8E-S123-AGB	8E-S123-TGB	8E-S123-UGB
Strata-X	8E-S100-AGB	8E-S100-TGB	8E-S100-UGB
Strata-X-C	8E-S029-AGB	8E-S029-TGB	8E-S029-UGB
Strata-X-CW	8E-S035-AGB	8E-S035-TGB	8E-S035-UGB
Strata-XL-AW	—	8E-S051-TGB	—
Strata-XL-A	—	8E-S053-TGB	—
Strata-XL	—	8E-S043-TGB	—
Strata-XL-C	—	8E-S044-TGB	—
Strata-XL-CW	—	8E-S052-TGB	—

## Strata-X Microelution SPE Plates

96-Well Plates (ea)	
Phase	2 mg / well
Strata-AW	8M-S038-4GA
Strata-A	8M-S123-4GA
Strata-X	8M-S100-4GA
Strata-X-C	8M-S029-4GA
Strata-X-CW	8M-S035-4GA

# Ordering Information

## Presston 100 Positive Pressure Manifold

Part No.	Description	Unit
AHO-9334	Presston 100 Positive Pressure Manifold, 96-Well Plate	1/Box
AHO-9342	Presston 100 Positive Pressure Manifold, 1 mL Tube Complete Assembly	1/Box
AHO-9347	Presston 100 Positive Pressure Manifold, 3 mL Tube Complete Assembly	1/Box
AHO-9343	Presston 100 Positive Pressure Manifold, 6 mL Tube Complete Assembly	1/Box



Phenomenex warrants that for a period of 12 months following delivery, the Presston 100 Positive Pressure Manifold you have purchased will perform in accordance with the published specifications and will be free from defects in materials or workmanship. In the event that the Presston 100 Positive Pressure Manifold does not meet this warranty, Phenomenex will repair or replace defective parts. Please visit [www.phenomenex.com/Presston](http://www.phenomenex.com/Presston) for complete warranty information.

## Kinetex Analytical LC Columns

5 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges <sup>‡</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJO-9298

for 2.1 mm ID



5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges <sup>‡</sup>
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJO-9296
F5	00B-4724-E0	00D-4724-E0	00F-4724-E0	00G-4724-E0	AJO-9320
Biphenyl	00B-4627-E0	00D-4627-E0	00F-4627-E0	00G-4627-E0	AJO-9207

for 4.6 mm ID

2.6 µm Minibore Columns (mm)						SecurityGuard ULTRA Cartridges <sup>‡</sup>
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJO-9298
Polar C18	00A-4759-AN	00B-4759-AN	—	00D-4759-AN	00F-4759-AN	AJO-9532
F5	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJO-9322
Biphenyl	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJO-9209
XB-C18	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
C18	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782
C8	00A-4497-AN	00B-4497-AN	00C-4497-AN	00D-4497-AN	00F-4497-AN	AJO-8784
HILIC	00A-4461-AN	00B-4461-AN	00C-4461-AN	00D-4461-AN	00F-4461-AN	AJO-8786
Phenyl-Hexyl	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJO-8788

for 2.1 mm ID

2.6 µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges <sup>‡</sup>
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	—	00B-4725-Y0	—	00D-4725-Y0	00F-4725-Y0	AJO-9297
Polar C18	—	00B-4759-Y0	—	00D-4759-Y0	00F-4759-Y0	AJO-9531
F5	—	00B-4723-Y0	—	00D-4723-Y0	00F-4723-Y0	AJO-9321
Biphenyl	—	00B-4622-Y0	—	00D-4622-Y0	00F-4622-Y0	AJO-9208
XB-C18	00A-4496-Y0	00B-4496-Y0	00C-4496-Y0	00D-4496-Y0	00F-4496-Y0	AJO-8775
C18	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJO-8775
C8	00A-4497-Y0	00B-4497-Y0	00C-4497-Y0	00D-4497-Y0	00F-4497-Y0	AJO-8777
HILIC	00A-4461-Y0	—	—	—	00F-4461-Y0	AJO-8779
Phenyl-Hexyl	—	00B-4495-Y0	—	00D-4495-Y0	00F-4495-Y0	AJO-8781

for 3.0 mm ID

## Luna Omega UHPLC Columns

1.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges <sup>‡</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
Polar C18	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJO-9505
PS C18	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJO-9508
C18	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJO-9502

for 2.1 mm ID



# CLINICAL SAMPLE PREPARATION

## SIMPLE | FAST | CLEAN

### Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auiinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

### China

t: +86 400-606-8099  
f: +86 (0)22 2532-1033  
phen@agela.com

### Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

### Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

### All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com



[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademark

Kinetex, Luna, and Strata are registered trademarks and MidBore, Novum, Phree,  $\beta$ -Gone, Presston, SecurityGuard, and Solvent Shielding Technology are trademarks of Phenomenex. Agilent is a registered trademark and Captiva is a trademark of Agilent Technologies, Inc. Supelco and HybridSPE are registered trademarks of Sigma-Aldrich Co., LLE. Cerilliant is a registered trademark of Cerilliant Corporation. API 4000 and API 5000 are trademarks and QTRAP is a registered trademark of AB SCIEX Pte, Ltd. AB SCIEX is being used under license.

### Disclaimer

Phenomenex is not affiliated with Agilent Technologies, Inc., Sigma-Aldrich Co., LLE, or Cerilliant Corporation.

Comparative separations may not be representative of all applications.

Novum is patent pending.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

Kinetex EVO is patented by Phenomenex. U.S. patent No. 7,563,367 and 8,658,038 and foreign counterparts.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

**CAUTION:** this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2017 Phenomenex, Inc. All rights reserved.