

Chiral Purity Testing of the Drug Substance Rivaroxaban Using Lux® 5 µm Cellulose-1 Chiral Column as per USP Monograph

Lauren Nakasone, Dr. Bryan Tackett, and Anania Assefa Amenu
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501-1430, USA

Overview


Rivaroxaban is a small molecule chiral drug used as an anticoagulant to treat blood clots. It is commercialized as the pure S-enantiomer. Separation of the enantiomers of chiral drugs has become of vital importance in analytical chemistry in recent years, because of differences in the biological activity and pharmacokinetic properties of drug enantiomers.

In this application note, we demonstrate the successful enantiomeric separation between Rivaroxaban and its enantiomer impurity, R-enantiomer, using the Lux Cellulose-1 chiral HPLC column according to the current USP monograph for Rivaroxaban.

System suitability per USP Monograph for enantiomeric purity of Rivaroxaban is for resolution of not less than (NLT) 1.5. The separation of the enantiomers was achieved and peaks corresponding to the R-enantiomer and Rivaroxaban had a resolution of 2.37, exceeding the system suitability requirement (**Figure 2**). The higher resolution indicates the powerful chiral recognition ability of Lux 5 µm Cellulose-1 chiral column.

The monograph for enantiomeric purity also requires a symmetry factor of not more than (NMT) 1.5 for the peak associated with Rivaroxaban. Lux 5 µm Cellulose-1 easily provided a symmetry factor of 1.24, meeting system suitability.

All reference solutions were prepared as indicated in the USP Monograph for Rivaroxaban. USP Rivaroxaban RS (Catalog No. 1604530), and USP Rivaroxaban R-Enantiomer RS (Catalog no. 1604541) were purchased from USP.



LC-UV Conditions

Column: Lux 5 µm Cellulose-1 (150 x 4.6 mm)

Part No.: [00F-4459-E0](#)

Mobile Phase: Ethanol and Heptane (30:70)

Pressure (bar): 69

Flow Rate: 1 mL/min

Injection: 15 µL

Temperature: 50 °C

Detector: UV @ 250 nm

System: Agilent® 1260 Binary UHPLC

Preparation of Solution

System Suitability Solution : 0.4 mg/mL of USP Rivaroxaban RS and 0.004 mg/mL of USP Rivaroxaban R-Enantiomer RS prepared as follows. To a suitable volumetric flask add a suitable amount of USP Rivaroxaban RS and USP Rivaroxaban R-Enantiomer RS. Add acetonitrile to about 50% of the flask volume, and dilute with ethanol to volume.

Figure 1. Rivaroxaban Chemical Structure






Figure 2. Resolution of Ricaroxaban R-Enantiomer (1) and Rivaroxaban (2)



No.	Peak	Time	Area	Height	Width	Symmetry Factor	Resolution
1	Rivaroxaban R-enantiomer	14.1	50.6	1.3	0.65		2.37
2	Rivaroxaban	17.0	6427.3	128.5	0.68	1.24	

Figure 3. Zoomed-in View of Chromatogram to Highlight Resolution



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Lux is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

