

APPLICATIONS

A Simple Quantitative Method for Monoclonal Antibody Coformulations

Ivan Lebedev and Brian Rivera

Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Monoclonal antibodies (mAbs) are well-established therapeutics for biopharmaceutical work. However, recent trends have investigated the administration of multiple mAbs for immunotherapy and passive immunization strategies. Although characterization may present unique challenges, quantitation of each individual mAb tend to rely on standard analytical workflows, but methods for drug substance of the mAb coformulation often consist of a more complex analysis.

One such consideration is a simple quantitative method for the mAb coformulation. Although this might seem like a straightforward analytical method, the similarities in physicochemical properties between any two mAbs may complicate the analysis. Many mAbs chosen for coformulation are standard format, i.e. IgG1 humanized or chimeric, and expressed from Chinese Hamster Ovary (CHO) cell lines. Consequently, coformulated mAbs will be very have similar molecular weight, glycosylation, and isoelectric point, limiting chromatographic method options.

As such, the analytical method must have the selectivity to separate these highly similar biomolecules. Intact reversed phase then would be an appropriate separation modality, as hydrophobic differences should be sufficient between most mAb samples. Here, we present a proof-of-concept, simple and quantitative method for mAb coformulations by intact protein reversed phase LC using the bioZen™ 2.6 µm WidePore C4 column.

Using a superficially porous, widepore particle morphology not only gives high efficiency separations but enables higher flow rates, e.g. 0.8 mL/min for a 2.1 mm ID. This ensures as shallow a gradient slope possible, which is critical for selectivity of intact proteins. The initial scouting method used was 30-40% B in 5 minutes; with a 150 mm length column this results in a gradient slope of approximately 0.6% B per column volume. Resulting chromatogram overlays for two mAbs- NIST and Trastuzumab- are shown in **Figure 1**. The retention time difference visually show a favorable result.

As seen in **Figure 2**, injecting a 1 mg/mL solution containing both mAbs shows baseline separation, with USP resolution being 3.2. Further optimization of the method by modifying the gradient from 40-50% B (**Figure 3**), using the same gradient slope, yields a marginal improvement of resolution of 4.

Finally, temperature was investigated to improve peak shape, with the intent of maximizing resolution. As observed in **Figure 4**, increasing temperature to 80°C yields an increase of resolution to 4.4.

In summary, a simple reversed phase method can be used for quantitation of coformulated mAbs. A scouting gradient can be used to confirm differences in retention time between two samples, with further optimization including gradient program optimization and modulation of temperature. Limitations of this method include discrepancies in concentration between two mAb samples. Further optimization then would have to be made to ensure applicable linearity for the method.

LC Conditions

- Column:** [bioZen 2.6 µm WidePore C4, 400 Å](#)
- Dimensions:** 150 x 2.1 mm
- Part No.:** [00F-4786-AN](#)
- Mobile Phase A:** 0.1 % TFA in Water
- Mobile Phase B:** 0.1 % TFA in Acetonitrile
- Gradient Program:** 30-40% B in 5 minutes (Fig. 1-2)
40-50% B in 5 minutes (Fig. 3-4)
- Flow Rate:** 0.8 mL/min
- Temperature:** 60°C (Fig. 1-3)
80°C (Fig. 4)
- Detection:** UV-Vis @ 280 nm
- Injection:** 1 µL
- Sample:** Trastuzumab (1 mg/mL)
NIST mAb (1 mg/mL)

Figure 1. Overlay of Trastuzumab and NIST mAb




Figure 2. Trastuzumab and NIST, 30-40% B Gradient, 60 °C




Figure 3. Trastuzumab and NIST, 40-50% B Gradient, 60 °C





Figure 4. Trastuzumab and NIST, 40-50% B Gradient, 80 °C



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

<p>Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com</p>	<p>India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com</p>	<p>Portugal t: +351 221 450 488 ptinfo@phenomenex.com</p>
<p>Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com</p>	<p>Ireland t: +353 (0)1 247 5405 eirinfo@phenomenex.com</p>	<p>Singapore t: +65 800-852-3944 sginfo@phenomenex.com</p>
<p>Belgium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com</p>	<p>Italy t: +39 051 6327511 italiainfo@phenomenex.com</p>	<p>Spain t: +34 91-413-8613 espinfo@phenomenex.com</p>
<p>Canada t: +1 (800) 543-3681 info@phenomenex.com</p>	<p>Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com</p>
<p>China t: +86 400-606-8099 cninfo@phenomenex.com</p>	<p>Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com</p>	<p>Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com</p>
<p>Denmark t: +45 4824 8048 nordicinfo@phenomenex.com</p>	<p>The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com</p>	<p>Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com</p>
<p>Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com</p>	<p>New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com</p>	<p>United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com</p>
<p>France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com</p>	<p>Norway t: +47 810 02 005 nordicinfo@phenomenex.com</p>	<p>USA t: +1 (310) 212-0555 info@phenomenex.com</p>
<p>Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com</p>	<p>Poland t: +48 22 104 21 72 pl-info@phenomenex.com</p>	<p>🌐 All other countries/regions Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com</p>

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen is a trademark of Phenomenex.

Comparisons may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.