

Meeting System Suitability Requirements for Atorvastatin Calcium per USP Monograph

Zandra Baja and Bryan Tackett, PhD
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview


Atorvastatin Calcium is a competitive inhibitor of HMG-CoA reductase, the rate-controlling enzyme of the metabolic pathway that produces cholesterol and other isoprenoids. It decreases the amount of LDL-cholesterol in the blood and reduces blood levels of triglycerides and slightly increases levels of HDL-cholesterol. The development of a quick and efficient analysis of Atorvastatin Calcium is of interest for generic drug manufacturers.

In this application note, we report meeting system suitability requirements for the USP Monograph for Atorvastatin Calcium Tablets Assay using a Luna 5 μm C18(2) column and a Luna Omega 5 μm C18 column, which provide similar selectivity as the original Ultramex 5 μm C18 column.

System suitability per USP Monograph for the Atorvastatin Calcium Tablets Assay is a percent relative standard deviation (%RSD) no more than (NMT) 1.0%, a symmetry factor NMT 1.5, and resolution not less than (NLT) 5.0 between Atorvastatin and Atorvastatin Related Compound H. The Luna 5 μm C18(2) column, Luna Omega 5 μm C18 column, and Ultramex 5 μm C18 column met all system suitability requirements.

All solutions were prepared as indicated in the USP Monograph for Atorvastatin Calcium Tablets. Atorvastatin Calcium RS (Catalog No. 1044516) and Atorvastatin Related Compound H RS (Catalog No. 1044582) were purchased from USP.

Figure 1. Atorvastatin Calcium



LC-UV Conditions

Column: Luna™ 5 μm C18(2)
Luna Omega 5 μm C18
Ultramex™ 5 μm C18

Dimension: 250 x 4.6 mm

Part No.: [00G-4252-E0](#)
[00G-4785-E0](#)
[00F-0048-E0](#)

Mobile Phase: Acetonitrile / stabilizer-free
Tetrahydrofuran / Buffer (27:20:53, v/v/v)

Buffer: 0.05 M Ammonium Citrate buffer,
pH 4.0

Flow Rate: 1.5 mL/min (Isocratic)

Injection Volume: 20 μL

Temperature: 30 °C

Detection: UV @ 244 nm


System: Agilent® 1260 Quaternary HPLC

Table 1. Preparation of Solutions


Solution	Composition
Solution A	Dissolve 9.62 g of Anhydrous Citric Acid in 900 mL of Water, adjust with Ammonium Hydroxide to a pH of 7.4, and dilute with Water to 1000 mL.
Diluent	Acetonitrile and Solution A (1:1)
Standard Solution	0.1 mg/mL of Atorvastatin Calcium RS in Diluent. Shake mechanically for 15 min or until dissolved.
System Suitability Solution	0.1 mg/mL of Atorvastatin Calcium RS and 0.01 mg/mL of Atorvastatin Related Compound H RS in Diluent. Shake mechanically for 30 min or until dissolved.




Figure 2. Standard Solution



Analyte	Retention Time (min)	Height	Area	Area %RSD
Atorvastatin Calcium	13.862	121.0	2550	0.2
N = 5 Injections				




Analyte	Retention Time (min)	Height	Area	Area %RSD
Atorvastatin Calcium	12.954	151.6	2594	0.1
N = 5 Injections				




Analyte	Retention Time (min)	Height	Area	Area %RSD
Atorvastatin Calcium	8.603	207.2	2550.4	0.2
N = 5 Injections				




Figure 3. System Suitability Solution



Peak No.	Analyte	Retention Time (min)	Symmetry Factor	Resolution
1	Atorvastatin Calcium	13.848	0.96	7.99
2	Related Compound H	18.961	0.95	



Peak No.	Analyte	Retention Time (min)	Symmetry Factor	Resolution
1	Atorvastatin Calcium	12.805	0.98	9.31
2	Related Compound H	17.500	0.95	



Peak No.	Analyte	Retention Time (min)	Symmetry Factor	Resolution
1	Atorvastatin Calcium	8.470	1.02	8.00
2	Related Compound H	11.276	1.00	



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Luna, Ultramex, and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

