

APPLICATION

The Use of Phenogel GPC Columns for Environmental and Biomonitoring Applications

Pierre-Luc Cloutier¹, Paule-Émilie Groleau¹, Mélanie Desrosiers¹, Brian Rivera², Michael Klein², and Ana Valdez²

¹Centre d'expertise en analyse environnementale du Québec, ministère du Développement durable, de l'Environnement et de la Lutte contre les changements climatiques, Québec City, QC, Canada

²Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Brian Rivera
Product Manager

In addition to chromatography, Brian also has a passion for ice cream-making, and enjoys experimenting with bold, new flavors.

This technical note provides examples and an application of preparative GPC (Gel Permeation Chromatography) columns for the automated purification of environmental samples in concordance to the 3640A EPA method. Results shown are the analysis and purification of PAHs by Phenogel™ GPC columns and the removal of biomaterial matrix interferences that are known to complicate the analysis of PAHs in biological tissues.

Introduction

Polycyclic aromatic hydrocarbons (PAHs) are contaminants resulting from incomplete combustion of organic materials. The Environmental Protection Agency has labeled 32 PAHs as priority pollutants because of their mutagenic and carcinogenic nature¹. As such, PAHs are monitored in sediments, soils, wildlife, water, and air samples².

Since GC is used for PAH analysis and determination, a general cleanup step is recommended to remove biomaterial matrix interferences, such as lipids, proteins, and steroids³. GPC is the primary method since Krahn and colleagues implemented the methodology, showing it's superior over gravity flow⁴. This technique is known to be robust, efficient, and enables lower maintenance of GC instrumentation.

The Phenogel 100 Å preparative GPC column and guard were used to remove lipids and macromolecules efficiently from mussels and salmon samples prior to PAH analysis⁵. Those matrices were selected due to their difficulty for PAH analysis and for sufficient purification to provide quality results at low levels (pg/g and ng/g).

Materials and Methods

A mixture of 49 native PAHs, ¹³C labeled and deuterated internal standards were purchased from Accustandard (New Haven, CT). Reference standard solution for GPC performance containing corn oil (250 mg/mL), bis(2-ethylhexyl)phthalate (5 mg/mL), methoxychlor (1 mg/mL), perylene (0.2 mg/mL), and sulfur (0.8 mg/mL) was purchased from Restek (Bellefonte, PA). The solution was diluted 50:1 in methylene chloride.

Fish and mussel tissue were prepared by homogenization in a blender and extracted by QuEChERS extraction with ethyl acetate to measure the efficiency of the matrix by the GPC preparative column.

Table 1 and **Table 2** list the internal standards and the 49 natural PAHs that were added in a blank to measure the recovery of the PAHs for the collected fraction.

Table 1.

Internal standards

D ₁₀ -2-Methylnaphthalene	D ₁₂ -Chrysene
D ₁₀ -Acenaphthene	¹³ C ₄ -Benzo[a]pyrene
¹³ C ₆ -Anthracene	D ₁₀ -Dibenz[a,h]anthracene
D ₁₀ -Pyrene	

Table 2.

Forty-nine natural PAHs

Naphthalene	Fluoranthene	Benzo[b]fluoranthene	Dibenz[a,h]anthracene
2-Methylnaphthalene	Pyrene	Benzo[k]fluoranthene	7H-Dibenzo[c,g]carbazole
1-Methylnaphthalene	2-Methylfluoranthene	Benzo[j]fluoranthene	Benzo[g,h,i]perylene
2-Chloronaphthalene	Benzo[c]phenanthrene	7,12-Dimethylbenz[a]anthracene	Anthanthrene
1-Chloronaphthalene	Benzo[c]acridine	Benzo[e]pyrene	Dibenzo[a,e]fluoranthene
1,3-Dimethylnaphthalene	Benz[a]anthracene	Benzo[a]pyrene	Dibenzo[a,i]pyrene
Acenaphthylene	Chrysene	Perylene	Coronene
Acenaphthene	3-Methylchrysene	3-Methylcholanthrene	Dibenzo[a,e]pyrene
2,3,5-Trimethylnaphthalene	2-Methylchrysene	Dibenzo[a,h]acridine	Dibenzo[a,i]pyrene
Fluorene	6-Methylchrysene	Dibenzo[a,i]acridine	Dibenzo[a,h]pyrene
Phenanthrene	5-Methylchrysene	Dibenz[a,i]anthracene	
Anthracene	4-Methylchrysene	Indeno[1,2,3-cd]pyrene	
Carbazole	1-Nitropyrene	Dibenz[a,c]anthracene	



HPLC analysis was performed using an Agilent[®] 1100 LC system (Agilent Technologies, Palo Alto, CA, USA) with an upper pressure limit of 400 bar, equipped with a binary pump, autosampler, UV-Vis detector and fraction collector. Phenogel[™] 5 μm 100 \AA 300 x 21.2 mm preparative column and 50 x 21.2 mm guard column (Phenomenex, USA) was used. The flow rate was 5 mL/min and the solvent used was methylene chloride. For both reference standards and samples, 1 mL was injected for each. Run times were 30 min per injection. Injections were monitored at 254 nm.

The eluate was then concentrated to 0.5 mL in isoctane and internal injection (D_8 -Naphthalene, D_8 -Acenaphthylene, D_{10} -Phenanthrene, D_{10} -Fluoranthene, D_{12} -Benz[a]anthracene, D_{12} -Benzo[e]pyrene, D_{12} -Benzo[g,h,i]perylene) standards were added prior to injection in GC-HRMS.

Results and Discussion

The system suitability was determined using the EPA standard solution. **Figure 1** shows 5 major compounds as the expected chromatogram by the US EPA Method 3640A, but with less solvent consumption and faster chromatography than previously reported methods. Compound retention times showed good repeatability (**Table 3**), with less than 0.2 % of RSD for the retention time for the calculated area that corresponds to the injected volume.

Figure 1.
Representative chromatogram for EPA 3640A test mix




Table 3.
Retention times for EPA 3640A test mix

Analyte	Retention Time (min)	Standard Deviation	%RSD
Corn Oil	11.779	0.019	0.16
Bis(2-ethylhexyl) phthalate	13.466	0.022	0.16
Methoxychlor	14.55	0.048	0.33
Perylene	19.362	0.083	0.43
Sulphur compounds	21.022	0.031	0.15

The lipid removal efficiency for fish and mussel samples was evaluated by a 2.25 mL injection of extracts in the GPC system and the recording of the UV absorbance at 254 nm. Fraction collection analysis was performed to selectively separate lipids from PAHs and the optimal fraction collection range was determined to be 15.5 to 20 minutes. Specificity was determined by the % recovery of deuterated standards by GC/MS (**Table 4**).

Table 4.
Recovery of ^{13}C -labeled and deuterated internal standards

Compounds	% Recovery
D_{10} -Methylnaphthalene	106.56
D_{10} -Acenaphthene	105.76
$^{13}\text{C}_6$ -Anthracene	106.54
D_{10} -Pyrene	99.15
D_{12} -Chrysene	94.02
$^{13}\text{C}_4$ -Benzo[a]pyrene	106.23
D_{10} -Dibenz[a,h]anthracene	120.40

Figure 2 and **Figure 3** show representative chromatograms from mussel and salmon, respectively. The major peaks observed between 10 and 15 min are attributed to lipids and high matrix content of the fish liver extracts. All the native PAHs were adequately recovered.

Figure 2.
Representative chromatogram for 5 g mussel extract





Figure 3.
Representative chromatogram for 2.5 g fish liver extract



Conclusion

This method showed that Phenogel™ 5µm 100 Å preparative GPC columns efficiently removed lipids from salmon and mussel samples with high recovery of the targeted 49 PAHs. The method showed good repeatability and could easily be automated for overnight injection.

Future studies will include an integrated method for the complete extraction and purification of different complex biological matrices. Additionally, PCBs, PBDEs, PCDD/Fs, and polychlorinated naphthalene have also shown interesting preliminary results and further development will be performed to analyze all persistent organic pollutants and PAHs simultaneously. Sediments and soils have also shown intriguing results and further validation will be performed for the chemical contaminants enumerated before.

References

1. Yan, Jian, Lei Wang, Peter P Fu, and Hongtao Yu. "Photomutagenicity of 16 Polycyclic Aromatic Hydrocarbons from the US EPA Priority Pollutant List." *Mutation Research/Genetic Toxicology and Environmental Mutagenesis (2004)*: 99-108.
2. Hunt CD, Slone E. Long-term monitoring using resident and caged mussels in Boston Harbor yield similar spatial and temporal trends in chemical contamination. *Marine Environ Res* 2010;70:343-57.
3. U.S. Environmental Protection Agency. "Extraction and Lipid Separation of Fish Samples for Contaminant Analysis and Lipid Determination." Standard Operating Procedure SOP No. HC521A
4. Krahn, M. M., C. A. Wigren, R. W. Pearce, L. K. Moore, R. G. Bogar, W. D. MacLeod, Jr., S-L. Chan, and D. W. Brown. 1988b. Standard analytical procedures of the NOAA National Analytical Facility, 1988: New HPLC cleanup and revised extraction procedures for organic contaminants. U.S. Dep. Commer. Tech. Memo. NMFS F/NWC-153, 52 p.
5. Cloutier, P.-L., Fortin, F., Fournier, M., Brousseau, P., Groleau, P.-E., Desrosiers, M. 2014. Method development for the determination of PAHs in biological samples. *Journal of Xenobiotics* 4 (4897) : 70-72.

Ordering Information

Phenogel SEC/GPC Columns

5 µm Analytical Columns (mm)			SecurityGuard™ Cartridges (mm)
300 x 7.8			4 x 3.0*
Pore Size	MW Range		
50 Å	100-3 K	00H-0441-K0	AJ0-9292
100 Å	500-6 K	00H-0442-K0	AJ0-9292
500 Å	1 K-15 K	00H-0443-K0	AJ0-9292
10 ³ Å	1 K-75 K	00H-0444-K0	AJ0-9292
10 ⁴ Å	5 K-500 K	00H-0445-K0	AJ0-9292
10 ⁵ Å	10 K-1,000 K	00H-0446-K0	AJ0-9292
10 ⁶ Å	60 K-10,000 K	00H-0447-K0	AJ0-9292
300 x 7.8			4 x 3.0*
Mixed Beds			
Linear(2)	100-10,000 K	00H-3259-K0	AJ0-9292

for 3.2-8.0 mm ID

5 µm Narrow Bore (NB) Columns (mm)			SecurityGuard™ Cartridges (mm)
300 x 4.6			4 x 3.0*
Pore Size	MW Range		
50 Å	100-3 K	00H-0441-E0	AJ0-9292
100 Å	500-6 K	00H-0442-E0	AJ0-9292
500 Å	1 K-15 K	00H-0443-E0	AJ0-9292
10 ³ Å	1 K-75 K	00H-0444-E0	AJ0-9292
10 ⁴ Å	5 K-500 K	00H-0445-E0	AJ0-9292
10 ⁵ Å	10 K-1,000 K	00H-0446-E0	AJ0-9292
10 ⁶ Å	60 K-10,000 K	00H-0447-E0	AJ0-9292
300 x 4.6			4 x 3.0*
Mixed Beds			
Linear(2)	100-10,000 K	00H-3259-E0	AJ0-9292

for 3.2-8.0 mm ID

10 µm Analytical Columns (mm)			SecurityGuard™ Cartridges (mm)
300 x 7.8			4 x 3.0*
Pore Size	MW Range		
50 Å	100-3 K	00H-0641-K0	AJ0-9292
100 Å	500-6 K	00H-0642-K0	AJ0-9292
500 Å	1 K-15 K	00H-0643-K0	AJ0-9292
10 ³ Å	1 K-75 K	00H-0644-K0	AJ0-9292
10 ⁴ Å	5 K-500 K	00H-0645-K0	AJ0-9292
10 ⁵ Å	10 K-1,000 K	00H-0646-K0	AJ0-9292
10 ⁶ Å	60 K-10,000 K	00H-0647-K0	AJ0-9292
300 x 7.8			4 x 3.0*
Mixed Beds			
Linear(2)	100-10,000 K	00H-3260-K0	AJ0-9292

for 3.2-8.0 mm ID

5 µm Preparative Columns (mm)		Guards	
300 x 21.2		50 x 21.2	
Pore Size	MW Range		
100 Å	500-6 K	00H-0442-P0	03B-0642-P0

10 µm Preparative Columns (mm)		Guards	
300 x 21.2		50 x 21.2	
Pore Size	MW Range		
100 Å	500-6 K	00H-0642-P0	03B-0642-P0

* SecurityGuard Analytical Cartridges require holder, Part No.: KJ0-4282

Guard Cartridge Holder

Part No.	Description
KJ0-4282	Reusable Holder (SecurityGuard Kit)

Note: SecurityGuard cartridges for Non-Aqueous Polymer GPC columns are not compatible with HFIP solvent.

Column Union

Part No.	Description
AQ0-8507	Zero Dead Volume Union, SS, with 10-32 fittings

Note: Additional union (AQ0-8507) may be necessary for SecurityGuard to fit in column oven with less than 30 cm length capacity.

Phenogel columns are routinely shipped in THF. Please contact your Phenomenex representative for information about other shipping solvents.



APPLICATION

Australia

t: +61 (0)2-9428-6444
 f: +61 (0)2-9428-6445
 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
 f: +43 (0)1-319-1300
 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
 t: +32 (0)2 511 8666 (Dutch)
 f: +31 (0)30-2383749
 beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
 f: +1 (310) 328-7768
 info@phenomenex.com

China

t: +86 (0)20 2282-6668
 f: +86 (0)20 2809-8130
 chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
 f: +45 4810 6265
 nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
 f: +33 (0)1 30 09 21 11
 franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
 f: +49 (0)6021-58830-11
 anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
 f: +91 (0)40-3012 2411
 indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
 f: +44 1625-501796
 eireinfo@phenomenex.com

Italy

t: +39 051 6327511
 f: +39 051 6327555
 italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
 f: +31 (0)30-2383749
 nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
 f: 001-310-328-7768
 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
 f: +31 (0)30-2383749
 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
 f: +64 (0)9-4780952
 nzinfo@phenomenex.com

Norway

t: +47 810 02 005
 f: +45 4810 6265
 nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
 f: +1 (310) 328-7768
 info@phenomenex.com

Spain

t: +34 91-413-8613
 f: +34 91-413-2290
 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
 f: +45 4810 6265
 nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
 f: +44 (0)1625-501796
 ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
 f: +1 (310) 328-7768
 info@phenomenex.com

**All other countries
Corporate Office USA** 

t: +1 (310) 212-0555
 f: +1 (310) 328-7768
 info@phenomenex.com


 guarantee

If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Phenogel and SecurityGuard are trademarks of Phenomenex. ChemStation is a trademark and Agilent is a registered trademark of Agilent Technologies, Inc.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162, 362

Caution: this patent applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.

© 2015 Phenomenex, Inc. All rights reserved.

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com