



## Lamivudine Related Substances Method per IP Monograph

Krishna Chaitanya Routhu<sup>1</sup>, Rajesh Babu Dandamudi, PhD<sup>1</sup>, and Bryan Tackett, PhD<sup>2</sup>

<sup>1</sup>India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India

<sup>2</sup>Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview

Lamivudine is an antiretroviral medication that belongs to a class of drugs known as nucleoside reverse transcriptase inhibitors (NRTIs) and plays a crucial role in combination therapy for managing the human immunodeficiency virus (HIV) infection. Additionally, Lamivudine is also utilized in the treatment of chronic hepatitis B virus (HBV) infection. The development of a quick and efficient analysis of Lamivudine and its related substances is of interest for India generic drug manufacturers. In this application note, we report the separation of Lamivudine and its related substances using a Luna Omega 5 µm C18 column and a Kinetex 5 µm C18 column according to the Indian Pharmacopoeia (IP) Monograph for Lamivudine, which references use of a 25 cm x 4.6 mm stainless-steel column packed with octadecylsilane bonded to porous silica (5 µm).


System suitability per IP Monograph for the Lamivudine Related Substances is resolution of not less than (NLT) 10.0 between Lamivudine and Salicylic Acid, theoretical plates NLT 5000, and a tailing factor not more than 1.5 for the peak due to Lamivudine.

The results clearly show that the system suitability criteria met and surpassed the minimum requirement. This illustrates that the Luna Omega 5 µm C18 and Kinetex 5 µm C18 columns meet all requirements for system suitability as outlined in the Indian Pharmacopoeia monograph for Lamivudine. The core-shell Kinetex 5 µm C18 column yielded a faster run time (>25 %) with comparable peak shape, while maintaining the resolution achieved with the fully porous Luna Omega C18 column.

All solutions were prepared as indicated in the IP Monograph for Lamivudine. Salicylic Acid was purchased from Sigma-Aldrich®. The following certified reference standards (CRS) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: Allee Kastner CS 30026 F - 67081 Strasbourg (France):

- Y0000425, Lamivudine CRS

Figure 1. Lamivudine



### LC-UV Conditions


|                          |   |
|--------------------------|---|
| <b>Column:</b>           | Luna™ Omega 5 µm C18<br>Kinetex™ 5 µm C18   |
| <b>Dimensions:</b>       | 250 x 4.6 mm  |
| <b>Part No.:</b>         | <a href="#">00G-4785-E0</a> (Luna Omega)<br><a href="#">00G-4601-E0</a> (Kinetex) |
| <b>Mobile Phase:</b>     | Mobile Phase Table 1  |
| <b>Flow Rate:</b>        | 1 mL/min (Isocratic)  |
| <b>Injection Volume:</b> | 10 µL   |
| <b>Temperature:</b>      | 35 °C   |
| <b>Detector:</b>         | UV @ 277 nm   |
| <b>System:</b>           | Waters® ACQUITY Arc® HPLC   |

Table 1. Preparation of Solutions


| Solution               | Composition  |
|------------------------|--|
| <b>Mobile Phase</b>    | Methanol / Buffer (5:95, v/v)<br><br>Buffer: 0.19 % w/v of Ammonium Acetate adjusted to pH 3.8 with Glacial Acetic Acid.   |
| <b>Test Solution</b>   | Dissolve 50 mg of <i>Lamivudine</i> CRS in 70 mL of <b>Mobile Phase</b> and dilute to 100 mL with <b>Mobile Phase</b> .  |
| Reference Solution (a) | Dilute 1 mL of the <b>Test Solution</b> to 100 mL with <b>Mobile Phase</b> . Dilute 1 mL of this solution to 10 mL with <b>Mobile Phase</b> .                            |
| Reference Solution (b) | Dissolve 5 mg of Salicylic Acid in <b>Mobile Phase</b> and dilute to 100 mL with <b>Mobile Phase</b> . Dilute 1 mL of this solution to 100 mL with <b>Mobile Phase</b> . |
| Reference Solution (c) | A solution containing 0.001 per cent w/v each of <i>Lamivudine</i> CRS and Salicylic Acid in <b>Mobile Phase</b> .   |



Figure 2. Reference Solution (c)



| Peak No.         | Analyte        | Retention Time (min) | Area   | Area %RSD | Tailing Factor | Theoretical Plates | Resolution |
|------------------|----------------|----------------------|--------|-----------|----------------|--------------------|------------|
| 1                | Lamivudine     | 13.14                | 292683 | 0.51      | 1.05           | 20536              | 30.81      |
| 2                | Salicylic Acid | 30.55                | 52125  | 0.65      | 0.99           | 25621              |            |
| N = 6 Injections |                |                      |        |           |                |                    |            |



| Peak No.         | Analyte        | Retention Time (min) | Area   | Area %RSD | Tailing Factor | Theoretical Plates | Resolution |
|------------------|----------------|----------------------|--------|-----------|----------------|--------------------|------------|
| 1                | Lamivudine     | 7.81                 | 302193 | 0.75      | 0.98           | 25787              | 31.83      |
| 2                | Salicylic Acid | 18.14                | 49676  | 0.72      | 1.10           | 25491              |            |
| N = 6 Injections |                |                      |        |           |                |                    |            |



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auiinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5019 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 6559 4364  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**🌐 All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/phx-terms-and-conditions-of-sales](http://www.phenomenex.com/phx-terms-and-conditions-of-sales).

**Trademarks**

Luna, Kinetex, and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY Arc are registered trademarks of Waters Technologies Corporation. Sigma-Aldrich is a registered trademark of Merck KGaA, Darmstadt, Germany.

**Disclaimer**

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Waters Technologies Corporation or Merck KGaA, Darmstadt, Germany.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

