

APPLICATIONS

Comparison of Polymeric and Silica-Based Columns for Oligonucleotide Analysis by Ion-Pair Reversed Phase

Ivan Lebedev, Brian Rivera
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Therapeutic oligonucleotides represent a recent breakthrough in the pharmaceutical industry. However, characterization of oligos, specifically by ion-pair reversed phase liquid chromatography (IP-RPLC), can be quite challenging. Oligos are manufactured by solid phase synthesis, where nucleotides are added in a step-wise manner. As such, process-related impurities such as n-1 and n+1 must be characterized. Further complicating methods are common mobile phase additives such as perfluorinated alcohols (HFIP) and an alkylamine ion-pair that facilitate the hydrophobic interaction of the negatively charged oligo. These mobile phase conditions combined with high temperature are often necessary to ensure good peak shapes but result in very harsh conditions that can encourage silica dissolution.

As such, the LC particle is an important consideration for column selection. The particle must be stable at neutral pH when run at moderately high temperature. It might then seem appropriate to run a polymer-based particle such as a polystyrene divinylbenzene to ensure method robustness and column stability. However, polymer-based particles typically do not have the efficiency necessary to separate the closely eluting impurities associated to oligos. Alternatively, superficially porous (core-shell) particles can be used to improve efficiency and resolution. Here, we investigate the chromatographic performance between two PS-DVB particles and the bioZen™ 2.6 µm Oligo column, which uses a pH stable organo-silica superficially porous particle packed in bio-inert hardware.

In comparing the 2.6 µm bioZen Oligo to a 5 µm PS-DVB of similar dimensions in the separation of a 2'-MOE Gapmer, there is a marked difference in chromatographic performance, with the 2.6 µm particle giving superior peak height (**Figure 1**). As such, the increase in chromatographic efficiency offered by the bioZen Oligo improves the separation of different process related impurities. This also may be due to the functionality of the C18 stationary phase that is bonded to the superficially porous particle, which has a better methylene selectivity. A similar performance for an RNA phosphorothioate is also observed, as seen in **Figure 2**. Notably, many short mer failure sequences are not detected using the 5 µm PS-DVB particle platform.

In comparing separation of a 5'-Amino C12 oligonucleotide on the 2.6 µm against a 4 µm PS-DVB column of similar dimensions (**Figure 3**), the superficially porous core-shell particle shows a slightly improved peak height and improved separation. A similar profile is observed with the BNA oligo (**Figure 4**), though the superficially porous particle yields slightly more peak areas and improved recoveries for some earlier eluting impurities. This improved recovery might be attributed to the bio-inert hardware of the bioZen Oligo, which tends to improve recovery and peak shape for oligos, especially using moderately high HFIP concentrations.

In summary, when selecting a particle for oligo analysis although a PS-DVB particle might seem like an appropriate choice, one may be making concessions for peak shape and efficiency, which might directly effect method sensitivity. Additionally, separation may not be optimal which can impact impurity analysis and characterization. The bioZen Oligo however, affords scientists a versatile combination of performance gains and robustness, along with improved recovery from bio-inert hardware.

LC Conditions

- Columns:** bioZen 2.6 µm Oligo ([00B-4790-AN](#))
5 µm PS-DVB
4 µm PS-DVB
- Dimension:** 50 x 2.1 mm
- Mobile Phase:** Figures 1, 2:
A: 12.5 mM HFIP, 0.5 mM Hexylamine in Water
B: 12.5 mM HFIP, 0.5 mM Hexylamine in Methanol
Figures 3, 4:
A: 100 mM HFIP, 4 mM TEA in Water
B: 100 mM HFIP, 4 mM TEA in Methanol
- Gradient:** 5-75 % B in 14 minutes (Figures 1,2)
5-30 % B in 14 minutes (Figures 3,4)
- Flow Rate:** 0.3 mL/min
- Injection:** 2 µL (1 µg/mL)
- Temperature:** 65 °C (Figure 1,2)
45 °C (Figure 3,4)
- Detection:** UV @ 260 nm
- Sample:** As indicated

Figure 1. Comparison of Impurity Profile, 2'-MOE Gapmer




Figure 2. Comparison of Impurity Profile, RNA Phosphorothioate




Figure 3. Comparison of Impurity Profile, 5'-Amino C12 Oligonucleotide





Figure 4. Comparison of Impurity Profile, BNA



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 All other countries/regions

Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen is a trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.