

Optimal Separation of Polar Anionic Pesticides From Fruits and Vegetables with Unique HPLC Column Selectivity

Luigi Margarucci, PhD¹, Aidan Harrison, PhD², Jorg Baute, PhD³, Richard Jack, PhD⁴, Ramkumar Dhandapani, PhD⁴, Sean Orlowicz⁴, Samuele Scurati, PhD⁵

¹Phenomenex Italy, Via M. Serenari, 15/D, 40013 Castel Maggiore (BO), Italy, ²Phenomenex Ltd., Queens Avenue, Hurdsfield Ind. Est., Macclesfield, Cheshire SK10 2BN, UK, ³Phenomenex Ltd., Deutshcland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany, ⁴Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA, ⁵AB SCIEX S.r.l. Via Montenapoleone, 8 21181 Milano CF e P.IVA: 06734220962 Italy

Introduction

Analysis of polar pesticides presents multiple challenges including adequate retention, separation of critical pairs, and reproducibility, to name a few. In addition, food matrices can add additional challenges due to the presence of complex matrix components including pigments, fats, and sugars that can interfere with the analyte of interest.

Often, polar, anionic analytes like Glyphosate will utilize QuEChERS or QuPPE sample preparation techniques, followed by HILIC LC-MS/MS methods for chromatographic retention and separation. Historically, these methods are not user friendly, and lack reproducibility necessary for a commercial application.

In this study, we are presenting a unique HPLC selectivity that provides optimal separation of various anionic polar pesticide classes including Glyphosate, Chlorate, Perchlorate, Ethephon, Phosphoric Acid-based pesticides, and N-Ac-Glu pesticides. The study demonstrates robust polar pesticide analysis from real sample matrix.

Materials and Methods

Sample Preparation: QuPPE-PO Method

Weigh Sample Homogenate into 50 mL Centrifuge Tube
Fresh fruits and vegetables (with high water content): 10 g ± 0.1 g
Previously re-hydrated dry fruit: e.g. 13.5 g ± 0.1 g (containing 5 g sample)
Dry commodities (e.g. herbs): 2 g ± 0.02 g

Adjust Water Content of Sample to 10 mL
(Mandatory for matrices w. <80 % water.
If no IL-IS used mandatory for ALL matrices)
e.g. +10 mL water to 2 g of dried mint; +2 mL water to 10 g potato;
+3.5 mL water to 10 g garlic

Add 100 µL Isotopically-labeled Internal Standard (IL-IS mix)

Add Extraction Solvent (10 mL Methanol containing 1 % Formic Acid)

Shake Thoroughly for 1 min to 15 min for Dry Commodities

Preferably Freeze-out Extract until Completely Frozen
e.g. >90 min at -18 °C or ca. 30 min at -80 °C

Centrifuge (5 min at >3,000 g but Preferably >10,000 g);
Preferably cryogenic centrifugation (e.g. at -10 °C)
(If centrifuge not refrigerated, proceed with centrifugation and following step to avoid redissolution of matrix)

dSPE to Remove Lipids for High Oil Content Samples;
(This step may be skipped if sample was centrifuged from
at -10 °C and <20 min)
Transfer 4 mL raw extract into a tube containing 200 mg
C18-sorbent, shake for 1 min and centrifuge (>3,000 g for 5 min)

Withdraw Supernatant and Filter in into a Plastic Autosampler Vial
(Use syringe filter of 0.2 µm pore size; e.g. H-PTFE)
(Plastic vials are recommended as some compounds tend to interact with glass)
(Withdraw vials supernatant quickly after centrifugation
to avoid that matrix components redissolve)

LC-MS/MS and IC-MS/MS Analysis

LC-MS/MS Conditions

Column: Venusil® 3 µm HILIC
Dimension: 100 x 2.1 mm
Part No.: VH931002-0
Injection: A: 0.2 % Formic Acid in Water
B: 0.2 % Formic Acid in Methanol

Gradient: Time (min)	%B
0	2
0.5	2
6	20
7	90
9	90
9.1	2

Flow Rate: 0.3 mL/min

Injection: 0.5 or 1 µL

Temperature: 40 °C

Detector: SCIEX® 7500


Sample Details

- Matrix Extracted following QuPPE Method:
Grain*, Kiwi*, Zucchini*, Rocket*, Soy**
- Raw Matrix, 10ppb Spike, 100ppb Spike
- Different Matrices Dilution depending on matrix composition (1:4* or 1:10**)
- Each sample injected in triplicate

Results


Figure 2. Real Matrix Analysis of Polar Pesticides.

Glyphosate




Matrix	RAW			10ppb			100ppb		
	Mean	CV		Mean	CV	Accuracy	Mean	CV	Accuracy
Grain DIL 1:4	3.27	14.7	14.46	2.32	111.9	96.67	2.56		93.4
Kiwi DIL 1:4	0	0	10.07	3.95	100.7	106.57	10.85		89.02
Zucchini DIL 1:4	0	0	10.75	3.99	107.5	99.02	1.86		99.02
Rocket DIL 1:10	0	0	10.17	4.17	101.7	86.41	4.66		86.41
Soy DIL 1:10	1.22	12.12	10.51	11.6	92.9	81.62	4.67		80.4

Chlorate




	RAW			10ppb			100ppb		
Matrix	Mean	CV		Mean	CV	Accuracy	Mean	CV	Accuracy
Grain DIL 1:4	0.15	24.12		9.41	3.25	92.6	87.14	0.91	86.99
Kiwi DIL 1:4	1.93	3.33		11.4	1.63	94.7	96.56	3.5	94.63
Zucchini DIL 1:4	1.92	4.85		11.3	1.28	93.8	92.35	3.14	90.43
Rocket DIL 1:10	32.65	4.88		41.2	0.91	85.5	103.25	6.28	70.6
Soy DIL 1:10	0.21	11.02		8.38	0.96	81.7	69.82	2.76	69.61

Perchlorate



Matrix	RAW			10ppb			100ppb		
	Mean	CV		Mean	CV	Accuracy	Mean	CV	Accuracy
Grain DIL 1:4	0.44	8.59		10.96	0.64	105.2	101.29	2.49	100.85
Kiwi DIL 1:4	0.47	5.75		11.01	1.98	105.4	104.04	1.03	103.57
Zucchini DIL 1:4	0.9	0.98		11.56	2.61	106.6	101.71	2	100.81
Rocket DIL 1:10	33.24	2.34		41.72	0.69	84.8	112.44	0.23	79.2

Ethephon




Matrix	RAW			10ppb			100ppb		
	Mean	CV		Mean	CV	Accuracy	Mean	CV	Accuracy
Grain DIL 1:4	0	0	9.43	2	94.3	89.12	5.8	89.12	
Kiwi DIL 1:4	0	0	9.09	4.88	90.9	89.28	9.23	89.28	
Zucchini DIL 1:4	0	0	10.23	2.35	102.3	91.6	3.73	91.6	
Rocket DIL 1:10	0	0	9.55	5.77	95.5	81.25	3.31	81.25	
Soy DIL 1:10	0.2	0	9.06	10.37	88.6	84.71	2.24	84.51	

Phosphonic Acid



Matrix	RAW			10ppb			100ppb		
	Mean	CV	Accuracy	Mean	CV	Accuracy	Mean	CV	Accuracy
Grain DIL 1:4	101.2	2.39	111.49	5.31	102.9	202.46	3.69	101.26	
Kiwi DIL 1:4	1.76	45.7	11.69	3.88	99.3	98.34	5.01	96.58	
Zucchini DIL 1:4	2.12	10.61	12.63	3.43	105.1	99.37	4.45	97.25	
Rocket DIL 1:10	24.67	4.46	35.72	1.19	110.5	105.95	1.38	81.28	

N Ac Glu




Matrix	RAW			10ppb			100ppb		
	Mean	CV		Mean	CV	Accuracy	Mean	CV	Accuracy
Grain DIL 1:4	0	0	15.53	3.31	155.3	143.63	6.33	143.63	
Kiwi DIL 1:4	1.76	19.2	13.26	1.73	115	123.77	7.52	122.01	
Zucchini DIL 1:4	0	0	9.81	7.13	98.1	100.03	2.64	100.03	
Rocket DIL 1:10	4.81	13.72	16.2	2.31	113.9	115.35	0.29	110.54	

Results

Figure 3. Retention Time Stability of Glyphosate.

RT Stability



Discussion

Column selectivity plays an important role in providing enhanced chromatographic resolution for critical pairs. In addition, retention of extremely polar analytes are very challenging. In this study, we present optimal separation of polar pesticides on a Venusil HILIC HPLC Column, which is a versatile selectivity with amide functionality that can be run in normal, reverse and HILIC mode. In this study, we have utilized the polar interactions in the Venusil HILIC stationary phase in reverse phase mode to obtain enhanced retention of polar pesticides. The chromatogram of standards on a SCIEX 7500 shows excellent retention and selectivity for polar pesticides (Figure 1).

Optimal concentration of 0.2% Formic Acid in the mobile phase provided a great balance of peak shape and retention. Traditional reverse phase columns do not retain analyte like Glyphosates, which can fall in the ion suppression zone in real samples and hence can show false positive or negative. With the Venusil HILIC column, enhanced polar selectivity from the un-encapped silica base and from the Amide ligand provides excellent retention which is evident from retention factor for polar pesticides that ranges from 0.7 to 6.6.

Real samples like Grains, Kiwi, Zucchini, Rocket, and Soy were analyzed with this method followed by spiking them with a known concentration of polar pesticides. The method proved to be precise, robust, and accurate for the polar pesticides (Figure 2). In addition, retention time stability of Glyphosate is presented as a representative in Figure 3. Consistent and robust retention of Glyphosate proves that the Venusil HILIC column is a robust stationary phase selectivity for the analysis of polar pesticides by LC-MS/MS.

Conclusion

Venusil HILIC is a versatile HPLC column selectivity that provides enhanced retention and selectivity for polar pesticides. Reproducible retention, optimal selectivity, and precise and accurate results prove that the Venusil HILIC column is the ideal choice for polar pesticide analysis.


In addition to providing consistent retention, the Venusil HILIC column offers short run time of less than 6 minutes for high throughput analysis.

Unlike traditional HILIC methods, the method demonstrated here provides stable retention in reverse phase by utilizing polar interactions from the stationary phase. Thus, the developed method is easy to adopt in labs running routine polar pesticide analysis.

Terms and Conditions
Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.
Trademarks
Venusil is a registered trademark of Bonna-Agela Technologies Inc. SCIEX is a registered trademark of AB SCIEX Pte. Ltd.
Disclaimer
Comparative separations may not be representative of all applications.
FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2022 Phenomenex, Inc. All rights reserved.

Results

Figure 1. Retention Profile of Polar Pesticides on a Venusil HILIC HPLC Column.



Analyte	RT (min)	Analyte	RT (min)
AMPA	0.95	Phosphonic Acid	2.9
MH	1.1	Ethephon	3
Glufosinate	1.4	Chlorate	3.7
MPPA	1.8	Fosetyl	3.9
Glyphosate	2	Perchlorate	4.5
N Acetyl Glu	2.2		