

# These Phases ROCK Your LC Laboratory



**Chromatographers in the following industries can benefit from the advantages of Core-Shell technology:**

-  **Agriculture**
-  **Forensics**
-  **Clinical**
-  **Life Science**
-  **Environmental**
-  **Pharmaceutical**
-  **Food and Beverage**
-  **Chemical/Industrial**
-  **Consumer Care**

## guarantee

If you are not completely satisfied with Kinetex core-shell columns, send in your comparative data to a similar product with the Kinetex column within 45 days for a FULL REFUND.

## Efficiency



## Table of Contents



The Technology .....	4
Core-Shell Advantage	
UHPLC .....	6
HPLC and UHPLC .....	7
HPLC .....	8
Prep LC .....	10
Significant Cost Savings .....	11
Core-Shell Fits Your System.....	12
Selectivity - The Phases.....	14
Method Development Flexibility .....	18
Simple Selection of the Suitable Column.....	19
Ordering Information .....	21


**Resolution**

**Productivity**

**Selectivity**




## Complete Scalable Solution from UHPLC to HPLC to PREP LC



## Better Performance than Fully Porous Particles

Using sol-gel processing techniques that incorporate nano structuring technology, a durable, homogeneous porous shell is grown on a solid silica core. This highly optimized process combined with industry leading column packing technology produces highly reproducible columns that generate extremely high plate counts.





Fully Porous		Kinetex Core-Shell		Average Efficiency Gain with Kinetex*
5 $\mu\text{m}$	vs	5 $\mu\text{m}$		90% Higher
3 $\mu\text{m}$	vs	2.6 $\mu\text{m}$		85% Higher
1.7 $\mu\text{m}$	vs	1.7 $\mu\text{m}$		20% Higher
1.7 $\mu\text{m}$	vs	1.3 $\mu\text{m}$		50% Higher

\* May not be representative of all applications

## Kinetex Core-Shell Technology

- Obtain higher throughput without sacrificing resolution
- Easy method transfer across LC system platforms
- Reduce solvent consumption with faster analysis
- Reach lower levels of detection and quantitation



2009

Small Molecules



2010

Synthetic  
Oligonucleotides



2011

Peptides / Proteins




2014

Kinetex EVO



## Kinetex EVO C18

Kinetex EVO C18 uses a patented organo-silica grafting process which incorporates uniform stabilizing ethane cross-linking to provide resistance to high pH attack while maintaining mechanical strength of the core-shell particle.



# Core-Shell Advantage for UHPLC





## Get the Most Performance Out of Your UHPLC System


Kinetex 1.3 and 1.7  $\mu\text{m}$  Core-Shell Technology produces increased efficiencies over traditional sub-2  $\mu\text{m}$  columns on the market, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity, so you can get the most out of every UHPLC analysis.



## Incredible Efficiency Gains



25% Higher Efficiency than the traditional fully porous sub-2  $\mu\text{m}$  column



68% Higher Efficiency than the traditional fully porous sub-2  $\mu\text{m}$  column

Easy Installation on UHPLC Systems:  
[www.phenomenex.com/kinetex](http://www.phenomenex.com/kinetex)

**Conditions for all columns:**

- Columns:** Kinetex 1.7  $\mu\text{m}$  C18  
Kinetex 1.3  $\mu\text{m}$  C18  
ACQUITY UPLC® BEH 1.7  $\mu\text{m}$  C18
- Dimensions:** 50 x 2.1 mm
- Mobile Phase:** Acetonitrile / Water (50:50)
- Flow Rate:** 0.5 mL/min
- Temperature:** Ambient
- Detection:** UV @ 254 nm
- Instrument:** Waters® ACQUITY® UPLC®
- Sample:** 1. Acetophenone  
2. Benzene  
3. Toluene  
4. Naphthalene



Efficiency calculated from peak 4 in each chromatogram. Waters, ACQUITY, and UPLC are registered trademarks, and BEH Technology is a trademark of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

## Column Protection with No Loss in Performance

Trap contaminants and microparticulates within the SecurityGuard™ ULTRA guard cartridge system and learn how to greatly extend UHPLC column life at:

[www.phenomenex.com/SecurityGuardULTRA](http://www.phenomenex.com/SecurityGuardULTRA)



# Core-Shell Advantage for HPLC and UHPLC


## A Versatile Upgrade for HPLC and UHPLC

On a low volume HPLC or UHPLC system Kinetex 2.6  $\mu\text{m}$  columns will perform like a fully porous sub-2  $\mu\text{m}$  column, providing up to 3x the efficiency of 5  $\mu\text{m}$  and double the efficiency of 3  $\mu\text{m}$  fully porous media. Dramatically improve the productivity and performance of your existing methods with the use of shorter Kinetex columns, all while decreasing your solvent usage!




## Sub-2 $\mu\text{m}$ Performance with Kinetex 2.6 $\mu\text{m}$ on HPLC and UHPLC Systems


### Kinetex 2.6 $\mu\text{m}$ C18



### Kinetex 2.6 $\mu\text{m}$ C18



### Agilent® ZORBAX® 1.8 $\mu\text{m}$ SB-C18



#### Conditions for all columns same except where noted:

**Columns:** Kinetex 2.6  $\mu\text{m}$  C18 (Agilent 1100\*)  
Kinetex 2.6  $\mu\text{m}$  C18 (Agilent 1290)  
ZORBAX 1.8  $\mu\text{m}$  SB-C18 (Agilent 1290)

**Dimensions:** 100 x 4.6 mm

**Mobile Phase:** A: Water with 0.1% TFA  
B: Acetonitrile with 0.1% TFA

Gradient:	Time (min)	% B
	0	10
	20	70

**Flow Rate:** 1.2 mL/min

**Temperature:** Ambient

**Detection:** UV @ 210 nm

**Sample:** Mupirocin degradants

\*Agilent 1100 was optimized with the Core-Shell Performance Enhancement Kit AQO-8892. Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc. Phenomenex is not affiliated with Agilent Technologies. Comparative separations may not be representative of all applications.

## Maximize Performance with Kinetex 2.6 $\mu\text{m}$

Decrease the system dwell volume of your HPLC instrument with the pre-cut tubing and column fittings found in the Core-Shell Performance Enhancement Kit and instantly improve the observed performance of your Kinetex core-shell 2.6  $\mu\text{m}$  column.

[www.phenomenex.com/enhancement](http://www.phenomenex.com/enhancement)




## Instantly Improve Your 3.5 µm Pharmacopoeia Methods

Immediately improve resolution, productivity, and sensitivity of your current 3.5 µm HPLC methods following the requirements of Ph. Eur. or USP with the new Kinetex 3.5 µm Core-Shell Technology. This core-shell particle was developed for use on standard or older model HPLC systems that may have low pressure limitations.




## Clarithromycin according to Ph. Eur. Monograph 1651

### Kromasil® 3.5 µm C18



Using Kinetex 3.5 µm XB-C18 for the Clarithromycin method provides the user with better resolution and narrower peaks for better quantification compared to a fully porous 3.5 µm C18 material.

### Kinetex 3.5 µm XB-C18



#### Conditions for both columns:

**Columns:** Kinetex 3.5 µm XB-C18  
Kromasil 3.5 µm C18

**Dimensions:** 100 x 4.6 mm

**Mobile Phase:** A: 4.76 g/L Potassium phosphate pH 4.4

B: Acetonitrile

Gradient	Time (min)	% B
	0	25
	32	65
	34	65

**Injection:** 10 µL

**Detection:** UV @ 205 nm

**Sample:** Clarithromycin

Kromasil is registered trademark of AkzoNobel Pulp And Performance Chemicals AB. Phenomenex is not affiliated with AkzoNobel. Comparative separations may not be representative of all applications.

## Allowable Adjustments to Pharmacopoeia Methods

Recent revisions of the Ph. Eur. and USP detail what “allowable adjustments” can be made to methods without resorting to revalidation. Find out what changes you are allowed to make by visiting:

[www.phenomenex.com/allowable](http://www.phenomenex.com/allowable)




# Core-Shell Advantage for HPLC

## Instantly Improve 5 $\mu$ m and 3 $\mu$ m Methods


Immediately improve resolution, productivity, and sensitivity of your current 3 $\mu$ m and 5 $\mu$ m HPLC methods with Kinetex 5 $\mu$ m Core-Shell Technology. This core-shell particle was specifically developed for use on standard or older model HPLC systems that may have low pressure limitations.



GL Sciences Inertsil<sup>®</sup> 5 $\mu$ m ODS-3



Phenomenex Kinetex 5 $\mu$ m C18



### Critical Advantages of Increased Resolution:

- Save time and money by adjusting method
- Separate critical pairs or difficult mixtures
- Less risk of failing system suitability requirements

#### Conditions for all columns:

**Columns:** Kinetex 5 $\mu$ m C18  
Inertsil ODS-3 5 $\mu$ m C18

**Dimensions:** 250 x 4.6 mm

**Mobile Phase:** A: Water with 0.1% TFA

B: Acetonitrile with 0.1% TFA

Gradient:	Time (min)	% B
	0	10
	40	70

**Flow Rate:** 1.0 mL/min

**Temperature:** Ambient

**Detection:** UV @ 210 nm

**Sample:** Pharmaceutical degradation sample

Inertsil is a registered trademark of GL Sciences Inc. Phenomenex is not affiliated with GL Sciences Inc. Comparative separations may not be representative of all applications.

### First and Only Core-Shell Material for Axia<sup>™</sup> Preparative Purifications

Axia packed Kinetex 5 $\mu$ m columns will provide incredible gains in efficiency and performance for any of your Prep LC methods. Learn more online at:

[www.phenomenex.com/kinetexprep](http://www.phenomenex.com/kinetexprep)



## Preparative Column Packing Technology



An advanced preparative column packing and hardware design, Axia incorporates patented Hydraulic Piston Compression technology that offers increased sorbent bed density and eliminates media bed collapse as a source of premature column failure in preparative HPLC columns. Unlike traditional column packing methods, the Axia packing method is completely automated and monitored by multiple sensors to allow for measurement and recording of all process parameters for every column. The result is a vastly improved packing process that offers the following benefits:

- Extended column lifetimes
- Improved reproducibility: Column-to-Column and Batch-to-Batch
- Efficiencies and peak symmetries on par with analytical separations
- Increased column stability under high flow rates





## Higher Efficiency!

Start with sharper peaks by taking advantage of the high efficiencies of Kinetex 5µm Axia preparative columns.



## Excellent Loadability!

With narrower peak widths than fully porous columns across every sample load, Axia packed Kinetex 5µm columns give you the capability of increased sample load and higher throughput for vastly improved purification performance and economics.




Waters and XBridge are registered trademarks of Waters Corp. OBD is a trademark of Waters Corp. Phenomenex is in no way affiliated with Waters Corp. Comparative separations may not be representative of all applications.

# Significant Cost Savings with Kinetex Core-Shell Columns



## Reduce Runtimes, Increase Productivity, and Cut Costs Today



**Column:** Kinetex 1.7 µm C18  
**Dimensions:** 100 x 2.1 mm  
**Part No.:** 00D-4475-AN  
**Mobile Phase:** A: 5 mM Ammonium formate pH 3.25 / Acetonitrile (95:5)  
 B: 5 mM Ammonium formate pH 3.25 / Acetonitrile (10:90)

Gradient	Time (min)	% B
	0	30
	1.5	50
	3	56.3
	5	95
	6	95
	6.1	30

**Flow Rate:** 0.4 mL/min  
**Temperature:** 50 °C  
**Detection:** PDA 210-300 nm, extracted channel 280 nm  
**Instrument:** Waters® ACQUITY® equipped with PDA

**Sample:**

1. Antidepressant drug (containing an HCl salt)
2. Hormone therapy #1 (containing a salt)
3. SERM drug (containing basic functional group)
4. CNS drug (containing basic functional group)
5. PPI drug (containing basic functional group)
6. CNS drug (containing basic functional group)
7. CNS drug (containing basic functional group)
8. Hormone therapy #2 (neutral)
9. Oral contraceptive hormone #1 (neutral)
10. Hormone therapy #3 (neutral)
11. Oral contraceptive hormone #2 (neutral)
12. Hormone therapy #4 (neutral)
13. Oral contraceptive hormone (neutral)
14. Hormone therapy #5 (neutral)
15. Hormone therapy #6 (acetate salt of 14)
16. Immunosuppressant drug (macromolecule, containing basic functional group)

It has been shown that the **1.7 µm Kinetex 100 x 2.1 mm column** was capable of resolving 16 different chemical entities with a 6 minute run time. This new analytical method will be used to replace 16 older methods thereby facilitating an **annualised cost saving for the site of €320,000 (\$346,000 USD).**

A. Charles, et. al., Pfizer Grange Castle,  
 Grange Castle Business Park, Clondalkin,  
 Dublin Republic of Ireland

# Core-Shell Fits Your System



## Typical UHPLC Instrumentation



### Kinetex on Agilent® 1290



### Kinetex on JASCO® X-LC



### Kinetex on Shimadzu® Nexera®



### Kinetex on Waters® ACQUITY® UPLC®



# Core-Shell Fits Your System




## Typical HPLC Instrumentation



### Kinetex on Agilent® 1100



### Kinetex on Shimadzu® Prominence® LC-20A




### Kinetex on Waters® Alliance®



### Kinetex on Gilson® HPLC




## Kinetex EVO C18




**Available in: 1.7  $\mu\text{m}$ , 2.6  $\mu\text{m}$  and 5  $\mu\text{m}$**   
Novel pH 1-12 stable C18 that delivers robust methods and improved peak shape for bases.

## Kinetex XB-C18




**Available in: 1.7  $\mu\text{m}$ , 2.6  $\mu\text{m}$ , 3.5  $\mu\text{m}$  and 5  $\mu\text{m}$**   
This unique C18 phase yields increased hydrogen bonding with hydrophobic selectivity, resulting in improved peak shape for basic compounds and increased retention of acidic compounds.

## Kinetex C18




**Available in: 1.3  $\mu\text{m}$ , 1.7  $\mu\text{m}$ , 2.6  $\mu\text{m}$  and 5  $\mu\text{m}$**   
Balanced C18 phase that provides the highest degree of hydrophobic selectivity relative to the other Kinetex phases.

## Kinetex Polar C18




**Available in: 2.6  $\mu\text{m}$**   
Combined C18 and polar modified surface that provide polar and non-polar retention alongside 100% aqueous stability

## Kinetex C8




**Available in: 1.7  $\mu\text{m}$ , 2.6  $\mu\text{m}$  and 5  $\mu\text{m}$**   
Moderate hydrophobic and steric selectivity is offered, bringing ultra-high performance to USP L7 and other octyl silane methods.

## Kinetex Biphenyl




**Available in: 1.7  $\mu\text{m}$ , 2.6  $\mu\text{m}$  and 5  $\mu\text{m}$**   
100% aqueous stable reversed phase chemistry with hydrophobic, aromatic, and enhanced polar selectivity.

## Kinetex F5




**Available in: 1.7  $\mu$ m, 2.6  $\mu$ m, and 5  $\mu$ m**  
Highly reproducible pentafluorophenyl phase exceptional for halogenated, conjugated, isomeric, or highly polar compounds.

## Kinetex Phenyl-Hexyl



**Available in: 1.7  $\mu$ m, 2.6  $\mu$ m and 5  $\mu$ m**  
Aromatic and moderate hydrophobic selectivity result in the great retention and separation of aromatic hydrocarbons.


## Kinetex HILIC





**Available in: 1.7  $\mu$ m, 2.6  $\mu$ m and 5  $\mu$ m**  
Used under HILIC running conditions, this phase provides the highest polar selectivity for retention and separation of hydrophilic compounds.

## pH Selectivity

With excellent performance across a large pH range, **Kinetex EVO C18** columns allow you to transform chromatograms and manipulate retention order, no matter the combination of compound functionalities. Now is your chance to break the mold and let your creative LC side flourish.



**Column:** Kinetex 5  $\mu$ m EVO C18  
**Dimensions:** 150 x 4.6 mm  
**Part No.:** 00F-4633-E0  
**Mobile Phase:** A: 20 mM Potassium phosphate  
B: Acetonitrile  
**Gradient:** 20-75% B in 10 minutes  
**Flow Rate:** 1.5 mL/min  
**Temperature:** 30  $^{\circ}$ C  
**Detection:** UV @ 254 nm  
**Sample:** 1. Ibuprofen  
2. Diphenhydramine  
3. Ethyl Benzene




# Selectivity – Polar C18




The Kinetex Polar C18 contains a C18 ligand alongside a polar modified surface that increases polar compound retention and improves resolution values. Additionally, the advanced proprietary bonding technology used with this phase ensures 100 % aqueous stability as well as balanced retention on non-polar compounds. This is an excellent all purpose phase for use with multi compound mixes that contain polar and nonpolar compounds, or even single class methods that have closely related compounds, impurities or metabolites.

## UHPLC Analysis of Cyclosporine and Impurities



**Column:** Kinetex 2.6  $\mu$ m Polar C18  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** 00B-4759-AN  
**Mobile Phase:** Acetonitrile/Tert-butyl methyl ether/Water/  
Phosphoric acid (430:50:520:1)  
**Flow Rate:** 0.30 mL/min  
**Temperature:** 80 °C  
**Detection:** UV @ 210 nm  
**Sample:** Cyclosporine

## Water Soluble Vitamins




**Column:** Kinetex 2.6  $\mu$ m Polar C18  
**Dimensions:** 100 x 4.6 mm  
**Part No.:** 00B-4759-E0  
**Mobile Phase:** A: 20 mM Potassium Phosphate  
B: Methanol  
**Gradient:**

Time (min)	% B
0	0
1	0
10	60

  
**Flow Rate:** 1.2 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 210 nm  
**Sample:** 1. Thiamine  
2. Nicotinamide  
3. Pyridoxal  
4. Pyridoxine  
5. Pantothenic Acid  
6. Riboflavin

## Multi-Class 206 Pesticide Panel Screen









**Column:** Kinetex 2.6  $\mu$ m Polar C18  
**Dimension:** 50 x 4.6 mm  
**Part No.:** 00B-4759-E0  
**Mobile Phase:** A: Water  
B: 0.1 % Formic Acid in Methanol  
**Gradient:** 5-100 % B in 5 min, hold 1 min  
**Flow Rate:** 0.7 mL/min  
**Temperature:** Ambient  
**Detection:** MS/MS (SCIEX API 4000™)  
**Sample:** 206 Pesticides.  
Find the full compound list online at  
[www.phenomenex.com/Application/Detail/24017](http://www.phenomenex.com/Application/Detail/24017)

# Selectivity – Combine Kinetex and Luna® Omega

Combine the available Kinetex core-shell and Luna Omega fully porous stationary phases to expand your chromatography reach no matter whether you are working with UHPLC, HPLC, or Preparative instrumentation. The Polar C18 selectivity offers balanced polar/non-polar retention, while the Biphenyl and F5 offer excellent orthogonal alternatives to help adjust retention times and peak order.

## Natural Cannabinoids



### Conditions for all columns:

**Columns:** Luna Omega 1.6µm Polar C18  
Kinetex 1.7µm Biphenyl  
Kinetex 1.7µm F5

**Dimension:** 100 x 2.1 mm

**Mobile Phase:** A: 20 mM Ammonium Formate pH 3.2  
B: Acetonitrile

Gradient:	Time (min)	% B
	0	60
	12	95
	13	95
	13.01	60
	15	60


**Flow Rate:** 0.4 mL/min

**Temperature:** 40 °C

**Detection:** UV @ 256 nm

- Sample:**
1. CBDV
  2. Cannabidiol
  3. CBG
  4. Cannabidiolic Acid
  5. CBG-A
  6. Cannabinol
  7. Delta 9 THC
  8. Delta 8 THC
  9. CBC
  10. THCA-A

## Recommended Selectivities If You're Working With:



## Upgrading Your Fully Porous Methods:

### Fully Porous 3 $\mu\text{m}$ - 5 $\mu\text{m}$

**Kinetex 5  $\mu\text{m}$**  – Drop-in for easy performance improvements with no backpressure increase

**Kinetex 3.5  $\mu\text{m}$**  – Drop-in for easy performance improvements of pharmacopoeia methods

**Kinetex 2.6  $\mu\text{m}$**  – Dramatically improve results with efficiency/peak capacity gains

### Fully Porous sub-2 $\mu\text{m}$

**Kinetex 2.6  $\mu\text{m}$**  – Get similar efficiencies at lower backpressure allowing for greater productivity gains

**Kinetex 1.7  $\mu\text{m}$**  – Up to 20% greater efficiencies resulting in drop-in improvements

**Kinetex 1.3  $\mu\text{m}$**  – Incredible efficiency gains on high end UHPLC systems

### Fully Porous Preparative LC

**Kinetex 5  $\mu\text{m}$**  – Drop-in for easy performance improvement with no backpressure increase

# Simple Selection of the Suitable Column

	5 $\mu$ m	3.5 $\mu$ m	2.6 $\mu$ m	1.7 $\mu$ m	1.3 $\mu$ m
UHPLC					
HPLC					
PREP LC					

Phase	Best Use	pH Stability	Available Particle Size(s)				
<b>Polar C18</b>	C18 provides all purpose non-polar interactions, while novel polar modified surface increases polar compound retention and provides 100% aqueous stability.	1.5 - 8.5*			2.6 $\mu$ m		
<b>EVO C18</b>	Robust reversed phase methods even in alkaline conditions with improved peak shape for polar basic compounds	1 - 12	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>C18</b>	All purpose phase that offers the hydrophobic retention and methylene selectivity chromatographers expect from a C18 column	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	1.3 $\mu$ m
<b>XB-C18</b>	C18 phase with protective butyl side chains for improved peak shape for basic compounds under neutral and acidic conditions	1.5 - 8.5*	5 $\mu$ m	3.5 $\mu$ m	2.6 $\mu$ m	1.7 $\mu$ m	
<b>C8</b>	USP L7 phase that provides less hydrophobic and methylene selectivity than a C18	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>F5</b>	Highly reproducible pentafluorophenyl propyl phase that offers a unique combination of polar, hydrophobic, aromatic, and shape selectivity	1.5 - 8.5	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>Biphenyl</b>	100% aqueous stable and allows for excellent reversed phase retention and enhanced polar and aromatic selectivity	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>Phenyl-Hexyl</b>	Reversed phase chemistry that allows for greater retention and separation of aromatic hydrocarbons	1.5 - 8.5*	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	
<b>HILIC</b>	Unbonded silica phase for HILIC conditions to provide selectivity for polar compounds	2.0 - 7.5	5 $\mu$ m		2.6 $\mu$ m	1.7 $\mu$ m	

\*pH stability under gradient conditions. pH stability is 1.5-10 under isocratic conditions.

Choose from an extensive selection of phases for greater flexibility in UHPLC/HPLC method development. Kinetex columns come in a variety of stationary phases to cover a full spectrum of applications ranging from acids and bases, to isomers and extremely polar compounds.

**Conditions for all columns:**

- Column:** Kinetex 2.6  $\mu$ m Biphenyl  
Kinetex 2.6  $\mu$ m C18  
Kinetex 2.6  $\mu$ m XB-C18  
Kinetex 2.6  $\mu$ m Phenyl-Hexyl

**Dimensions:** 50 x 4.6 mm

**Mobile Phase:** A:Water  
B: Acetonitrile





**Gradient:** 20 to 60 % B in 6 minutes

**Flow Rate:** 1.85 mL/min

**Temperature:** 30 °C

**Detection:** UV @ 220 nm


- Samples:**
1. Estriol
  2. Hydrocortisone
  3. Prednisone
  4. Cortisone
  5. Corticosterone
  6.  $\beta$ -Estradiol
  7. Cortisone Acetate
  8. 17-Hydroxyprogesterone
  9. 21-Hydroxyprogesterone
  10. Deoxycorticosterone



# Protect Any UHPLC Column

Protect your UHPLC column, including Kinetex® core-shell columns, from damaging contaminants and microparticulates with the SecurityGuard ULTRA guard cartridge system!

- Simple to use
- Extend column lifetime
- Pressure rated to 20,000 psi (1,378 bar)
- Fits virtually all manufacturers' columns 2.1 to 4.6 mm ID



See it in action:  
[www.phenomenex.com/SecurityGuardULTRA](http://www.phenomenex.com/SecurityGuardULTRA)

# Ordering Information

5 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges <sup>‡</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJO-9298
<b>F5</b>	00A-4724-AN	00B-4724-AN	00D-4724-AN	00F-4724-AN	AJO-9322
<b>Biphenyl</b>	00A-4627-AN	00B-4627-AN	00D-4627-AN	—	AJO-9209
<b>XB-C18</b>	00A-4605-AN	00B-4605-AN	00D-4605-AN	—	AJO-8782
<b>C18</b>	00A-4601-AN	00B-4601-AN	00D-4601-AN	00F-4601-AN	AJO-8782
<b>C8</b>	—	00B-4608-AN	00D-4608-AN	—	AJO-8784
<b>Phenyl-Hexyl</b>	—	00B-4603-AN	00D-4603-AN	—	AJO-8788

for 2.1 mm ID

5 µm MidBore™ Columns (mm)				SecurityGuard™ ULTRA Cartridges <sup>‡</sup>
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
<b>EVO C18</b>	00B-4633-YO	00D-4633-YO	00F-4633-YO	AJO-9297
<b>F5</b>	00B-4724-YO	00D-4724-YO	00F-4724-YO	AJO-9321
<b>Biphenyl</b>	00B-4627-YO	00D-4627-YO	00F-4627-YO	AJO-9208
<b>XB-C18</b>	00B-4605-YO	00D-4605-YO	00F-4605-YO	AJO-8775
<b>C18</b>	00B-4601-YO	00D-4601-YO	00F-4601-YO	AJO-8775
<b>C8</b>	00B-4608-YO	00D-4608-YO	—	AJO-8777
<b>Phenyl-Hexyl</b>	00B-4603-YO	00D-4603-YO	—	AJO-8781

for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard™ ULTRA Cartridges <sup>‡</sup>
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
<b>EVO C18</b>	00B-4633-EO	00D-4633-EO	00F-4633-EO	00G-4633-EO	AJO-9296
<b>F5</b>	00B-4724-EO	00D-4724-EO	00F-4724-EO	00G-4724-EO	AJO-9320
<b>Biphenyl</b>	00B-4627-EO	00D-4627-EO	00F-4627-EO	00G-4627-EO	AJO-9207
<b>XB-C18</b>	00B-4605-EO	00D-4605-EO	00F-4605-EO	00G-4605-EO	AJO-8768
<b>C18</b>	00B-4601-EO	00D-4601-EO	00F-4601-EO	00G-4601-EO	AJO-8768
<b>C8</b>	00B-4608-EO	00D-4608-EO	00F-4608-EO	00G-4608-EO	AJO-8770
<b>Phenyl-Hexyl</b>	00B-4603-EO	00D-4603-EO	00F-4603-EO	00G-4603-EO	AJO-8774

for 4.6 mm ID

5 µm Semi-Preparative Columns (mm)			SecurityGuard™ SemiPrep Cartridges <sup>***</sup>
Phases	150 x 10	250 x 10	10 x 10
<b>EVO C18</b>	00F-4633-NO	00G-4633-NO	AJO-9306 /3pk
<b>F5</b>	—	00G-4724-NO	AJO-9323
<b>C18</b>	00F-4601-NO	00G-4601-NO	AJO-9278
<b>Biphenyl</b>	00F-4627-NO	00G-4627-NO	AJO-9280

for 10 mm ID

5 µm Axia™ Packed Preparative Columns (mm)					SecurityGuard™ PREP Cartridges*
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	15 x 21.2
<b>EVO C18</b>	00B-4633-PO-AX	00D-4633-PO-AX	00F-4633-PO-AX	00G-4633-PO-AX	AJO-9304 /ea
<b>F5</b>	—	—	00F-4724-PO-AX	00G-4724-PO-AX	AJO-9324
<b>Biphenyl</b>	00B-4627-PO-AX	00D-4627-PO-AX	00F-4627-PO-AX	00G-4627-PO-AX	AJO-9272
<b>XB-C18</b>	00B-4605-PO-AX	00D-4605-PO-AX	00F-4605-PO-AX	00G-4605-PO-AX	AJO-9145
<b>C18</b>	00B-4601-PO-AX	00D-4601-PO-AX	00F-4601-PO-AX	00G-4601-PO-AX	AJO-9145
<b>C8</b>	00B-4608-PO-AX	00D-4608-PO-AX	00F-4608-PO-AX	00G-4608-PO-AX	AJO-9205
<b>Phenyl-Hexyl</b>	00B-4603-PO-AX	00D-4603-PO-AX	00F-4603-PO-AX	00G-4603-PO-AX	AJO-9147
<b>HILIC</b>	—	00D-4606-PO-AX	00F-4606-PO-AX	00G-4606-PO-AX	AJO-9277

for 21.2 mm ID

<sup>‡</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

\* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223

\*\* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277

\*\*\* SemiPrep SecurityGuard Cartridges require holder, Part No.: AJO-9281

5 µm Axia Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges**
Phases	50 x 30	100 x 30	150 x 30	250 x 30	15 x 30
					/ea
<b>EVO C18</b>	00B-4633-U0-AX	00D-4633-U0-AX	00F-4633-U0-AX	00G-4633-U0-AX	AJO-9305
					/ea
<b>F5</b>	00B-4724-U0-AX	00D-4724-U0-AX	00F-4724-U0-AX	00G-4724-U0-AX	AJO-9325
<b>Biphenyl</b>	—	—	00F-4627-U0-AX	—	AJO-9273
<b>XB-C18</b>	00B-4605-U0-AX	00D-4605-U0-AX	00F-4605-U0-AX	00G-4605-U0-AX	AJO-9204
<b>C18</b>	00B-4601-U0-AX	00D-4601-U0-AX	00F-4601-U0-AX	00G-4601-U0-AX	AJO-9204
<b>C8</b>	00B-4608-U0-AX	00D-4608-U0-AX	00F-4608-U0-AX	00G-4608-U0-AX	AJO-9217
<b>Phenyl-Hexyl</b>	00B-4603-U0-AX	00D-4603-U0-AX	00F-4603-U0-AX	00G-4603-U0-AX	AJO-9216

for 30 mm ID

3.5 µm Analytical Columns (mm)			SecurityGuard ULTRA Cartridges†
Phases	100 x 4.6	150 x 4.6	3/pk
<b>XB-C18</b>	00D-4744-E0	00F-4744-E0	AJO-8768

for 4.6 mm ID

2.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
<b>XB-C18</b>	00B-4496-A0	00D-4496-A0	00F-4496-A0

2.6 µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges†
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJO-9298
<b>Polar C18</b>	00A-4759-AN	00B-4759-AN	—	00D-4759-AN	00F-4759-AN	AJO-9530
<b>F5</b>	00A-4723-AN	00B-4723-AN	—	00D-4723-AN	00F-4723-AN	AJO-9322
<b>Biphenyl</b>	00A-4622-AN	00B-4622-AN	—	00D-4622-AN	00F-4622-AN	AJO-9209
<b>XB-C18</b>	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
<b>C18</b>	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782
<b>C8</b>	00A-4497-AN	00B-4497-AN	00C-4497-AN	00D-4497-AN	00F-4497-AN	AJO-8784
<b>HILIC</b>	00A-4461-AN	00B-4461-AN	00C-4461-AN	00D-4461-AN	00F-4461-AN	AJO-8786
<b>Phenyl-Hexyl</b>	00A-4495-AN	00B-4495-AN	00C-4495-AN	00D-4495-AN	00F-4495-AN	AJO-8788

for 2.1 mm ID

2.6 µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges†
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
<b>EVO C18</b>	—	00B-4725-Y0	—	00D-4725-Y0	00F-4725-Y0	AJO-9297
<b>Polar C18</b>	—	00B-4759-Y0	—	00D-4759-Y0	00F-4759-Y0	AJO-9531
<b>F5</b>	—	00B-4723-Y0	—	00D-4723-Y0	00F-4723-Y0	AJO-9321
<b>Biphenyl</b>	—	00B-4622-Y0	—	00D-4622-Y0	00F-4622-Y0	AJO-9208
<b>XB-C18</b>	00A-4496-Y0	00B-4496-Y0	00C-4496-Y0	00D-4496-Y0	00F-4496-Y0	AJO-8775
<b>C18</b>	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJO-8775
<b>C8</b>	00A-4497-Y0	00B-4497-Y0	00C-4497-Y0	00D-4497-Y0	00F-4497-Y0	AJO-8777
<b>HILIC</b>	00A-4461-Y0	—	—	—	00F-4461-Y0	AJO-8779
<b>Phenyl-Hexyl</b>	—	00B-4495-Y0	—	00D-4495-Y0	00F-4495-Y0	AJO-8781

for 3.0 mm ID

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

\* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223

\*\* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277

\*\*\* SemiPrep SecurityGuard Cartridges require holder, Part No.: AJO-9281

2.6 µm Analytical Columns (mm)						SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
<b>EVO C18</b>	—	00B-4725-E0	—	00D-4725-E0	00F-4725-E0	AJO-9296
<b>Polar C18</b>	—	00B-4759-E0	—	00D-4759-E0	00F-4759-E0	AJO-9532
<b>F5</b>	—	00B-4723-E0	—	00D-4723-E0	00F-4723-E0	AJO-9320
<b>Biphenyl</b>	—	00B-4622-E0	—	00D-4622-E0	00F-4622-E0	AJO-9207
<b>XB-C18</b>	—	00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	AJO-8768
<b>C18</b>	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJO-8768
<b>C8</b>	—	00B-4497-E0	00C-4497-E0	00D-4497-E0	00F-4497-E0	AJO-8770
<b>HILIC</b>	—	00B-4461-E0	00C-4461-E0	00D-4461-E0	00F-4461-E0	AJO-8772
<b>Phenyl-Hexyl</b>	—	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0	AJO-8774

for 4.6 mm ID

1.7 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	—	00B-4726-AN	00D-4726-AN	00F-4726-AN	AJO-9298
<b>F5</b>	—	00B-4722-AN	00D-4722-AN	00F-4722-AN	AJO-9322
<b>Biphenyl</b>	—	00B-4628-AN	00D-4628-AN	00F-4628-AN	AJO-9209
<b>XB-C18</b>	00A-4498-AN	00B-4498-AN	00D-4498-AN	00F-4498-AN	AJO-8782
<b>C18</b>	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJO-8782
<b>C8</b>	00A-4499-AN	00B-4499-AN	00D-4499-AN	00F-4499-AN	AJO-8784
<b>HILIC</b>	00A-4474-AN	00B-4474-AN	00D-4474-AN	—	AJO-8786
<b>Phenyl-Hexyl</b>	—	00B-4500-AN	00D-4500-AN	00F-4500-AN	AJO-8788

for 2.1 mm ID

1.7 µm MidBore Columns (mm)				SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 3.0	50 x 3.0	100 x 3.0	3/pk
<b>XB-C18</b>	00A-4498-Y0	00B-4498-Y0	00D-4498-Y0	AJO-8775
<b>C18</b>	—	00B-4475-Y0	00D-4475-Y0	AJO-8775
<b>C8</b>	00A-4499-Y0	00B-4499-Y0	00D-4499-Y0	AJO-8777
<b>HILIC</b>	—	00B-4474-Y0	—	AJO-8779

for 3.0 mm ID

1.3 µm Minibore Columns (mm)		
Phases	30 x 2.1	50 x 2.1
<b>C18</b>	00A-4515-AN	00B-4515-AN

<sup>†</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000.

#### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Phenomenex, Aeris, Clarity, Luna and Kinetex are registered trademarks of Phenomenex. Axia, Oligo-MS, Novum, Phenex, MidBore, and SecurityGuard are trademarks of Phenomenex. Waters, Alliance, ACQUITY, UPLC, and XBridge are registered trademarks, OBD and BEH Technology are trademarks of Waters Corporation. Agilent, Poroshell, and Zorbax are registered trademarks of Agilent Technologies, Inc. Ascentis and Supelco are registered trademarks of Sigma-Aldrich Co. ACE is a registered trademark and UltraCore is a trademark of Advanced Chromatography Technologies Limited. API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license. Kromasil is a registered trademark of AkzoNobel Pulp And Performance Chemicals AB. Gilson is a registered trademark of Gilson, Inc. Shimadzu and Prominence are registered trademarks of Shimadzu Corporation. Inertsil is a registered trademark of GL Sciences, Inc.

#### Disclaimer

Phenomenex is not affiliated with Advanced Chromatography Technologies Limited, Agilent Technologies, AkzoNobel, Gilson, GL Sciences, Sigma-Aldrich Co., or Waters Corporation.

Comparative separations may not be representative of all applications.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

© 2017 Phenomenex, Inc. All rights reserved.

**Protect Your  
LC Investment!**

**novum**  
simplified liquid extraction



Remove interferences such as proteins, phospholipids, salts, and more in a quick and simple procedure.

[www.phenomenex.com/Novum](http://www.phenomenex.com/Novum)

**PHENEX**  
SYRINGE FILTERS



Convenient and inexpensive way to remove particulates from HPLC samples.

[www.phenomenex.com/Phenex](http://www.phenomenex.com/Phenex)

**guarantee**

If you are not completely satisfied with Kinetex core-shell columns, send in your comparative data to a similar product with the Kinetex column within 45 days for a FULL REFUND.

# These Phases ROCK Your LC Laboratory



## Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auiinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

## China

t: +86 400-606-8099  
f: +86 (0)22 2532-1033  
phen@agela.com

## Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

## India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

## Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

## All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com