

APPLICATION

Choosing the Right UHPLC Column for Analyzing Strongly Basic Molecules

Zeshan Aqeel and Jeff Layne
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA




Overview

One of the biggest challenges in reversed phase chromatography has always been maintaining acceptable peak shape for strongly basic analytes. With typical silica-based reversed phase media, it is impossible to completely eliminate unbonded silanol groups, even with extensive end-capping. These residual silanols can interact via ion-exchange mechanisms with the charged basic functional groups on analyte molecules, leading to peak tailing and, in some cases, adsorption that can confound accurate quantitation. These unwanted secondary interactions can be overcome through the use of ion-pairing or ion-masking agents, such as alkyl sulfonic acid additives (e.g. hexane sulfonic acid) or triethylamine (TEA), but these sorts of additives are notoriously finicky and can lead to method robustness problems as subtle alterations in the concentrations of these additives can have dramatic effects on retention, as well as peak shape. Peak tailing can be especially problematic when using typical LC-MS compatible mobile phases, such as 0.1 % formic acid or ammonium formate, as these volatile buffers do not seem to be as effective as UV compatible mobile phases like potassium phosphate or trifluoroacetic acid (TFA) at blocking or minimizing silanol interactions.

To overcome this challenge, we have designed the Luna[®] Omega 1.6 μ m PS C18 phase which contains a proprietary positively charged functional group on the surface, in addition to the standard C18 bonding and trimethylsilane (TMS) end-capping. Although the exact mechanisms are not fully defined, the presence of this positive charge on the surface results in a dramatic improvement in peak shape for basic analytes, even when using weak buffer systems such as 0.1 % formic acid. In the example below, we are attempting to isolate and quantify a small impurity eluting immediately before the peak for the basic antihistamine diphenhydramine (pKa ~8.98).

Using an LC-MS compatible mobile phase of 0.1 % formic acid in water and acetonitrile under a simple linear gradient, a hybrid particle C18 UHPLC column (Waters[®] ACQUITY[®] BEH 1.7 μ m C18; **Figure 1**) displays an exceptionally broad and tailing peak shape, and the early eluting impurity is not observed (either absorbed or buried under the main peak). However, when analyzed using the Luna[®] Omega 1.6 μ m PS C18 column under identical running conditions (**Figure 2**), you can see the vast improvement in peak width and reduced peak tailing, allowing full visualization and separation of the early eluting impurity, marked with the asterisk (*).




Luna Omega PS C18

Unique, 100 % aqueous stable mixed-mode phase that provides both polar and non-polar retention. The surface contains a positive charged ligand which aids in the retention of acidic compounds through ionic interactions, while the C18 ligand promotes general reversed phase retention. The positively charged surface also improves basic compound peaks shape through ionic repulsion.



Figure 1.
Waters[®] ACQUITY[®] BEH 1.7 μm C18




HPLC Conditions

Column: As specified
Dimension: 50 x 2.1 mm
Mobile Phase: A: 0.1 % Formic acid in Water
 B: 0.1 % Formic acid in Acetonitrile
Gradient:

Time (min)	% B
0	5
5	95

Flow Rate: 0.5 mL/min
Temperature: 40 °C
Detection: UV @ 254 nm
Sample: 1. Diphenhydramine

Figure 2.
Luna[®] Omega PS C18



Phenomenex is not affiliated with Waters Technologies Corporation. Comparative separations may not be representative of all applications.

Luna® Omega Ordering Information

1.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
Polar C18	00B-4748-A0	00D-4748-A0	00F-4748-A0
C18	00B-4742-A0	00D-4742-A0	00F-4742-A0

1.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges†
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
Polar C18	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJO-9505
PS C18	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJO-9508
C18	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJO-9502

for 2.1 mm ID

3 µm Minibore Columns (mm)					SecurityGuard™ Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0*
Polar C18	00A-4760-AN	00B-4760-AN	00D-4760-AN	00F-4760-AN	AJO-7600
PS C18	00A-4758-AN	00B-4758-AN	00D-4758-AN	00F-4758-AN	AJO-7605

for ID: 2.0 - 3.0 mm

3 µm MidBore™ Columns (mm)				SecurityGuard™ Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
Polar C18	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJO-7600
PS C18	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJO-7605

for ID: 2.0 - 3.0 mm

3 µm Analytical Columns (mm)					SecurityGuard™ Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
Polar C18	00B-4760-E0	00D-4760-E0	00F-4760-E0	00G-4760-E0	AJO-7601
PS C18	00B-4758-E0	00D-4758-E0	00F-4758-E0	00G-4758-E0	AJO-7606

for ID: 3.2-8.0 mm

5 µm Minibore Columns (mm)					SecurityGuard™ Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0*
Polar C18	00A-4754-AN	00B-4754-AN	00D-4754-AN	00F-4754-AN	AJO-7600
PS C18	00A-4753-AN	00B-4753-AN	00D-4753-AN	00F-4753-AN	AJO-7605

for ID: 2.0 - 3.0 mm

5 µm MidBore™ Columns (mm)				SecurityGuard™ Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
Polar C18	00B-4754-Y0	00D-4754-Y0	00F-4754-Y0	AJO-7600
PS C18	00B-4753-Y0	00D-4753-Y0	00F-4753-Y0	AJO-7605

for ID: 2.0 - 3.0 mm

5 µm Analytical Columns (mm)					SecurityGuard™ Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
Polar C18	00B-4754-E0	00D-4754-E0	00F-4754-E0	00G-4754-E0	AJO-7601
PS C18	00B-4753-E0	00D-4753-E0	00F-4753-E0	00G-4753-E0	AJO-7606

for ID: 3.2-8.0 mm

5 µm Axia™ Packed Preparative Columns (mm)			SecurityGuard™ Cartridges (mm)
Phases	150 x 21.2	250 x 21.2	15 x 21.2**
Polar C18	00F-4754-P0-AX	00G-4754-P0-AX	AJO-7603
PS C18	00F-4753-P0-AX	00G-4753-P0-AX	AJO-7608

for ID: 18-29 mm

5 µm Axia™ Packed Preparative Columns (mm)				SecurityGuard™ Cartridges (mm)
Phases	150 x 30	250 x 30	250 x 50	15 x 30.0*
Polar C18	00F-4754-U0-AX	00G-4754-U0-AX	00G-4754-V0-AX	AJO-7604
PS C18	00F-4753-U0-AX	00G-4753-U0-AX	00G-4753-V0-AX	AJO-7609

for ID: 30-49 mm

† SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000
 * SecurityGuard Analytical Cartridges require holder, Part No.: KJO-4282
 ** PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223
 †† PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277



If Luna analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.



APPLICATION

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

**All other countries
Corporate Office USA **

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Luna is a registered trademark and SecurityGuard, Axia, and MidBore are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Waters Technologies Corporation. Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6, 162, 362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

© 2017 Phenomenex, Inc. All rights reserved.